
10x Genomics® Sample Preparation Demonstrated Protocol

Moss Protoplast Suspension for Single Cell RNA Sequencing





10xGenomics.com

Notices

Manual Part Number

CG000090 Rev C

Legal Notices

© 2017 10x Genomics, Inc. All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, Inc., is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. "10x", "10x Genomics", "Changing the Definition of Sequencing", "Chromium", "GemCode", "Loupe", "Long Ranger", "Cell Ranger" and "Supernova" are trademarks of 10x Genomics, Inc. All other trademarks are the property of their respective owners. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of 10x Product(s) in practicing the methods set forth herein has not been validated by 10x, and such non-validated use is NOT COVERED BY 10X STANDARD WARRANTY, AND 10X HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE.

Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics, Inc., terms and conditions of sale for the Chromium[™] Controller or the Chromium Single Cell Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics, Inc that it currently or will at any time in the future offer or in any way support any application set forth herein.

Customer Information and Feedback

For technical information or advice, please contact our Customer Technical Support Division online at any time.

Email: <u>support@10xgenomics.com</u> 10x Genomics 7068 Koll Center Parkway Suite 401 Pleasanton, CA 94566 USA

Table of Contents

Demonstrated Protocol			1
	1.	Overview	2
	2.	Getting Started	2
	2.1.	Tips & Safety	2
	2.2.	Materials	2
	2.3.	Literature	3
	2.4.	Preparation – Wash Solution & Dissociation Solution	3
	3.	Preparation of Protoplast Suspension	3
	4.	Results	4

Demonstrated Protocol

Moss Protoplast Suspension for Single Cell RNA Sequencing

1. Overview

This Demonstrated Protocol outlines how to obtain a single protoplast suspension from spreading earthmoss, *Physcomitrella patens* (Grandsen or Villesexel ecotypes), in preparation for use in 10x Genomics[®] Single Cell Protocols. Protoplasts are plant cells without a cell wall but with an intact plasma membrane.

2. Getting Started

2.1. Tips & Safety

Best practices for handling any protoplast or cell sample include using sterile technique, nuclease-free reagents and consumables, minimizing pipetting steps, and using wide-bore pipette tips where possible to minimize cell or protoplast damage.

To determine accurate protoplast or cell counts, best practices include sampling the protoplast or cell suspension at least twice and at least two counts on each sample (*i.e.* a minimum of four counts in total, based on two independent draws from the suspension). Consult Technical Note *Guidelines on Accurate Target Cell Counts* (Document CG000091) for more information. A hemacytometer was used for protoplast counting.

CRITICAL!

Follow tissue and material supplier recommendations and local laboratory procedures and regulations for the safe handling, storage and disposal of biological materials.

2.2. Materials

Supplier	Description	Part Number (US)
Rainin	Tips LTS 1ML Filter RT-L1000FLR	17007954
	Pipet-Lite LTS Pipette L-1000XLS+	17014382
Corning/Falcon	Cell Strainer 40 µm	352340
Thermo Fisher	INCYTO™ C-Chip™ Disposable Hemacytometer	22-600-101
Scientific	Trypan Blue Stain (0.4%)	15250061
	Cell Scrapers	08-100-241
Eppendorf	DNA LoBind Tubes 2.0 ml*	022431048
	Model 5430 Microcentrifuge	022620584
	Model 5430 Rotor 6X15/50ML with lid	05-401-512
	Model 5430 Rotor 30X1.5ML with lid	05-401-503
VWR	Sterile Polypropylene Centrifuge Tubes with Flat Caps, 50 ml	82018-050
	Sterile Polypropylene Centrifuge Tubes with Flat Caps, 15 ml	21008-103
	Standard Orbital Shaker, Model 1000	89032-088
Sigma-Aldrich	D-Mannitol	M4125
	Driselase	D8037

*No substitutions are allowed. Items have been validated by 10x Genomics and are required for Single Cell workflow, training and system operations.

2.3. Literature

This Protocol was developed in accordance with Liu, Y.C., Vidali, L. *Efficient Polyethylene Glycol (PEG) Mediated Transformation of the Moss Physcomitrella patens*. J. Vis. Exp. (50), e2560, doi:10.3791/2560 (2011). Both Grandsen or Villesexel ecotypes were used and produced similar results.

2.4. Preparation – Wash Solution & Dissociation Solution

- a) Prepare 50 ml 8% Mannitol solution in distilled water (Wash Solution).
- b) Dissolve 2% Driselase solution in **10 ml** 8% Mannitol solution (Dissociation Solution).

3. Preparation of Protoplast Suspension

		-
- N	[]]	

This Protocol was demonstrated using moss from a total of four 10 cm petri dishes, but can be scaled down based on sample availability.

- a) Using a cell scraper, gently transfer the moss to a 50 ml centrifuge tube. Use Wash Solution to help transfer the moss into the tube.
- b) Add Wash Solution up to 10 ml.
- c) Gently mix by inverting the tube 5 times.
- d) Remove as much Wash Solution as possible by pouring; leaving moss solids in the tube.
- e) Add **10 ml** Dissociation Solution to the 50 ml centrifuge tube containing the moss.
- f) Incubate on the orbital shaker (speed setting 3) at **room temperature** for **1 h**.
- g) Place a 40 μm cell strainer on top of a 50 ml centrifuge tube and filter the protoplast suspension into the tube.
- h) Centrifuge the filtered suspension at **300 rcf** for **5 min**.
- i) Remove the supernatant without disturbing the protoplast pellet. Leave only enough solution to cover the protoplast pellet.
- j) Add 10 ml Wash Solution into the tube. Using a 10 ml serological pipette, resuspend the protoplast pellet by gently pipetting 10 times or until protoplasts are completely suspended.
- k) Centrifuge at **300 rcf** for **5 min**.
 - l) Repeat steps i k for a total of two washes.
 - m) Resuspend the protoplasts in approximately **1 ml** Wash Solution.
 - n) Determine the protoplast concentration using a hemocytometer. The target concentration is **2000 protoplasts/µl** (2 x 10⁶ protoplasts/ml).

NOTE

Repeat


The protoplasts are concentrated further than in cell-based 10x Genomics[®] Demonstrated Protocols to minimize the volume of protoplast suspension added to the Single Cell Master Mix, thereby minimizing the final concentration of Mannitol in the Single Cell Master Mix.

o) Determine the protoplast viability using a trypan blue staining method, for example.

- p) If necessary, pass the protoplast suspension through a 40 µm cell strainer to remove large debris particles and cell clumps.
- q) If necessary, dilute the protoplasts with additional Wash Solution to achieve the target protoplast concentration.
- r) Once the target protoplast concentration is obtained, place the protoplasts on ice.
- s) Proceed immediately with the 10x Genomics[®] Single Cell Protocol to avoid prolonged (>30 min) incubation of the prepared protoplasts on ice.

4. Results

Ideally, the protoplasts should be well singulated with little or no large debris particles. Very small particles of debris may be seen on the hemacytometer but do not appear to adversely affect the performance of 10x Genomics Single Cell Protocols. The percent viability of protoplasts is expected to be ~90% based on trypan blue staining.



Protoplasts are well singulated (white arrows). A few large and small debris particles can be seen in the sample but do not appear to adversely affect the performance of the 10x Genomics Single Cell Protocols.