

## DEMONSTRATED PROTOCOL

# Tumor Dissociation for Single Cell RNA Sequencing

## Overview

This protocol outlines best practices for dissociation of mouse melanoma, color tumor, and breast tumor for use in 10x Genomics Single Cell RNA protocols. Modifications to this protocol (e.g. dissociation time, resuspension buffer, enzyme concentration, centrifugation speed and time) may be required when working with other tumor types.

## Additional Guidance


Consult Demonstrated Protocol Cell Preparation Guide (Document CG00053) for Tips & Best Practices and Technical Note Guidelines for Accurate Target Cell Counts using 10x Genomics Single Cell Solutions (Document CG00091) for more information on determining accurate cell counts.

**Cells carry potentially hazardous pathogens. Follow material supplier recommendations and local laboratory procedures and regulations for the safe handling, storage, and disposal of biological materials.**

## Cell Sourcing

Tumor/Species	Description	Supplier	Part Number
CT26 colon/mouse	Whole colon tumor in Miltenyi tissue storage solution Tumor size: 300-500 mm <sup>3</sup>	Charles River Laboratories	CT26
4T1-e274 breast/mouse	Whole breast tumor in Miltenyi tissue storage solution Tumor size: 300-500 mm <sup>3</sup>	Charles River Laboratories	-
B16F10 melanoma/mouse	Whole colon tumor in Miltenyi tissue storage solution Tumor size: 300-500 mm <sup>3</sup>	Charles River Laboratories	B16F10

## Preparation - Buffers

Tumor Dissociation	
Reagents	Composition
Enzyme D*	1 vial lyophilized Enzyme D in 3 ml RPMI 1640/DMEM
Enzyme R*	1 vial lyophilized Enzyme R in 2.7 ml RPMI 1640/DMEM Thoroughly mix reconstituted Enzyme R immediately before withdrawing the required reaction volume.
Enzyme A*	1 vial lyophilized Enzyme A in 1 ml RPMI 1640/DMEM
Red Blood Cell Lysis (maintain at 4°C)	
Buffers	Composition
Wash Buffer	0.04% BSA in 1X DPBS
Red Blood Cell Lysis Solution	1X Red Blood Cell Lysis Solution in double distilled water  DO NOT use deionized water. Discard unused solution at the end of the day.

\*Part of Tumor Dissociation kit. Aliquot reconstituted enzymes and store at -20°C. Reconstituted enzyme is stable for 6 months.


## Specific Reagents & Consumables

Vendor	Item	Part Number
Thermo Fisher Scientific	Dulbecco's Phosphate-Buffered Saline (DPBS), No Calcium, No Magnesium	14190144
	UltraPure BSA (50 mg/ml)	AM2616
	Nuclease-free Water (not DEPC-treated)	AM9937
	Trypan Blue Stain, 0.4%	T10282
	Countess II FL Automated Cell Counter	AMQAF1000
Miltenyi Biotec	Countess II FL Automated Cell Counting Chamber Slides	C10228
	MACS Multistand	130-042-303
	Tumor Dissociation Kit, mouse	130-096-730
	Red Blood Cell Lysis Solution (10X)	130-094-183
	MACS SmartStrainers (70 µm)	130-098-462
Corning	gentleMACS Octo Dissociator with Heaters	130-096-427
	gentleMACS C Tubes	130-093-237
	RPMI 1640	10-040-CV
ATCC	Dulbecco's Modified Eagle Medium (DMEM)	30-2002
Eppendorf	DNA LoBind Tubes, 2.0 ml	022431048




## Protocol Overview

### 1. Tumor Dissociation



### 2. Red Blood Cell Lysis



## Protocol

This protocol was demonstrated using commercially-sourced fresh mouse tumors (300-500 mm<sup>3</sup>). The tissues were not frozen.

### 1. Tumor Dissociation

- Wash the tumor tissue in a 50-ml centrifuge tube by adding **10 ml** chilled 1X DPBS.
- Place in a petri dish and cut the tumor tissue to small pieces of ~2-4 mm<sup>3</sup>.
- Prepare enzyme mix in a gentleMACS C Tube by adding:
  - 2.35 ml RPMI 1640 or DMEM
  - 100 µl Enzyme D
  - 50 µl Enzyme R
  - 12.5 µl Enzyme A
- Transfer the tumor tissue pieces to the C Tube containing the enzyme mix.
- Tightly close the C Tube and attach upside down to a sleeve of a gentleMACS Octo Dissociator with Heaters. The tissue in the C Tube should be close to the stator of the dissociator.
- Select the dissociator program based on tumor tissue texture (soft, medium or hard tissue).
- Run the program 37C\_m\_TDK\_1 for melanoma and colon tumor (soft and medium tissue, respectively). Run the program 37C\_m\_TDK\_2 for breast tumor (hard tissue).
- At the end of the run, detach the C Tube from the dissociator.
- Centrifuge at **300 rcf** for **30 sec** at **room temperature**.
- Remove the supernatant without disturbing the cell pellet.
- Add **10 ml** RPMI 1640 or DMEM and gently pipette mix to resuspend the cell pellet.
- Filter the cell suspension through a prewetted 70-µm MACS SmartStrainer placed on a 50-ml centrifuge tube.
- Wash the strainer with **10 ml** RPMI1640 or DMEM and collect the wash in the tube with the cell suspension.
- Centrifuge the cell suspension at **300 rcf** for **7 min** at **room temperature**.
- Remove supernatant without disturbing the cell pellet. Proceed **immediately** to Red Blood Cell Lysis.

### 2. Red Blood Cell Lysis

The red blood cells present in the dissociated tumor cell suspension are lysed by the Red Blood Cell Removal Solution. Prolonged exposure to this solution could lyse the tumor cells.

- Add **1 ml** chilled 1X Red Blood Cell Removal Solution to the cell pellet from step 1o and gently pipette mix to resuspend the cells.



DO NOT vortex.

- Incubate for **10 min** at **4°C**.
- Add **10 ml** chilled Wash Buffer.
- Centrifuge at **300 rcf** for **10 min** at **4°C**.
- Remove supernatant without disturbing the cell pellet.
- Add **5 ml** chilled Wash Buffer and gently pipette mix to resuspend the cell pellet.
- Determine the cell concentration and viability using a Countess II FL Automated Cell Counter or hemocytometer. If the percentage of viable cells <70%, follow Demonstrated Protocol Removal of Dead Cells from Single Cell Suspensions for Single Cell RNA Sequencing (Document CG000093). Multiple rounds of dead cell removal may be performed to obtain a high viability sample.
- Add appropriate volume chilled Wash Buffer to the cell suspension and gently pipette mix to achieve the target cell concentration of 700-1,200 cells/µl.
- Proceed **immediately** with the 10x Genomics Single Cell protocol.

## Results

After tumor dissociation, trypan blue staining of the sample showed a high percentage (68%) of dead cells. To reduce the number of dead cells after tumor dissociation, the protocol outlined in Demonstrated Protocol Removal of Dead Cells for Single Cell RNA Sequencing (Document CG000093) was followed. After dead cell removal, the dead cell fraction was reduced from 68% to 12%, thus increasing the percentage of viable cells in the sample. Results may vary between different tumor types.

#### Before Dead Cell Removal



#### After Dead Cell Removal



## Troubleshooting

Problem	Possible Solution
High fraction of non-viable cells after tissue dissociation	Optimize dissociation protocol specific to tissue type to enhance sample quality
	Reduce fraction of dead cells. Refer to Demonstrated Protocol Removal of Dead Cells for Single Cell RNA Sequencing (Document CG000093)
	Gently handle cell suspensions by following best practices and reduce cell processing times
High fraction of cellular debris in final solution	Add 1-2 additional wash steps
	Filter cell suspension with the appropriate strainer
	Sort cells using flow cytometry
	Clean sample by density centrifugation
<700 cells/ $\mu$ l after final resuspension	Concentrate cell suspension to achieve target concentration of 700-1,200 cells/ $\mu$ l

© 2020 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: [www.10xgenomics.com/](http://www.10xgenomics.com/) trademarks. 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: [www.10xgenomics.com/](http://www.10xgenomics.com/) patents. The use of products described herein is subject to 10x Genomics Terms and Conditions of Sale, available at <http://www.10xgenomics.com/legal-notices>, or such other terms that have been agreed to in writing between 10x Genomics and user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of 10x Genomics products in practicing the methods set forth herein has not been validated by 10x Genomics, and such non-validated use is NOT COVERED BY 10x GENOMICS STANDARD WARRANTY, AND 10x GENOMICS HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics terms and conditions of sale for the Chromium Controller or the Chromium Single Cell Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics that it currently or will at any time in the future offer or in any way support any application set forth herein.

**Contact:**  
[support@10xgenomics.com](mailto:support@10xgenomics.com)  
 10x Genomics  
 6230 Stoneridge Mall Road  
 Pleasanton, CA 94588 USA

