

TECHNICAL NOTE

Sequencing Metrics & Base Composition of Chromium Single Cell ATAC Libraries

Introduction

The Chromium Single Cell ATAC Reagent Kits workflow produces Single Cell ATAC libraries for detecting accessible chromatin regions for epigenetic profiling. This Technical Note presents a comparison of sequencing metrics for Single Cell ATAC v1.1 and v2 libraries across Illumina® platforms. The expected base percentage profiles and Phred quality scores based on a control library is described to provide general guidance on the expected range of sequencing metrics on Illumina® platforms. Individual results may vary depending on the specific sequencing instrument and/or particular sample and loading characteristics.

Method

Single Cell ATAC v1.1 and v2 libraries were generated from a 1:1 mixture of human GM12878 and mouse EL4 cells following the Chromium Next GEM Single Cell ATAC Reagent Kits v1.1 (Document CG000209) and Chromium Next GEM Single Cell ATAC Reagent Kits v2 (Document CG000496) User Guides. Approximately 500 nuclei were targeted using the Chromium Next GEM Chip H Single Cell to generate libraries, and two libraries were pooled on each flowcell for each kit version.

The library fragments contain an accessible chromatin genomic insert, a 10x Barcode, a Chromium sample index N, standard Illumina® paired-end P5 and P7 sequences, and Read 1 (Read 1N) and Read 2 (Read 2N) sequences (Figure 1). The libraries were quantified using a KAPA DNA Quantification Kit and sequenced with 1% PhiX on various Illumina® sequencers to evaluate performance. The following platforms were tested using the run layout specified in Table 1:

- MiSeq™ - one flowcell
- NextSeq™ 500/550 - one flowcell
- NextSeq™ 1000/2000 - one flowcell
- HiSeq™ 2500 (RR) - one flowcell (ATAC v1.1 only)
- HiSeq™ 4000 - one flowcell (ATAC v1.1 only)
- NovaSeq™ 6000 - one S1 flowcell




Figure 1. Schematic of final library construct of a Chromium Single Cell ATAC library. The 10x Barcode (16 bp) is sequenced in the i5 index read and the Sample Index N is sequenced in the i7 index read. The accessible chromatin insert is sequenced in Read 1N and Read 2N.

Table 1. Sequencing configuration.

Sequencing Read	Number of Cycles
Read 1N	50 cycles
i7 Index	8 cycles
i5 Index	16 cycles
Read 2N	49* cycles

*50 cycles is also acceptable

Data were analyzed using the Cell Ranger ATAC pipeline (2.1), which processes the sequencing data to align reads, detect cells, and identify and annotate open chromatin regions in clusters of cells.

Results

To assess performance across Illumina® platforms, metrics for the same Single Cell ATAC library run on different sequencers are shown in Table 2. The loading concentration and cluster density for each sequencer is specified. The sequencing yield for each run is in line with the expected yield for each sequencing platform.

Q30 quality scores: In general, the Q30 quality scores for read 1 (Read 1N) were higher than those for the sample index read (i7), followed by the quality of the 10x barcode read (i5) and read 2 (Read 2N) (Table 2).

Data by Cycle: Two features available from the 'data by cycle' panel in the Illumina® Sequencing Analysis Viewer (SAV) software were assessed between the sequencing platforms. All sequencers tested had similar profiles. Representative examples for the % Bases by cycle (Figure 2A) and quality score of ≥ 30 ($\% \geq Q30$; Figure 2B) distribution are shown for the Illumina® MiSeq™ run data.

The fluctuation in the % Bases across the cycles arises from the different sequence features contained within the Single Cell ATAC libraries, including the 10x Barcode and sequence bias of the Tn5 transposase. Table 3 provides an interpretation of the expected profile for the Single Cell ATAC libraries.

The % of Phred quality scores $\geq Q30$ across the cycles for the Illumina® MiSeq™ platform are shown in Figure 2B. The Phred quality score reflects the base calling accuracy and is used to determine how much of the data from a given sequencing run can be used. Sequencing runs with lower quality scores can result in a significant portion of reads being unmappable.

In summary, % Bases by cycle and $\% \geq Q30$ quality score distribution showed highly consistent profiles for all 6 sequence platforms.

Conclusion

The sequencing metrics and base composition of sequencing reads generated for Chromium Single Cell ATAC libraries described in this Technical Note serve as guidelines for assessing the quality of Single Cell ATAC library sequencing. Additional factors that may contribute to overall success of a sequencing run and impact downstream application performance metrics include:

- Sample preparation to obtain a high quality single nuclei suspension.
- Final libraries with fragment length of 200–2,000 bp and a significant number of inserts between 200–600 bp for optimal cluster formation on Illumina® flowcells.
- Reliable and accurate library quantification using the KAPA DNA Quantification Kit using the average insert size determined by Agilent Bioanalyzer QC.
- Sequencing platform loading concentration.




Figure 2. 'Data by Cycle' plot from the Illumina® SAV software displaying the percentage of base calls (A) and Q30 quality scores (B) from representative Single Cell ATAC v1.1 and v2 libraries sequenced on Illumina® MiSeq™ using the MiSeq™ Reagent Kit v3.

Table 2. Sequencing metrics obtained when sequencing the same Single Cell ATAC libraries across different sequencers.

Instrument	Assay	Loading Concentration (pM)	Cluster Density or % PF	Yield per Lane (Gb)		% ≥Q30			
				Read 1N	Read 2N	Read 1N	i7	i5	Read 2N
MiSeq™	ATAC v1.1	11	792 K/mm ²	0.9	0.8	97.4	97.1	97.0	95.2
	ATAC v2		735 K/mm ²	0.8	0.8	97.7	96.7	97.7	94.7
NextSeq™ 500	ATAC v1.1	1.7	91.9%*	16.3	16.0	95.7	95.4	93.6	94.6
	ATAC v2		91.2%*	16.3	15.9	95.8	94.6	92.5	94.2
NextSeq™ 2000	ATAC v1.1	650	82.1%*	26.8	26.2	95.4	92.7	90.5	94.6
	ATAC v2		82.8%*	27.0	26.5	94.8	94.6	90.9	95.0
HiSeq™ 2500RR	ATAC v1.1	11	1159 K/mm ²	9.8	9.6	96.0	92.9	91.0	87.1
HiSeq™ 4000	ATAC v1.1	180	80.6%*	19.1	18.7	98.9	95.0	90.5	88.9
NovaSeq™ 6000	ATAC v1.1	300	79.9%*	50.1	49.0	95.4	92.6	90.0	93.9
	ATAC v2		82.0%*	51.4	50.3	95.4	93.0	89.9	94.1

*Percent PF is reported for NextSeq™, 2000 HiSeq™ 4000, and NovaSeq™ instead of cluster density due to the patterned flowcell.

Table 3. Interpretation of base percentages from Illumina® SAV 'Data by Cycle' plot.

Read Number	Cycle Number	Expected Profile
Read 1N	1-16	Base percentages fluctuate due to the transposase sequence bias.
	17-50	Base percentages reflect the expected base composition of the accessible chromatin in the genome.
i7	51-58	Base percentages fluctuate due to the 8 bp sample index.
i5	59-74	Base percentages fluctuate due to sequences from the 16 bp 10x Barcode. Each base is represented in roughly equal proportions.
Read 2N	75-90	Base percentages fluctuate due to the transposase sequence bias.
	91-124	Base percentages reflect the expected base composition of the accessible chromatin in the genome.

References

- Chromium Single Cell ATAC Reagent Kits v1.1 User Guide (Document CG000209)
- Chromium Single Cell ATAC Reagent Kits v2 User Guide (Document CG000496)

© 2022 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: www.10xgenomics.com/trademarks. 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: www.10xgenomics.com/patents. The use of products described herein is subject to 10x Genomics Terms and Conditions of Sale, available at www.10xgenomics.com/legal-notices, or such other terms that have been agreed to in writing between 10x Genomics and user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of 10x Genomics products in practicing the methods set forth herein has not been validated by 10x Genomics, and such non-validated use is NOT COVERED BY 10x GENOMICS STANDARD WARRANTY, AND 10x GENOMICS HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics terms and conditions of sale for the Chromium Controller or the Chromium Single Cell Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics that it currently or will at any time in the future offer or in any way support any application set forth herein.

Contact:
support@10xgenomics.com
 10x Genomics
 7068 Koll Center Parkway
 Suite 401
 Pleasanton, CA 94566 USA

