

## TECHNICAL NOTE

# Chromium Next GEM Single Cell ATAC v1.1: Reagent & Workflow Updates

## Introduction

The Chromium Single Cell ATAC Solution provides a comprehensive, scalable approach to determine the regulatory landscape of chromatin in hundreds to thousands of cells in a single sample. This document highlights the key changes in Chromium Next GEM Single Cell ATAC v1.1 reagents along with workflow updates for generating and analyzing Illumina-ready sequencing libraries.

Refer to the Chromium Next GEM Single Cell ATAC Reagent Kits v1.1 User Guide for the complete protocol.

## Chromium Next GEM Single Cell ATAC v1.1 Workflow




Figure 1. Chromium Next GEM Single Cell ATAC v1.1 workflow.

## Chromium Next GEM Chip H

The Chromium Next GEM Chip H is assembled in the Chromium Next GEM Secondary Holder. The location of the rows labeled 1-3 and the NO FILL (not used) wells are unique to the Chromium Next GEM Chip H (Figure 2A). GEMs are generated by combining a Master Mix containing transposed nuclei and enzymes, barcoded Single Cell ATAC Gel Beads v1.1, and Partitioning Oil onto the Chromium Next GEM Chip H (Figure 2B).




Figure 2. Chromium Next GEM Chip H assembled in a Chromium Next GEM Secondary Holder (A). GEM generation in Chromium Next GEM Chip H (B).






Next GEM reagents are specific to Next GEM products and should not be used interchangeably with non-Next GEM reagents.

## Sample Prep




20X Nuclei Buffer (name and PN change only) for nuclei isolation included in the Chromium Next GEM ATAC Library Kit v1.1.

## Library Prep

Protocol Steps	Single Cell ATAC	Next GEM Single Cell ATAC v1.1
<b>Reagent Kits</b>		
<b>10x Genomics Reagents</b>	Chromium Single Cell ATAC Reagent Kits  Chromium i7 Multiplex Kit N, Set A	Chromium Next GEM Single Cell ATAC Reagent Kits v1.1  Single Index Kit N, Set A
<b>Transposition</b>		
<b>Nuclei Buffer Name</b>	Nuclei Buffer	20X Nuclei Buffer (name change only)
<b>Transposition Mix</b>	ATAC Buffer	ATAC Buffer B
<b>GEM Generation &amp; Barcoding</b>		
<b>Reagents in Master Mix</b>	Barcoding Reagent	<b>Barcoding Reagent B</b> (updated reagent volumes in Master Mix; 60 µl Master Mix volume added to Transposed Nuclei)
<b>Gel Bead</b>	Single Cell ATAC Gel Beads	<b>Single Cell ATAC Gel Bead v1.1</b> (updated handling instructions; vortex 30 sec, centrifuge 5 sec before use)
<b>Chip Loading (Chips cannot be used interchangeably)</b>		
	 <p><b>Master Mix + Sample Gel Beads Partitioning Oil</b></p> <p>75 µl - row labeled 1 40 µl - row labeled 2 240 µl - row labeled 3</p>	 <p>● 70 µl - row labeled 1 ● 50 µl - row labeled 2* ● 40 µl - row labeled 3**</p> <p>* After loading the Gel Beads, wait 30 sec before loading Partitioning Oil. **After loading Partitioning Oil, attach the gasket and run the chip in the Chromium Controller <b>immediately</b></p>
<b>Chromium Controller</b>		
<b>Firmware Version</b>	3.16 or higher	4.0 or higher
<b>Run Time</b>	~8.5 min	~18 min
<b>GEM Recovery</b>		
	From top row labeled "NO FILL"	From top row labeled "3"
		
<b>Library Construction</b>		
<b>SI-PCR Cycles</b>	-	Updated cycle numbers

## Sequencing



Recommendations for sequencing libraries generated using the Chromium Next GEM Single Cell ATAC Reagent Kits v1.1 protocol are same as sequencing recommendations for libraries generated using the Chromium Single Cell ATAC Reagent Kit protocol.

## Software



Libraries generated using the Chromium Next GEM Single Cell ATAC Reagent Kits v1.1 can be analyzed and visualized using the current versions of Cell Ranger ATAC and Loupe Cell Browser.

## Chromium Next GEM Single Cell ATAC v1.1 – Product List & Documents

Product list for generating Chromium Single Cell ATAC Libraries:

REAGENT KITS	REACTIONS	PART NUMBER (PN)
Chromium Next GEM Single Cell ATAC Library & Gel Bead Kit v1.1	16 rxns 4 rxns	1000175 1000176
Chromium Next GEM Chip H Single Cell Kit	48 rxns 16 rxns	1000161 1000162
Single Index Kit N, Set A	96 rxns	1000212
INSTRUMENT		
Chromium Controller & Next GEM Accessory Kit	-	120223 (12 month warranty) 120246 (24 month warranty)
SOFTWARE		
Cell Ranger ATAC ( <a href="#">DOWNLOAD</a> )		
Loupe Cell Browser ( <a href="#">DOWNLOAD</a> )		
DOCUMENTS		
User Guide : Chromium Next GEM Single Cell ATAC Reagent Kits v1.1 ( <a href="#">CG000209</a> )		
Demonstrated Protocol : Nuclei Isolation for Single Cell ATAC Sequencing ( <a href="#">CG000169</a> )		
Demonstrated Protocol : Nuclei Isolation from Mouse Brain Tissue for Single Cell ATAC Sequencing ( <a href="#">CG000212</a> )		

© 2020 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: [www.10xgenomics.com/trademarks](http://www.10xgenomics.com/trademarks). 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: [www.10xgenomics.com/patents](http://www.10xgenomics.com/patents). The use of products described herein is subject to 10x Genomics Terms and Conditions of Sale, available at [www.10xgenomics.com/legal-notices](http://www.10xgenomics.com/legal-notices), or such other terms that have been agreed to in writing between 10x Genomics and user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of 10x Genomics products in practicing the methods set forth herein has not been validated by 10x Genomics, and such non-validated use is NOT COVERED BY 10x GENOMICS STANDARD WARRANTY, AND 10x GENOMICS HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics terms and conditions of sale for the Chromium Controller or the Chromium Single Cell Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics that it currently or will at any time in the future offer or in any way support any application set forth herein.

**Contact:**  
[support@10xgenomics.com](mailto:support@10xgenomics.com)  
10x Genomics  
6230 Stoneridge Mall Road  
Pleasanton, CA 94588 USA

