

DEMONSTRATED PROTOCOL

Thawing Dissociated Tumor Cells for Single Cell RNA Sequencing

Overview

This protocol outlines thawing of frozen dissociated tumor cells for use with 10x Genomics Single Cell protocols. The protocol was demonstrated using dissociated tumor cells from glioblastoma, renal cell carcinoma, and endometrial cancer patients.

While this Demonstrated Protocol is specific to dissociated tumor cells, the protocol may be used as a basis for handling other primary/fragile cells in preparation for use in the 10x Genomics Single Cell protocols.

Additional Guidance

Consult Demonstrated Protocol Cell Preparation Guide (Document CG00053) for Tips & Best Practices on handling cells and Technical Note Guidelines for Accurate Target Cell Counts using 10x Genomics Single Cell Solutions (Document CG000091) for determining accurate cell counts.

Cells carry potentially hazardous pathogens. Follow material supplier recommendations and local laboratory procedures and regulations for the safe handling, storage, and disposal of biological materials.

Cell Sourcing

Cryopreserved dissociated tissue cells from glioblastoma, renal cell carcinoma, and endometrial cancer patients were acquired from Discovery Life Sciences. Refer to www.dls.com for specifics.

Preparation-Buffers


Thawing
Maintain at 37°C

Media/Buffers	Composition
Complete Growth Medium	10% FBS in cell culture media (e.g., IMDM/RPMI)


Specific Reagents & Consumables

Vendor	Item	Part Number
Thermo Fisher Scientific	Gibco IMDM	12440-053
	Trypan Blue Stain (0.4%)	T10282
	Countess II FL Automated Cell Counter	AMQAF1000
	Countess II FL Automated Cell Counting Chamber Slides	C10228
Miltenyi	MACS SmartStrainers, 30 µm	130-098-458
Bel-Art	Flowmi Cell Strainer, 40 µm (alternative to Miltenyi product)	H13680-0040
Corning	Corning RPMI 1640	10-040-CM
VWR	Seradigm Premium Grade Fetal Bovine Serum (FBS)	97068-085
	Sterile Polypropylene Centrifuge Tubes with Flat Caps, 50 ml	82018-050
Eppendorf	DNA LoBind Tubes, 2.0 ml	022431048

Protocol Overview



Resuspension




Protocol

Thawing

Set up a water bath to 37°C before starting cell thawing. All cell washes are performed at room temperature.

- a. Remove cryovials from storage and **immediately** thaw in the water bath at 37°C for 2-3 min. Gently shake the cryovial while in water bath.

 **DO NOT** submerge the entire vial in the water bath. Remove from the water bath when a tiny ice crystal remains.

- b. In a biosafety hood, gently transfer thawed cells to a 50-ml centrifuge tube using a **wide-bore** pipette tip. Rinse the cryovial with 1 ml warm complete growth medium and add the rinse dropwise (1 drop per 5 sec) to the 50-ml tube while gently shaking the tube.
- c. Dilute cells in the 50-ml tube by adding 20 ml warm complete growth media. Add media dropwise. Slowly increase addition speed as the volume in the tube increases. Gently mix by inverting the tube 2x.
- d. Centrifuge at 150 rcf for 10 min. A swinging bucket rotor can be used to increase cell recovery.
- e. Remove most of the supernatant, leaving ~200 µl and resuspend cell pellet in this volume using a **regular-bore** pipette tip.
- f. Based on starting cell concentration and viability, add an appropriate volume complete growth medium to obtain a concentration of ~2,500 cells/µl.

For example, 5 x 10⁶ cells with a 70% viability can be resuspended in ~2.0 ml of medium.
- g. Determine cell concentration and viability using a Countess II Automated Cell Counter or a hemocytometer. The viability should be ~70% or more.
- h. Visually inspect cells for the presence of cellular clumps or debris.

If no visible cell clumps or debris are present:

- i. Centrifuge at 150 rcf for 10 min.
- ii. Remove the supernatant without disrupting the cell pellet.
- iii. Proceed to Resuspension.

If large cell clumps or debris are present:

- i. Add 20 ml warm media to the cells.
- ii. Pass the sample through a 30 µm MACS SmartStrainer.
- iii. Centrifuge at 150 rcf for 10 min.
- iv. Remove the supernatant without disrupting the cell pellet.
- v. Proceed to Resuspension.

Exceptionally large clumps can be removed by mechanical disruption using a 10-ml serological pipette.

If small cell debris is present:

- i. Add 20 ml warm media to the cells.
- ii. Pass the sample through a 30 µm MACS SmartStrainer.
- iii. Centrifuge at 150 rcf for 10 min.
- iv. Remove the supernatant without disrupting the cell pellet.
- v. **Repeat i-iv.**
- vi. Proceed to Resuspension.

Resuspension

- a. Based on the starting cell concentration and assuming ~50% cell loss, add an appropriate volume complete growth medium to obtain a concentration of 700-1,200 cells/µl.

Gently pipette mix using a **regular-bore** pipette tip until a single cell suspension is achieved.
- b. Determine cell concentration and viability using a Countess II Automated Cell Counter or a hemocytometer. The targeted final concentration is 700-1,200 cells/µl.
- c. Once the final cell concentration is achieved, place cells on ice.
- d. Proceed **immediately** to the 10x Genomics Single Cell protocols.

© 2019 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: www.10xgenomics.com/trademarks. 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: www.10xgenomics.com/patents. The use of products described herein is subject to 10x Genomics Terms and Conditions of Sale, available at <http://www.10xgenomics.com/legal-notices>, or such other terms that have been agreed to in writing between 10x Genomics and user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of 10x Genomics products in practicing the methods set forth herein has not been validated by 10x Genomics, and such non-validated use is NOT COVERED BY 10x GENOMICS STANDARD WARRANTY, AND 10x GENOMICS HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics terms and conditions of sale for the Chromium Controller or the Chromium Single Cell Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics that it currently or will at any time in the future offer or in any way support any application set forth herein.

Contact:
support@10xgenomics.com
 10x Genomics
 6230 Stoneridge Mall Road
 Pleasanton, CA 94588 USA

