

Performance of Blood Collection Tubes and PBMC Isolation Methods for Single Cell Gene Expression

Introduction

Blood products such as Peripheral Blood Mononuclear Cells (PBMCs), Bone Marrow Aspirates (BMAs) and leukocytes provide insightful information for assessing disease and immune states across patients. They are easily accessible, minimally invasive, and provide readouts reflective of other organ states. Although existing methods for isolating these cell types are broadly used across traditional bulk assays, the impact of anticoagulants and isolation methods on single cell assay performance is not fully understood. This Technical Note compares different blood collection tubes and PBMC isolation methods, demonstrating how they can be implemented into single cell workflows.

To assess the performance of different anticoagulants and separation methods, whole blood was collected in either BD Vacutainer tubes (ACD-A, Sodium Citrate, Heparin, or EDTA anticoagulants) or CPT tubes (Heparin or Citrate anticoagulants). PBMCs were isolated from the whole blood samples using CPT or SepMate density centrifugation methods as described in the Isolation of Leukocytes, Bone Marrow, Peripheral Blood Mononuclear Cells for Single Cell RNA Sequencing Demonstrated Protocol (CG000392). Two replicates were performed for each experimental condition.

Figure 1 delineates the sample preparation for these experiments. Samples were processed using the Chromium Single Cell 3' v3.1 (Dual Index) Reagent Kits User Guide (CG000315) and Single Cell 5' v2 Reagent Kits User Guide (CG000331).




Figure 1. Sample preparation involving different collection tubes and isolation methods. Analysis may be performed with Cell Ranger and Loupe Browser Software

Blood Anticoagulant Comparison

Cell Viability and Single Cell Behavior

PBMCs isolated from whole blood vacutainer tubes showed high cell viability across anticoagulants and isolation methods. Single cell behavior can be assessed by the barcode rank plot and fraction of reads in cells. In the barcode rank plots, a distinct knee and cliff are expected, showing good separation between cells and background. This also correlates to the high fraction of reads in cells, demonstrating that most reads are coming from called cells rather than background. Barcode rank plots overlap well across all anticoagulants and show consistent high fraction of reads in cells.




Figure 2. Cell Viability and single cell behavior data for Single Cell 3' and 5' libraries. Rank plots contain data from both Single Cell 3' and 5' libraries (B). Fraction Reads in Cells plots are for Single Cell 3' libraries (C) and Single Cell 5' libraries (D).

Sensitivity and Sequencing Saturation

To compare sensitivity, median genes and UMIs are plotted against sequenced raw read pairs. Median genes and UMIs captured across different anticoagulants and across replicates are comparable for both Single Cell 3' v3.1 and Single Cell 5' v2 libraries. At 20,000 read pairs per cell, both assays reach ~50% sequence saturation.




Figure 3. Sensitivity and sequencing saturation data for Single Cell 3' v3.1 (A) and 5' v2 (B) libraries. Sensitivity and sequencing comparisons for Multiome and Single Cell 3' v3.1 are shown in (C). Sequencing saturation curves for Single Cell 3' v3.1 (D) and Single Cell 5' v2 (E) assays.

Gene Expression Profiles Are Consistent Across Anticoagulants

The single cell barcode UMAP plots below show gene expression data analyzed and visualized using Cell Ranger aggr and Loupe Cell Browser, respectively. These UMAP plots for PBMCs collected with different anticoagulants show significant overlap as evidenced by the consistent cluster structure across both Single Cell 3' v3.1 and Single Cell 5' v2 assays.




Figure 4. UMAP plots for Single Cell 3' v3.1 (A) and 5' v2 (B) libraries.

Cell Populations Are Consistent Across Anticoagulants

Major cell populations in each cluster were identified with known PBMC markers. Comparable cell type proportions amongst the different anticoagulants are shown in the figure below, demonstrating that most anticoagulants conserve main PBMC cell populations. Slight variation in red blood cell (RBC) retention is observed across anticoagulant conditions.




Figure 5. Cell type proportions for Single Cell 3' v3.1 (A) and 5' v2 (B) libraries.

Global Gene Expression Profiles Are Consistent Across Anticoagulants

Global gene expression across different anticoagulants is highly correlated. Differentially expressed genes due to differential RBC contamination are observed.




Figure 6. Global gene expression profiles for Single Cell 3' v3.1 (A) and 5' v2 (B) libraries.

Isolation Method Comparison

Gene Expression Profiles Are Consistent Across Isolation Methods

Whole blood was isolated using CPT and SepMate tubes. Major cell populations in each cluster were identified with known PBMC markers. Comparable cell type proportions amongst the different isolation methods are observed.




Figure 7. Global gene expression clustering and cell population frequencies for Single Cell 3' v3.1 (A) and 5' v2 (B) libraries .




Figure 7. Global gene expression for Single Cell 3' v3.1 (A) and 5' v2 (B) libraries.

Global Gene Expression Profiles Are Consistent Across Isolation Methods

Conclusion

In summary, blood collection tubes of different, widely available types as well as different isolation methods are compatible with 10x Genomics Single Cell Gene Expression solutions. As described in this Technical Note, utilizing different collection tubes and isolation methods yielded comparable gene expression profiles and cell population clustering between samples. Users who are already using these blood collection tubes and isolation methods can continue to do so knowing they are fully compatible with 10x Genomics Single Cell Gene Expression assays.

References

- Garvin, Tyler, Robert Aboukhalil, Jude Kendall, Timour Baslan, *et. al.*, Interactive analysis and assessment of single-cell copy-number variations. *Nature methods* 12, no. 11 (2015): 1058
- Chromium Single Cell 3' v3.1 (Dual Index) Reagent Kits User Guide (CG000315)
- Chromium Single Cell 5'v2 (Dual Index) Reagent Kits User Guide (CG000331)

© 2021 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: www.10xgenomics.com/trademarks. 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: www.10xgenomics.com/patents. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of 10x Genomics products in practicing the methods set forth herein has not been validated by 10x Genomics, and such non-validated use is NOT COVERED BY 10X GENOMICS STANDARD WARRANTY, AND 10X GENOMICS HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics terms and conditions of sale for the Chromium Controller or the Chromium Single Cell Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics that it currently or will at any time in the future offer or in any way support any application set forth herein.

Contact:

support@10xgenomics.com

10x Genomics

6230 Stoneridge Mall Road

Pleasanton, CA 94588 USA

