Chromium Next GEM Single Cell 5' HT v2: Reagents, Workflow & Data Overview

Introduction

Chromium Next GEM Single Cell 5' HT v2 is a high throughput solution on the Chromium X series instrument for analyzing hundreds of thousands of cells per run, with up to 16 samples in a single chip. In combination with Feature Barcode technology, the high throughput assay also enables simultaneous cell surface protein detection in single cells. This Technical Note highlights sample preparation, reagents, and workflow specifics for Single Cell 5' HT v2, along with information about data analysis. A comparison of representative data derived from the Single Cell 5' HT v2 assay versus the standard Single Cell 5' V2 assay is also presented.

Chromium Next GEM Single Cell 5' HT Workflow

Chromium Next GEM Single Cell 5' HT v2 workflow (referred to as high throughput or HT) is similar to the Chromium Next GEM Single Cell 5' v2 workflow (referred to as standard), with specific updates that are indicated by an "HT" icon adjacent to the updated steps in the Single Cell 5' HT v2 User Guides (see Product List & Documents section for link to user guides).

Figure 1 provides a high level overview of the Single Cell 5' HT v2 workflow that includes a Chromium Next GEM chip designed to run up to 16 reactions on the Chromium X instrument.




Figure 1. Chromium Next GEM Single Cell 5' HT v2 workflow.



The Chromium Next GEM Single Cell 5' HT v2 assay enables partitioning of 2,000–20,000 cells per channel of the Chromium Next GEM Chip N. Up to 16 samples loaded in two rows (labeled 2A & 2B) of the chip can be processed on the Chromium X per run. Nanoliter-scale Gel Beads-in-emulsion (GEMs) are generated by combining a Master Mix containing cells and enzymes, 10x Barcoded Single Cell 5' HT Gel Beads (loaded on row labeled 1), and Partitioning Oil (loaded on rows labeled 3A & 3B) onto the chip (Figure 2) and running the chip on Chromium X. DNA molecules generated in a GEM share a common 10x Barcode. Libraries are generated and sequenced from the DNA molecules and 10x Barcodes are used to associate the individual reads back to the individual partitions. Similar to the standard assay, Single Cell 5' Gene Expression and V(D)J libraries can be generated alone or in combination with Cell Surface Protein libraries using the high throughput assay.




Figure 2. Chromium Next GEM Chip N is used for the Single Cell 5' v2 HT (high throughput) assay (A). Up to 16 samples loaded in two rows (labeled 2A & 2B) of the chip can be processed on the Chromium X per run (B). During the run (C), thousands of cells are partitioned into nanoliter-scale Gel Beads-in-emulsion (GEMs) by combining a Master Mix containing cells and enzymes (rows labeled 2A & 2B) and 10x Barcoded Single Cell 5' HT Gel Beads (row labeled 1) in two microfluidic channels (for each sample) that connect with corresponding Partitioning Oil well (rows labeled 3A & 3B). GEMs are retrieved from rows 3A and 3B to generate sequencing-ready single cell libraries.

The key differences between the Single Cell 5' HT v2 and the standard Single Cell 5' v2 assay are presented in the following section.

	Single Cell 5' v2 Standard	Single Cell 5' HT v2 High throughput	
Sample Prep			
	Recommendations for preparing single cell suspensions are unchanged between the Single Cell 5' HT v2 and the standard Single Cell 5' v2 assays. Visit the 10x Genomics Support website for specifics. Optimal cell stock concentration for chip loading is same for both assays. • 700-1,200 cells/µl		
10x Genomics Reagents			
	Chromium Next GEM Single Cell 5' Reagent Kits v2	Chromium Next GEM Single Cell 5' HT Reagent Kits v2	
	16 rxn & 4 rxn kits	48 rxn & 8 rxn kits (See <u>Product List & Documents</u> for details)	
GEM Generation & Barco	oding		
Master Mix Volume	36.3 µl/sample	72.6 µl /sample (no change in Master Mix reagents)	
Cell Suspension Vol. Table	-	Updated volume	
Gel Beads	Single Cell 5' Gel Beads v2	Single Cell 5' HT Gel Bead v2 (in updated gel bead holder)	
Chip Loading			
	Chromium Next GEM Chip K is assembled in Next GEM Secondary Holder Chip Gasket attached after chip loading	Chromium Next GEM Chip N is assembled in Chromium X Secondary Holder Chip Gasket HT attached before chip loading	
	Load up to 8 samples/chip	Load up to 16 samples/chip Only an even number of reactions can be run on this chip	
	 3 μ μ Master Mix + Sample - row labeled 1 50 μl Gel Beads - row labeled 2 45 μl Partitioning Oil - row labeled 3 Bottom row is NO FILL 	Iso µl Gel Beads - row labeled 1 Ho µl Master Mix + Sample - rows labeled 2A & 2B Ho µl Partitioning Oil - rows labeled 3A & 3B	

Technical Note | Chromium Next GEM Single Cell 5' HT v2: Reagents, Workflow & Data Overview

	Single Cell 5' v2 Standard			Single Cell 5' HT v2 High throughput		
Instrument						
	Chromium Controller			-		
	Now also compa	tible with Chrom	ium X/iX	Chromium X		
	Run time ~18 mir	n		Run time ~ 18 min		
GEM Transfer						
	Transfer 100 µl GEMs from row labeled 3			Transfer 90 μI GEMs from row 3A twice (total 180 μl) Transfer 90 μI GEMs from row 3B twice (total 180 μl)		
	100 µl GEMs/sample			180 µl GEMs /sa	ample	
	1 chip well » tran	sfer GEMs to 1 tu	ibe	1 chip well » transfer GEMs to 2 tubes		
cDNA Amplification						
cDNA Amp. PCR Cycles	Targeted Cell Recovery	Total Cycles Low RNA Cells	Total Cycles High RNA Cells	Targeted Cell Recovery	Total Cycles Low RNA Cells	Total Cycles High RNA Cells
	500-2,000	16	14	≤4,000	16	14
	2,001-6,000	14	12	4,000-12,000	14	12
	6,001-10,000	13	11	>12,000	13	11
cDNA Cleanup	After cDNA Cleanup, for each sample the cleanup product is present 1 tube.			After cDNA Cleanup, for each sample recombine the cleanup product present in 2 tubes to 1 tube.		
Library Construction						
cDNA Carry Forward	Carry forward 50 ng amplified cDNA for 5' GEX Library Construction			Carry forward only 20 µl Sample for 5' GEX Library Construction		
Recommended SI PCR	cDNA Input	(ng) T	otal Cycles	cDNA Input	(ng) T	otal Cycles
Cycles	1-25		16	0.25-25		14-16
	26-50		14	25-150		12-14
				150-500)	10-12
				500-1,00	0	8-10
				1,000-1,50	00	6-8
				>1,500		5
Single Cell 5' GEX Library	assay have the s	ame configuration	umi tso	ries generated us	Sample Index (17: Index (17: Inde	9

	Single Cell 5' v2 Standard	Single Cell 5' I High throug		
V(D)J Amplification and	Library Construction			
V(D)J Amplified Product Carry Forward	Carry forward 50 ng V(D)J amplified sample for V(D)J Library Construction	Carry forward only 10 µl V(D). for V(D)J Library Construction		
Sample Index PCR Cycles	Total 8 cycles	Cycle number depends upon t forward V(D)J amplified produ		
		V(D)J Amplified Product Input	Total Cycles	
		<25	9	
		25-150	8	
		151-500	7	
		501-1,000	6	
		>1,000	5	
V(D)J Library	Sample Index (15:10)	Sample Index (i7:10)		
	P5 TruSeq 10x UMI TSO Read 1 Barcode	V D J C TruSeq Read 2	P7	
Sequencing				
	Paired-end, dual indexing	Paired-end, dual indexing		
	Same sequencing parameters	Same sequencing parameters		
Software				
	Sequencing data derived from libraries generated visualized using the latest versions of Cell Ranger Support website.			

Multiplet Rate in Single Cell 5' HT v2 Assay

Multiplets are defined as two or more cells that have the same cell-associated barcode sequence. The multiplet rate in a single cell assay is dependent on the loading of cells in GEMs according to Poisson statistics and barcode collisions.

The table below shows empirically derived multiplet rate comparison for the standard Single Cell 5' v2 and the Single Cell 5' HT v2 assays. The multiplet rate data derived from human HEK293T and mouse NIH/3T3 cells that were mixed (1:1) and profiled using both the Single Cell 5' HT v2 and standard Single Cell 5' v2 assays is shown in Figure 3.

The multiplet rate when normalized to the same cell load is approximately half for Single Cell 5' HT v2 (1,598 multiplets in 19,632 cells detected, 0.4% multiplets per 1,000 cells, Fig. 3B) compared to standard Single Cell 5' v2 (861 multiplets in 10,140 cells detected, 0.8% multiplets per 1,000 cells, Fig. 3C).

	# of Cells	Recovered	B High Throughput
Multiplet Rate	Single Cell 5' HT v2 High throughput	Single Cell 5' v2 Standard	60k • Mouse • Multiple • Multiple
0.8%	2,000	1,000	N AUK
1.6%	4,000	2,000	S 20k
2.4%	6,000	3,000	0 50k 100k
3.2%	8,000	4,000	Human UMI Counts
4.0%	10,000	5,000	C Standard Kit
4.8%	12,000	6,000	တ္ င မ Mouse စ Multiple
5.6%	14,000	7,000	
6.4%	16,000	8,000	• Mouse • Multiple • Multiple
7.2%	18,000	9,000	0 20k
8.0%	20,000	10,000	0 0 50k 100k Human UMI Counts

Figure 3. Multiplet rates based on cell recovery in Single Cell 5' HT v2 and standard assays (A). Scatter plot of human and mouse UMI counts detected in a mixture of HEK293T and NIH/3T3 cells. Cell barcodes mapping to human (green), mouse (blue) or both, multiplets (gray), are shown for (B) Single Cell 5' HT v2 assay (1,598 multiplets in 19,632 cells detected, 0.4% multiplets per 1,000 cells) compared to the standard Single Cell 5' v2 (C) assay (861 multiplets in 10,140 cells detected, 0.8% multiplets per 1,000 cells).


Results

The representative Data Highlight provides a methods overview along with comparison of key results derived from bone marrow mononuclear cells. The libraries generated using the specified Single Cell 5' HT v2 and standard Single Cell 5' v2 reagents and protocols, were sequenced, and the data were analyzed and visualized using Cell Ranger 6.1 and Loupe Browser.

The results shown in Figures 4-10 clearly demonstrate that the high throughput and standard Single Cell 5' v2 assays yield comparable data in terms of library complexity, mapping rates, gene expression, and cell surface protein detection. Additionally, the scale of the Single Cell 5' HT v2 assay further enables rare V(D)J clonotypes to be identified in these samples.

Data Highlight

Methods Overview



Mononuclear cells were isolated from fresh bone marrow aspirates following the Demonstrated Protocol Isolation of Leukocytes, Bone Marrow and Peripheral Blood Mononuclear Cells for Single Cell RNA Sequencing (CG000392) and cryopreserved. Thawed cells were labeled with antibodies (TotalSeq C Human Universal Cocktail) for cell surface protein detection. The labeled cells were sorted to remove dead cells (7AAD+) and enrich for viable cells that were loaded separately onto a Single Cell 5' v2 Next GEM Chip N (20,000 cells) and Chip K (10,000 cells) and were processed on Chromium X followed by library preparation, sequencing, and data analysis.




Figure 4. Comparable library complexity and chemistry correlation was observed between the data derived from Single Cell 5' HT v2 and Single Cell 5' v2 assays run on Chromium X. Comparable median genes per cell (A), median counts per cell (B), and the UMI correlation (C) plots are shown.

Data Highlight Contd.

Comparable Read Mapping Rates




Figure 5. Comparable read mapping rates between the Single Cell 5' HT v2 and the standard Single Cell 5' v2 data.

Gene Expression Clustering




Figure 6. Single Cell 5' Gene Expression clustering data derived from the Single Cell 5' HT v2 (A) and standard Single Cell 5' v2 assays (B) is shown. Comparable t-SNE clusters and cell types were observed between the Single Cell 5' HT v2 and the standard Single Cell 5' v2 data.

Cell Type Annotation




Figure 7. Cell type annotation using both 5' gene expression and cell surface protein labeling data for Cell 5' HT v2 (A) and standard Single Cell 5' v2 assays (B) is shown.

Data Highlight Contd.

Comparable Cell Types Abundance




Figure 8. Comparable abundance of cell types was observed between the Single Cell 5' HT v2 and the standard Single Cell 5' v2 data.

Comparable V(D)J Clonotypes




Figure 9. Shown are honeycomb plots in which each T-cell receptor clonotype is represented by a cluster of hexagons colored according to the gene expression assay (standard or HT) used. Each hexagon represents an individual T cell with a TCR specific sequence. T-cell repertoires are comparable between the data derived from Single Cell 5' HT v2 (blue hexagons) and standard Single Cell 5' v2 (red hexagons) assays. The scale of HT enables detection of rare TCR clonotypes.






Figure 10. Shown are honeycomb plots in which memory and effector T-cell receptor clonotypes are represented by a cluster of hexagons colored according to the gene expression assay (standard or HT) used. Each hexagon represents an individual T cell with a TCR specific sequence. Enhanced TCR clonality was observed in data derived from Single Cell 5' HT v2 (blue hexagons) compared to the standard Single Cell 5' v2 assays (green hexagons).

Conclusions

The data generated using the Single Cell 5' HT v2 and the standard Single Cell 5' v2 assays are comparable in terms of library complexity, mapping rates, gene expression, and cell surface protein based cell detection. Additionally, the scale of the Single Cell 5' HT v2 assay enables detection of rare V(D)J clonotypes.

Chromium Next GEM Single Cell 5' HT v2 - Product List & Documents

Product list for generating Chromium Single Cell 5' Gene Expression Libraries using the high throughput Single Cell 5' HT v2 assay :			
Reagent Kits	Reactions	Part Number (PN)	
	8 rxns	1000374	
Chromium Next GEM Single Cell 5' HT Kit v2	48 rxns	1000356	
Chromeium Neut CEM Chin N Single Cell Vit	16 rxns	1000375	
Chromium Next GEM Chip N Single Cell Kit	80 rxns	1000357	
Chromium Single Cell Human TCR Amplification Kit	16 rxns	1000252	
Chromium Single Cell Human BCR Amplification Kit	16 rxns	1000253	
Chromium Single Cell Mouse TCR Amplification Kit	16 rxns	1000254	
Chromium Single Cell Mouse BCR Amplification Kit	16 rxns	1000255	
Dual Index Kit TT Set A	96 rxns	1000215	
Additional Kits for Feature Barcode technology Protocols	Reactions	Part Number (PN)	
5' Feature Barcode Kit	16 rxns	1000256	
Dual Index Kit TN Set A	96 rxns	1000250	
Instrument			
		1000331 (12 month warranty)	
Chromium X Upgrade Kit	-	1000332 (24 month warranty)	
User Guides			
Chromium Next GEM Single Cell 5' HT Reagent Kits v2 (Dual Index)	CG000423		
Chromium Next GEM Single Cell 5' HT Reagent Kits v2 (Dual Index) with Feature Barcode technology for Cell Surface Protein	CG000424		
Chromium X series (X/iX) with Readiness Test	CG000396		
Software			
Cell Ranger Analysis Pipeline (DOWNLOAD)			
Loupe Cell Browser (DOWNLOAD)			

Document Revision Summary

Document Number	CG000425
Title	Chromium Next GEM Single Cell 5' HT v2: Reagents, Workflow & Data Overview
Revision	Rev A
Revision Date	August 2021

© 2021 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of forducts described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: www10xgenomics.com/trademarks. 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: www.10xgenomics.com/patents. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of 10x Genomics products in practicing the methods set forth herein has not been validated by 10x Genomics, and such non-validated use is NOT COVERED BY 10X GENOMICS STANDARD WARRANTY, AND 10X GENOMICS HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics terms and conditions of sale for the Chromium Controller or the Chromium Single Cell Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics that it currently or will at any time in the future offer or in any way support any application set forth herein.

Contact:

support@10xgenomics.com 10x Genomics 6230 Stoneridge Mall Road Pleasanton. CA 94588 USA

