

## **TECHNICAL NOTE**

# Importance of Cell Stock Concentration for Accurate Target Cell Recovery

#### INTRODUCTION

10x Genomics® Single Cell Protocols require suspensions of viable, single cells (Single Cell Protocols – Cell Preparation Guide - CG00053). Accurate counting and assessment of cell viability are required as part of the workflow outlined in the Single Cell Protocols (Chromium™ Single Cell 3' Reagent Kits v2 User Guide – CG00052, Chromium™ Single Cell V(D)J Reagent Kits User Guide – CG000086) and are essential to determine the number of cells loaded onto the microfluidic chip for single cell partitioning. Manual (e.g. hemocytometer) and automated (e.g. Countess® II Automated Cell Counter) counting methods are common techniques to determine cell concentrations. A number of factors can affect the final cell count (Technical Note Guidelines for Accurate Target Cell Counts Using 10x Genomics Single Cell Solutions – CG000091). This Technical Note highlights cell stock concentration as one factor that can impact the accuracy of the final cell count and ultimately, the agreement between targeted and calculated cell recovery.

#### **METHODS**

We prepared four different single cell suspensions from HEK293T cells, each at the following cell concentrations:

- Cell Suspension #1: 340 cells/μl
- Cell Suspension #2: 500 cells/μl
- Cell Suspension #3: 780 cells/μl
- Cell Suspension #4: 1180 cells/μl

Samples were counted in four replicates with the Countess II Automated Cell Counter. For each single cell suspension, we prepared three Chromium Single Cell 3' v2 libraries with different target cell counts (1000, 3000 and 5000 cells) following the protocol outlined in the *Chromium Single Cell 3' Reagent Kits v2 User Guide* – CG00052 (see Table 1).

#### **RESULTS**

To assess the impact of different cell stock concentrations on cell recovery rates, we compared the *Estimated Number of Cells* reported by Cell Ranger with the targeted cell counts for each sample (see Figure 1). The difference between the targeted and calculated number of recovered cells is listed for each sample. Barcode Rank plots are also presented which illustrate the distribution of barcodes, ranked according to the number of UMIs that are associated with each barcode. Note that each cell partition is represented by a unique barcode. The y-axis represents the number of UMI counts that is associated with each barcode that is presented at the x-axis. Each barcodes are ranked from left to right based on the decreasing number of UMI counts associated for a given barcode. Cell-containing partitions are shown in green. Background partitions are shown in grey.

Cells Targeted	Volume Added to RT Master Mix	Cell Suspension #1 (340 cells/µl)	Cell Suspension #2 (500 cells/μl)	Cell Suspension #3 (780 cells/µl)	Cell Suspension #4 (1180 cells/µl)
1000	Cells	5.1 µl	3.5 μl	2.2 μl	1.5 μl
	Water	28.7 µl	30.3 μl	31.6 μl	32.3 µl
3000	Cells	15.3 µl	10.4 µl	6.7 μl	4.4 µl
	Water	18.5 µl	23.4 μl	27.1 μl	29.4 μl
5000	Cells	25.5 μl	17.3 µl	11.2 μl	7.4 µl
	Water	8.3 µl	16.5 µl	22.6 μl	26.4 μl

Table 1. Cell suspension and water volumes used for different cell suspensions to achieve three different cell target recovery counts.




Fig. 1. Ranked Barcode plots for all 12 libraries were prepared from four cell suspensions with different target cell counts.

Libraries were sequenced using paired-end sequencing (26bp Read 1 and 98bp Read 2) with a single sample index (8bp i7) on an Illumina® HiSeq® 4000 platform with approximately 5000 reads per cell. The resulting sequencing data was processed with Cell Ranger™ 2.0.

We determined the relative difference in the number of recovered cells calculated by Cell Ranger<sup>TM</sup> vs. the targeted cell counts for each library (Figure 2). Cell Suspensions #1 (340 cells/ $\mu$ l) and #2 (500 cells/ $\mu$ l) showed the greatest deviation from the targeted cell counts, 16 – 23% and 8 – 10% respectively. Target accuracy was improved with cell suspensions that were prepared at higher stock concentrations, Cell Suspensions #3 (780 cells/ $\mu$ l) and #4 (1180 cells/ $\mu$ l).




Fig 2. Relative difference of recovered cells from targeted cell counts.

The total number of cells required as input for 10x Genomics® Single Cell Solutions is determined by the User and their target for number of recovered cells. It is recommended that, for each assay, the input cell suspension be prepared to a concentration that allows pipetting of the desired number of cells in  $2.5 \,\mu$ l to  $15 \,\mu$ l. Transferring cell suspension volumes less than  $2.5 \,\mu$ l increases variance due to pipetting inaccuracy, while using volumes larger than  $15 \,\mu$ l increases the risk of introducing unwanted debris or inhibitors into the reactions.

We examined the relationship between the volume of single cell suspension added to the Master Mix for each Cell Suspension (#1-4) with the relative difference in the calculated recovered cells from the targeted number of cells (Figure 3). Expected cell recovery appeared positively correlated with both the input cell concentration and the volume of cells transferred into each reaction. Note that the lowest, targeted cell load (1000 cells) for Cell Suspension #4 required pipetting 1.5  $\mu$ l into the Master Mix, which resulted in a 5% difference between the number of expected and recovered cells.




Fig. 3. Relative difference of recovered cells from expected cells as a function of the volume of single cell suspension added to the Master Mix.

#### **DISCUSSION**

Libraries were made from HEK293T cells at different cell stock concentrations. The impact of initial cell concentration on cell recovery targets was assessed in the analysis of sequencing data. To achieve the targeted number of recovered cells, the optimal input cell stock concentration was found to be between 700 and 1200 cells/µl. Cell suspensions that are outside this optimal concentration range may result in unreliable cell counts. If samples are outside this range, we recommend adjusting the cell stock concentrations accordingly. If the concentration is below this range, the number of cells counted may not accurately represent the total cell count. For example, if only 100 cells are counted, the variation can be as high as 10%. A 10% or higher variation in recovered cells will contribute to an unexpected deviation seen in the estimated number of cells reported by Cell Ranger™. If the concentration is above the range specified here the cells may be difficult to count.

#### CONCLUSION

Targeting the optimal cell number for analysis in 10x Genomics® Single Cell applications is dependent on multiple factors including the cell stock concentration. To minimize variability in cell counts loaded onto the microfluidic chip and to obtain an accurate recovered cell count, cell stock concentrations should be within the optimal range specified in this Technical Note. Adhering to the best practices described here will improve overall agreement of targeted cell counts and cell counts calculated by Cell Ranger.

#### REFERENCES

- Chromium™ Single Cell 3' Reagent Kits v2 User Guide (CG00052)
- Single Cell Protocols Cell Preparation Guide CG00053
- Chromium™ Single Cell V(D)J Reagent Kits User Guide (CG000086)
- Guidelines for Accurate Target Cell Counts Using 10x Genomics® Single Cell Solutions (CG000091)
- 4 10x Genomics® | CG000127 Rev A Technical Note Importance of Cell Stock Concentration for Accurate Target Cell Recovery

### **Notices**

#### **Document Number**

CG00127 Rev A Technical Note

#### **Legal Notices**

© 2017 10x Genomics, Inc. All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, Inc., is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. "10x", "10x Genomics", "Changing the Definition of Sequencing", "Chromium", "GemCode", "Loupe", "Long Ranger", "Cell Ranger" and "Supernova" are trademarks of 10x Genomics, Inc. All other trademarks are the property of their respective owners. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of 10x Product(s) in practicing the methods set forth herein has not been validated by 10x, and such non-validated use is NOT COVERED BY 10X STANDARD WARRANTY, AND 10X HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE.

Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics, Inc., terms and conditions of sale for the Chromium™ Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics, Inc that it currently or will at any time in the future offer or in any way support any application set forth herein.

#### **Customer Information and Feedback**

For technical information or advice, please contact our Customer Technical Support Division online at any time.

Email: support@10xgenomics.com

10x Genomics 7068 Koll Center Parkway

Suite 401

Pleasanton, CA 94566 USA