# Chromium Next GEM Single Cell 3' HT v3.1: Reagents, Workflow & Data Overview

### Introduction

Chromium Next GEM Single Cell 3' HT v3.1 is a high throughput solution on the Chromium X series instrument for analyzing hundreds of thousands of cells per run, with up to 16 samples in a single chip. In combination with Feature Barcode technology, the high throughput assay also enables simultaneous cell surface protein detection, CRISPR screening, and cell multiplexing in single cells. This Technical Note highlights sample preparation, reagents, and workflow specifics for Single Cell 3' HT v3.1, along with information about data analysis. A comparison of representative data derived from the Single Cell 3' HT v3.1 assay versus the standard Single Cell 3' v3.1 assay is also presented.

### **Chromium Next GEM Single Cell 3' HT Workflow**

Chromium Next GEM Single Cell 3' HT v3.1 workflow (referred to as high throughput or HT) is similar to the Chromium Next GEM Single Cell 3' v3.1 workflow (referred to as standard), with specific updates that are indicated by an "HT" icon adjacent to the updated steps in the Single Cell 3' HT v3.1 User Guides (see Product List & Documents section for link to user guides).

Figure 1 provides a high level overview of the Single Cell 3' HT v3.1 workflow that includes a Chromium Next GEM chip designed to run up to 16 reactions on the Chromium X instrument.





Figure 1. Chromium Next GEM Single Cell 3' HT v3.1 workflow.



The Chromium Next GEM Single Cell 3' HT v3.1 assay enables partitioning of 2,000–20,000 cells per channel of the Chromium Next GEM Chip M (2,000– 60,000 cells per channel with CellPlex). Up to 16 samples loaded in two rows (labeled 2A & 2B) of the chip can be processed on the Chromium X per run. Nanoliter-scale Gel Beads-in-emulsion (GEMs) are generated by combining a Master Mix containing cells and enzymes, 10x Barcoded (~3,500,000 barcodes) Single Cell 3' HT Gel Beads (loaded in row labeled 1), and Partitioning Oil (loaded in rows labeled 3A & 3B) onto the chip (Figure 2) and running the chip on Chromium X. DNA molecules that are generated in a GEM share a common 10x Barcode. Libraries are generated and sequenced from the DNA molecules and 10x Barcodes are used to associate the individual reads back to the individual partitions. Similar to the standard assay, Single Cell 3' Gene Expression libraries can be generated alone or in combination with Cell Surface Protein, CRISPR Screening, and Cell Multiplexing libraries using the high throughput assay.



**Figure 2.** Chromium Next GEM Chip M is used for the Single Cell 3' v3.1 HT (high throughput) assay (A). Up to 16 samples loaded in two rows (labeled 2A & 2B) of the chip can be processed on the Chromium X per run (B). During the run (C), thousands of cells are partitioned into nanoliter-scale Gel Beads-in-emulsion (GEMs) by combining a Master Mix containing cells and enzymes (rows labeled 2A & 2B) and 10x Barcoded Single Cell 3' HT Gel Beads (row labeled 1) in two microfluidic channels (for each sample) that connect with corresponding Partitioning Oil well (rows labeled 3A & 3B). GEMs are retrieved from rows 3A and 3B to generate sequencing-ready single cell libraries.

The key differences between the Single Cell 3' HT v3.1 and the standard Single Cell 3' v3.1 assays are presented in the table below. Refer to the relevant user guides for complete information.

Sample Prep         Sample Prep       Recommendations for preparing single cell sub-sub-sub are unchanged between the Single v3.1 and the standard Single Cell 3' v3.1 assays. Visit the 10x Genomics Support website for Optimal cell stock concentration for chip loading is same for both assays. 	
v3.1 and the standard Single Cell 3' v3.1 assays. Visit the 10x Genomics Support website for Optimal cell stock concentration for chip loading is same for both assays.         v700-1,200 cells/µl         • 1,300-1,600 cells/µl (with cell multiplexing)         10x Genomics Reagents         Chromium Next GEM Single Cell 3'         Reagent Kits v3.1         16 rxn & 4 rxn kits         CBEM Generation & Barcoting         GEM Generation & Barcoting         Master Mix Volume         31.9 µl/sample         Gel Beads         Single Cell 3' Gel Bead v3.1         Gel Beads         Single Cell 3' Gel Bead v3.1         Chromium Next GEM Chip G is assembled in Next GEM Chip G is assembled in Next GEM Scondary Holder         Chip Gasket attached after chip loading	
• 1,300-1,600 cells/µl (with cell multiplexing) <b>10x Genomics Reagents</b> Chromium Next GEM Single Cell 3' Reagent Kits v3.1         Reagent Kits v3.1       Chromium Next GEM Single Cell 3' Reagent Kits v3.1         16 rxn & 4 rxn kits       48 rxn & 8 rxn kits (See Product List & Documents for det (See Product List & Documents for det (See Product List & Documents for det (See Product List & Documents for det (no change in Master Mix reagents)         Master Mix Volume       31.9 µl/sample <b>63.6</b> µl/sample (no change in Master Mix reagents)         Cell Suspension Vol. Table       -       Updated volume         Gel Beads       Single Cell 3' Gel Bead v3.1 (in updated gel bead holder)       Image: Single Cell 3' HT Gel Bead v3.1 (in updated gel bead holder)         Chip Loading       Chromium Next GEM Chip G is assembled in Next GEM Secondary Holder Chip Gasket attached after chip loading       Chromium Next GEM Chip M is assemily chip Gasket HT attached before chip for the condary Holder	
10x Genomics Reagents         Chromium Next GEM Single Cell 3' Reagent Kits v3.1       Chromium Next GEM Single Cell 3' Reagent Kits v3.1         16 rxn & 4 rxn kits       48 rxn & 8 rxn kits (See Product List & Documents for detain (See Product List	
Chromium Next GEM Single Cell 3' Reagent Kits v3.1Chromium Next GEM Single Cell 3' HT Reagent Kits v3.116 rxn & 4 rxn kits48 rxn & 8 rxn kits (See Product List & Documents for det (See Product List & Documents for detGEM Generation & Barcours63.6 µl/sample (no change in Master Mix reagents)Master Mix Volume31.9 µl/sample63.6 µl/sample (no change in Master Mix reagents)Cell Suspension Vol. Table-Updated volumeGel BeadsSingle Cell 3' Gel Bead v3.1Single Cell 3' HT Gel Bead v3.1 (in updated gel bead holder)Chrip LoadingChromium Next GEM Chip G is assembled in Next GEM Secondary Holder Chip Gasket attached after chip loadingChromium Next GEM Chip M is assembled in Chromium X Secondary Holder Chip Gasket HT attached before chip I	
Reagent Kits v3.1       Reagent Kits v3.1         16 rxn & 4 rxn kits       48 rxn & 8 rxn kits         CEM Generation & Barcotion         Master Mix Volume       31.9 µl/sample         Master Mix Volume       31.9 µl/sample         Cell Suspension Vol.       -         Table       Updated volume         Gel Beads       Single Cell 3' Gel Bead v3.1         Chip Loading       Chromium Next GEM Chip G is assembled in Next GEM Secondary Holder         Chip Gasket attached after chip loading       Chromium Next GEM Chip M is assembled in Chip Gasket HT attached before chip I	
GEM Generation & Barcoding       G3.6 µl/sample (no change in Master Mix reagents)         Master Mix Volume       31.9 µl/sample       G3.6 µl/sample (no change in Master Mix reagents)         Cell Suspension Vol. Table       -       Updated volume         Gel Beads       Single Cell 3' Gel Bead v3.1       Single Cell 3' HT Gel Bead v3.1 (in updated gel bead holder)         Chip Loading       Chromium Next GEM Chip G is assembled in Next GEM Secondary Holder Chip Gasket attached after chip loading       Chromium Next GEM Chip M is assembled in Chip Gasket HT attached before chip I	
Master Mix Volume       31.9 μl/sample       63.6 μl/sample (no change in Master Mix reagents)         Cell Suspension Vol. Table       -       Updated volume         Gel Beads       Single Cell 3' Gel Bead v3.1       Single Cell 3' HT Gel Bead v3.1 (in updated gel bead holder)         Chip Loading       Chromium Next GEM Chip G is assembled in Next GEM Chip M is assembled in Chip Gasket attached after chip loading       Chromium X Secondary Holder Chip G is assembled in Chip Gasket HT attached before chip H	ails)
Cell Suspension Vol. Table       -       Updated volume         Gel Beads       Single Cell 3' Gel Bead v3.1       Single Cell 3' HT Gel Bead v3.1         Chip Loading       Chromium Next GEM Chip G is assembled in Next GEM Secondary Holder Chip Gasket attached after chip loading       Chromium Next GEM Chip M is assembled in Chip Gasket HT attached before chip H	
Table       Single Cell 3' Gel Bead v3.1       Single Cell 3' HT Gel Bead v3.1         Gel Beads       Single Cell 3' Gel Bead v3.1       Single Cell 3' HT Gel Bead v3.1         Chip Loading       Chromium Next GEM Chip G is assembled in Next GEM Secondary Holder Chip Gasket attached after chip loading       Chromium Next GEM Chip M is assembled in Chromium X Secondary Holder Chip Gasket HT attached before chip I	
Chip Loading       Chromium Next GEM Chip G is assembled in Next GEM Secondary Holder       Chromium Next GEM Chip G is assembled in Chromium X Secondary Holder         Chip Gasket attached after chip loading       Chip Gasket HT attached before chip I	
Chromium Next GEM Chip G is assembled in Next GEM Secondary HolderChromium Next GEM Chip M is assembled Chromium X Secondary HolderChip Gasket attached after chip loadingChip Gasket HT attached before chip I	
Next GEM Secondary HolderChromium X Secondary HolderChip Gasket attached after chip loadingChip Gasket HT attached before chip	
	bled in
	oading
Load up to 8 samples/chip Load up to 16 samples/chip	
Only an even number of reactions can b this chip	be run on
70 μl Master Mix + Sample - row labeled 1130 μl Gel Beads - row labeled 1	
50 µl Gel Beads - row labeled 2 <b>140 µl</b> Master Mix + Sample - rows labeled	ed 2A & 2B
45 μl Partitioning Oil - row labeled 3 <b>140 μl</b> Partitioning Oil - rows labeled 3	A & 3B
Bottom row is NO FILL -	

	Single Cel Stand		Single Cell 3 High throu		
Instrument					
	Chromium Controller		-		
	Now also compatible with	Chromium X/iX	Chromium X		
	Run time ~18 min		Run time <b>~18 min</b>		
GEM Transfer					
	Transfer 100 µl GEMs from row labeled 3		Transfer <b>90 µl</b> GEMs from row 3A twice (total 180 µl) Transfer <b>90 µl</b> GEMs from row 3B twice (total 180 µl)		
	100 µl GEMs/sample		180 µl GEMs /sample		
	1 chip well » transfer GEMs	s to 1 tube	1 chip well » transfer GEMs	to 2 tubes	
cDNA Amplification					
cDNA Amp. PCR Cycles	Targeted Cell Recovery/Well	PCR Cycles	Targeted Cell Recovery/Well	PCR Cycles	
	<500 cells	13	<12,000 cells	12	
	500-6,000 cells	12	>12,000 cells	11	
	>6,000 cells	11			
cDNA Cleanup	After cDNA Cleanup, for e cleanup product is presen		After cDNA Cleanup, for ea the cleanup product prese		
Library Construction					
Single Cell 3' Gene Expression Library	configuration. Sample Index (i5:10	-	Index	assay have the same ((i7:10) P7	
Sequencing					
	Paired-end, dual indexing		Paired-end, dual indexing		
	-		Same sequencing paramete	rs	
	Single Cell 3' v3.1 is also available as a single index kit. Sequencing type for these libraries is paired-end, single indexing		Not available as single index kit		
Software					
			using either protocol can be a and Loupe Browser available (		

# Multiplet Rate in Single Cell 3' HT v3.1 Assay

Multiplets are defined as two or more cells that have the same cell-associated barcode sequence. The multiplet rate in a single cell assay is dependent on the loading of cells in GEMs according to Poisson statistics and barcode collisions. The table below shows empirically derived multiplet rate comparison for the standard Single Cell 3' v3.1 and the Single Cell 3' HT v3.1 assays. The multiplet rate data derived from human HEK293T and mouse NIH/3T3 cells that were mixed (1:1) and profiled using both the Single Cell 3' HT v3.1 and standard Single Cell 3' v3.1 assays is shown in Figure 3. The multiplet rate when normalized to the same cell load is approximately half for Single Cell 3' HT v3.1 (1,558 multiplets in 18,293 cells detected, ~0.4% multiplets per 1,000 cells, Fig. 3B) compared to standard Single Cell 3' v3.1 (931 multiplets in 9,513 cells detected, ~0.8% multiplets per 1,000 cells, Fig. 3C).

	e Cell 3' HT v3.1		В	Doowarad	# of Cells	
<ul><li>Human</li><li>Mouse</li><li>Multiple</li></ul>		100k 80k			Single Cell 3' HT v3.1 High throughput	Multiplet Rate
		60k 40k		1,000	2,000	0.8%
		20k		2,000	4,000	1.6%
	lk 100k	00		3,000	6,000	2.4%
	n UMI Counts	Hu		4,000	8,000	3.2%
	le Cell 3' v3.1		С	5,000	10,000	4.0%
<ul> <li>Human</li> <li>Mouse</li> <li>Multiple</li> </ul>		100k 80k		6,000	12,000	4.8%
• Multiple		60k		7,000	14,000	5.6%
		40k		8,000	16,000	6.4%
		20k		9,000	18,000	7.2%
	lk 100k n UMI Counts	00		10,000	20,000	8.0%


**Figure 3.** Multiplet rates based on cell recovery in Single Cell 3' HT v3.1 and standard assays (A). Scatter plot of human and mouse UMI counts detected in a mixture of HEK293T and NIH/3T3 cells. Cell barcodes mapping to human (green), mouse (blue) or both, multiplets (gray), are shown for (B) Single Cell 3' HT v3.1 assay (1,558 multiplets in 18,293 cells detected, 0.4% multiplets per 1,000 cells) compared to the standard Single Cell 3' v3.1 (C) assay (931 multiplets in 9,513 cells detected, 0.8% multiplets per 1,000 cells).

### Results

The representative Data Highlight provides a Methods Overview along with comparison of key results derived from 7 non-small-cell lung cancer (NSCLC) patients. The libraries were generated using the specified Single Cell 3' HT v3.1 and standard Single Cell 3' v3.1 reagents and protocols, were sequenced, and the data were analyzed and visualized using Cell Ranger 6.1 and Loupe Browser. The results shown in Figures 4-12 clearly demonstrate that the high throughput and standard Single Cell 3' v3.1 assays yield comparable data in terms of library complexity, mapping rates, gene expression, cell multiplexing, and cell surface protein detection. Additionally, the scale of the Single Cell 3' HT v3.1 assay further enhances the ability to detect rare cell types in these samples.


### **Data Highlight**

### **Methods Overview**



Dissociated tumor cells (DTCs) from 7 patients with non-small-cell lung cancer (NSCLC) were thawed and each sample was labeled with antibodies (TotalSeq B TBNK panel) for cell surface protein detection. Individual samples were then labeled with Cell Multiplexing Oligos (CMOs). All labeled samples were pooled, followed by cell sorting to remove dead cells (7-AAD+) and enrich for viable cells that were loaded separately onto a Single Cell 3' v3.1 Next GEM Chip M and Chip G targeting 40,000 and 20,000 cell recovery respectively. The chips were run on Chromium X followed by library preparation, sequencing, and data analysis as described in the respective user guides.

#### **Comparable Library Complexity & Correlation**



**Figure 4.** Comparable library complexity and chemistry correlation was observed between the data derived from Single Cell 3' HT v3.1 and Single Cell 3' v3.1 assays run on Chromium X. For the 36,293 cells recovered in the Single Cell 3' HT v3.1 assay and the 16,442 cells recovered in the standard assay, comparable median genes per cell (A), median counts per cell (B), and UMI correlation (C) were observed.

## Data Highlight Contd.

#### **Comparable Read Mapping Rates**





Figure 5. Comparable read mapping rates between the Single Cell 3' HT v3.1 and the standard Single Cell 3' v3.1 data.

### **Cell Multiplexing Oligo based Projection**



**Figure 6.** Chromium Single Cell 3' Cell Multiplexing libraries derived from multiplexed dissociated tumor cell samples from each of the seven donors showed similar cell multiplex tag clustering in both the Single Cell 3' HT v3.1 (A) and standard Single Cell 3' v3.1 (B) data.




#### **Overlapping Gene Expression Clustering & Cell Multiplexing Tag**

**Figure 7.** Overlapping Single Cell 3' Gene Expression clustering data derived from the Single Cell 3' HT v3.1 and standard Single Cell 3' v3.1 assays is shown in the aggregated t-SNE plot (A). The same plot with the superimposed data from 7 donors (based on cell multiplexing tags) is shown in panel B. Donor-specific clusters are observed, highlighting the diversity within 7 cancer patients.


### Data Highlight Contd.

#### Gene Expression Based Cell Clustering & Cell Type Identification




**Figure 8**. Similar cell clusters and cellular populations were detected in the tumor cells profiled using aggregated gene expression data from the Single Cell 3' HT v3.1 (A) and the standard Single Cell 3' v3.1 assays (B). The t-SNE plots show overlapping gene expression based cell clustering along with highly concordant cell type identification based on manual annotation.

### **Cell Surface Protein Based Cell Identification**



**Figure 9.** Chromium Single Cell 3' Gene Expression and Cell Surface Protein libraries generated using the Single Cell 3' HT v3.1 (A) and the standard Single Cell 3' v3.1 (B) assays show comparable cell type identification. The cell types were manually annotated based on single cell gene expression and data derived from cells labeled with TotalSeq B TBNK antibody panel.


#### **Tumor Gene Expression Clusters**



**Figure 10.** Seven comparable tumor gene clusters were observed in manually annotated data derived using the Single Cell 3' HT v3.1 (A) and the standard Single Cell 3' v3.1 (B) assays.


# Data Highlight Contd.

### **Cell Type Identification**



**Figure 11.** Tumor cells derived from 7 donors were manually annotated based on cell surface protein data derived using the Single Cell 3' HT v3.1 and the standard Single Cell 3' v3.1 assays. For each of the 7 donors (D1-D7), the graph shows absolute number of identified cell types for both assays, including CD45+ immune cells, normal lung cells, and various other tumor cells. Each donor displays a unique profile of tumor cell populations with specific genes upregulated. Additional rare cell types were detected in the Single Cell 3' HT v3.1 data

### **Rare Tumor Cell Type Detection**



**Figure 12.** Data derived from Donor 6 sample shows overlapping Single Cell 3' Gene Expression for 6 out of 7 tumor clusters between the Single Cell 3' HT v3.1 and standard Single Cell 3' v3.1 assays as shown the aggregated t-SNE plot (A) The individual t-SNE plots derived from the two assays (B, C), show detection of the cluster (ELF3/EPCAM+) only in the Single Cell 3' HT v3.1 data. These results indicate that the scale of the Single Cell 3' HT v3.1 assay enables identification of rare tumor clusters to further enhance the understanding of tumor physiology.

### Conclusions

The data generated using the Single Cell 3' HT v3.1 and the standard Single Cell 3' v3.1 assays are comparable in terms of library complexity, mapping rates, gene expression, and cell surface protein based cell detection. The ability to annotate cell types is further enhanced in Single Cell 3' HT v3.1 assay due to the ability to multiplex up to 60,000 cells per chip channel. The Single Cell 3' HT v3.1 data also highlights different tumor markers present in samples from various donors along with detection of rare cell types. These can be powerful tools for understanding tumor physiology and enabling discovery of better drug targets.

# Chromium Next GEM Single Cell 3' HT v3.1 – Product List & Documents

Product list for generating Chromium Single Cell 3' Gene Expression Libraries using the high throughput Single Cell 3' HT v3.1 assay :

Reagent Kits	Reactions	Part Number (PN)
	48 rxns	1000348
Chromium Next GEM Single Cell 3' HT Kit v3.1	8 rxns	1000348
Chromium Next GEM Chip M Single Cell Kit	80 rxns	1000349 (orderable only with 1000348)
	16 rxns	<b>1000371</b> (orderable only with 1000370)
Dual Index Kit TT Set A	96 rxns	1000215
Additional kits for Feature Barcode technology protocols		
3' Feature Barcode Kit	16 rxns	1000262
3' CellPlex Kit	48 rxns	1000261
Dual Index Kit NT Set A	96 rxns	1000242
Dual Index Kit NN Set A	96 rxns	1000243
Instrument		
Chromium X Upgrade Kit		1000331 (12 month warranty) 1000332 (24 month warranty)
User Guides		
Chromium Next GEM Single Cell 3' HT Reagent Kits v3.1 (Dual Index)	CG00416	
Chromium Next GEM Single Cell 3' HT Reagent Kits v3.1 (Dual Index) with Feature Barcode technology for Cell Surface Protein	CG00417	
Chromium Next GEM Single Cell 3' HT Reagent Kits v3.1 (Dual Index) with Feature Barcode technology for CRISPR Screening	CG00418	
Chromium Next GEM Single Cell 3' HT Reagent Kits v3.1 (Dual Index) with Feature Barcode technology for Cell Multiplexing	CG00419	
Chromium Next GEM Single Cell 3' HT Reagent Kits v3.1 (Dual Index) with Feature Barcode technology for Cell Surface Protein & Cell Multiplexing	CG00420	
Chromium Next GEM Single Cell 3' HT Reagent Kits v3.1 (Dual Index) with Feature Barcode technology for CRISPR Screening & Cell Multiplexing	CG00421	
Chromium X Series (X/iX) with Readiness Test	CG00396	
Software		
Cell Ranger Analysis Pipeline (DOWNLOAD)		
Loupe Browser (DOWNLOAD)		

### **Document Revision Summary**

Document Number	CG000422
Title	Chromium Next GEM Single Cell 3' HT v3.1: Reagents, Workflow & Data Overview
Revision	Rev A
<b>Revision</b> Date	August 2021

© 2021 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: www.10xgenomics.com/trademarks. 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: www.10xgenomics.com/patents. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of 10x Genomics products in practicing the methods set forth herein has not been validated by 10x Genomics, and such non-validated use is NOT COVERED BY 10X GENOMICS STANDARD WARRANTY, AND 10X GENOMICS HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics terms and conditions of sale for the Chromium Controller or the Chromium Single Cell Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics that it currently or will at any time in the future offer or in any way support any application set forth herein.

#### **Contact:**

support@10xgenomics.com 10x Genomics 6230 Stoneridge Mall Road Pleasanton. CA 94588 USA

