


Xenium In Situ Gene Expression – Protocol Planner

Introduction

Xenium In Situ measures gene expression in tissue sections derived from either formalin fixed and paraffin embedded (FFPE) or fresh frozen (FF) tissue samples placed on Xenium Slides. This Protocol Planner provides an overview of the workflow along with the Xenium Analyzer overview. To enable efficient planning, a breakdown of key protocol steps and times, list of user-acquired reagents and consumables, and information about supporting documentation that will be available for executing the Xenium Gene Expression workflow is also provided.

10x Genomics Xenium Reagent Kits are not listed in this document.



Contents

3	Introduction
4	Xenium Analyzer Overview
4	Key Specifications
5	Installation & Training
6	Site Preparation
9	FFPE Samples
9	Workflow Documents
10	Key Protocol Steps & Timing
11	Reagents & Consumables - Tissue Sectioning
12	Reagents & Consumables - Deparaffinization
13	Fresh Frozen Samples
13	Workflow Documents
14	Key Protocol Steps & Timing
15	Reagents & Consumables - Tissue Sectioning
16	Reagents & Consumables - Fixation & Permeabilization
17	Probe Hybridization, Ligation & Amplification
17	Key Protocol Steps & Timing
18	Reagents & Consumables
19	Xenium Analyzer
19	Key Protocol Steps & Timing
20	Reagents & Consumables
21	Appendix
21	Reagents & Consumables for Quencher Removal and H&E Staining <i>(optional)</i>

Introduction

The protocol planner provides resources and guidelines to prepare a laboratory for seamless planning and execution of the Xenium In Situ Gene Expression workflow. The key topics covered in this document are highlighted below.

Xenium Analyzer Overview

The document provides a high level overview of key instrument specifications, delivery, installation, training along with guidelines for site preparation.

Workflow Documents

A list of documents to support various steps of the in situ gene expression workflow is provided for formalin-fixed & paraffin-embedded (FFPE) and fresh frozen (FF) and tissue samples. These documents, along with many additional resources, will be available on the 10x Genomics Support website once the Xenium Analyzer is installed and ready to use.

Key Protocol Steps & Timing

A breakdown of the off-instrument and on-instrument workflow steps, the time required to perform each step, and safe stopping points are provided.

Reagents & Consumables (not supplied by 10x Genomics)

The reagents & consumables for various steps of the Xenium In Situ Gene Expression workflow are listed in this document. The Appendix also includes a list of items for optional H&E staining along with post-run quencher removal (only if staining slides after the instrument run). The listed items have been tested by 10x and perform optimally with the assay. These items will not be supplied by 10x Genomics and should be acquired from the indicated vendors. Refer to the manufacturer's website for regional part numbers. For items with multiple options, choose one based on availability and preference. **Substituting materials may adversely affect system performance.** This list may not include some standard laboratory equipment.



Note that some reagents and consumables, such as PBS, Tween, ice buckets etc., are common across multiple steps of the workflow and need not be bought individually for each step.

Gene Panel Selection

Prior to executing the Xenium In Situ Gene Expression workflow, ensure that a compatible gene panel has been selected. 10x Genomics provides the option of using pre-designed gene panels. Additionally, the pre-designed panel may be customized by adding genes of interest.

- 10x Genomics pre-designed panels: Xenium Mouse Brain Gene Expression Panel and Xenium Human Breast Gene Expression Panel
- Custom gene panels: Contact 10x Genomics via email at customerservice@10xgenomics.com for information about designing custom gene panels that are compatible with pre-designed panels. The lead time for acquiring custom panels is ~4 weeks (~1 week for gene selection, 3 weeks for ordering and shipping).

Visit the 10x Genomics website for additional information.

1.0 Xenium Analyzer Overview

Xenium is an end-to-end platform from 10x Genomics that provides highly sensitive, targeted gene expression information at sub-cellular resolution. This platform is powered by the Xenium Analyzer, a versatile instrument for fully automated high-throughput in situ analysis.

1.1 Key Specifications

Parameter	Xenium Analyzer Specifications		
Weight			
Xenium Analyzer	~425 lb/192.7 kg	<i>Total weight of system: ~932 lb (422.5 kg)</i>	
Xenium Analysis Computer	~57 lb/25.8 kg		
Vibration Isolation Table	~450 lb/204 kg		
Dimensions			
	L	W	H
Xenium Analyzer	52.5"/133.3 cm	27"/68.5 cm	31"/ 78.7 cm <i>59"/149.8 cm - door open</i>
Xenium Analysis Computer	7"/17.8 cm	26.5"/67.3 cm	18"/ 45.7 cm
Vibration Isolation Table	53"/134.6 cm	30"/76.2 cm	29"/73.6 cm
UPS (<i>APC SRT3000XLT* or similar; not provided by 10x Genomics</i>)	3.4"/8.5 cm	25"/63.5 cm	17"/43.2 cm

**Use equivalent regional models*



For detailed specifications, consult the Xenium Analyzer Site Preparation Survey (CG000587). Specifications will also be available in the Xenium Analyzer User Guide (CG000584).

1.2 Xenium Analyzer - Installation & Training

An overview of the delivery and installation process is provided below.

1. Site Prep Survey (CG000587)

After an introductory call with 10x Genomics, fill & return the Xenium Analyzer Site Preparation Survey to 10x Genomics

2. Site Readiness Visit

On-site visit by a 10x Genomics Rep. to verify that the site is ready to receive the shipment

3. Shipment

Shipment is received on-site & stored without unboxing until installation

*(*shipped items listed below)*

4. Installation (~5-7 d)

Instrument is installed and verified by a 10x Service Engineer

5. On-site Training (~5 d)

On-site training by a 10x FAS

(remote workflow & data analysis trainings will also be provided by 10x and may happen prior to the on-site training; ~2 d)

Ready for Xenium In Situ Gene Expression!


1.3 Xenium Analyzer - Site Preparation

Space Requirements

It is critical to install the instrument in a location away from any vibration sources, such as equipment with compressors (refrigerators, freezers, etc.), motors (centrifuges, shakers, etc.), doors, and busy walkways. Additionally, the installation space should not have soft floor types, such as linoleum or carpet. The instrument should not be placed in direct sunlight or next to other heat generating sources.

The illustration below provides the dimensions and configuration of the space required for installing the Xenium Analyzer. Note that the clearances specified are required for instrument installation, operation, service access, and ventilation. The Xenium Analysis Computer, user-acquired UPS, and compressor (optional) may be placed under the Vibration Isolation Table or adjacent to it with the indicated clearances.

To float the instrument on the Vibration Isolation Table, on-site CDA (compressed dry air) is highly recommended. Alternatively, a low noise/vibration compressor may be used (such as, Air Compressor, Low Noise, 3.5 Liter Capacity, 110 VAC, Model ACGP from Newport or equivalent).



Fill and share the Xenium Analyzer Site Preparation Survey (CG000587) with 10x Genomics which will be followed by more in-depth discussion with a 10x Representative.

1.3 Xenium Analyzer - Site Preparation *contd.*

Power Supply



The Xenium Analyzer and the Xenium Analysis Computer require uninterrupted power supply for a successful run (~2-4 days/run). Standard emergency generator-backed power is often not uninterruptible and a brief power outage is typical before power resumes. Any interruption in the power supply will terminate the run, resulting in the potential loss of samples, reagents, and data that cannot be replaced/recovered by 10x Genomics.

A user-supplied uninterruptible power supply (UPS) is highly recommended but not required during installation. It is recommended that the instrument should be connected to a UPS during runs (provides ~5 min backup power for 2,000 W). Additionally, connecting the UPS to an emergency generator-backed power supply is recommended.

UPS

The recommended UPS specifications are provided below.

- **Minimum UPS power rating (VA):** 3,000 VA
- **Minimum UPS power rating (W):** 2,700 W
- **Nominal Input Voltage:** 208 V
- **UPS design topology type:** On Line
- **Minimum backup run time for 2,000 W:** ~5 min

APC SRT3000XLT or similar UPS may be used.


Follow manufacturer's instructions for UPS setup and ensure that the battery is connected.

1.3 Xenium Analyzer - Site Preparation *contd.*

The configurations of the Xenium Analyzer and the Xenium Analysis Computer with a UPS (recommended) and without a UPS are illustrated below.



UPS should be plugged into an independent circuit. Refer to the user-acquired UPS installation documentation for UPS input circuit requirements.



Power Cords

10x Genomics will ship two region-specific power cords that are compatible with the regional wall sockets and the Xenium Analyzer IEC C20 input and the Xenium Analysis Computer IEC C14 input. If the instrument is connected to a UPS, region-specific compatible power cords in compliance with the local standards need to be acquired by the user.


Network Connectivity

Networking capabilities allow for egress of output files to shared network drives and enable faster troubleshooting via remote support of the Xenium Analyzer. Users will have the ability to enable and disable remote access to their instrument directly. The user needs to inform the IT department of their institution regarding the network/Internet access.

Contact support@10xgenomics.com for additional information.

2.0 FFPE Samples

2.1 Workflow Documents



Documents will be available on the 10x Genomics Support website.

2.2 FFPE Samples - Key Protocol Steps & Timing

»Tissue Sectioning & Section Placement (off-instrument; ~3 h)

Demonstrated Protocol CG000578




Sections placed on the Xenium slide can be stored at room temperature in a desiccator for up to 4 weeks.

»Tissue Section Deparaffinization & Decrosslinking (off-instrument; ~4.5 h)

Demonstrated Protocol CG000580

Steps	Timing
1.1 Buffer Preparation	30 min
1.2 Deparaffinization	3 h (includes 2 h baking step at 60°C)
1.3 Cassette Assembly	10 min
1.4 Decrosslinking	45 min

 Proceed **immediately** to Probe Hybridization, Ligation & Amplification

»Probe Hybridization, Ligation & Amplification (off-instrument; ~2 days)

Refer to the [Probe Hybridization, Ligation & Amplification](#) section for details.

User Guide CG000582

»Xenium Analyzer (on-instrument; ~2-4 days)

Refer to the [Xenium Analyzer](#) section for details.

User Guide CG000584

2.3 FFPE Samples - Reagents & Consumables (not supplied by 10x Genomics)

FFPE Tissue Sectioning & Section Placement

For items with multiple options listed, choose option based on availability and preference. Refer to the manufacturer's website for regional part numbers.

For FFPE Tissue Sectioning & Section Placement			
Item	Description	Vendor	Part Number
<input type="checkbox"/>	Microtome	Epredia HM 355S Automatic Microtome <i>Or any standard histology grade microtome</i>	Fisher Scientific 23-900-672
<input type="checkbox"/>	Microtome blade	Epredia MX35 Premier Disposable Microtome Blades, Low Profile	Fisher Scientific 3052835
<input type="checkbox"/>	Cool-Cut, <i>Optional</i>	Thermo Scientific Cool-Cut	Fisher Scientific 77-112-0
<input type="checkbox"/>	Section transfer system (STS) <i>Optional</i>	Thermo Scientific Section Transfer System (STS),	Fisher Scientific 771200
<input type="checkbox"/>	Probes	Fisherbrand Fine Precision Probe	Fisher Scientific 12-000-153
<input type="checkbox"/>	Forceps	Fisherbrand Curved Medium Point General Purpose Forceps	Fisher Scientific 16-100-110
<input type="checkbox"/>	Blank Slides <i>Optional, for sectioning practice</i>	Superfrost Plus Microscope Slides	Geyer 194242
<input type="checkbox"/>	Water bath	Tissue Floating Bath, Lighted <i>Or equivalent</i>	Fisher Scientific A84600061
		Epredia Digital Round Tissue Section Water bath <i>If using optional Section Transfer System</i>	Fisher Scientific A84600061
<input type="checkbox"/>	Section dryer oven <i>Optional, but recommended</i>	Epredia High Capacity Section Dryer <i>Or equivalent. Thermal cycler may also be used for section drying</i>	Fisher Scientific A84600051
<input type="checkbox"/>	Brushes	Camel Hair Brushes <i>Or equivalent paintbrush</i>	Ted Pella 11859
<input type="checkbox"/>	Fan <i>For drying slides</i>	Personal Rechargeable Fan <i>Or equivalent</i>	Holmes 085-01-0117
<input type="checkbox"/>	Cutting Mat	WellTech Cutting Mat	WellTech Precision Lab -
<input type="checkbox"/>	Wax Trimmer <i>Optional</i>	Electronic Microscopy Sciences Paraffin Block Trimmer Wax Trimmer, 115 VAC	Fisher Scientific NC0310844
Additional Materials			
<input type="checkbox"/>	Razor blades		
<input type="checkbox"/>	Ice bucket (4-5 L)		
<input type="checkbox"/>	Ultrapure/Milli-Q Water for Water Bath, <i>from Milli-Q Integral Ultrapure Water System or equivalent</i>		

This list may not include some standard laboratory equipment.

FFPE Tissue Sections: Deparaffinization & Decrosslinking

For items with multiple options listed, choose option based on availability and preference. Refer to the manufacturer's website for regional part numbers.

For FFPE Tissue Sections: Deparaffinization & Decrosslinking			
Item	Description	Vendor	Part Number
Xylene	Xylene, Reagent Grade	Millipore Sigma	214736
<input type="checkbox"/> or	Xylene, Histological Grade	Millipore Sigma	534056
Neo-clear	Neo-clear Xylene Alternative Substitute	Millipore Sigma	1098435000
<input type="checkbox"/> Ethanol	Ethyl Alcohol, 200 Proof, anhydrous	Millipore Sigma	E7023
	Ethanol absolute ≥99.5%, TechniSolv, pure (Europe)	VWR	83813.360DP
<input type="checkbox"/> Nuclease-free Water	Nuclease-free Water (not DEPC-treated)	Thermo Fisher Scientific	AM9932/ AM9937
<input type="checkbox"/> PBS	PBS - Phosphate Buffered Saline (10X) pH 7.4, RNase-free	Thermo Fisher Scientific	AM9624
<input type="checkbox"/> Urea	Urea Solution, 8M	Millipore Sigma	51457
<input type="checkbox"/> Forceps	Tweezers, 4" Wafer Handling	Excelta Corp	491P-SA-PI
<input type="checkbox"/> Staining jar/dishes	Coplin Jar	VWR	100500-232
	Staining Dishes	VWR	25608-906
<input type="checkbox"/> Section dryer oven	Epredia High Capacity Section Dryer <i>Or equivalent. Thermal cyclers may also be used for section drying</i>	Fisher Scientific	A84600051

Additional Materials


- Water Bath
Alternatively, Thermomixer with SmartBlock -1.5 mL, 2.0 mL from Eppendorf (5382000023/5360000038/5362000035) or equivalent may be used
- Slide drying rack
- Fume Hood
- Vortex
- Ultrapure/Milli-Q Water for Water Bath,
from Milli-Q Integral Ultrapure Water System or equivalent

This list may not include some standard laboratory equipment.

Refer to the [Probe Hybridization, Ligation & Amplification](#) section and the [Xenium Analyzer](#) section for reagents & consumables required. The information in these two sections applies to both FFPE and FF samples.

3.0 Fresh Frozen Samples

3.1 Workflow Documents



FF (Fresh Frozen) Samples

Documents will be available on the 10x Genomics Support website.

3.2 Fresh Frozen Samples - Key Protocol Steps & Timing

»Tissue Sectioning & Section Placement (off-instrument; ~3 h)

Demonstrated Protocol CG000579



Sections placed on the Xenium slide can be stored at room temperature in a desiccator for up to 4 weeks.

»Tissue Section Fixation + Permeabilization (off-instrument; ~2.5 h)

Demonstrated Protocol CG000581

Steps	Timing
1.1 Buffer Preparation	30 min
1.2 Slide Preparation	5 min
1.3 Fixation	30 min
1.4 Permeabilization	65 min
1.5 Cassette Assembly	10 min



Proceed **immediately** to Probe Hybridization, Ligation & Amplification

FF (Fresh Frozen)
Samples

»Probe Hybridization, Ligation & Amplification (off-instrument; ~2 days)

Refer to the [Probe Hybridization, Ligation & Amplification](#) section for details.

User Guide CG000582

»Xenium Analyzer (on-instrument; ~2-4 days)

Refer to the [Xenium Analyzer](#) section for details.

User Guide CG000584

3.3 Fresh Frozen Samples - Reagents & Consumables

(not supplied by 10x Genomics)

Fresh Frozen (FF) Tissue Sectioning & Section Placement

For items with multiple options listed, choose option based on availability and preference. Refer to the manufacturer's website for regional part numbers.

For FF Tissue Sectioning & Section Placement			
Item	Description	Vendor	Part Number
Tissue Freezing			
<input type="checkbox"/>	Isopentane	Isopentane (2-Methylbutane)	Millipore Sigma 270342
<input type="checkbox"/>	Forceps	Specimen Forceps, Straight, 203 mm (8")	VWR 82027-436
		Specimen Forceps, Straight, 152 mm (6")	VWR 82027-438
Frozen Tissue Embedding			
<input type="checkbox"/>	Embedding Compound	TissueTek O.C.T. Compound	VWR 25608-930
<input type="checkbox"/>	Embedding Molds	EpreDia Peel-A-Way Disposable Embedding Molds	Fisher Scientific 12-20
Frozen Tissue Sectioning			
<input type="checkbox"/>	Blank Slides <i>Optional, for sectioning practice</i>	Superfrost Plus Microscope Slides	Fisher Scientific 12-550-15
<input type="checkbox"/>	Cryostat	CryoStar NX70 Cryostat	Fisher Scientific 957020
<input type="checkbox"/>	Brushes	Flat cryostat brush, 10 mm <i>Or equivalent</i>	Fisher Scientific 14-071-00
<input type="checkbox"/>	Specimen Chuck	Thermo Scientific CryoStar NX70 Specimen Chuck	Fisher Scientific 14-071-413
<input type="checkbox"/>	Microtome Blade	MX35 Ultra Microtome Blade, Low Profile	Fisher Scientific 3051835
<input type="checkbox"/>	Slide Mailer	Simport Scientific LockMailer Tamper Evident Slide Mailer	Fisher Scientific 22-038-399
<input type="checkbox"/>	Anti-Roll Plate <i>Optional</i>	Glass Anti-Roll Plate	Fisher Scientific A78930200
Additional Materials			
<input type="checkbox"/>	Dry Ice		
<input type="checkbox"/>	Razor blades		
<input type="checkbox"/>	Ice bucket		
<input type="checkbox"/>	Aluminum Foil		

This list may not include some standard laboratory equipment.

Fresh Frozen (FF) Tissue Sections: Fixation & Permeabilization

For items with multiple options listed, choose option based on availability and preference. Refer to the manufacturer's website for regional part numbers.

For FF Tissue Sections: Fixation & Permeabilization

Item	Description	Vendor	Part Number
<input type="checkbox"/> PBS	PBS - Phosphate Buffered Saline (10X) pH 7.4, RNase-free	Thermo Fisher Scientific	AM9624
<input type="checkbox"/> Nuclease-free Water	Nuclease-free water (not-DEPC treated)	Thermo Fisher Scientific	AM9932/ AM9937
<input type="checkbox"/> Formaldehyde or Paraformaldehyde	Formaldehyde (37% by Weight/Molecular Biology) Paraformaldehyde 16% Aqueous Solution, EM Grade	Thermo Fisher Scientific Electron Microscopy Sciences	BP531-500 15710
<input type="checkbox"/> Ethanol	Ethyl Alcohol, 200 Proof, anhydrous Ethanol absolute ≥99.5%, TechniSolv, pure (Europe Only)	Millipore Sigma VWR	E7023 83813.360DP
<input type="checkbox"/> Tween-20	Tween 20 Surfact-Amps Detergent Solution (10% solution)	Thermo Fisher Scientific	28320
<input type="checkbox"/> Methanol	Methanol, for HPLC	Millipore Sigma	34860
<input type="checkbox"/> SDS	Sodium dodecyl sulfate solution (for molecular biology, 10% in H2O)	Millipore Sigma	71736
<input type="checkbox"/> Forceps	Tweezers, 4" Wafer Handling	Excelta Corp	491P-SA-PI
<input type="checkbox"/> Slide Mailers	Sim port Scientific LockMailer Tamper Evident Slide Mailer	Fisher Scientific	22-038-399

Additional Materials

- Dry Ice
- Thermal Cycler (C1000 Touch Thermal Cycler with 96-Deep Well Reaction Module, Bio-Rad, 1851197)
Currently no alternate recommendations are available; ONLY supported thermal cycler
- Slide drying rack
- Serological pipettes
- Fume Hood
- Vortex
- Ice bucket
- Ultrapure/Milli-Q Water for Water Bath,
from Milli-Q Integral Ultrapure Water System or equivalent

This list may not include some standard laboratory equipment.

FF (Fresh Frozen) Samples

4.0 Probe Hybridization, Ligation & Amplification

4.1 Key Protocol Steps & Timing

(off-instrument; for both FFPE & FF samples)

»Probe Hybridization, Ligation & Amplification (off-instrument; ~2 days)

User Guide CG000582

Steps		Timing	Stop & Store
Day 1			
Step 1: Probe Hybridization			
1.1	Buffer Preparation	20 min	
1.2	Probe Hybridization	16-24 h (overnight)	
Day 2			
Step 2: Post Hybridization Wash			
2.1	Post Hybridization Wash	35 min	
Step 3: Ligation			
3.1	Ligation	~2 h	
Step 4: Amplification			
4.1	Amplification	~2 h	
4.2	Post-Amplification wash	15 min	 4°C overnight or ≤4 days (in the dark)
Step 5: Autofluorescence Quenching			
5.1	Autofluorescence Quenching	45 min	 4°C overnight or ≤4 days (in the dark)
5.2	Nuclei Staining	20 min	 4°C overnight or ≤4 days (in the dark)

Probe Hyb., Lig. & Amplification

4.2 Probe Hybridization, Ligation & Amplification - Reagents & Consumables

(not supplied by 10x Genomics)

For items with multiple options listed, choose option based on availability and preference. Refer to the manufacturer's website for regional part numbers.

For Probe Hybridization, Ligation & Amplification			
Item	Description	Vendor	Part Number
<input type="checkbox"/>	Nuclease-free water	Nuclease-free Water (not DEPC-Treated)	Thermo Fisher Scientific AM9932/ AM9937
<input type="checkbox"/>	TE Buffer	TE Buffer, TRIS-EDTA, 1X Solution, pH 8.0	Thermo Fisher Scientific BP24731
<input type="checkbox"/>	PBS	PBS - Phosphate Buffered Saline (10X) pH 7.4, RNase-free	Thermo Fisher Scientific AM9624
<input type="checkbox"/>	Tween 20	Tween 20 Surfact-Amps Detergent Solution (10% solution)	Thermo Fisher Scientific 28320
<input type="checkbox"/>	Ethanol	Ethyl Alcohol, 200 Proof, anhydrous	Millipore Sigma E7023-500ML
		Ethanol absolute ≥99.5%, TechniSolv, pure (Europe)	VWR 83813.360DP
<input type="checkbox"/>	1.5 ml tubes	DNA LoBind Tubes, 1.5 ml	Eppendorf 022431021
		Low DNA Binding Tubes, 1.5 ml	Sarstedt 72.706.700
<input type="checkbox"/>	2.0 ml tubes	DNA LoBind Tubes, 2.0 ml	Eppendorf 022431048
		Low DNA Binding Tubes, 2.0 ml	Sarstedt 72.695.700
<input type="checkbox"/>	15 ml tubes	15 ml PP Centrifuge Tubes	Corning 730791
<input type="checkbox"/>	50 ml tubes	Self-Standing Polypropylene Centrifuge Tubes (50 ml), sterile	Corning 430921
<input type="checkbox"/>	Pipette tips	Tips LTS 200UL Filter RT-L200 FLR	Rainin 30389240
		Tips LTS 1ML Filter RT-L1000 FLR	Rainin 30389213
		Tips LTS 20UL Filter RT-L20 FLR	Rainin 30389226
<input type="checkbox"/>	Pipettes	Pipet-Lite LTS Pipette L-20XLS+	Rainin 17014392
		Pipet-Lite LTS Pipette L-100XLS+	Rainin 17014384
		Pipet-Lite LTS Pipette L-200XLS+	Rainin 17014391
		Pipet-Lite LTS Pipette L-1000XLS+	Rainin 17014382

Additional Materials

- Water Bath**
Alternatively, Thermomixer with SmartBlock -1.5 mL, 2.0 mL from Eppendorf (5382000023/5360000038/5362000035) or equivalent may be used
- Mini centrifuge**
- Vortex**
- Ice Bucket**
- Ultrapure/Milli-Q Water for Water Bath,**
from Milli-Q Integral Ultrapure Water System or equivalent

This list may not include some standard laboratory equipment.

5.0 Xenium Analyzer

5.1 Key Protocol Steps & Timing

(on-instrument; for both FFPE & FF samples)

»Xenium Analyzer (on-instrument; ~2-4 days)

User Guide CG000584

Steps	Timing	
	Hands-on Time	Total Time
Day 1		
Thaw Decoding Reagents	5 min	16-24 h (overnight)
Day 2		
Prepare Buffers	1 h	1 h
Initialize Instrument	-	5-10 min
Input Experimental Details	5-10 min	5-10 min
Load Instrument	~5 min	~5 min
Overview Scan	-	1 h
Select Region & Initiate Run	~10 min	~10 min
Day 4-6		
Run Time	-	2-4 days
Post-Run Cleanup	5 min	10 min

Xenium Analyzer
(Instrument workflow)

5.2 Xenium Analyzer - Reagents & Consumables (not supplied by 10x Genomics)

For items with multiple options listed, choose option based on availability and preference. Refer to the manufacturer's website for regional part numbers.

For Reagent Bottle Buffer Preparation			
Item	Description	Vendor	Part Number
<input type="checkbox"/>	Nuclease-free Water	Nuclease-free Water (not DEPC-treated)	Thermo Fisher Scientific AM9932/ AM9937
<input type="checkbox"/>	PBS-T	Phosphate Buffered Saline with 0.05% Tween 20, pH 7.4 Phosphate Buffered Saline with 0.05% Tween 20, pH 7.4 (select one based on availability)	Millipore Sigma P3563-10PAK Millipore Sigma PPB005-20PAK
<input type="checkbox"/>	PBS <i>Alternate for making PBS-T</i>	PBS - Phosphate Buffered Saline (10X) pH 7.4, RNase-free	Thermo Fisher Scientific AM9624
<input type="checkbox"/>	Tween 20	Tween 20 Surfact-Amps Detergent Solution (10% solution) (use one ampule per use)	Thermo Fisher Scientific 28320
<input type="checkbox"/>	100% DMSO	Dimethyl sulfoxide (molecular biology grade) Dimethyl sulfoxide (molecular biology grade) Dimethyl sulfoxide, Fisher BioReagents (>99.7%) (select one based on availability)	Millipore Sigma D8418-250ML Millipore Sigma D8418-1L Fisher Scientific BP231-1
<input type="checkbox"/>	KCl	Potassium Chloride (KCl, sterile), 500 ml Potassium Chloride (KCl, sterile), 1L KCl (2 M), RNase-free (conc. in working solution will be 50 mM; select one based on availability)	Teknova P0330 Teknova P0335 Invitrogen AM9640G
Additional Materials			
<input type="checkbox"/>	Centrifuge (with microplate pkg.)	Allegra X-14 Series Benchtop Centrifuge 120 V Or equivalent; fits deep-well 96 well plates (~2 ml vol.)	Beckman Coulter Coulter -
<input type="checkbox"/>	Serological Pipettes	10 ml, 25 ml, 50 ml, 100 ml	
<input type="checkbox"/>	Serological Pipette Controller	Compatible with 10, 25, 50 & 100 ml serological pipettes	
<input type="checkbox"/>	Graduated Cylinders	100 ml and other volumes as needed	
<input type="checkbox"/>	Pipette Tips	Tips LTS 1ML Filter RT-L1000FLR Or equivalent	Rainin 30389213
<input type="checkbox"/>	Pipettes	Pipet-Lite LTS Pipette L-1000XLS+ Or equivalent	Rainin 17014382
<input type="checkbox"/>	Glass Bottles with Cap	Pyrex Reusable Media Storage Bottles (500 ml and 1 L) Or equivalent	
<input type="checkbox"/>	Compressed Canned Air for cleaning		
<input type="checkbox"/>	Lens-cleaning Paper or Lint-free Laboratory Wipes		
<input type="checkbox"/>	70% Isopropanol		
<input type="checkbox"/>	Laboratory Balance		
<input type="checkbox"/>	Ultrapure water	Ultrapure/Milli-Q water, from Milli-Q Integral Ultrapure Water System or equivalent	

This list may not include some standard laboratory equipment.

Appendix

Quencher Removal and H&E Staining - Reagents & Consumables (not supplied by 10x Genomics) Optional; Only if following 10x Genomics H&E protocol

For items with multiple options listed, choose option based on availability and preference. Refer to the manufacturer's website for regional part numbers.

For Quencher Removal				
<i>Post-run quencher removal is required only if staining slides after the instrument run</i>				
Item	Description	Vendor	Part Number	
<input type="checkbox"/>	Sodium hydrosulfite	Sodium hydrosulfite, technical grade	Sigma Aldrich	157953
<input type="checkbox"/>	Forceps	Tweezers, 4' Water Handling	Excelta Corp	491P-SA-PI
<input type="checkbox"/>	PBS <i>(optional)</i>	PBS - Phosphate Buffered Saline (10X) pH 7.4, RNase-free	Thermo Fisher Scientific	AM9624
<input type="checkbox"/>	Slide Mailers	Sim port Scientific LockMailer Tamper Evident Slide Mailer	Fisher Scientific	22-038-399
For H&E Staining				
Item	Description	Vendor	Part Number	
<input type="checkbox"/>	Hematoxylin	Hematoxylin Solution, Mayer's	Sigma Aldrich	MHS16
<input type="checkbox"/>	Eosin	Eosin Y Solution, Alcoholic	Leica	3801615
<input type="checkbox"/>	Bluing Reagent	Bluing Solution	Dako	CS702
<input type="checkbox"/>	Mounting Media	Surgipath SUB-X Mounting Media	Leica	3801741
<input type="checkbox"/>	Ethanol	Ethyl Alcohol, 200 Proof, anhydrous	Millipore Sigma	E7023
		Ethanol absolute ≥99.5%, TechniSolv, pure (Europe)	VWR	83813.360DP
<input type="checkbox"/>	Xylene	Xylene, Reagent Grade	Millipore Sigma	214736
		Xylene, Histological Grade	Millipore Sigma	534056
<input type="checkbox"/>	Forceps	Tweezers, 4' Water Handling	Excelta Corp	491P-SA-PI
<input type="checkbox"/>	Filter Paper	Fisherbrand Qualitative Grade Plain Filter Paper Circles	Fisher Scientific	09-795-H
<input type="checkbox"/>	Coverslips	Fisherbrand Cover Glasses: Rectangles	Fisher Scientific	12-544-EP
		Cover Glasses, Rectangles	VWR	16004-322
Additional Materials				
<input type="checkbox"/>	Vortex			
<input type="checkbox"/>	Staining jar/dishes			
<input type="checkbox"/>	Wide-bore pipette tips			
<input type="checkbox"/>	Ultrapure water	Ultrapure/Milli-Q water, from Milli-Q Integral Ultrapure Water System or equivalent		

This list may not include some standard laboratory equipment.

Document Revision Summary

Document Number	CG000601
Title	Xenium In Situ Gene Expression - Protocol Planner
Revision	Rev A to Rev B
Revision Date	December 2022

Specific Changes

- Minor updates in workflow overview (page 1)
- Updated to include stop and store guidance after tissue sectioning (pages 10, 14)
- Included additional reagents & materials (Isopropanol, alternate vendor for purchasing KCl, laboratory balance - page 20)
- Included reagents and consumables for Quencher Removal and H&E staining (page 21)

General Changes

Updated for general minor consistency of language and terms throughout

© 2022 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document, without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: www.10xgenomics.com/trademarks. 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: www.10xgenomics.com/patents. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of 10x Genomics products in practicing the methods set forth herein has not been validated by 10x Genomics, and such non-validated use is NOT COVERED BY 10X GENOMICS STANDARD WARRANTY, AND 10X GENOMICS HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics terms and conditions of sale for the Chromium Controller or the Chromium Single Cell Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics that it currently or will at any time in the future offer or in any way support any application set forth herein.

Contact:

support@10xgenomics.com

10x Genomics
6230 Stoneridge Mall Road
Pleasanton, CA 94588 USA

