DEMONSTRATED PROTOCOL CG000553 | Rev B

# Tissue Fixation & Dissociation for Chromium Fixed RNA Profiling

## Introduction

This protocol outlines how to perform tissue fixation followed by dissociation of the fixed tissue for use with Chromium Fixed RNA Profiling workflow. Storage recommendations for the fixed tissue and post-storage processing are listed in the Appendix.

## **Additional Guidance**

This protocol was demonstrated using 25 mg of various mouse and human tissue types (flash-frozen and fresh), including brain, heart, liver, spleen, testes, ileum, duodenum, pancreas, kidney, lung, colon, tonsil, lymph, etc. See Appendix for a complete list of tissues tested along with corresponding cell yields. Optimize the protocol based on tissue type, tissue density, and dissociation efficiency.

Tissue and cells carry potentially hazardous pathogens. Follow material supplier recommendations and local laboratory procedures and regulations for the safe handling, storage and disposal of biological materials.

# **Specific Reagents & Consumables**

Vendor	Item Part Number			
For Tissue M	For Tissue Mincing & Fixation			
10x Genomics	Chromium Next GEM Single Cell Fixed RNA Sample Preparation Kit	1000414		
Millipore Sigma	Protector RNase Inhibitor	3335399001		
	Albumin, Bovine Serum, 10% Aqueous Solution, Nuclease- Free	126615		
Thermo Fisher Scientific	UltraPure Bovine Serum Albumin (BSA, 50 mg/ml) Alternative to Millipore Sigma product	AM2616		
	Formaldehyde (37% by Weight/Molecular Biology), Fisher BioReagents	BP531-25		
	Nuclease-free Water (not DEPC-Treated)	AM9937		
	Pyrex Petri Dish*	08-748D		
	General-Purpose Forceps*	10-270		
	Standard Dissecting Scissors*	08-951-20		
	Wide-Bore Pipette Tips RT LTS 1000 uL*	FLW 768A/8 (30389218)		
	*May be procured from alternat	ive vendors		
Corning	Phosphate-Buffered Saline, 1X without Calcium and Magnesium	21-040-CV		
VWR	Tris Buffer, 1M sterile solution, E199-100ML pH 8.0			
Additional Materials				
Blade; Lab Weighing Scale				

This list may not include some standard laboratory equipments.



For Tissue Dissociation			
Millipore Sigma	Liberase TL Liberase TH (optional)	5401020001 5401135001	
Miltenyi Biotec	gentleMACS Octo Dissociator with Heaters gentleMACS C Tubes	130-096-427 130-093-237	
Corning	Corning RPMI 1640	10-040-CV	
Thermo Fisher Scientific	Gibco RPMI 1640	11875093	

Choose either Corning or Gibco RPMI

For Sample Filtration			
Sysmex	Sterile Single-Pack CellTrics Filters (30 µm)	04-004-2326	
Miltenyi Biotec	MACS SmartStrainers (30 μm)	130-098-458	

Choose either Sysmex or Miltenyi Biotec filter.

For Cell Counting			
Nexcelom Biosciences	ViaStain AOPI Staining Solution	CS2-0106-5mL	
	Cellaca MX High-throughput <sup>†</sup> Automated Cell Counter	MX-112-0127	
Thermo Fisher Scientific	Countess II FL Automated Cell Counter <sup>†</sup>	AMAQAF1000	
	Countess II FL Automated Cell Counting Chamber Slides	C10228	
	Ethidium Homodimer-1	E1169	

<sup>&</sup>lt;sup>†</sup>Choose either Countess or Cellaca.

For Storage & Post-Storage Processing			
Acros Organics	Glycerol, 99.5%, for molecular 327255000 biology, DNAse, RNAse and Protease free Alternative to Millipore Sigma product		
Millipore Sigma	Glycerol for molecular biology, ≥99.0% Alternative to Acros Organics product	G5516-100ML	
	Protector RNase Inhibitor	3335402001	
VWR	Vacuum Filter/Storage System (0.2 µm)	29442-936	
Additional Materials			
Eppendorf	DNA LoBind Tubes 2.0 ml	022431048	
	ThermoMixer C	5382000023	

This list may not include some standard laboratory equipments.

# **Preparation - Buffers**

All buffer preparations should be fresh.

Buffers for Fixation - Prepare fresh (1 ml)			
Fixation Buffer Maintain at room temperature	Stock	Final	Per 25 mg tissue (µl)
Nuclease-free Water	-	-	791.9
Conc. Fix & Perm Buffer* (10x Genomics PN 2000517)	10X	1X	100.0
Formaldehyde	37%	4%	108.1
Tissue Resuspension Buffer Maintain at 4°C	Stock	Final	Per sample (µI)
PBS	1X	0.496X	496.0
Tris buffer (pH 8.0; mM)	1000	50	50.0
BSA (RNase free)	10%	0.02%	2.0
RNase Inhibitor (U/µl)	40	0.24	6.0
Nuclease-free Water	-		446.0
Quenching Buffer Maintain at 4°C	Stock	Final	Per 25 mg tissue (µl)
Nuclease-free Water	-	-	875.0
Conc. Quench Buffer* (10x Genomics PN 2000516)	8X	1X	125.0
Dissociation Colution Dynnaus 2 ml			

#### **Dissociation Solution** Prepare 2 ml

- Reconstitute 5mg Liberase by adding 1 ml sterile water.
   Agitate at 2-8 °C until dissolved. Store stock solution in single-use aliquots at -15 to -25°C
- Prepare RPMI + 0.2 mg/ml Liberase (Add 80 µl Liberase stock solution into 1,920 µl of RPMI, mix, maintain at 4°C)
- Warm Dissociation Solution for 10 min at 37°C before use.

#### **Additional Buffers**

PBS + 0.04% BSA (maintain at 4°C)

# **Buffers for Storage of Fixed Samples - Prepare fresh**

#### **50% Glycerol Solution**

(needed for long-term storage of fixed samples)

- Mix an equal volume of water and 99% Glycerol, Molecular Biology Grade
- Filter through a 0.2 µm filter
- Store at room temperature in 2-ml LoBind tubes

#### **Additional Buffers**

0.5X PBS + 0.02% BSA (maintain at  $4^{\circ}$ C) For post-storage processing; use RNase-free BSA

\*Included in the 10x Genomics Chromium Next GEM Single Cell Fixed RNA Sample Preparation Kit, 16 rxns (PN-1000414).

Formaldehyde should always be used with adequate ventilation, preferably in a fume hood. Follow appropriate regulations.

# **Tips & Best Practices**

The following recommendations are critical for optimal performance of the Chromium Fixed RNA Profiling assay.

## **Tissue Quality & Processing**

- Perform pilot experiments to determine if the tissue type is suitable for fixation and dissociation for preparing single cell suspension(s).
- Ensure that variation in cell yields from different tissue types is considered when performing the protocol.
- The tissue should be minced only on a glass surface. DO NOT use plastic petri dish for mincing tissue.

## **Centrifugation & Pellet Resuspension**

- Use a swinging-bucket rotor for higher cell/nuclei recovery.
- Centrifugation speed and time may need optimization depending upon the sample type.
- When working with samples with low cell numbers (i.e. <500,000 cells), complete removal of the supernatant is not required. Up to 30  $\mu$ l supernatant may be left behind to optimize cell recovery without significantly impacting assay performance.
- After each buffer addition step, gently mix cells/ nuclei 5x, or until the pellet is completely resuspended, without introducing bubbles.




## **Fixed Cell Counting**


- Samples should have minimal debris for best results; debris can have associated RNA that can contribute to non-cell background.
- Accurate sample counting is critical for optimal assay performance.
- It is strongly recommended that the sample be stained with a fluorescent nucleic acid dye such as Ethidium Homodimer-1 or AO/PI Staining Solution and counted using an automated fluorescent cell counter (Countess II Automated Cell Counter or a Cellaca counter).
- See Appendix for details on fixed cell counting.

## **Fixed Sample Storage**

- Fixed tissue pieces and cell suspension derived from fixed dissociated tissues can be stored at 4°C for up to 1 week or at -80°C for up to 6 months after resuspending in appropriate reagents.
- Sample storage guidelines are provided in the Appendix.

## **Protocol Overview: Tissue Fixation & Dissociation**





## **Tissue Fixation & Dissociation Protocol**

This protocol was demonstrated using 25 mg fresh or flash frozen tissue. For tissues with high cell density, less than 25 mg may also yield sufficient cell numbers. See Appendix for a complete list of tissues tested along with corresponding cell yields from 25 mg of tissue. Chromium Next GEM RNA Profiling Sample Fixation Kit. 16 rxns (PN-1000414) was used for tissue fixation.

#### 1. Fix Tissue



If using frozen tissue, maintain the tissue on dry ice before and after weighing. Use a pre-chilled petri dish placed on ice while mincing the tissue.

- **a.** Weigh tissue to determine Fixation Buffer volume. 1 ml Fixation Buffer/25 mg tissue will be used at step 1c.
- **b.** Place the tissue on a pre-chilled glass petri dish maintained on ice and using a blade, mince tissue finely (enables passing through a 1 ml wide-bore pipette tip).
- c. Using a wide-bore (1.5 mm) 1-ml pipette, add 1 ml Fixation Buffer/25 mg tissue (see table below). Pass the tissue up and down. Aspirate the solution containing the minced tissue and transfer to a 2 ml or larger centrifuge tube. Pass the tissue up and down. If larger tissue pieces fail to pass through the pipette, use dissection scissors to cut the tissue pieces further until they pass through the pipette. Maintain on ice.

Tissue Mass (mg)	Fixation Buffer (ml)
12-37 mg	1 ml
38-62 mg	2 ml
63-87 mg	3 ml
88-112 mg	4 ml

d. Incubate for 16-24 h at 4°C.



DO NOT agitate or mix the sample during incubation.

Fixation time and temperature should be consistent across all samples in an experiment.

- e. Centrifuge at **850 rcf** for **5 min** at **room** temperature.
- **f.** Remove the supernatant without disturbing the tissue pellet.

- **g.** Add **2 ml** chilled PBS and resuspend the tissue pellet.
- **h.** Centrifuge at **850 rcf** for **5 min** at **room temperature**.
- i. Remove the supernatant without disturbing the tissue pellet.
- **j.** Add **1 ml** Tissue Resuspension Buffer or Quenching Buffer, resuspend the tissue pellet, and maintain on ice.



Fixed tissue pieces can be stored at 4°C for up to 1 week or at -80°C for up to 6 months. See Appendix for guidance on short-term storage and post-storage processing of fixed tissue pieces.

- **k.** Centrifuge at **850 rcf** for **5 min** at **room temperature**.
- **l.** Remove the supernatant without disturbing the fixed tissue pellet.
- m. Proceed to tissue dissociation.

#### 2. Dissociate Fixed Tissue

- **a.** Warm Dissociation Solution for **10 min** at **37°C** before use.
- **b.** Add **2 ml** pre-warmed Dissociation Solution to the sample.



Use 2 ml Dissociation Solution for up to 100 mg tissue. Use only 2 ml Dissociation Solution per C tube.

**c.** Dissociate tissue. Using an Octo Dissociator is highly recommended. Alternatively, tissue may be dissociated manually.

#### • Octo Dissociator:

Transfer to Miltenyi C tubes. Use multiple Miltenyi C tubes if tissue exceeds 100 mg and use 2 ml Dissociation Solution/C tube. Ensure that no fixed tissue pieces are stuck on the walls of the C tube. Place the C tube in the Octo Dissociator, apply Heating Units, and run the following program:

gentleMACS Program			
1	temp ON		
2	spin 50 rpm, 20' 0"		
3	spin 1000 rpm, 30"		
4	spin -1000 rpm, 30"		
5	end		



Dissociation can be monitored by taking a 10  $\mu$ l aliquot for counting. If the cell concentration is lower than expected, additional dissociation may be required. See Troubleshooting guidance.

Detach the C tubes and proceed to step 2d.
 OPTIONAL: Centrifuge the C tubes at 300 rcf briefly to collect all the cells at the bottom of the tube. Resuspend the pellet in the supernatant and proceed to step 2d.

#### OR

#### • Manual Dissociation:

Incubate for **20 min** at **37°C** shaking the tube intermittently. Using a P1000 pipette tip or a silanized glass pipette, triturate the tissue pieces 15-20X (until solution begins to turn cloudy) to obtain a single cell suspension. Proceed to step 2d.



If tissue is tough to triturate and/or solution is not cloudy after trituration, see Troubleshooting guidance

- **d.** Pass the dissociated tissue through a 30 μm filter to remove debris and undissociated tissue pieces.
- **e.** Perform an additional wash of the 30 μm filter by adding **2 ml** PBS to the filter. Collect the filtrate in the same tube as step 2d.



Use of swinging-bucket rotor is recommended for higher cell recovery.

- f. Centrifuge at 850 rcf for 5 min.
- **g.** Remove the supernatant without disturbing the pellet.
- **h.** Resuspend pellet in **1 ml** chilled Quenching buffer.
- Determine the cell concentration using a Countess II FL Automated Cell Counter or hemocytometer. See Appendix for Fixed Cell Counting.



For accurate cell counting, it is strongly recommended that the cell/nuclei suspension be stained with a fluorescent nucleic acid dye such as Ethidium Homodimer-1 or AO/PI Staining Solution and counted using an automated fluorescent cell counter.

If the cell suspension contains aggregates or debris larger than 30  $\mu$ m, pass through a 30  $\mu$ m filter and recount.



For lower than expected cell yields, see Troubleshooting guidance.

j. Proceed immediately to appropriate Chromium Fixed RNA Profiling protocols (see References) or store the sample after resuspending in appropriate reagents.



Samples can be stored at 4°C for up to 1 week or at -80°C for up to 6 months, depending upon the buffers used. See Appendix for guidance on short and long-term storage and post-storage processing of fixed tissue pieces and cells.

# **Appendix**

## **Fixed Sample Storage Guidance**

Fixed samples can be stored for short or long-term.

# Short-term Storage at 4°C (for fixed tissue pieces & cells)

a. Thaw Enhancer for 10 min at 65°C. Vortex and centrifuge briefly. Keep warm and verify no precipitate before use.



DO NOT keep the thawed reagent on ice, or the solution will precipitate. Once thawed, Enhancer can be kept at 42°C for up to 10 min.

- **b.** Add **0.1 volume** pre-warmed Enhancer to fixed tissue pieces or cells in Tissue Resuspension Buffer or fixed cells in Quenching Buffer. For example, add 100 µl Enhancer to 1,000 µl fixed cells in Quenching Buffer. Pipette mix.
- **c.** Store sample at **4°C** for up to **1 week**.

# Long-term Storage at −80°C (for fixed tissue pieces & cells)

a. Thaw Enhancer for 10 min at 65°C. Vortex and centrifuge briefly. Keep warm and verify no precipitate before use.



DO NOT keep the thawed reagent on ice, or the solution will precipitate. Once thawed, Enhancer can be kept at 42°C for up to 10 min.

- **b.** Add **0.1 volume** of pre-warmed Enhancer to fixed sample in Quenching Buffer. For example, add 100 µl Enhancer to 1,000 µl fixed sample in Quenching Buffer. Pipette mix.
- **c.** Add 50% Glycerol for a final concentration of 10%. For example: add 275 µl 50% Glycerol to 1,100 µl fixed sample in Quenching Buffer and Enhancer. Pipette mix.
- **d.** Store at **-80°C** for up to **6 months**.



Storing fixed cells at -80°C is recommended for best results.

## Fixed Tissue Pieces - Post-Storage Processing

Samples may undergo a color change during storage (e.g. black, light gray, or green), however this will not impact assay performance.

If samples were stored at -80°C, thaw at room temperature until no ice is present.

- a. Centrifuge sample at **850 rcf** for **5 min** at **room** temperature.
- **b.** Remove the supernatant without disturbing the tissue pellet.
- **c.** Proceed to Dissociate Fixed Tissue protocol (step 2a).

## Fixed Cells - Post-Storage Processing

Samples may undergo a color change during storage (e.g. black, light gray, or green), however this will not impact assay performance.

If samples were stored at -80°C, thaw at room temperature until no ice is present.

- a. Centrifuge sample at 850 rcf for 5 min at room temperature.
- **b.** Remove the supernatant without disturbing the pellet.
- **c.** Resuspend cell pellet in **1 ml** Quenchinng Buffer or 0.5X PBS + 0.02% BSA (optionally supplemented with 0.2 U/µl RNase Inhibitor) and maintain on ice.



Use RNase-free BSA at this step. See Specific Reagents & Consumables for details.

- **d.** Determine cell concentration of the fixed sample using an Automated Cell Counter (Countess II/Cellaca MX) or hemocytometer. See Appendix for Fixed Cell Counting.
- e. Proceed immediately to appropriate Chromium Fixed RNA Profiling protocols (see References).

# **Fixed Sample Shipping Guidance**

- **a.** Fixed single cells resuspended in Quenching Buffer or fixed tissue pieces resuspended in Tissue Resuspension Buffer and supplemented with Enhancer can be shipped with a cold pack. See Short-term Storage for details.
- **b.** Fixed single cells or fixed tissue pieces resuspended in Quenching Buffer supplemented with Enhancer and Glycerol can be shipped on dry ice. See Long-term Storage for details.

# **Fixed Cell/Nuclei Counting**

- Accurate sample counting is critical for optimal assay performance.
- It is strongly recommended that the fixed sample be stained with a fluorescent nucleic acid dye such as Ethidium Homodimer-1 or AO/PI staining solution and counted using an automated fluorescent cell counter (Countess II Automated Cell Counter or a Cellaca counter).

# **Counting Using Ethidium Homodimer-1**

This protocol provides instructions for counting samples using Ethidium Homodimer-1 and the Countess II FL Automated Cell Counter (with RFP light cube) to enable accurate quantification even in the presence of sub-cellular debris. The optimal cell concentration for the Countess is 1,000-4,000 cells/ $\mu$ l. Refer to manufacturer's instructions for details on operations.

- Vortex Ethidium Homodimer-1, centrifuge briefly, and dilute the concentrated stock as per manufacturer's instructions (~1:100 dilution).
- Aliquot **10 µl** diluted Ethidium Homodimer-1 in a tube.
- Gently mix the sample. Immediately add 10 μl sample to 10 μl diluted Ethidium Homodimer-1.
   Gently pipette mix 10x.
- Transfer **10 µl** sample to a Countess II Cell Counting Slide chamber.
- Insert the slide into the Countess II FL Cell Counter. Image the sample using the RFP setting for fluorescent illumination and filtering.



Optimize focus and exposure settings.

 Confirm the absence of large clumps using the bright-field mode. Make sure the cell counter is circling RFP positive cells. Note the RFP-positive concentration. Multiply by dilution factor 2 to determine cell concentration.

Samples stained with Ethidium Homodimer-1 can also be counted using Cellaca counter. Refer to manufacturer's instructions for details.

# **Counting using AO/PI Staining Solution**

This protocol provides instructions for counting sample using AO/PI staining solution and the Cellaca Counter to enable accurate quantification even in the presence of sub-cellular debris. The optimal cell concentration for the Cellaca Counter is 100-10,000 cells/µl. Refer to manufacturer's instructions for details on operations.

- Add 25 μl AO/PI Staining Solution into Mixing Row of Cellaca plate
- Gently mix the sample. If the sample is too concentrated, a 1:1 dilution in PBS can also be prepared. For example, add 15  $\mu$ l fixed cell suspension to 15  $\mu$ l PBS.
- Add 25 μl sample to Mixing Row of plate containing AO/PI Staining Solution. Gently pipette mix 8x.
- Transfer stained sample to Loading Row of Cellaca plate.

## **Cell Yields from Fixed Dissociated Tissues**

Listed below are the cell yields from 25 mg of various mouse and human tissue types (flash-frozen or fresh) that were processed as described in this protocol. The final cell yields are impacted by the cell density and morphology of the tissue, and by the overall tissue quality prior to fixation. Less than 25 mg of some of the listed tissues may also yield sufficient cells.

	Tissue Type	Cell Yields/25 mg tissue	
rissue rype		Flash-frozen tissue	Fresh tissue
	Adipose	Not recommended Nuclei isolation & fixation recommended	
	Bladder	0.7 x 10 <sup>6</sup>	Not tested
	Brain (Forebrain)	1 x 10 <sup>6</sup>	0.79 x 10 <sup>6</sup>
	Brain (Cerebellum)	Not tested	3.1 x 10 <sup>6</sup>
	Brain (Whole)	1.6 x 10 <sup>6</sup>	0.84 x 10 <sup>6</sup>
	Colon	2.2 x 10 <sup>6</sup>	Not tested
	Eye	1.6 x 10 <sup>6</sup>	Not tested
	Heart	0.99 x 10 <sup>6</sup>	2 x 10 <sup>6</sup>
	Intestine	3.1 x 10 <sup>6</sup>	1.6 x 10 <sup>6</sup>
	Kidney	$2.7 \times 10^6$	2.5 x 10 <sup>6</sup>
Mouse	Liver	1.3 x 10 <sup>6</sup>	0.57 x 10 <sup>6</sup>
Ĭ	Lung	2 x 10 <sup>6</sup>	3.1 x 10 <sup>6</sup>
	Muscle	0.41 x 10 <sup>6</sup>	Not tested
	Ovary	3.2 x 10 <sup>6</sup>	Not tested
		1.8 x 10 <sup>6</sup>	0.92 x 10 <sup>6</sup>
	Pancreas	Digestion using Liberase TH (1 mg/ml) recommended	
	Spinal Cord	1.5 x 10 <sup>6</sup>	Not tested
	Spinal Cord	High debris	s observed
	Spleen	14 x 10 <sup>6</sup>	6.3 x 10 <sup>6</sup>
	Stomach	2.5 x 10 <sup>6</sup>	Not tested
	Testis	4.2 x 10 <sup>6</sup>	Not tested
	Thymus	15 x 10 <sup>6</sup>	26 x 10 <sup>6</sup>
	Tongue	1.4 x 10 <sup>6</sup>	Not tested

	Tissue Type	<b>Cell Yields/25 mg tissue</b> Flash-frozen tissue
	Jejunum	1.37 x 10 <sup>6</sup>
	Colon*	0.04 x 10 <sup>6</sup>
	Lung	3.4 x 106
	Kidney	2 x 10 <sup>6</sup>
	Liver	1.8 x 10 <sup>6</sup>
	Spleen	9.6 x 10 <sup>6</sup>
⊑	lleum	0.14 x 10 <sup>6</sup>
Human	Testes	2.05 x 10 <sup>6</sup>
	Duodenum	0.32 x 10 <sup>6</sup>
	Breast Cancer	0.96 x 10 <sup>6</sup>
	Melanoma	3.18 x 10 <sup>6</sup>
	Ovarian Cancer	0.93 x 10 <sup>6</sup>
	Pancreatic Cancer	1.78 x 10° cells

<sup>\*</sup>Lower than expected yield from 25 mg colon tissue; using >25 mg colon tissue recommended for this assay.

# **Troubleshooting**

# Problem

## Solution

Lower than expected cell yield after dissociation

- Increase the tissue dissociation time.
- If using an Octo Dissociator, after transferring the dissociated cells from the C tube, perform an additional PBS rinse of the C tube and pass the rinse through the 30 µm strainer in step 2e to collect additional cells.
- If using gentleMACS Octo Dissociator, intact/large tissue pieces at the end of run or run fails midway. Run a "spin only" program on the Octo Dissociator with steps 3-4 (from the protocol in step 2C).
- If manually dissociating tough tissue, use the back of a 1-ml syringe plunger to push any undissociated tissue pieces through the 30 µm strainer.
- If cell yield is still low, replace Liberase TL with Liberase TH in Dissociation solution, and increase final Liberase concentration from 0.2 mg/ml to 1 mg/ml. If cell yield is still low, increase the dissociation time.

## References

- 1. Chromium Fixed RNA Profiling Reagent Kit for Singleplexed Samples with Feature Barcode technology for Cell Surface Protein User Guide (CG000477)
- 2. Chromium Fixed RNA Profiling Reagent Kit for Multiplexed Samples User Guide (CG000527)

# **Document Revision Summary**

**Document Number** CG000553

Title Tissue Fixation & Dissociation for Chromium Fixed RNA Profiling

Revision Rev A to Rev B

**Revision Date** February 2023

#### **Specific Changes**

- Updated vendor recommendations for buying RPMI and Liberase TH (optional)
- · Updated step 1c to include tissue weights and corresponding fixation buffer volumes
- Included Quenching Buffer as an option in step 1j (Fix Tissue) and step c (Fixed Cells Post-Storage Processing)
- Updated long term fixed sample storage recommendation to -80°C
- · Updated to include mouse tissues tested
- · Updated Troubleshooting section to include additional guidance

#### **General Changes**

• Updated for general minor consistency of language and terms throughout

© 2023 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: www.10xgenomics.com/trademarks. 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: www.10xgenomics.com/patents. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of 10x Genomics products in practicing the methods set forth herein has not been validated by 10x Genomics, and such non-validated use is NOT COVERED BY 10X GENOMICS STANDARD WARRANTY, AND 10X GENOMICS HERBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics terms and conditions of sale for the Chromium Controller or the Chromium Single Cell Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics that it currently or will at any time in the future offer or in any way support any application set forth herein.

#### **Contact:**

support@10xgenomics.com

10x Genomics 6230 Stoneridge Mall Road Pleasanton, CA 94588 USA

