

Visium CytAssist Spatial Gene Expression for Fixed Frozen – Rehydration, H&E Staining, Imaging & Decrosslinking

Introduction

The Visium CytAssist Spatial Gene Expression for Fixed Frozen assay is designed to analyze mRNA in tissue sections derived from fixed frozen tissue samples. The Visium CytAssist instrument requires a glass slide with intact tissue sections as input. This protocol outlines rehydration, Hematoxylin & Eosin (H&E) staining, imaging, and decrosslinking of tissue for use with the 10x Genomics Visium CytAssist Spatial Gene Expression for Fixed Frozen assay. Rehydrated, stained, and decrosslinked tissue sections are inputs for the downstream Visium CytAssist Spatial Gene Expression workflow.

Additional Guidance

Consult the Visium CytAssist Spatial Gene Expression for Fixed Frozen - Tissue Preparation Guide (CG000663) for complete information on fixing, cryopreserving, embedding, and sectioning of tissue blocks and placing sections on slides. Consult the Visium CytAssist Spatial Gene Expression Imaging Guidelines (CG000521) to verify imaging settings prior to starting this Demonstrated Protocol. After completing this Demonstrated Protocol (CG000662), proceed with the Visium CytAssist Spatial Gene Expression - User Guide (CG000495), Rev D or later.

Contents

2	Reagent Kits	21	2. Tissue Imaging
4	Recommended Thermal Cyclers	21	2.1 Imaging System Recommendations
5	Workflow Overview	22	2.2 Imaging
6	Specific Reagents & Consumables	22	2.3 Coverslip Removal
7	Tips & Best Practices	24	3. Decrosslinking
15	1. Rehydration & H&E Staining	25	3.1 Preparation - Buffers
15	1.1 Preparation - Buffers	26	Protocol Overview
16	Protocol Overview	27	3.2 Destaining
18	1.2 Rehydration	29	3.3 Decrosslinking
19	1.3 H&E Staining	30	Troubleshooting
20	1.4 Coverslipping		

Reagent Kits

Visium Spatial Gene Expression for FFPE Reagent Kits

Refer to SDS for handling and disposal information

Visium CytAssist Slide and Cassettes, 6.5 mm, 2 rxns PN-1000519

Visium CytAssist Slide and Cassettes, 6.5 mm 2 rxns

PN-1000519

Store at ambient temperature

	#	PN
Visium Cassette, 8 port	1	3000811
Visium Tissue Slide Cassette*		
Visium CytAssist moveable gasket small (pre-assembled with translator)	2	3000814
Visium CytAssist moveable translator (pre-assembled with gasket)	2	3000816
Visium CytAssist moveable Cassette, frame	2	3000813
Visium CytAssist Slide Seals, 40 pack*	1	2000284
Visium CytAssist Spatial Gene Expression Slide v2, 6.5 mm	1	2000549

10x
GENOMICS

*Only these items are used in this protocol.

Visium CytAssist Slide and Cassettes, 11 mm, 2 rxns PN-1000518

Visium CytAssist Slide and Cassettes, 11 mm 2 rxns

PN-1000518

Store at ambient temperature

	#	PN
Visium Cassette, 2 port	1	3000812
Visium Tissue Slide Cassette*		
Visium CytAssist moveable gasket large	2	3000815
Visium CytAssist moveable Cassette, frame	2	3000813
Visium CytAssist Slide Seals, 40 pack*	1	2000284
Visium CytAssist Spatial Gene Expression Slide v2, 11 mm	1	2000701

10x
GENOMICS

*Only these items are used in this protocol.

Reagent Kits

Visium Spatial Gene Expression for FFPE Reagent Kits

Refer to SDS for handling and disposal information

Visium FFPE Reagent Kit v2 - Small PN-1000436

Enough reagent is provided for processing two 6.5 mm slides or one 11 mm slide.

Visium FFPE Reagent Kit – Small		
PN-1000436		
Store at -20°C		
	#	PN
○ Amp Mix B	1	2000567
● Extension Enzyme	1	2000389
● Extension Buffer	1	2000409
● RNase Enzyme	1	3000593
● RNase Buffer B	1	2000551
● Tissue Removal Enzyme	1	3000387
● Tissue Removal Buffer B**	1	2000543
● Tissue Removal Buffer Enhancer**	1	2000557
● Decrosslinking Buffer*	1	2000566
● TS Primer Mix B	1	2000537
● Block and Stain Buffer	2	2000554

*Only this reagent is used in this protocol.

**These tubes may not be included in the kit. They are not used in this assay.

10x Genomics Accessories

Product	#	Kit and Part Number	Part Number (Item)
Low Profile Plate Insert	2	Visium CytAssist Reagent Accessory Kit: 1000499	3000823
10x Magnetic Separator	1		120250

Recommended Thermal Cyclers

Supplier	Description	Part Number
Bio-Rad	C1000 Touch Thermal Cycler with 96-Deep Well Reaction Module	1851197
Eppendorf	MasterCycler Pro (discontinued)	North America 950030010 International 6321 000.019
	MasterCycler X50s	North America 6311000010
Thermo Fisher Scientific	Veriti 96-Well Thermal Cycler (discontinued)	4375786
Analytik Jena	Biometra TAdvanced 96 SG with 96-well block (silver, 0.2 mL) and gradient function	846-x-070-241




If using thermal cyclers other than the BioRad C1000, ramp rates should be adjusted for all the steps as described below:


- Eppendorf MasterCycler X50s: 3°C/sec heating and 2°C/sec cooling
- Analytik Jena Biometra TAdvanced: 2°C/sec heating and cooling

Thermal cycler must be able to accommodate the Low Profile Plate Insert (also referred to as the Low Profile Thermocycler Adapter):

- Well depth: 4.5 mm
- Distance between block and heated lid: 12 mm
- Reaction block dimensions 95.5 x 73 mm



Workflow Overview



Visit the [10x Genomics Support website](https://support.10xgenomics.com) for the most current documentation.

Specific Reagents & Consumables

For each item, a number of vendor options are listed. Choose item based on availability and preference.

Substituting materials may adversely affect system performance.

Item	Alternatives/Options	Vendor	Part Number
Ethanol	Ethyl Alcohol, 200 Proof	Millipore Sigma	E7023
	Ethanol absolute ≥99.5%, TechniSolv, pure (Europe Only)	VWR	83813.360DP
Eosin	Eosin Y-solution, Alcoholic	Millipore Sigma	HT110116
Hematoxylin	Hematoxylin Solution, Mayer's	Millipore Sigma	MHS16
Bluing Reagent	Bluing Reagent, Dako	Agilent	CS70230-2
1X PBS	Phosphate-Buffered Saline, 1X without calcium and magnesium, PH 7.4	Corning	21-040-CV
Glycerol	Glycerol Solution	Millipore Sigma	49781
	Glycerol	Acros Organics	327255000
0.1 N HCl	Hydrochloric Acid Solution, 0.1 N <i>Or any equivalent HCl</i>	Fisher Chemical	SA54-1
Coplin Jar/ Staining Dishes	Coplin Jar	VWR	100500-232
	Staining Dishes	VWR	25608-906
Green Marker, <i>Optional, if annotating slide</i>	Sharpie Argyle Green Permanent Marker	Sharpie	1785396
Slide Holders	Slide Holders, 24-place	VWR	25608-868
Coverslips	Fisherbrand Cover Glasses: Rectangles	Fisher Scientific	12-544-EP
	Cover Glasses, Rectangular	VWR	16004-322
Pipettes	Pipet-Lite Multi Pipette L8-200XLS+	Rainin	17013805
	Pipet-Lite LTS Pipette L-2XLS+	Rainin	17014393
	Pipet-Lite LTS Pipette L-10XLS+	Rainin	17014388
	Pipet-Lite LTS Pipette L-20XLS+	Rainin	17014392
	Pipet-Lite LTS Pipette L-100XLS+	Rainin	17014384
	Pipet-Lite LTS Pipette L-200XLS+	Rainin	17014391
	Pipet-Lite LTS Pipette L-1000XLS+	Rainin	17014382
Wide Bore Pipette Tips	Tips RT LTS 200UL FLW	Rainin	30389241
	Tips RT LTS 1000UL FLW	Rainin	30389218
Pipette Tips	Tips LTS 200UL Filter RT-L200 FLR	Rainin	30389240
	Tips LTS 1ML Filter RT-L1000 FLR	Rainin	30389213
	Tips LTS 20UL Filter RT-L20 FLR	Rainin	30389226
Additional Materials			
1000-ml Beakers (6)	-	-	-
Ultrapure/Milli-Q Water, <i>from Milli-Q Integral Ultrapure Water System or equivalent</i>	-	-	-

Tips & Best Practices



Tips & Best Practices section includes additional guidance



Signifies critical step requiring accurate execution



Troubleshooting section includes additional guidance

General Reagent Handling


- Thoroughly mix reagents before use.

Pipette Calibration

- Follow manufacturer’s calibration and maintenance schedules.

Slide Handling

- Always wear gloves when handling slides.
- DO NOT touch the tissue sections on slides.
- Keep tissue slides in a container on dry ice.
- Minimize exposure of slides to sources of particles and fibers.
- When immersing slides in prepared solutions and water, ensure that the tissue sections are completely submerged.
- Keep slides flat on the bench when adding reagents to the tissue.
- Ensure that no absorbent surface is in contact with the reagents on slides during incubation.



Slide Incubation Guidance

Incubation at a specified temperature

Incubation using a Thermal Cycler:

- Position a Low Profile Plate Insert (also referred to as Low Profile Thermocycler Adapter) on a thermal cycler that is set at the incubation temperature. Move Low Profile Thermocycler Adapter back and forth to ensure that it is seated properly.
- Ensure that the Low Profile Thermocycler Adapter is in contact with the thermal cycler surface uniformly.
- When incubating a slide, position the slide on the Low Profile Thermocycler Adapter with the tissue surface facing up.
- Ensure that the entire bottom surface of the slide is in contact with Low Profile Thermocycler Adapter.
- When incubating a slide encased in a cassette, place the assembled unit on the Low Profile Thermocycler Adapter with the wells facing up. Cassettes should always be sealed when on the Low Profile Thermocycler Adapter.
- Allow Low Profile Thermocycler Adapter to cool before removing it from the thermal cycler.



Incubation at room temperature



- Place the Visium CytAssist Tissue Slide Cassette on a flat, clean, non-absorbent work surface.

Visium CytAssist Tissue Slide Cassette

- The Visium CytAssist Tissue Slide Cassette encases the slide and creates a leakproof well for adding reagents on tissue slides.
- The cassette is a single use item.
- Gaskets are adjusted by the user to ensure that the tissue section or area of interest is encased in a well.
- Refer to Visium CytAssist Tissue Slide Cassette Assembly & Removal instructions for details.



- Prior to use, inspect the moveable gasket to ensure that the gasket perimeter and corners are free of excess silicone prior to assembly.
- Excess silicone should be safely removed with forceps or a pipette tip prior to assembly.
- Assembly should occur against a white background for easy tissue visualization during alignment.
- Practice assembly with a blank slide (75 x 25 x 1 mm).
- Place slides in the cassette only when specified.



Visium CytAssist Tissue Slide Cassette Assembly




Wear fresh gloves while assembling Tissue Slide Cassette



Exercise caution when handling slide edges to prevent injury.

Break cassette into two halves by bending each half at the hinge until they snap apart




Place tissue slide into lower half of cassette with tissue facing up




The 6.5 mm gasket can be adjusted horizontally and rotated in two directions (180°) while 11 mm gasket can be rotated in four directions (90°). Determine the appropriate configuration that allows the gasket to encompass the tissue area of interest.




Securely combine gasket with top half of cassette until the gasket snaps into place.




Gently place top half of cassette over bottom half. DO NOT assemble together until Step 7.




Adjust gasket such that gasket is over the tissue area of interest. The 6.5 mm gasket can be adjusted horizontally as well as vertically.



Apply even pressure on top of cassette until it clicks shut. Verify that clip is completely secured over hinges.




Turn cassette over and verify tissue area of interest is within gasket. DO NOT move gasket once cassette is closed. If necessary, open cassette and recenter gasket.



Visium CytAssist Tissue Slide Cassette Removal

Pull clip up to detach upper and lower halves of cassette



Open cassette by continuing to lift clip upward. If slide sticks to gasket, continue to apply even upward pressure to separate slide from gasket




Hold slide by the label and lift slide out from lower half




Reagent Addition to Wells

- Place the assembled slide in the Visium CytAssist Tissue Slide Cassette flat on a clean work surface.
- Dispense reagents along the side of the wells without touching the tissue sections and without introducing bubbles.
- Always cover the tissue section completely when adding reagents to the well. A gentle tap may help spread the reagent more evenly.



Reagent Removal from Wells

- Slightly tilt the cassette while removing the reagent.
- Place the pipette tip to the bottom edge of the wells.
- Remove reagents along the side of the wells without touching the tissue sections and without introducing bubbles.
- Remove all liquid from the wells in each step. To ensure complete removal, check the bottom of the well by tilting the cassette slightly. A meniscus at the bottom of the well will indicate the presence of liquid in the well. Repeat removal steps until no reagent remains.



Visium CytAssist Slide Seal Application & Removal


Application

- If applying a Visium Slide Seal to a Tissue Slide Cassette, the seal must be cut in half lengthwise.
- Four pre-cut seals per tissue section are needed for the entire Visium CytAssist Spatial Gene Expression assay.
- Cut the seal as shown in the image below. Ensure scissors are cleaned with a RNase decontamination solution.
- Place the CytAssist Tissue Slide Cassette flat on a clean work surface.
- Remove the back of the adhesive Visium Slide Seal.
- Align the Visium Slide Seal with the surface of the cassette and apply while firmly holding the cassette with one hand.
- Press on the Visium Slide Seal to ensure uniform adhesion.


Removal

- Place the CytAssist Tissue Slide Cassette flat on a clean work surface.
- Carefully pull Visium Slide Seal up and over from the edge while firmly holding the cassette.
- Ensure that no liquid splashes out of the wells.


Cut Slide Seal in Half Lengthwise



Visium Slide Seal Application



Visium Slide Seal Removal



Visium CytAssist Tested Slides

The following slides have been tested for use with the Visium CytAssist Tissue Slide Cassette and instrument.

Item	Length (mm)	Width (mm)	Thickness (mm)
Epredia Shandon ColorFrost Plus	75.0	25.0	1.0
Fisherbrand SuperFrost Plus	75.0	25.0	1.0
Sigma-Aldrich Poly Prep Slides	75.0	25.0	1.0
VWR SuperFrost Plus Micro Slide, Premium	75.0	25.0	1.0


If unsure of slide part number, refer to the "generic" slide diagram below for general guidance (images not to scale). Diagrams for verifying that tissue sections are placed in the allowable area can also be found in the Visium CytAssist Quick Reference Cards - Accessory Kit (CG000548). The diagrams demonstrate allowable areas that are far enough away from frosted sections to not interfere with gasket closure during the CytAssist assay. Frosted sections include the opaque area of the slide as well as any etching or writing on the slide.



While slides are specified as being 25 mm x 75 mm, manufacturing tolerances may lead to dimensions that are too small or large to be compatible with 10x Genomics products. Tissue slide dimensions must be within 24.8 mm - 25.3 mm in width and 74.4 mm - 76.2 mm in length to fit the Visium CytAssist Tissue Slide Cassette.

Minimum slide dimensions: 24.8 x 74.4 mm


Maximum slide dimensions: 25.3 x 76.2 mm



1. Rehydration & H&E Staining

1.0 Overview

This chapter provides guidance on rehydration and H&E staining of tissue slides containing fixed frozen tissue sections. Ensure that microscope settings have been verified and imaging programs have been created prior to starting this protocol. Consult the Visium Spatial Gene Expression Imaging Guidelines Technical Note (CG000521) for more information.



1.1 Preparation - Buffers

For Rehydration	
Prepare fresh. If using 50-ml centrifuge tubes, two slides can be faced back to back (tissue facing out) per tube.	
Items	Preparation & Handling
<input type="checkbox"/> 100% Ethanol	Label one 50-ml centrifuge tube as 100% Ethanol Tube. Dispense 30 ml 100% ethanol. <i>Alternatively, use a coplin jar.</i>
<input type="checkbox"/> 70% Ethanol	Label one 50-ml centrifuge tube as 70% Ethanol Tube. Dispense 30 ml 70% ethanol. <i>Alternatively, use a coplin jar.</i>
<input type="checkbox"/> Milli-Q or UltraPure Water	Label two 50-ml centrifuge tubes as Water Tube 1 and 2. Dispense 30 ml water in each. <i>Alternatively, use a coplin jar.</i>
<input type="checkbox"/> 1X PBS	Label one 50-ml centrifuge tube as 1X PBS Tube. Dispense 30 ml 1X PBS. <i>Alternatively, use a coplin jar.</i>



Reagent volume of 30 ml is sufficient for both 50-ml centrifuge tubes and coplin jars.

For H&E Staining

Items	Preparation & Handling
<input type="checkbox"/> Milli-Q or UltraPure Water	Label six 1000-ml beakers as Water Beakers 1 – 6. Dispense 800 ml of water into each beaker. Dispensed volumes in each beaker can be used for two slides.

For Coverslipping


Items	Preparation & Handling
<input type="checkbox"/> Mounting Medium	The dilution below is not necessary if stock glycerol is already at 85%. Invert or pipette gently with a wide-bore pipette tip to mix. Briefly centrifuge to remove bubbles.

Mounting Medium	Stock	Final	1X (µl)	2X+ 15% (µl)
Glycerol	100%	85%	85.0	195.5
Nuclease-free Water	100%	15%	15.0	34.5
Total	-	-	100.0	230.0




DO NOT let the attached coverslip dry.
DO NOT use Cytoseal or nail polish for securing the coverslip.

Protocol Overview



H&E Staining & Coverslipping ⌚ ~15 min




1.2 Rehydration

- a. Place a Low Profile Thermocycler Adapter on a thermal cycler and preheat thermal cycler to **37°C**.
- b. Retrieve slide containing tissue from **-80°C** and place on dry ice.
- c. Place slide on the Low Profile Thermocycler Adapter with the tissue side facing up and incubate **10 min** at **37°C**. DO NOT close the thermal cycler lid.



Ensure tissue is dry before proceeding to next step. Damp tissue may lead to detachment. DO NOT exceed 20 min incubation time.




- d. Remove slide from thermal cycler. Gently immerse slide in the 1X PBS Tube and incubate for **5 min**.



When immersing slides in solutions, ensure that the tissue sections are completely submerged.


- e. Gently immerse slide in the Milli-Q Water Tube 1 and incubate for **3 min**.
- f. Gently immerse slide in the 100% Ethanol Tube for **3 min**.
- g. Gently immerse slide in the 70% Ethanol Tube for **3 min**.
- h. Gently immerse slide in the Milli-Q Water Tube 2 for **20 sec**.
- i. Proceed **immediately** to H&E Staining & Coverslipping.



To ensure even staining, DO NOT let the slides dry.

1.3 H&E Staining

- a. Place slide on a flat, clean, non-absorbent work surface. Some residual droplets may remain.
- b. Add **1 ml** Hematoxylin per slide to uniformly cover the entire tissue section.
- c. Incubate **3 min** at **room temperature**.
- d. Discard reagent by draining and/or holding slide at an angle with the bottom edge in contact with a laboratory wipe.



- e. Immerse slide 5x in Water Beaker 1.
- f. Immerse slide 15x in Water Beaker 2.
- g. Immerse slide 15x in Water Beaker 3.
- h. Wipe excess liquid from the back of the slide without touching the tissue section.
- i. Place slide on a flat, clean, non-absorbent work surface. Some droplets may remain.
- j. Add **1 ml** Bluing Buffer per slide to uniformly cover the entire tissue section. Incubate **1 min** at **room temperature**.
- k. Discard reagent by draining and/or holding slide at an angle with the bottom edge in contact with a laboratory wipe.
- l. Immerse slide 10x in Water Beaker 4.
- m. Wipe excess liquid from the back of the slide without touching the tissue section. Place slide on a flat, clean, non-absorbent work surface. Some droplets may remain.
- n. Add **1 ml** Eosin per slide to uniformly cover the entire tissue section. Incubate **1 min** at **room temperature**. DO NOT use diluted Eosin.
- o. Discard reagent by draining and/or holding slide at an angle with the bottom edge in contact with a laboratory wipe.
- p. Immerse slide for **30 sec** in Water Beaker 5.
- q. Immerse slide 10x in Water Beaker 6.


- r. Wipe excess liquid from the back of the slide without touching the tissue section.




DO NOT air dry the slides.

1.4 Coverslipping

- Place slide on a flat, clean, non-absorbent work surface. Some residual droplets may remain.
- Using a **wide-bore** pipette tip, add **100 µl** Mounting Medium to uniformly cover the entire tissue section.
- Apply the coverslip at an angle on one end of the slide. Slowly lower the coverslip, without introducing bubbles. Allow Mounting Medium to spread and settle.
- If needed, remove any large excess of Mounting Medium by carefully wicking away from the edge of the coverslip with a laboratory wipe. Be careful not to move the coverslip and disturb the tissue.




- e. Once coverslipping is complete, **immediately** proceed with imaging or store slide laying flat in a squared slide mailer or tissue slide book at **4°C** in the dark for up to **two weeks**. Ensure multiple slides do not come in contact with one another. Storage can be done before or after imaging. **DO NOT** exceed two weeks of storage time.



2. Tissue Imaging

2.0 Overview

This chapter provides guidance on imaging tissue slides containing H&E stained fixed frozen tissue sections and coverslip removal.



2.1 Imaging System Recommendations

The following table shows imaging systems used by 10x Genomics in the development of this protocol. Any equivalent imaging system can be used as an alternative.

Supplier	Model	Configuration
Thermo Fisher Scientific	EVOS M7000	Inverted
Leica	Aperio Versa 8	Upright
	Leica DMI8	Inverted
MetaSystems	Metafer	Upright
Nikon	Nikon Eclipse Ti2	Inverted
BioTek	Cytation 7	Inverted or Upright
Keyence	Keyence BZX800	Inverted
Olympus	VS200	Upright
Zeiss	Imager.Z2	Upright

Brightfield Recommended Configuration
Color camera (3 x 8 bit, 2,424 x 2,424 pixel resolution)
White balancing functionality
2.18 µm/pixel minimum capture resolution
Exposure times 2-10 milli sec


2.2 Imaging

- Image each tissue section individually at the desired magnification using brightfield imaging settings.
- Consult the Visium CytAssist Spatial Gene Expression Imaging Guidelines Technical Note (CG000521) for additional information.
- After imaging, proceed **immediately** to Coverslip Removal or store slide laying flat in a squared slide mailer or tissue slide book at **4°C** in the dark for up to **two weeks**. Storage can be done before or after imaging. **DO NOT** exceed two weeks of storage time.



2.3 Coverslip Removal

- Dispense **800 ml** Milli-Q Water in a beaker. Up to 10 slides may be processed using this beaker.
- Immerse slide sideways/horizontal in the beaker containing **800 ml** water with the coverslipped surface fully sideways.
- Hold slide in water until the coverslip slowly separates away from the slide.



To avoid damaging the tissue sections or causing tissue detachment, **DO NOT** move the slide up and down, shake forcibly, or manually move the coverslip.

- Gently immerse slide 30x in water to ensure all Mounting Medium is removed.
- Wipe the back of the slide with a laboratory wipe. Place on a flat, clean, non-absorbent work surface and air dry for a minimum of **5 min** until the tissue is mostly dry. **DO NOT** exceed **20 min**.

- f. Incubate slide on the Low Profile Thermocycler Adapter with the thermal cycler lid open for **3 min** at **37°C** to dry the slide.




- g. Proceed **immediately** to Decrosslinking.

3. Decrosslinking

3.0 Overview


This chapter provides guidance on performing decrosslinking to release RNA that was sequestered by formalin fixation from rehydrated and H&E stained tissues. Ensure that the coverslip is removed and tissue section is dry before destaining and decrosslinking.

Once sections have been decrosslinked, step 1 (Probe Hybridization) of the Visium CytAssist Spatial Gene Expression User Guide (CG000495) should be immediately performed.




3.1 Preparation - Buffers

For Decrosslinking						
Items	Preparation & Handling					
<input type="checkbox"/> 0.1 N HCl	If necessary, prepare 0.1 N HCl using nuclease-free water.					
<input type="checkbox"/> Diluted Decrosslinking Buffer	Thaw Decrosslinking Buffer at room temperature. Vortex and centrifuge briefly after preparing Diluted Decrosslinking Buffer. Store excess stock buffer at 4°C.					




6.5 mm Gaskets						
Diluted Decrosslinking Buffer	Stock	Final	1X (µl)	2X+ 10% (µl)	4X+ 10% (µl)	
Decrosslinking Buffer	10X	1X	25	55	110	
Nuclease-free Water	-		225	495	990	
Total	-	-	250	550	1,100	



11 mm Gaskets						
Diluted Decrosslinking Buffer	Stock	Final	1X (µl)	2X+ 10% (µl)	4X+ 10% (µl)	
Decrosslinking Buffer	10X	1X	50	110	220	
Nuclease-free Water	-		450	990	1,980	
Total	-	-	500	1,100	2,200	

Protocol Overview

Destaining and Decrosslinking ~1 h



Once the sections have been decrosslinked, proceed directly to step 1 (Probe Hybridization) of the Visium CytAssist Spatial Gene Expression User Guide (CG000495)

3.2 Destaining

■ denotes volumes for 6.5 mm gaskets and ▲ denotes volumes for 11 mm gaskets

- a. Place a Low Profile Thermocycler Adapter in the thermal cycler. Prepare the thermal cycler with the following incubation protocol and start the program.

Lid Temperature	Reaction Volume	Run Time
42°C (lid may be set to lowest setting if instrument does not enable 42°C)	100 µl	15 min

Step	Temperature	Time
Pre-equilibrate	42°C	Hold
Destaining	42°C	00:15:00
Hold	22°C	Hold

- b. Place the slide in the Visium CytAssist Tissue Slide Cassette.




See *Tips & Best Practices* for assembly instructions. Practice assembly with a blank slide.

- c. Add ■150 µl or ▲300 µl 0.1 N HCl along the side of the wells to uniformly cover the tissue sections, without introducing bubbles. Tap cassette gently to ensure uniform coverage.
- d. Remove HCl from the wells.



- e. Add ■100 µl or ▲200 µl 0.1 N HCl along the side of the wells to uniformly cover the tissue sections, without introducing bubbles. Tap cassette gently to ensure uniform coverage.
- f. Apply pre-cut slide seal on cassette and place the cassette on the Low Profile Thermocycler Adapter at 42°C.



- g. Close the thermal cycler lid. Skip Pre-equilibrate step to initiate Destaining.

- h.** Remove the cassette from the Low Profile Thermocycler Adapter and place on a flat, clean work surface. Some color remaining in the tissue after Destaining is normal.




3.3 Decrosslinking

- a. Prepare the thermal cycler with the following incubation protocol and start the program.

Lid Temperature	Reaction Volume	Run Time
70°C	100 µl	30 min
Step	Temperature	Time
Pre-equilibrate	70°C	Hold
Decrosslinking	70°C	00:30:00
Cooling	22°C	00:10:00
Hold	22°C	Hold

- b. Peel back slide seal and using a pipette, remove all the HCl from the well corners.



- c. Add **150 µl or 300 µl** Diluted Decrosslinking Buffer along the side of the wells.
- d. Remove Diluted Decrosslinking Buffer from the wells.
- e. Add **100 µl or 200 µl** Diluted Decrosslinking Buffer along the side of the wells.
- f. Re-apply slide seal on the cassette and place the cassette on the Low Profile Thermocycler Adapter at **70°C**.
- g. Close the thermal cycler lid. Skip Pre-equilibrate step and initiate Decrosslinking.
-  h. Proceed **immediately** to Visium CytAssist Spatial Gene Expression User Guide (CG000495).

Troubleshooting

Notes

Tissue Detachment

- Ensure compatible glass slides are used to minimize tissue detachment. Refer to Visium CytAssist Tested Slides.

Bubbles

- Avoid bubble formation during coverslipping. Introduction of bubbles can be mitigated by applying the coverslip at an angle and slowly lowering it onto the slide, allowing air to escape. Briefly centrifuge Mounting Medium to remove bubbles before use. Avoid introducing bubbles when pipetting Mounting Medium onto slide.




Bubbles may cause blackening of tissue

Slide in image is representative

Uneven Staining


- Ensure fresh staining reagents are applied to the tissue uniformly and adequate washes are performed. A gentle tap may help reagents spread uniformly. Consider filtering hematoxylin.



Slide in image is representative

Incorrect Staining Protocol

- Ensure that the correct staining protocol with fresh reagents was followed.



Incorrect staining protocol (right image) may result in poor staining performance.

Slide in image is representative

Document Revision Summary

Document Number	CG000662
Title	Visium CytAssist Spatial Gene Expression for Fixed Frozen – Rehydration, H&E Staining, Imaging & Decrosslinking Demonstrated Protocol
Revision	Rev A
Revision Date	March 2023

LEGAL NOTICE © 2023 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: www.10xgenomics.com/trademarks. 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: www.10xgenomics.com/patents. The use of products described herein is subject to 10x Genomics Terms and Conditions of Sale, available at www.10xgenomics.com/legal-notices, or such other terms that have been agreed to in writing between 10x Genomics and user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of 10x Genomics products in practicing the methods set forth herein has not been validated by 10x Genomics, and such non-validated use is NOT COVERED BY 10X GENOMICS STANDARD WARRANTY, AND 10X GENOMICS HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics terms and conditions of sale for the Chromium Controller or the Chromium Single Cell Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics that it currently or will at any time in the future offer or in any way support any application set forth herein.

Contact:

support@10xgenomics.com

10x Genomics

6230 Stoneridge Mall Road

Pleasanton, CA 94588 USA

