

TECHNICAL NOTE

Sequencing Metrics & Base Composition of Single Cell 5' v2 Dual Index Libraries

Introduction

The Chromium Next GEM Single Cell 5' v2 (dual index), standard, and High Throughput (HT) workflow produces sequencing-ready Gene Expression, V(D)J, Cell Surface Protein, CRISPR, and/or Barcode Enabled Antigen Mapping (BEAM) libraries from the same single cells. This enables simultaneous profiling of cellular features along with gene expression profiling. This Technical Note presents a comparison of sequencing metrics for various Single Cell 5' v2 Dual Index library types across Illumina platforms. Results may vary depending on the sequencing instrument and/or sample and loading characteristics.

Single Cell 5' v2 Dual Index Libraries

The dual index library types that can be generated using the Chromium Next GEM Single Cell 5' v2 (dual index) or the Chromium Next GEM Single Cell 5' HT v2 (dual index) reagents and protocols are shown in the schematics below.

The libraries include cDNA/V(D)J insert or Feature Barcode constructs which begin with P5 and end with P7 sequences necessary for binding to the Illumina flow cell. Read 1 is used to sequence 16 bp 10x Barcodes and 10 bp UMI and Read 2 is used for priming and sequencing the cDNA insert or the Feature Barcode. The two 10 bp sample indexes are sequenced in the i5 and i7 reads.

Chromium Single Cell 5' v2
Gene Expression
Dual Index Library



Chromium Single Cell 5' v2
V(D)J Amplified
Dual Index Library



Chromium Single Cell 5' v2
Cell Surface Protein
Dual Index Library



Chromium Single Cell 5' v2
CRISPR Screening
Dual Index Library



Chromium Single Cell 5' v2
BEAM Dual Index Library



Methods Overview

Single Cell 5' v2 dual index Gene Expression libraries alone or in combination with V(D)J, Cell Surface Protein, CRISPR Screening, and/or BEAM libraries were generated from both standard and HT workflows as described in the respective User Guides (see References). The libraries were sequenced in the following combinations:

- Gene Expression, V(D)J, and Cell Surface Protein libraries
- Gene Expression and V(D)J libraries
- V(D)J and Cell Surface Protein libraries
- Gene Expression and CRISPR Screening libraries
- Gene Expression, V(D)J, and BEAM libraries
- Gene Expression, V(D)J, Cell Surface Protein, and BEAM libraries

Libraries were generated from a target of 1,000 human peripheral blood mononuclear cells (PBMC) using the standard Single Cell 5' v2 assay. Libraries were generated from a target of 20,000 human bone marrow cells using the Single cell 5' v2 HT assay.

These combinations were selected as the most representative library pooling strategies based on common experimental designs. The libraries were quantified and sequenced as indicated in the results (Tables 1-6).

Results Overview

Tables 1-4 show representative sequencing metrics and base composition data derived from the indicated libraries. The Q30 quality scores, representative Data by Cycle plots, and other metrics for each sequencer/workflow is shown for the standard assay. Figure 1 shows representative "Data by Cycle" plots from the Illumina SAV software displaying the percentage of base calls and Q30 quality scores of

Single Cell 5' v2 standard and the HT derived Gene Expression, V(D)J, and Cell Surface Protein libraries run on the NovaSeq. Individual results may vary depending on the specific sequencing instrument and/or particular sample and loading characteristics.

Conclusions

In summary, % Bases by cycle and % \geq Q30 Quality Score distribution showed highly consistent profiles for all sequencing platforms and workflows tested. Furthermore, sequencing performance between the library types generated using the Single Cell 5' v2 standard and the HT assay are comparable. These data serve as guidelines for assessing the quality of Single Cell 5' v2 Dual Index library sequencing. Additional factors that may contribute to overall success of a sequencing run and impact downstream application performance metrics include:

- Sample preparation to obtain a high quality single cell suspension.
- Final libraries with fragment lengths within the expected size range for each library type, for optimal cluster formation on Illumina flow cells.
- Reliable and accurate library quantification using the KAPA DNA Quantification Kit using the average insert size determined by Agilent Bioanalyzer or Agilent TapeStation QC.
- Sequencing platform loading concentration.

Gene Expression, V(D)J, & Cell Surface Protein Dual Index Libraries

Four Chromium Single Cell 5' Gene Expression, eight V(D)J (four TCR amplified, four BCR amplified) and four Cell Surface Protein (dual index) libraries were pooled and sequenced on indicated Illumina sequencers. Table 1 shows 'Data by Cycle' plots from the Illumina SAV software displaying the percentage of base calls and Q30 quality scores along with additional metrics.

Sequencing configuration & run parameters

Minimum sequencing depth:

Gene Expression 20,000 read pairs/cell; V(D)J amplified 5,000 read pairs/targeted cell*; Cell Surface Protein 5,000 read pairs/cell



Paired-end, dual indexing




- Read 1: 26 cycles
- i7 Index: 10 cycles
- i5 Index: 10 cycles
- Read 2: 90 cycles

Table 1: Representative Plots and Sequencing Data

Plots shown are from a pool of four Gene Expression libraries, eight V(D)J libraries, and four Cell Surface protein libraries sequenced on a NovaSeq S2 flow cell.

NovaSeq S2



		% ≥ Q30				Yield per Lane (Gb)		Reads Mapped to Reference (%)		
		R1	i7	i5	R2	R1	R2	Gene Expression	V(D)J	Cell Surface Protein
HiSeq 2500 RR										
	Loading Conc. (pM): 10 Cluster density: 973 K/mm ² Phix (%): 1	97.5	96.7	95.3	95.8	4.1	14.6	98.3±0.2	78.1±11.0	96.7±0.2
NovaSeq 6000										
	Loading Conc. (pM): 300 % PF**: 79.5 Phix (%) 1	93.2	94.1	89.6	90.5	13.7	45.3	86.7±1.8	79.4±2.0	96.3±0.2
NextSeq 2000										
	Loading Conc. (pM): 650 % PF**: 84.0 Phix (%) 1	92.5	92.2	92.7	93.1	15.2	50.1	89.1±1.5	75.8±14.0	96.2±0.2

* Targeted cell refers to the number of V(D)J expressing cells in the sample. For TCR amplified libraries, targeted cell refers to the number of T cells. For BCR amplified libraries, targeted cell refers to the number of B cells

**Percent Pass Filter (% PF) is reported for NovaSeq and NextSeq 2000 instead of cluster density due to the patterned flow cell

Gene Expression & V(D)J Dual Index Libraries


Four Chromium Single Cell 5' Gene Expression and eight V(D)J (four TCR and four BCR amplified) libraries were pooled and sequenced on indicated Illumina sequencers. 'Data by Cycle' plots from the Illumina SAV software displaying the percentage of base calls and Q30 quality scores along with additional metrics are shown in Table 2.





Sequencing configuration & run parameters:

- Minimum sequencing depth:
Gene Expression 20,000 read pairs/cell; V(D)J amplified 5,000 read pairs/targeted cell*
- Paired-end, dual indexing
- Read 1: 26 cycles
 - i7 Index: 10 cycles
 - i5 Index: 10 cycles
 - Read 2: 90 cycles

Table 2: Representative Plots and Sequencing Data

Plots shown are from a pool of four Gene Expression libraries and eight V(D)J libraries, sequenced on a MiSeq or NextSeq 500. Percentage of bases >Q30 were slightly higher in MiSeq, resulting in slightly higher usable reads for both library types.



		% ≥Q30				Yield per Lane (Gb)		Reads Mapped to Reference (%)	
		R1	i7	i5	R2	R1	R2	Gene Expression	V(D)J
MiSeq									
	Loading Conc. (pM): 10 Cluster density: 1,040 K/mm ² Phix (%): 1	98.2	97.9	97.7	95.6	0.6	2.0	96.7±0.2	80.3±9.2
NextSeq 500									
	Loading Conc. (pM): 1.5 Cluster density: 236 K/mm ² Phix (%): 1	95.3	93.5	93.8	89.1	4.4	15.8	88.0±1.6	77.0±11.6
NovaSeq 6000									
	Loading Conc. (pM): 300 % PF**: 80 Phix (%): 1	93.8	94.5	90.8	90.5	13.9	45.5	86.7±2.0	79.5±10.6
NextSeq 2000									
	Loading Conc. (pM): 650 % PF**: 83.3 Phix (%): 1	92.0	91.5	92.2	91.8	15.0	49.4	88.6±1.4	76.5±13.0

* Targeted cell refers to the number of V(D)J expressing cells in the sample. For TCR amplified libraries, targeted cell refers to the number of T cells. For BCR amplified libraries, targeted cell refers to the number of B cells

**Percent Pass Filter (% PF) is reported for NovaSeq and NextSeq 2000 instead of cluster density due to the patterned flow cell

V(D)J & Cell Surface Protein Dual Index Libraries

Eight Chromium Single Cell V(D)J (four TCR amplified and four BCR amplified), and four Cell Surface Protein libraries were pooled and sequenced on indicated Illumina sequencers. 'Data by Cycle' plots from the Illumina SAV software displaying the percentage of base calls and Q30 quality scores along with additional metrics are shown in Table 3.

Sequencing configuration & run parameters:

Minimum sequencing depth: V(D)J amplified 5,000 read pairs/targeted cell*; Cell Surface Protein 5,000 read pairs/cell



Paired-end, dual indexing






- Read 1: 26 cycles
- i7 Index: 10 cycles
- i5 Index: 10 cycles
- Read 2: 90 cycles

Table 3: Representative Plots and Sequencing Data

Plots shown are from a pool of eight V(D)J libraries and four Cell Surface Protein libraries, sequenced on a HiSeq 2500 in Rapid Run mode

HiSeq 2500



		% ≥Q30				Yield per Lane (Gb)		Reads Mapped to Reference (%)	
		R1	i7	i5	R2	R1	R2	V(D)J	Cell Surface Protein
MiSeq									
	Loading Conc. (pM): 10 Cluster density: 987 K/mm ² Phix (%): 1	98.3	97.7	97.5	95.9	0.6	2.0	80.1±9.2	96.7±0.1
NextSeq 550									
	Loading Conc. (pM): 1.5 Cluster density: 211 K/mm ² Phix (%): 1	95.9	94.2	94.0	91.4	3.9	14.0	75.8±13.1	96.2±0.1
HiSeq 2500 RR									
	Loading Conc. (pM): 10 Cluster density: 969 K/mm ² Phix (%): 1	98.0	96.7	95.7	96.0	4.3	15.2	76.6±12.6	96.7±0.2
NovaSeq 6000									
	Loading Conc. (pM): 300 % PF**: 75 Phix (%): 1	92.0	90.8	85.6	87.9	77.4	255.0	78.7±10.6	96.3±0.2
NextSeq 2000									
	Loading Conc. (pM): 650 % PF**: 83 Phix (%): 1	92.8	91.5	92.1	93.7	15.0	49.3	81.3±8.9	96.1±0.2

*Targeted cell refers to the number of V(D)J expressing cells in the sample. For TCR amplified libraries, targeted cell refers to the number of T cells. For BCR amplified libraries, targeted cell refers to the number of B cells

**Percent Pass Filter (% PF) is reported for NovaSeq and NextSeq 2000 instead of cluster density due to the patterned flow cell

Gene Expression & CRISPR Screening Dual Index Libraries

Twelve Chromium Single Cell Gene Expression, and ten CRISPR Screening libraries were pooled and sequenced on indicated Illumina sequencers. 'Data by Cycle' plots from the Illumina SAV software displaying the percentage of base calls and Q30 quality scores along with additional metrics are shown in Table 4.

Sequencing configuration & run parameters:

Minimum sequencing depth:

Gene Expression 20,000 read pairs/targeted cell;

CRISPR Screening 5,000 read pairs/cell


Paired-end, dual indexing




- Read 1: 26 cycles
- i7 Index: 10 cycles
- i5 Index: 10 cycles
- Read 2: 90 cycles

Table 4: Representative Plots and Sequencing Data

Plots shown are from a pool of twelve Gene Expression and ten CRISPR screening libraries, sequenced on a NextSeq 2000.

NextSeq 2000



		% ≥ Q30				Yield per Lane (Gb)		Reads Mapped to Reference (%)	
		R1	i7	i5	R2	R1	R2	GEX	CRISPR Screening
MiSeq									
	Loading Conc. (pM): 10 Cluster density: 94 K/mm ² Phix (%): 1	98.6	98.5	98.3	95.7	0.4	1.5	98.6±0.2	88.6±6.0
NextSeq 500									
	Loading Conc. (pM): 1.2 Cluster density: 95 K/mm ² Phix (%): 1	98.1	97.0	96.7	95.1	1.7	5.9	95.1±0.4	90.3±6.0
NextSeq 2000									
	Loading Conc. (pM): 650 % PF**: 72 Phix (%): 1	95.8	90.0	92.1	95.7	12.0	42.7	95.0±0.5	93.8±5.8

**Percent Pass Filter (% PF) is reported for NovaSeq and NextSeq 2000 instead of cluster density due to the patterned flow cell

Gene Expression, V(D)J, & BEAM Dual Index Libraries

Four Chromium Single Cell 5' Gene Expression, four V(D)J (two TCR amplified, two BCR amplified) and four BEAM (two BEAM-T and two BEAM-Ab) (dual index) libraries were pooled and sequenced on indicated Illumina sequencers. Table 5 shows 'Data by Cycle' plots from the Illumina SAV software displaying the percentage of base calls and Q30 quality scores with additional metrics.

Sequencing configuration & run parameters:

Minimum sequencing depth: Gene Expression 20,000 read pairs/targeted cell*; V(D)J amplified 5,000 read pairs/cell; BEAM 5,000 read pairs/cell;

* If in-depth analysis of Gene Expression data is not required in a BEAM experiment, the Gene Expression library may be sequenced to 5,000 read pairs/cell, for a pooling ratio of 1:1:1 for Gene Expression:V(D)J:BEAM libraries. 10% PhiX is recommended for 1:1:1 pooling ratios



Paired-end, dual indexing

- Read 1: 26 cycles
- i7 Index: 10 cycles
- i5 Index: 10 cycles
- Read 2: 90 cycles



Table 5: Representative Plots and Sequencing Data

Plots shown are from a pool of four Gene Expression, four V(D)J (two TCR amplified and two BCR amplified), and four BEAM (two BEAM-T and two BEAM-Ab) libraries, sequenced on a NovaSeq SP flow cell. Libraries were pooled at a ratio of 4:1:1 or 1:1:1.




NovaSeq SP 4:1:1



NovaSeq SP 1:1:1, with 10% PhiX





Pooling ratio 4:1:1 for GEX:V(D)J:BEAM libraries

		% ≥Q30				Yield per Lane (Gb)		Reads Mapped to Reference (%)		
		R1	i7	i5	R2	R1	R2	Gene Expression	V(D)J	BEAM
NextSeq 550										
	Loading Conc. (pM): 1.6 Cluster density: 153 K/mm ² Phix (%): 1	97.9	96.9	97.0	94.7	2.4	10.8	78.7±0.02	88.0±0.02	97.9±0.00
NovaSeq (SP flow cell)										
	Loading Conc. (pM) 300 %PF**: 74.8 Phix (%): 1	92.1	91.0	87.8	86.6	12.0	42.6	78.2±0.02	87.8±0.02	98.3±0.00
NextSeq 2000										
	Loading Conc. (pM): 650 %PF**: 83.2 Phix (%): 1	96.1	93.2	94.0	94.1	13.8	49.3	78.2±0.02	87.9±0.02	98.23±0.00

**Percent Pass Filter (% PF) is reported for NovaSeq and NextSeq 2000 instead of cluster density due to the patterned flow cell

Pooling ratio 1:1:1 for GEX:V(D)J:BEAM libraries

		% ≥Q30				Yield per Lane (Gb)		Reads Mapped to Reference (%)		
		R1	i7	i5	R2	R1	R2	Gene Expression	V(D)J	BEAM
NovaSeq (SP flow cell)										
	Loading Conc. (pM) 300 %PF** 68.4 Phix (%): 10***	91.7	88.5	83.3	87.2	10.9	39.0	78.7±0.02	87.8±0.04	98.3±0.00
MiSeq										
	Loading Conc. (pM) 10 Cluster density: 1,200 K/mm² Phix (%): 1	97.3	96.0	95.9	92.7	0.6	2.3	82.3±0.01	87.7± 0.02	96.9±0.00

*** 10% PhiX is recommended for 1:1:1 pooling ratios on NovaSeq

Gene Expression,V(D)J, Cell Surface Protein, & BEAM Dual Index Libraries

Four Chromium Single Cell 5' Gene Expression, four V(D)J (two TCR amplified, two BCR amplified) four Cell Surface Protein, and four BEAM (two BEAM-T and two BEAM-Ab) (dual index) libraries were pooled and sequenced on indicated Illumina sequencers. Table 6 shows 'Data by Cycle' plots from the Illumina SAV software displaying the percentage of base calls and Q30 quality scores along with additional metrics.

Sequencing configuration & run parameters

Minimum sequencing depth: Gene Expression 20,000 read pairs/cell*; V(D)J amplified 5,000 read pairs/cell; Cell Surface Protein 5,000 read pairs/cell, BEAM 5,000 read pairs/cell

* If in-depth analysis of Gene Expression data is not required in a BEAM experiment, the Gene Expression library can be sequenced to 5,000 read pairs/cell, for a pooling ratio of 1:1:1:1 for Gene Expression:V(D)J:BEAM:Cell Surface Protein libraries. 10% PhiX is recommended for 1:1:1:1 pooling ratios.


Paired-end, dual indexing



- Read 1: 26 cycles
- i7 Index: 10 cycles
- i5 Index: 10 cycles
- Read 2: 90 cycles

Table 6: Representative Plots and Sequencing Data

Plots shown are from a pool of four Gene Expression libraries, four V(D)J libraries, four Cell Surface Protein libraries, and four BEAM libraries sequenced on a NovaSeq SP flow cell. Libraries were pooled at a ratio of 4:1:1:1

NovaSeq S2



		% ≥Q30				Yield per Lane (Gb)		Reads Mapped to Reference (%)			
		R1	i7	i5	R2	R1	R2	Gene Expression	V(D)J	Cell Surface Protein	
NextSeq 550											
	Loading Conc. (pM) 1.4 Cluster density: 196 K/mm ² Phix (%): 1	97.1	95.8	96.1	93.4	2.9	10.5	78.5±0.02	87.9±0.04	90.4±0.04	97.8±0.00
NovaSeq (SP flow cell)											
	Loading Conc. (pM): 300 % PF**: 79.5 Phix (%) 1	91.1	88.7	85.5	85.7	11.0	39.0	77.6±0.02	87.3±0.02	89.3±0.05	98.2±0.00
NextSeq 2000											
	Loading Conc. (pM): 650 %PF**: 78.9 Phix (%): 1	94.8	91.8	92.2	93.3	13.1	46.7	78.0±0.02	87.9±0.02	88.0±0.06	98.2±0.00

**Percent Pass Filter (% PF) is reported for NovaSeq and NextSeq 2000 instead of cluster density due to the patterned flow cell

Gene Expression, V(D)J, and Cell Surface Protein Dual Index Libraries from Single Cell 5' v2 standard & HT assays

A pool (1:1) of two Chromium Single Cell 5' v2 Gene Expression, two BCR, and two Cell Surface Protein (standard assay) libraries were compared to a pool (1:1) of two Chromium Single Cell 5' v2 Gene Expression, two BCR, and two Cell Surface Protein (HT assay) libraries. "Data by Cycle" plot from the Illumina SAV software displaying the percentage of base calls and Q30 quality scores are shown in Figure 1.

Sequencing configuration & run parameters:

Minimum sequencing depth:

Gene Expression 20,000 read pairs/targeted cell,

V(D)J amplified 5,000 read pairs/targeted cell*, Cell Surface Protein 5,000 read pairs/cell*

Paired-end, dual indexing


- Read 1: 26 cycles
- i7 Index: 10 cycles
- i5 Index: 10 cycles
- Read 2: 90 cycles

Figure 1: Representative Plots

Plots shown are from a pool (1:1) of two Gene Expression libraries, two BCR, and two Cell Surface protein libraries, sequenced together on a NovaSeq flow cell.

NovaSeq

Data from Single Cell 5' v2 standard assay



Data from Single Cell 5' v2 HT assay




Figure 1. Representative plots derived from Single Cell 5' v2 standard assay (top panel) and HT assay (bottom panel) libraries sequenced on a Novaseq flow cell. Metrics were comparable across all other sequencers tested (data not shown).

* Targeted cell refers to the number of V(D)J expressing cells in the sample. For TCR amplified libraries, targeted cell refers to the number of T cells. For BCR amplified libraries, targeted cell refers to the number of B cells

Note

Depending on the size of the antibody panel, Cell Surface Protein libraries may vary in library complexity. Follow Illumina best practices for working with low complexity libraries, such as optimizing the percentage of Phi-X spike-in.

References

- Chromium Next GEM Single Cell 5' Reagent Kits v2 (Dual Index) User Guide (CG000331)
- Chromium Next GEM Single Cell 5' Reagent Kits v2 (Dual Index) with Feature Barcode technology for Cell Surface Protein & Immune Receptor Mapping User Guide (CG000330)
- Chromium Next GEM Single Cell 5' Reagent Kits v2 (Dual Index) with Feature Barcode technology for CRISPR Screening User Guide (CG000510)
- Chromium Next GEM Single Cell 5' Reagent Kits v2 (Dual Index) with Feature Barcode technology for CRISPR Screening and Cell Surface Protein (CG000511)
- Chromium Next GEM Single Cell 5' HT Reagent Kits v2 (Dual Index) (CG000423)
- Chromium Next GEM Single Cell 5' HT Reagent Kits v2 (Dual Index) with Feature Barcode technology for Cell Surface Protein & Immune Receptor Mapping (CG000424)
- Chromium Next GEM Single Cell 5' HT Reagent Kits v2 (Dual Index) with Feature Barcode technology for CRISPR Screening (CG000512)
- Chromium Next GEM Single Cell 5' HT Reagent Kits v2 (Dual Index) with Feature Barcode technology for CRISPR Screening and Cell Surface Protein (CG000513)
- Chromium Next GEM Single Cell 5' Reagent Kits v2 (Dual Index) with Feature Barcode technology for Barcode Enabled Antigen Mapping (BEAM) (CG000591)
- Chromium Next GEM Single Cell 5' Reagent Kits v2 (Dual Index) with Feature Barcode technology for Barcode Enabled Antigen Mapping (BEAM) and Cell Surface Protein (CG000592)
- Chromium Next GEM Single Cell 5' HT Reagent Kits v2 (Dual Index) with Feature Barcode technology for Barcode Enabled Antigen Mapping (BEAM) (CG000593)
- Chromium Next GEM Single Cell 5' HT Reagent Kits v2 (Dual Index) with Feature Barcode technology for Barcode Enabled Antigen Mapping (BEAM) and Cell Surface Protein (CG000594)

Document Revision Summary

Document Number	CG000401
Title	Sequencing Metrics & Base Composition of Single Cell 5' v2 Dual Index Libraries
Revision	Rev C to Rev D
Revision Date	March 2023
Specific Changes	1. Updated sequencing data for Barcode Enabled Antigen Mapping (BEAM) libraries
General Changes	2. Updated for general minor consistency of language and terms throughout

© 2022 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: www.10xgenomics.com/trademarks. 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: www.10xgenomics.com/patents. The use of products described herein is subject to 10x Genomics Terms and Conditions of Sale, available at www.10xgenomics.com/legal-notices, or such other terms that have been agreed to in writing between 10x Genomics and user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of 10x Genomics products in practicing the methods set forth herein has not been validated by 10x Genomics, and such non-validated use is NOT COVERED BY 10x GENOMICS STANDARD WARRANTY, AND 10x GENOMICS HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics terms and conditions of sale for the Chromium Controller or the Chromium Single Cell Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics that it currently or will at any time in the future offer or in any way support any application set forth herein.

Contact:
support@10xgenomics.com
10x Genomics
6230 Stoneridge Mall Road
Pleasanton, CA 94588 USA

