Visium CytAssist Spatial Applications Imaging Guidelines

Introduction

The Visium CytAssist platform is designed to analyze tissue sections placed on blank slides. Before running a Visium CytAssist assay, a high resolution image is captured with a microscope that is combined with the image taken by the Visium CytAssist instrument to overlay gene expression data. Successful NGS data and protein visualization are highly dependent on good imaging practices. This Technical Note provides hardware recommendations, general image acquisition and analysis guidelines, and examples of images that are suitable for downstream analysis with Space Ranger. Individual results may vary depending on the specific imaging system, and/or sample characteristics.

This Technical Note applies to the following CytAssist assays and sample types:

- Visium CytAssist Spatial Gene Expression
 - Hematoxylin & Eosin (H&E) stained formalin fixed & paraffin embedded (FFPE), fresh frozen & OCT embedded (FF), and fixed frozen & OCT embedded (FxF) tissue sections
 - Immunofluorescence (IF) stained FFPE tissue sections
- Visium CytAssist Spatial Gene and Protein Expression
 - H&E and IF stained FFPE tissue sections

Contents

- 2 General Imaging Guidelines
- 2 1.1 Slide Handling
- 2 1.2 Coverslip Use
- 3 1.3 Imaging System Requirements and Recommendations
- 5 1.4 Image Format Requirements
- 6 Brightfield, Fluorescence, and CytAssist Imaging Guidelines
- 6 2.0 Brightfield Imaging Guidelines
- 6 2.1 Fluorescence Imaging Guidelines
- 7 2.2 CytAssist Imaging Guidelines
- 8 Imaging Examples
- **10** References

Hematoxylin and Eosin (H&E) and Immunofluorescence (IF) Stained FFPE Tissue Sections



Figure 1. H&E stained mouse embryo (left) and IF stained mouse brain (right) tissue.



General Imaging Guidelines

1.0 Imaging Guidelines Overview

This section provides guidelines on slide handling, coverslip use, imaging system requirements, and output image requirements to ensure good image quality and compatibility with the Space Ranger analysis pipeline.

1.1 Slide Handling

- Wear a clean pair of gloves when handling slides.
- DO NOT touch or wipe the side of the slide with tissue.
- Ensure slides are clean. If necessary, use a laboratory wipe or compressed air to clean the bottom (side without tissue) of the slide before imaging.
- Place slides gently and evenly on the imaging stage.
- If using an inverted system, slides can be imaged without inversion when using a 4X or 5X objective. These slides may also be imaged if using an objective with a correction collar or high working distance that accommodates imaging through 1 mm slides.

1.2 Coverslip Use

- Coverslip application is optional for FF samples.
- Imaging through a coverslip will ensure optimal image quality for objectives with a magnification of 10X or higher. Use a coverslip appropriate for the chosen objective. Most objectives are corrected for e.g. a #1.5 (0.17 mm) coverslip. Some specialized objectives are corrected such that a coverslip is not necessary, while objectives with a correction collar (20X or higher) can accommodate other types of coverslips

- If using a coverslip, samples should be loosely mounted with a nonpermanent mounting medium. The coverslip will be removed before the assay workflow.
- It is critical that all excess mounting media be removed before imaging, as mounting media may damage the imaging system. Inverted systems are especially prone to this type of damage. The amount of mounting media may be reduced depending on tissue size.
- If a coverslip is not used for an H&E stained sample, blackening of hematoxylin may occur.
- Coverslip application and removal guidance are provided in documents found in the References section.

1.3 Imaging System Requirements and Recommendations

A good quality image depends on even illumination, an objective and camera combination that provides acceptable resolution and image size, correct imaging settings, and automation (i.e. stitching and shading correction) that provides successful image processing. Example image systems tested by 10x Genomics are provided in Table 3.

Objectives and Cameras

- The level of detail, clarity, and sharpness of an image depends on the objective's magnification, numerical aperture (NA), and the pixel size of the detection camera.
- Air objectives with magnifications of 4X, 5X, or 10X are recommended for image acquisition. DO NOT use immersion media on air objectives. Choose an objective based on imaging needs. See Table 1 for more information. Objectives with higher magnification and NA may be used, but are not required.
- The resolution, i.e. the smallest discernible detail in the image, is primarily dependent on the objective's NA. Objectives with NA 0.1 can distinguish tissue morphology and subcompartments (resolution of $3-4 \mu m$), while single cell resolution is achieved at NA 0.5 or higher (resolution of 0.7 μm).

- The image pixel size should be equal or 2-3 times lower than the smallest discernible detail in the image. 10x Genomics recommends an image pixel size in the range of 0.8 μ m 4 μ m. The pixel size is typically displayed in the imaging software, saved in the image information, or can be calculated by dividing the camera pixel size by the objective's magnification.
- A color camera is required for H&E stained tissue and a monochromatic capable camera/detector for IF stained tissue sections.
- Researchers who wish to perform image-based cell segmentation should use a higher image resolution, e.g. with a 20X or 40X objective.

Filters for Fluorescence Imaging

• Commercially available fluorescently conjugated primary and/or secondary antibodies may require specific filter sets. Ensure filter sets are compatible with imaging system and fluorophore choice.

Objective Magnification / NA	Image Size (pixels) for 8 x 8 mm area	lmage Pixel Size (µm)	Image Resolution
4X / 0.1	2,000–11,500	0.7–4	Tissue morphology and compartments with > 10 cells (50–200 μm)
5X/ 0.16			
4X / 0.2	4,000–11,500	0.7–2	Tissue compartments with 3 – 10 cells (20–100 $\mu m)$
5X / 0.25			
10X / 0.25	5,000-20,000	0.4–1.6	
10X / 0.45 10X / 0.5	10,000–20,000	0.4-0.8	Tissue compartments with 1 – 10 cells (3–10 μm)
20X / 0.8	15,000–45,000	0.2 – 0.5	Subcellular resolution, nuclear features (0.5–2.5 μm)

Table 1. Image resolution obtained with example objectives and image size settings

Image Acquisition and Processing

- H&E stained tissue sections are imaged with a color camera. Slide scanners and some micro-scopes are automated for optimal imaging. For microscopes that require manual setting adjustments, verify the following:
 - When imaging with a 10X or higher objective, adjust the microscope for Köhler illumination: This adjustment of the condenser focus and position provides uniform and strong illumination of the tissue and ensures good contrast and clarity. Refer to the manufacturer directions for performing this adjustment.
 - Adjust the illumination power and exposure time to obtain a bright image without saturation. Poorly exposed images cannot be corrected after acquisition.
 - Adjust the color balance on the camera by performing a "white balance" on an empty area.
 - Stitch the individual tiles with shading correction, if available in the software.

- Fluorescence imaging is performed with a monochromatic camera/detector. For microscopes that require manual setting adjustments, verify the following:
 - Adjust the illumination power, detector gain (if available), and exposure time to obtain a clear and nonsaturated image. Poorly exposed images cannot be corrected after acquisition.
 - 10x Genomics recommends acquiring images with correction for the field of view flatness. Without this correction, the fluorescence intensity is typically lower at the edge of the field of view compared to its center.
 - Stitch the individual tiles without shading correction that could alter the acquired fluorescence levels.

1.4 Image Format Requirements

- A computer system with sufficient resources to handle large images (0.5–10 GB) and the resulting data should be used for image processing. The compute system should be able to stitch images via the microscope's native software or third party software such as FIJI/ImageJ. Preinstall the computer system with the most recent version of Loupe Browser to support manual alignment and tissue selection for IF and H&E stained tissues. Space Ranger is a linux pipeline and will be run on a separate linux server.
- Space Ranger and Loupe Browser accept TIFF/ bigTIFF (file extension .tif or TIFF) or JPEG (file extension .jpg or JPEG) image formats. Space Ranger and Loupe Browser support pyramid TIFF or multiple single TIFF/bigTIFF images that

contain multiple resolutions. Table 2 provides more information. Images in a different format must be converted to TIFF/bigTIFF or JPEG in the native software or other software (i.e. FIJI/ ImageJ). Consult the 10x Genomics Support site for recommendations.

• Images have a minimum size requirement of 2,000 pixels (Visium CytAssist Spatial Gene Expression slide v2, 6.5 mm) and 4,000 pixels (Visium CytAssist Spatial Gene Expression slide v2, 11 mm) in at least one direction to be accepted in Space Ranger. This minimum size does not guarantee images with good resolution. Refer to Section 1.2. (Imaging system requirements and recommendations) for imaging best practices.

Image Type	Image Format		
Brightfield Image	3 x 8-bit color TIFF/bigTIFF or JPEG		
	16-bit grayscale TIFF/bigTIFF or JPEG		
	3 x 8-bit color pyramid TIFF and multiple 8-bit TIFF/bigTIFFs at different resolution (Space Ranger 1.2.2 or newer and Loupe Browser 5.0.1 or newer)		
Fluorescence Image	8 or 16-bit grayscale single, multi-page TIFF/bigTIFF		
	8 or 16-bit grayscale multiple, single-page TIFF/bigTIFF or JPEG		
	3 x 8-bit colored image TIFF/bigTIFF or JPEG		

 Table 2. Acceptable Image Formats for Brightfield and Fluorescence Images

Brightfield, Fluorescence, and CytAssist Imaging Guidelines

2.0 Brightfield Imaging Guidelines

- If performing Hematoxylin and Eosin (H&E) staining and imaging, do so before the Visium CytAssist workflow. Consult documentation listed in the References section.
- If the tissue is over or underexposed, adjust the illumination power and exposure such that the tissue morphology appears clearly. Poorly exposed images cannot be corrected after acquisition.
- Images intended for automatic processing in Space Ranger should be roughly axis-aligned, although rotations of less than 15 degrees are acceptable.
- After image acquisition, stitch image tiles together with the microscope's native software or third party software such as FIJI/ImageJ.
- Export stitched images as a 24-bit color or 8-bit/16-bit monochrome TIFF/BigTIFF or JPEG images. Pyramid TIFF or TIFF files that contain multiple images at different resolutions are compatible with the latest versions of Space Ranger and Loupe Browser.
- For information on image analysis and tissue alignment, refer to the Visium section of the 10x Genomics Software Support Website.

2.1 Fluorescence Imaging Guidelines

- Immunofluorescence staining and imaging should be performed before the Visium CytAssist workflow and is only required if alignment of specific protein staining with gene expression is desired.
- Consult the IF staining Demonstrated Protocols listed in the References section for the full staining protocol, inclusion of autofluorescence quenching, and information on antibody optimization.

Image Export

- Fluorescence images can be provided to Space Ranger as individual grayscale images or as a single, combined, color image. Save images in one of the following ways:
 - A multi-page TIFF/bigTIFF file that contains one or more 8-bit or 16-bit grayscale images. The file can contain up to six total pages/images with the same bit depth.
- Individual channel images as grayscale, 8-bit or 16-bit TIFF/bigTIFF or JPEG files with up to six total files from the same tissue section. Ensure that each monochrome image has been acquired at the same magnification, has the same bit depth, dimension, orientation, and file format.
- A color image with merged channels, where individual channels are colored and combined into one image. This format can be a 24-bit RGB color TIFF/bigTIFF or JPEG file. Saving and exporting images in this format will result in the inability to independently adjust each channel after merging.
- Table 3 shows imaging systems used or tested by 10x Genomics. Any equivalent imaging system can be used as an alternative and is dependent on user needs and applications. For additional information, contact imaging system manufacturer.

2.2 CytAssist Imaging Guidelines

- CytAssist captures a single field of view image of the tissue section, which can be automatically registered to a high resolution H&E or IF image of the same tissue section (if provided) in the Space Ranger analysis pipeline or can be registered using Loupe Browser 6.2 or later.
- If the tissue section is larger than the Visium CytAssist slide Capture Area, the area of interest may be optionally annotated after high resolution imaging using the recommended marker provided in the Demonstrated Protocols.
- If the tissue slide is annotated, annotation marks will appear in the CytAssist image, but they will not be transferred on to the data visualization image. Annotation marks interfere with the auto-image alignment in Space Ranger, but most cases can be resolved by performing manual image alignment in Loupe Browser.

Supplier	Description	Configuration	Image Acquisition Method
Thermo Fisher Scientific	EVOS M7000	Inverted	Semiautomated
Leica	Aperio Versa 8 Leica DMi8	Upright Inverted	Automated Automated
MetaSystems	Metafer	Upright	Automated
Nikon	Nikon Eclipse Ti2	Inverted	Manual
BioTek	Cytation 7	Inverted or Upright	Semiautomated
Keyence	Keyence BZX800	Inverted	Semiautomated
Olympus	VS200	Upright	Automated
Zeiss	Imager.Z2	Upright	Manual

Table 3. Tested imaging systems.

Imaging Examples

The following imaging artifacts may cause poor alignment of the microscope image with NGS data, reduced image resolution, image analysis failure, and/or difficulty in data interpretation. Ensure that optimal imaging settings are verified with prestained tissue before beginning the Visium CytAssist workflow. The examples below are from representative tissue sections. Refer to assay-specific documentation for compatible species and tissue types.



Technical Note | Visium CytAssist Spatial Applications Imaging Guidelines

Setting	Correct	Incorrect
Focus	A	C
White Balancing	R	R
Fluorescence Signal		Undersaturated
Shade Correction		

References

Visium CytAssist Spatial Gene Expression:

- Visium CytAssist Spatial Gene Expression for FFPE - Deparaffinization, H&E Staining, Imaging, & Decrosslinking Demonstrated Protocol (CG000520)
- Visium CytAssist Spatial Gene Expression for Fresh Frozen – Methanol Fixation, H&E Staining, Imaging & Destaining Demonstrated Protocol (CG000614)
- Visium CytAssist Spatial Gene Expression for Fixed Frozen – Rehydration, H&E Staining, Imaging & Decrosslinking Demonstrated Protocol (CG000662)
- Visium CytAssist Spatial Gene Expression for FFPE - Deparaffinization, Decrosslinking, IF Staining & Imaging Demonstrated Protocol (CG000519)

Visium CytAssist Spatial Gene and Protein Expression:

- Visium CytAssist Spatial Gene and Protein Expression for FFPE - Deparaffinization, H&E Staining, Imaging & Decrosslinking Demonstrated Protocol (Document CG000658)
- Visium CytAssist Spatial Gene and Protein Expression for FFPE - Deparaffinization, Decrosslinking, Immunofluorescence Staining & Imaging Demonstrated Protocol (Document CG000659)

Document Revision Summary

Document Number	CG000521		
Title	Visium CytAssist Spatial Applications Imaging Guidelines		
Revision	Rev D		
Revision Date	May 2023		
	Added documentation references for Gene and Protein Expression		
Specific Changes	 Updated document title to Visium CytAssist Spatial Applications Imaging Guidelines 		

© 2023 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: www10xgenomics.com/trademarks. 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: www.10xgenomics.com/patents. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PRCEDURES.

The use of 10x Genomics products in practicing the methods set forth herein has not been validated by 10x Genomics, and such non-validated use is NOT COVERED BY 10X GENOMICS STANDARD WARRANTY, AND 10X GENOMICS HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics terms and conditions of sale for the Chromium Controller or the Chromium Single Cell Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics that it currently or will at any time in the future offer or in any way support any application set forth herein.

Contact:

support@10xgenomics.com 10x Genomics

6230 Stoneridge Mall Road Pleasanton, CA 94588 USA

