


CG000656 Rev A



USER GUIDE

Automated Fixed RNA Profiling Library Construction

FOR USE WITH

Chromium Automated Library Construction Kit-B, 24 rxns PN-1000581

Chromium Automated Library Construction Kit-B, 4 rxns PN-1000597

Dual Index Kit TS Set A, 96 rxns PN-1000251

Dual Index Kit TT Set A, 96 rxns PN-1000215

Notices

Document Number

CG000656 • Rev A

Legal Notices

© 2023 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: www.10xgenomics.com/trademarks. 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: www.10xgenomics.com/patents. The use of products described herein is subject to 10x Genomics Terms and Conditions of Sale, available at www.10xgenomics.com/legal-notices, or such other terms that have been agreed to in writing between 10x Genomics and user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

Instrument & Licensed Software Updates Warranties

Updates to existing Instruments and Licensed Software may be required to enable customers to use new or existing products. In the event of an Instrument failure resulting from an update, such failed Instrument will be replaced or repaired in accordance with the 10x Limited Warranty, Assurance Plan or service agreement, only if such Instrument is covered by any of the foregoing at the time of such failure. Instruments not covered under a current 10x Limited Warranty, Assurance Plan or service agreement will not be replaced or repaired.

Support

Email: support@10xgenomics.com

10x Genomics

6230 Stoneridge Mall Road

Pleasanton, CA 94588 USA

Document Revision Summary

| | |
|------------------------|---|
| Document Number | CG000656 |
| Title | Automated Fixed RNA Profiling Library Construction User Guide |
| Revision | Rev A |
| Revision Date | May 2023 |

Table of Contents

| | |
|--|-----------|
| Introduction | 6 |
| Workflow Overview | 7 |
| Automated Fixed RNA Profiling Library Construction Kit | 8 |
| Additional Kits, Reagents & Equipment | 11 |
| Protocol Steps & Timing | 13 |
| Chromium Connect | 14 |
| Instrument Orientation | 15 |
| Deck Orientation | 16 |
| CSV Setup | 18 |
| Items & Reagents | 20 |
| Tips & Best Practices | 21 |
| Gather Items & Reagents for Fixed RNA Profiling Library Construction | 21 |
| Thaw & Prep Reagents for Fixed RNA Profiling Library Construction | 22 |
| Quick Planner Card | 23 |
| Carrier Loading Guidelines | 24 |
| Carrier Loading Guidelines | 25 |
| Fixed RNA Profiling Library Construction Guidelines | 26 |
| DNA Input | 27 |
| Sample Index PCR | 28 |
| Post Library Construction QC | 29 |
| Post Library Construction Quantification & Pooling | 30 |
| Deck Orientation – Library Quantification | 31 |
| Gather Items & Reagents | 32 |
| Post Library Construction Quantification | 33 |
| Deck Orientation – Library Pooling | 35 |
| Gather Items & Reagents | 36 |
| Library Pooling | 36 |
| Sequencing | 37 |
| Sequencing Libraries | 38 |
| Illumina Sequencer Compatibility | 38 |
| Sample Indices | 38 |
| Fixed RNA – Gene Expression Library Sequencing Parameters | 39 |
| Fixed RNA – Cell Surface Protein Library Sequencing Parameters | 39 |

| | |
|--|--------------------|
| Library Loading | 39 |
| Library Pooling | 40 |
| Data Analysis and Visualization | 40 |
| Appendix | 41 |
| Agilent TapeStation Traces | 42 |
| LabChip Traces | 43 |
| Post Library Construction Cleanup – SPRIselect | 44 |

Introduction

Workflow Overview

Automated Library Construction Kit

Additional Kits, Reagents & Equipment

Recommended Thermal Cyclers


Protocol Steps & Timing

Workflow Overview

This User Guide provides an overview of the Automated Fixed RNA Profiling Library Construction workflow to generate Gene Expression (singleplexed/multiplexed) or Cell Surface Protein (Feature Barcode; singleplexed) libraries. These libraries can be generated on Chromium Connect from manually prepared Amplified Ligated and Extended Product (singleplexed/multiplexed) or Amplified DNA from Cell Surface Protein Feature Barcode (singleplexed) respectively. Consult the Chromium Fixed RNA Profiling Reagent Kits User Guide for Multiplexed Samples (CG000527) or the Chromium Fixed RNA Profiling Reagent Kits User Guide for Singleplexed Samples with Feature Barcode technology for Cell Surface Protein (CG000477) for details of the manual steps.

SPRIselect cleaned Pre-Amplification PCR product (20 µl), prepared as per the manual workflow, is the starting material for the automated Sample Index PCR and final Library Construction steps on Chromium Connect. During Automated Fixed RNA Profiling Library Construction, the Sample Index PCR cycle for a specific library type can be selected through the Chromium Connect User Interface as per recommendations provided in this User Guide (also specified in the relevant manual user guides).

For instrument operation, refer to the Chromium Connect Instrument User Guide (CG000180).



Automated Fixed RNA Profiling Library Construction Kit

Refer to SDS for handling and disposal information

Chromium Automated Library Construction Kit-B, 24 rxns PN-1000581

Reagent volumes and colors are different in each of the module types, not all module tubes contain reagents.

Automated Library Construction Kit-B, Module 1, 24 rxns PN-1000579 (store at 4°C)

Automated Library Construction Kit-B
Module 1, 24 rxns

| | # |
|---|----------------|
| <input checked="" type="radio"/> Module 1 | 24 tube strips |

10xGenomics.com

10x
GENOMICS

Automated Library Construction Kit-B,
Module 1



Automated Library Construction Kit-B, Module 2, 24 rxns PN-1000580 (store at -20°C)

Automated Library Construction Kit-B
Module 2, 24 rxns

| | # |
|--------------------------------|----------------|
| <input type="radio"/> Module 2 | 24 tube strips |

10xGenomics.com

10x
GENOMICS

Automated Library Construction Kit-B,
Module 2



Chromium Automated Library Construction Kit-B, 4 rxns PN-1000597

Reagent volumes and colors are different in each of the module types, not all module tubes contain reagents.

Automated Library Construction Kit-B, Module 1, 4 rxns PN-1000595 (store at 4°C)

Automated Library Construction Kit-B
Module 1, 4 rxns

| | # |
|--|---------------|
| <input checked="" type="checkbox"/> Module 1 | 4 tube strips |

10xGenomics.com **10x**
GENOMICS

Automated Library Construction Kit-B,
Module 1



Automated Library Construction Kit-B, Module 2, 4 rxns PN-1000596 (store at -20°C)

Automated Library Construction Kit-B
Module 2, 4 rxns

| | # |
|-----------------------------------|---------------|
| <input type="checkbox"/> Module 2 | 4 tube strips |

10xGenomics.com **10x**
GENOMICS

Automated Library Construction Kit-B,
Module 2



Dual Index Kit TS Set A, 96 rxns PN-1000251 (store at -20°C)
for Fixed RNA Profiling Gene Expression library construction

Dual Index Kit TS Set A

| | # | PN |
|---------------------------|---|---------|
| Dual Index Plate TS Set A | 1 | 3000511 |

Dual Index Kit TT Set A, 96 rxns PN-1000215 (store at -20°C)
for Fixed RNA Profiling Cell Surface Protein library construction

Dual Index Kit TT Set A

| | # | PN |
|---------------------------|---|---------|
| Dual Index Plate TT Set A | 1 | 3000431 |

Additional Kits, Reagents & Equipment

The items in the table below have been validated by 10x Genomics and are required for the Chromium Connect Automated Fixed RNA Library Construction protocol. DO NOT substitute any of the listed materials.

| Supplier | Description | Part Number (US) |
|---|--|------------------------|
| Plastics | | |
| Hamilton | CO-RE/CO-RE II Tips 50 µl Filtered Tips* | 235948 |
| | CO-RE/CO-RE II Tips 300 µl Filtered Tips* | 235903 |
| | 60 ml Reagent Reservoir Self-Standing | 194051 |
| | Hamilton PCR ComfortLid | 814300 |
| <i>*CO-RE pipette tips will be phased out and replaced by new CO-RE II pipette tips (same part number as CO-RE tips) in 2022/2023. CO-RE II tips include a new sealing surface to interface with the CO-RE II stop disk. Geometry that interfaces with the current CO-RE stop disk is identical between the two tip designs and performance remains unaffected.</i> | | |
| Eppendorf | 96-well Full-Skirted Plate | 951020460 |
| | 96-well Semi-Skirted Plate (Blue color listed; other colors are acceptable) | 951020362 |
| Thermo Fisher Scientific | MicroAmp 8-Tube Strip, 0.2 ml | N8010580 |
| | MicroAmp 8-Cap Strip, clear | 4323032 |
| Kits & Reagents | | |
| Thermo Fisher Scientific | Nuclease-free Water | AM9937 |
| Millipore Sigma | Ethanol, Pure (200 Proof, anhydrous) | E7023-500ML |
| Qiagen | Qiagen Buffer EB | 19086 |
| Equipment | | |
| 10x Genomics | 10x Vortex Adapter | 330002 |
| | Benchtop Vortex | standard lab equipment |
| | Benchtop Centrifuge | standard lab equipment |
| | Plate Centrifuge | standard lab equipment |
| | Benchtop Thermal Cycler (optional) | standard lab equipment |
| Additional materials ONLY for optional assays – qPCR and pooling | | |
| Bio-Rad | 10% Tween 20 | 1662404 |
| | 96-well PCR Plates | HSP9665 |
| Thermo Fisher Scientific | 2 ml-Screw-cap Tubes, NonKnurl | 3488NK |
| | 0.5 ml-Screw-cap Tubes, NonKnurl | 3472NK |
| KAPA Biosystems | KAPA Library Quantification Kit for Illumina Platforms | KK4824 |
| Qiagen | Qiagen Buffer EB | 19086 |
| Additional materials for Chromium Connect maintenance | | |
| Use only indicated cleaning agents. DO NOT use bleach or organic oxidizers | | |
| Thor Labs | Lens tissues | MC-5 |
| VWR | Microcide SQ Broad Spectrum Disinfectant | 25099 |
| Contec | 70% Isopropanol (alternative to VWR disinfectant) | SB167030IR |

Additional Kits, Reagents & Equipment

| Supplier | Description | Part Number (US) |
|--|--|----------------------------|
| Quantification & Quality Control | | |
| Agilent | 2100 Bioanalyzer Laptop Bundle (<i>discontinued</i>) (Replacement 2100 Bioanalyzer Instrument/ 2100 Expert Laptop Bundle) | G2943CA G2939BA/G2953CA |
| Choose Bioanalyzer, TapeStation, LabChip, or Qubit based on availability & preference. | High Sensitivity DNA Kit | 5067-4626 |
| | 4200 TapeStation | G2991AA |
| | High Sensitivity D1000 ScreenTape/Reagents | 5067-5592/ 5067-5593 |
| | High Sensitivity D5000 ScreenTape/Reagents | 5067-5584/ 5067-5585 |
| Thermo Fisher Scientific | Qubit 4.0 Fluorometer Qubit dsDNA HS Assay Kit | Q33238 Q32854 |
| PerkinElmer | LabChip GX Touch HT Nucleic Acid Analyzer DNA High Sensitivity Reagent Kit | CLS137031 CLS760672 |

Protocol Steps & Timing

| Steps | Timing | |
|---|---|---------------------|
| MANUAL | Sample Preparation Prepare single cell or nuclei suspension for Fixed RNA Profiling workflows (multiplex workflow or singleplex workflow with optional Cell Surface Protein Labeling prior to fixation) | Dependent on sample |
| | Sample Fixation | 2-3 h |
| | Probe Hybridization | 16-24 h |
| | GEM Generation & Barcoding | ~4 h |
| | Post-hybridization wash | |
| | Prepare GEM Master Mix + Sample Dilution | |
| | Load Chip Q | |
| | Run Chromium X/iX | |
| | Transfer GEMs | |
| | GEM Incubation | |
| GEM Recovery and Pre-Amplification | 1.5 h | |
| Post-GEM Incubation - Recovery | | |
| Pre-Amplification PCR | | |
| DNA Cleanup - SPRIselect | | |



Use only 20 µl pre-amplified, cleaned-up DNA product (derived from the multiplex workflow or the singleplex workflow with optional Cell Surface Protein Labeling) for Automated Fixed RNA Profiling Library Construction.

| | | | |
|------------------|--|------------------------------|---------|
| AUTOMATED | Fixed RNA Profiling - Library Construction <ul style="list-style-type: none"> • Sample Index PCR • SPRIselect Library Cleanup • Optional qPCR set-up & Pooling | ~2.5 h Walk-away time | |
| | MANUAL | Post Library Construction QC | ~60 min |


Chromium Connect

Instrument Orientation

Deck Orientation



CSV Setup

Instrument Orientation






Refer to the Chromium Connect Instrument User Guide (CG000180) and Quick Reference Cards (CG000256) for more information.

Deck Orientation



Refer to the Chromium Connect Instrument User Guide (CG000180) and Quick Reference Cards (CG000256) for more information.

| Deck Layout Reagents/Consumables | | Automated Single Cell Fixed RNA Profiling Library Construction | |
|----------------------------------|----------------|--|--|
| Carrier | Zone | Item | |
| 1 Stationary | Zone 1 (Black) | 22°C Block, Reagent Strips, Module 1 | |
| | Zone 3 (White) | 10°C Block, Reagent Strips, Module 2 | |
| | - | Magnetic Plate | |
| 2 Stationary | - | Thermal Cycler | |
| | - | Full Skirted PCR Plate (within Thermal Cycler) | |
| | Zone 4 | ComfortLids (2) | |
| 3* Sliding | Position W | Waste Reservoirs | |
| | Position TS | Tube Strip (position 4) | |
| | Position CP | CO-RE Paddles | |
| | Position E | Ethanol Reservoir | |
| | Position Lid | Lid for Ethanol Reservoir | |
| 4 Sliding | Zone 5 | Sample Index Plate TS or TT Set A <i>Verify name & PN</i> | |
| | Zone 6 | Pipette Tips 50 µl | |
| | Zone 7 | Semi Skirted Plate | |
| 5 Sliding | Zone 8 | Pipette Tips 300 µl | |



CSV Setup

Sample information can also be uploaded using a CSV file at the run set-up screen. Use the folder icon to search a network file system or USB drive. Navigate to the appropriate folder with CSV file and click "SELECT".

The Chromium Connect Input Sample Template (CG000309) is available on the 10x Genomics support website.

Run Set-up Screen

The image displays two screenshots of the Chromium Connect Run Set-up Screen. The left screenshot shows the main setup form with the following fields and options:

- Progress bar: FRNA LIBR > **Setup** > Load > Run > Complete
- Experiment Name: Enter Experiment Name (with folder, search, and refresh icons)
- Instruction Level: Standard (dropdown)
- Sample Input Type: GEX Multiplex (dropdown)
- qPCR Setup?: No (dropdown)
- Pooling?: No (dropdown)
- Buttons: CANCEL, NEXT

The right screenshot shows a file selection dialog titled "Select File" with the following elements:

- Progress bar: FRNA LIBR > **Setup** > Load > Run > Complete
- File list: A table with columns "Name" and "Z:".
- File Name: [] CSV files (*.csv) (dropdown)
- Buttons: SELECT, CANCEL
- Buttons: CANCEL, NEXT

Sample CSV File

A sample CSV file is shown below. The columns highlighted in blue are mandatory to start a run. Any missing fields/corrections can be added during sample information setup. Final selections will be recorded in the final run report CSV file.

| SAMPLEPARAMETERS | SAMPLENAME | SIINDEX | VOLUME | CellCount | Expression | CYCLES | USERDEFINED1 | USERDEFINED2 | USERDEFINED3 | USERDEFINED4 |
|--------------------|----------------|---------|--------|-----------|------------|--------|--------------|--------------|--------------|--------------|
| ID1 | aaa | A1 | | | | | | | | |
| ID2 | bbb | B1 | | | | | | | | |
| ID3 | ccc | C1 | | | | | | | | |
| ID4 | ddd | D1 | | | | | | | | |
| ID5 | eee | E1 | | | | | | | | |
| ID6 | fff | F1 | | | | | | | | |
| ID7 | ggg | G1 | | | | | | | | |
| ID8 | hhh | H1 | | | | | | | | |
| RUNPARAMETERS | SELECTION | | | | | | | | | |
| runName | Sample Run | | | | | | | | | |
| Instruction Level | Standard | | | | | | | | | |
| Sample Input Type? | GEX Singleplex | | | | | | | | | |
| qPCR Setup? | No | | | | | | | | | |
| Pooling? | No | | | | | | | | | |
| SI Cycles | 14 | | | | | | | | | |
| Notes | | | | | | | | | | |

Uploading Sample Information Using a CSV File

The following tables provide specific guidelines on sample entry in the CSV file.

| Sample Parameters | Information |
|-------------------|---|
| Sample Name | Alphanumeric and up to 32 characters |
| Sample Index | Location on dual index plate to be used for each sample during SI PCR |

Up to four user-defined fields (LIMS data) can be passed through the instrument for additional sample tracking.

| Run Parameters | Information |
|-------------------|--|
| Run Name | Alphanumeric and up to 32 characters |
| Instruction Level | Standard, Advanced, Expert Refer to the Chromium Connect Instrument User Guide (CG000180) for details |
| Sample Input Type | GEX Multiplex, GEX Singleplex, FB Singleplex |
| qPCR Setup | Opt-in for optional assay step: Yes/No |
| Pooling | Opt-in for optional assay step: Yes/No |
| SI Cycles | User defined field. Refer to appropriate section in this User Guide for guidance on optimal cycles. |



Items & Reagents

Tips & Best Practices

Reagent Handling

- Fully thaw and thoroughly mix reagents before use.
- Ensure that there are no bubbles at the bottoms of modules and plate wells.
- Follow the prompts on the instrument touchscreen for handling modules during setup and use.
- Prepare and dispense 80% ethanol off-deck to avoid spilling on consumables.

Gather Items & Reagents for Fixed RNA Profiling Library Construction

Follow prompts on the Chromium Connect touchscreen to gather the listed items and reagents for loading the Deck Carriers. Gather the quantities specified for each of the items and reagents.

| Item | Qty |
|---|---------------------|
| Nuclease-free Water | 10 ml |
| Ethanol, Pure (200 Proof, anhydrous) | 40 ml |
| Hamilton | |
| ComfortLids | 2 |
| 50 µl CO-RE/CO-RE II Pipette Tips, with filter (Black, Conductive) | 1 rack |
| 300 µl CO-RE/CO-RE II Pipette Tips, with filter (Black, Conductive) | 1 rack |
| Reagent Reservoir, 60 ml | 3 |
| Eppendorf | |
| 96-well Semi Skirted Plate | 1 |
| 96-well Full Skirted Plate | 1 |
| Thermo Fisher Scientific | |
| MicroAmp 8-Tube Strip, 0.2 ml | 1 |
| 10x Genomics | |
| Automated Fixed RNA Profiling Library Construction Kit | |
| Automated Library Construction Kit-B, Module 1 (stored at 4°C) <i>Black tube strip</i> | 1 tube strip/sample |
| Automated Library Construction Kit-B, Module 2 (stored at -20°C) <i>White tube strip</i> | 1 tube strip/sample |
| Dual Index Plate TS or TT Set A (stored at -20°C) <i>Verify name & PN</i> | 1 plate |

See [Additional Kits, Reagents & Equipment list](#) for performing optional assays and/or QC.

Thaw & Prep Reagents for Fixed RNA Profiling Library Construction

Follow prompts on the touchscreen to thaw and prepare reagents. Some important guidelines are highlighted below.

| ACTION | GUIDELINES <i>Follow touchscreen prompts for specifics and timing</i> |
|--|---|
| Thaw Reagents | <ul style="list-style-type: none"> • Thaw reagents as indicated on the touchscreen. Verify no precipitate is present. • Ensure that the correct thawing locations and temperatures are used. • During reagent thaw load the consumables following touchscreen prompts. |
| Prepare Ethanol | <ul style="list-style-type: none"> • Prepare 50 ml 80% Ethanol in Nuclease-free water and dispense in Ethanol Reservoir when prompted. |
| Automated Fixed RNA Profiling Library Construction Modules | <ul style="list-style-type: none"> • Thaw Modules as prompted on the touchscreen. • Confirm that there are no bubbles at the bottoms of any module tubes. |
| Dual Index Plate TS or TT, Set A | <ul style="list-style-type: none"> • Use the indicated plate. Verify name & PN. Vortex Dual Index Plate for 15 sec at maximum speed and centrifuge at 300 rcf for 1 min. |



Confirm that there are no bubbles at the bottoms of any module tubes and Dual Index Plate wells.

Quick Planner Card

Gather the listed items & reagents before running the assay. Follow the touchscreen prompts for detailed information.

Gather indicated items prior to running the assay

| | |
|--|--|
| <input type="checkbox"/> Nuclease free water – 10 ml | <input type="checkbox"/> Semi skirted plate, 96 well – 1 per run |
| <input type="checkbox"/> Ethanol, Pure (200 Proof, anhydrous) – 40 ml <input type="checkbox"/> Combine 40 ml EtOH and 10 ml nuclease free water to prepare 80% EtOH | <input type="checkbox"/> Full skirted plate, 96 well – 1 per run |
| <input type="checkbox"/> ComfortLids – 2 per run | <input type="checkbox"/> 50 µl Black CO-RE/CO-RE II Pipette Tips, with filter • 1 rack |
| <input type="checkbox"/> Micro Amp 8-tube strips, 0.2 ml – 1 per run | |
| <input type="checkbox"/> Reagent reservoirs, 60 ml – 3 per run | <input type="checkbox"/> 300 µl Black CO-RE/CO-RE II Pipette Tips, with filter • 1 rack |

| 10x Reagents | Storage | Preparation & Handling |
|---|---------|---|
| <input type="checkbox"/> Reagent Module 1 (black tube strip) • 1 tube strip per sample | 4°C | Maintain at room temperature for 30 min. Vortex, centrifuge at 300 rcf for 1 min. |
| <input type="checkbox"/> Reagent Module 2 (white tube strip) • 1 tube strip per sample | -20°C | Thaw at 4°C or on ice. Maintain chilled until ready to load. Before loading, invert mix (DO NOT vortex), centrifuge at 300 rcf for 1 min. |
| <input type="checkbox"/> Dual Index Plate TS or TT, Set A • 1 tube strip per sample | -20°C | Vortex for 15 sec at maximum speed, centrifuge at 300 rcf for 1 min. |



Carrier Loading Guidelines


Carrier Loading Guidelines

Follow the instructions on the touchscreen to load the carriers.

Carriers


- Handle the carriers as prompted.
- Ensure that Carriers 3, 4, and 5 are completely slid out and placed on an off-deck workspace before loading.

Align the carriers to the corresponding Deck Rails when sliding them in or out of the deck.



Modules


- Load one tube strip/sample of each of the indicated modules in the corresponding positions on the Carrier, starting from back to front row.
- DO NOT skip any rows when loading.
- Use pinhole alignment to place module tube strips in the correct orientation (as shown on the touchscreen).



Label Tube Strip Orientation

- The tube strip to receive the final libraries should be placed in Position 4 of the Tube Strip Holder.
- Label tube strip orientation for collecting final libraries.

Tube Strip (TS) Holder



Consult the Chromium Connect User Guide (CG000180) for more information.

Fixed RNA Profiling Library Construction Guidelines

Sample Index PCR

Post Sample Index PCR & Size Selection – SPRIselect

Post Library Construction QC

DNA Input

Follow the guidelines in the relevant manual workflow user guides (listed below) to prepare and cleanup pre-amplification DNA product as input for Automated Fixed RNA Profiling Library Construction.

- Chromium Fixed RNA Profiling Reagent Kits for Singleplexed Samples with Feature Barcode technology for Cell Surface Protein (CG000477).
- Chromium Fixed RNA Profiling Reagent Kits for Multiplexed Samples (CG000527).

Transfer DNA Input

- a. Centrifuge **100 μ l** cleaned-up pre-amplification DNA (prepared as per the manual workflow) at **500 rcf for 3 min** at room temperature.

Representative pre-amplification DNA

Example 1 (before centrifuging)



Example 1 (after centrifuging)



Example 2 (before centrifuging)



Example 2 (after centrifuging)



The pellet formed after centrifugation may vary between different samples. Ensure only the supernatant is transferred without disturbing the pellet.

- b. Transfer **20 μ l** supernatant (without disturbing the pellet) to column 1 of a full skirted plate.
After transferring 20 μ l, the remaining volume of pre-amplified DNA may be stored at 4°C for ≤ 72 h or at -20°C for ≤ 4 weeks.



Always use only **20 μ l** cleaned-up pre-amplification DNA product as input for Automated Fixed RNA Profiling Library Construction.

- c. Load onto the Chromium Connect deck, following the prompts on the instrument touchscreen.

Sample Information: Enter the sample information (sample name, singleplex or multiplex information, targeted cell recovery, etc.) in the "Sample ID" field on the instrument touchscreen.

Sample Index PCR

- Enter the applicable SI-PCR cycle number through the Chromium Connect User Interface as per recommendations provided below for specific library types (also specified in the manual workflow user guides).
- Cycle number selected will apply to all the samples in the run.

Chromium Fixed RNA Profiling – Gene Expression Library (multiplexed)

| Targeted Cell Recovery | Total Cycles* | | | |
|------------------------|----------------|--------------------|--|-------------------------------------|
| | for Cell Lines | for PBMCs & Nuclei | for Cells from Fixed & Dissociated Tissues** | for Cells from FFPE Tissue Sections |
| 500-2,000 | 12 | 16 | 15-16 | 17 |
| 2,000-4,000 | 11 | 15 | 14-15 | 16 |
| 4,000-7,000 | 10 | 14 | 13-14 | 15 |
| 7,000-12,000 | 9 | 13 | 12-13 | 14 |
| 12,000-25,000 | 8 | 12 | 11-12 | 13 |
| 25,000-50,000 | 7 | 11 | 10-11 | 12 |
| 50,000-128,000 | 6 | 10 | 9-10 | 11 |

Chromium Fixed RNA Profiling – Gene Expression Library (singleplexed)

| Targeted Cell Recovery | Total Cycles* | | | |
|------------------------|----------------|--------------------|--|-------------------------------------|
| | for Cell Lines | for PBMCs & Nuclei | for Cells from Fixed & Dissociated Tissues** | for Cells from FFPE Tissue Sections |
| 500-2,000 | 12 | 16 | 15-16 | 17 |
| 2,000-4,000 | 11 | 15 | 14-15 | 16 |
| 4,000-7,000 | 10 | 14 | 13-14 | 15 |
| 7,000-10,000 | 9 | 13 | 12-13 | 14 |

;

Chromium Fixed RNA Profiling – Cell Surface Protein Library (singleplexed)

| Targeted Cell Recovery | Total Cycles |
|------------------------|--------------|
| 500-2,000 | 14-15 |
| 2,000-4,000 | 13-14 |
| 4,000-10,000 | 12-13 |


*Optimization of cycle number may be needed based on the total RNA content of the sample. The ideal target library concentration is 50 - 200 nM. However, if the concentration is between 10-50 nM or between 200-500 nM and if the libraries do not contain low or high molecular weight peaks, sequencing can still be performed. For dissociated tumor cells, cycle numbers for cell lines can be used as a starting point. For dissociated primary cells, cycle numbers for PBMCs can be used as a starting point.

**For cells derived from the fixed and dissociated tissue samples, the cycle number will depend on the RNA expression level of the tissue and on overall quality of the tissue prior to fixation. Additional optimization may be required.

Post Library Construction QC


Representative Trace - Gene Expression Library

Run 1 μ l sample at 1:40 dilution on an Agilent Bioanalyzer High Sensitivity chip.



Representative Traces - Cell Surface Protein Library

Run 1 μ l sample at 1:20 dilution on an Agilent Bioanalyzer High Sensitivity chip.



The number of distinct peaks present may vary. An additional peak (~115 bp) may be present due to carryover adapters/primer dimers and can negatively impact sequencing performance. To reduce the additional peak, performing a SPRIselect (0.9X) cleanup before sequencing the library is recommended. See [Appendix](#) for representative traces and [Post Library Construction Cleanup - SPRIselect](#) steps. Note that ~40% of material may be lost when performing the SPRIselect cleanup.

Manually select the analysis region of ~150-300 bp and determine the average fragment size from the Bioanalyzer trace. This will be used as the insert size for library quantification.

Alternate QC Methods

- Agilent TapeStation
- PerkinElmer LabChip

See [Appendix](#) for representative traces

Post Library Construction Quantification & Pooling

Deck Orientation – Library Quantification

Post Library Construction Quantification


Deck Orientation – Library Pooling

Library Pooling

Deck Orientation – Library Quantification

Library quantification using qPCR is recommended for accurate pooling and loading on sequencers. If the option is selected during run-setup, automated qPCR plate set-up can be run directly on Chromium Connect after library generation and final library QC. Alternatively, the option can be selected from the instrument home screen, at the user's convenience. Up to 8 samples can be quantified on a 96 well reaction plate, including duplicates for each sample. The minimum sample volume required is **25 µl**. Only **6 µl** of the sample will be used for qPCR plate setup.

The Chromium Connect deck layout for Library Quantification setup is shown below.



Gather Items & Reagents

Follow prompts on the Chromium Connect touchscreen to gather the listed items and reagents for loading the Deck Carriers for Library Quantification.

Gather the quantities specified for each of the items and reagents.

| Item | Qty |
|---|-------------------|
| Hamilton | |
| ComfortLid | 1 |
| 50 µl CO-RE/CO-RE II Pipette Tips, with filter (Black, Conductive) | 1-2 racks |
| 300 µl CO-RE/CO-RE II Pipette Tips, with filter (Black, Conductive) | 1-2 racks |
| 60-ml Reagent Reservoir | 1 |
| Eppendorf | |
| 96-well Semi Skirted Plate | 1 |
| Thermo Fisher Scientific | |
| 2-ml Tube with Screw Cap | 6 |
| Bio-Rad | |
| 96-well Hard-Shell Full Skirted Plate | 1 |
| Reagent | Qty |
| Qiagen Buffer EB | 50 ml |
| Nuclease-free Water | 1 ml |
| 10% Tween-20 | 250 µl |
| Libraries (in up to four 8-tube strips) | up to 8 libraries |
| KAPA Library Quantification Kit, thawed | |
| SYBR FAST Master Mix | 5 ml |
| Primer Mix | 1 ml |
| Standards | 6 |

Post Library Construction Quantification

- Prepare reagents as prompted on the touchscreen.
- Vortex and centrifuge KAPA standards and libraries before use.
- Retrieve previously prepared Master Mix + Primer Mix
OR
Add **1 ml** Primer Mix to **5 ml** SYBR FAST Master Mix.
- Prepare specified Quantification Master Mix in the 2-ml tube using the guidance below.

| # Sample | Master Mix + Primer Mix (μ l) | Water (μ l) | Total Vol (μ l) |
|----------|------------------------------------|------------------|----------------------|
| 8 | 1305 | 435 | 1740 |
| 7 | 1200 | 400 | 1600 |
| 6 | 1095 | 365 | 1460 |
| 5 | 990 | 330 | 1320 |
| 4 | 885 | 295 | 1180 |
| 3 | 780 | 260 | 1040 |
| 2 | 675 | 225 | 900 |
| 1 | 570 | 190 | 760 |

Volumes listed take into account volume for 6 standards

See Tips & Best Practices for handling instructions. Load the chip in the order listed below. See dispensing instructions and illustration for each row. Always dispense slowly without introducing bubbles. Raising and depressing the pipette plunger should each take ~5 sec. Raise the pipette tips at the same rate as the liquid is rising, keeping the tips slightly submerged.


Post Library Construction Quantification

- Follow the touchscreen prompts for loading, scanning, and executing the run.
- During the run, the following steps will be executed by the instrument:
 - KAPA Master Mix transfer to the 96-well Hard Shell Full Skirted Plate (layout below)
 - Diluent transfer to dilution plate
 - Serial dilutions of libraries
 - Addition of library dilutions, KAPA Standards, and negative controls to the plate

Total reaction volume (20 μ l)=
 16 μ l Master Mix
 +
 4 μ l Library Dilution/
 KAPA Standard/
 Negative Control (NTC)

Dilutions:

1:12,500
 1:62,500
 1:312,500
 1:1,562,500



- After the run is completed, follow the unloading instructions on the touchscreen.
- Cap and store libraries at $4^{\circ}\text{C} \leq 72 \text{ h}$ or $-20^{\circ}\text{C} \leq 4 \text{ weeks}$.
- Remove Full Skirted Plate. Seal plate and centrifuge at **300 rcf** for **1 min** at **22°C**.
- Follow the manufacturer's recommendations for qPCR-based quantification. For library quantification for sequencer clustering, determine the concentration based on average size (bp) derived from the Bioanalyzer/TapeStation trace.

| Step | Temperature | Run Time |
|------|-------------------------------------|----------|
| 1 | 95°C | 00:05:00 |
| 2 | 95°C | 00:00:30 |
| 3 | 60°C | 00:00:45 |
| 4 | Go to Step 2, 29X (Total 30 cycles) | |


Resource:

Use the Chromium Connect Library Quantification Worksheet (CG000157) provided on the 10x Genomics Support website for calculating library concentrations.

Deck Orientation – Library Pooling

The libraries may be pooled on the Chromium Connect instrument and used for sequencing, taking into account the preferred cell numbers and per-cell read depth requirements for each library. Samples utilizing the same sample index should not be pooled together, or run on the same flow cell lane, as this would prevent correct sample demultiplexing.

The Chromium Connect deck layout for Library Pooling is shown below.



Gather Items & Reagents

Follow prompts on the Chromium Connect touchscreen to gather the listed items and reagents for loading the Deck Carriers for Library Pooling.

Gather the quantities specified for each of the items and reagents.

| Item | Qty |
|---|--------------------|
| Hamilton | |
| 50 µl CO-RE/CO-RE II Pipette Tips, with filter (Black, Conductive) | 1-2 racks |
| 300 µl CO-RE/CO-RE II Pipette Tips, with filter (Black, Conductive) | 1-2 racks |
| Reagent Reservoir, 60 ml | 1 |
| Eppendorf | |
| 96-well Semi Skirted Plate | 1 |
| Thermo Fisher Scientific | |
| 0.5-ml Tube with Screw Cap | 6 |
| Reagent | Qty |
| Qiagen Buffer EB | 50 ml |
| Libraries (in up to four 8-tube strips) | up to 32 libraries |

Library Pooling

- Follow the touchscreen prompts for loading, scanning, and executing the run.
- Briefly vortex and centrifuge libraries in the 8-tube strips.
- Ensure a minimum **25 µl** library volume is available in the tubes.
- After run is complete, follow touchscreen prompts to unload and store the libraries.
- Unload remaining items and clean as prompted on the touchscreen.

Resource:

Use the Chromium Connect Library Pooling Worksheet (CG000466) provided on the 10x Genomics Support website to calculate volumes to be pooled. The calculated volumes can be input into the instrument either manually, or via the CSV file generated from this workbook.




Sequencing


Sequencing Libraries

Chromium Fixed RNA Profiling – Gene Expression and Cell Surface Protein libraries comprise standard Illumina paired-end constructs which begin with P5 and end with P7. These libraries include 16 bp 10x GEM Barcodes (10x Barcode) encoded at the start of TruSeq Read 1 (Read 1T). Sample index sequences are incorporated as the i5 and i7 index reads. TruSeq Read 1 (Read 1T) and Small RNA Read 2 (Read 2S) are used in paired-end sequencing of Fixed RNA – Gene Expression libraries. TruSeq Read 1 (Read 1T) and TruSeq Read 2 (Read 2T) are used for paired-end sequencing of Fixed RNA – Cell Surface Protein library. Sequencing these libraries produces a standard Illumina BCL data output folder.

Chromium Fixed RNA Profiling – Gene Expression Library



Chromium Fixed RNA Profiling – Cell Surface Protein Library



Illumina Sequencer Compatibility

The compatibility of the listed sequencers has been verified by 10x Genomics. Some variation in assay performance is expected based on sequencer choice. For more information about performance variation, visit the 10x Genomics Support website.

- iSeq
- MiSeq
- NextSeq 500/550
- NextSeq 1000/2000
- NovaSeq

Sample Indices

Each sample index in the Dual Index Kit TT Set A (PN-1000215) or Dual Index Kit TS Set A (PN-1000251) is a mix of one unique i7 and one unique i5 sample index. If multiple samples are pooled in a flow cell lane, the sample index name (i.e. the Dual Index TS Set A plate well ID, SI-TS-) is needed in the sample sheet used for generating FASTQs with “cellranger mkfastq.” Samples utilizing the same sample index should not be pooled together, or run on the same flow cell lane, as this would not enable correct sample demultiplexing.

Fixed RNA – Gene Expression Library Sequencing Parameters

| | |
|-------------------------|---|
| Sequencing Depth | Minimum 10,000 read pairs per cell |
| Sequencing Type | Paired-end, dual indexing |
| Sequencing Read | Read 1: 28 cycles i7 Index: 10 cycles i5 Index: 10 cycles Read 2: 90 cycles (<i>minimum required Read 2 length is 50 bp</i>) |

Fixed RNA – Cell Surface Protein Library Sequencing Parameters

Pooling Fixed RNA – Gene Expression & Cell Surface Protein libraries is recommended for sequencing to maintain nucleotide diversity.

| | |
|-------------------------|---|
| Sequencing Depth | Minimum 5,000 read pairs per cell |
| Sequencing Type | Paired-end, dual indexing |
| Sequencing Read | Read 1: 28 cycles i7 Index: 10 cycles i5 Index: 10 cycles Read 2: 90 cycles (<i>minimum required Read 2 length is 25 bp</i>) |

Library Loading

Once quantified and normalized, the libraries should be denatured and diluted as recommended for Illumina sequencing platforms. Refer to Illumina documentation for denaturing and diluting libraries. Refer to the 10x Genomics Support website for more information.

Fixed RNA – Gene Expression alone or in combination with Cell Surface Protein libraries

Library Loading:

| Instrument | Loading Concentration (pM) | PhiX (%) |
|------------------|----------------------------|----------|
| MiSeq | 12 | 5 |
| NextSeq 500/550 | 1.6 | 5 |
| NexSeq 1000/2000 | 650 | 5 |
| NovaSeq | 150*/300 | 10 |

* Use 150 pM loading concentration for Illumina XP workflow.

Library Pooling

Fixed RNA – Gene Expression and Cell Surface Protein libraries may be pooled for sequencing, taking into account the differences in cell number and per-cell read depth requirements between each library. Samples utilizing the same sample index should not be pooled together or run on the same flow cell lane, as this would not enable correct sample demultiplexing.

Library Pooling Examples:

| Libraries | Sequencing Depth (read pairs per cell) | Library Pooling Ratio |
|--|---|-----------------------|
| Example | | |
| Fixed RNA – Gene Expression Library | 10,000 | 2 |
| Fixed RNA – Cell Surface Protein Library | 5,000 | 1 |

Data Analysis and Visualization

Sequencing data may be analyzed using Cell Ranger or 10x Genomics Cloud Analysis and visualized using Loupe Browser. Key features for these tools are listed below. For detailed product-specific information, visit the 10x Genomics Support website.

Cell Ranger

Cell Ranger is a set of analysis pipelines that processes Chromium Single Gene Expression data to align reads, generate Feature Barcode matrices and perform clustering and gene expression analysis.

- Input: Base call (BCL) and FASTQ
- Output: BAM, MEX, CSV, HDF5, Web Summary, .cloupe/.loupe
- Operating System: Linux

Cloud Analysis

Cloud Analysis is currently only available for US customers.

Cloud Analysis allows users to run Cell Ranger analysis pipelines from a web browser while computation is handled in the cloud.

- Key features: scalable, highly secure, simple to set up and run
- Input: FASTQ
- Output: BAM, MEX, CSV, HDF5, Web Summary, .cloupe/.loupe

Loupe Browser

Loupe Browser is an interactive data visualization tool that requires no prior programming knowledge.

- Input: .cloupe
- Output: Data visualization, including t-SNE and UMAP projections, custom clusters, differentially expressed genes
- Operating System: MacOS, Windows

Appendix


Agilent TapeStation Traces

LabChip Traces


Agilent TapeStation Traces

Agilent TapeStation High Sensitivity D5000 ScreenTape was used.

Chromium Fixed RNA Profiling – Gene Expression Library



Chromium Fixed RNA Profiling – Cell Surface Protein Library




All traces are representative. Samples were run at 1:40 dilution.


LabChip Traces

DNA High Sensitivity Reagent Kit was used.

Chromium Fixed RNA Profiling – Gene Expression Library



Chromium Fixed RNA Profiling – Cell Surface Protein Library




All traces are representative. Samples were run at 1:20 dilution.

Post Library Construction Cleanup – SPRIselect

| Item | 10x PN | Preparation & Handling | Storage |
|--|--------|---------------------------------|--------------|
| Equilibrate to Room Temperature | | | |
| <input type="checkbox"/> Beckman Coulter SPRIselect Reagent | - | Manufacturer's recommendations. | - |
| Obtain | | | |
| <input type="checkbox"/> Qiagen Buffer EB | - | - | Ambient |
| <input type="checkbox"/> 10x Magnetic Separator | 230003 | - | Ambient |
| <input type="checkbox"/> Prepare 80% Ethanol Prepare 15 ml for 8 reactions | - | Prepare fresh. | Ambient |
| Retrieve | | | |
| <input type="checkbox"/> Chromium Fixed RNA Profiling – Cell Surface Protein Library QC trace shows additional peak (~115 bp) | | | 4°C or -20°C |


- a. Retrieve **35 µl** undiluted library that shows an additional peak (~115 bp) during QC (representative trace below).



- b. Vortex to resuspend SPRIselect Reagent. Add **31 µl** SPRIselect Reagent (**0.9X**) to each sample/library. Pipette mix 15x (pipette set to 55 µl).
- c. Incubate **5 min** at **room temperature**.
- d. Place on the magnet•**High** until the solution clears.
- e. Remove the supernatant.
- f. Add **200 µl** 80% ethanol to the pellet. Wait **30 sec**.
- g. Remove the ethanol.
- h. Repeat steps f and g for a total of 2 washes.

Post Library Construction Cleanup – SPRIselect

- i. Centrifuge briefly. Place on the magnet•**Low**.
- j. Remove any remaining ethanol. Air dry for **2 min**.
- k. Remove from the magnet. Add **35.5 µl** Buffer EB. Pipette mix 15x. If beads still appear clumpy, continue pipette mixing until fully resuspended.
- l. Incubate **2 min** at **room temperature**.
- m. Place on the magnet•**Low** until the solution clears.
- n. Transfer **35 µl** sample to a new tube. Perform **Post Library Construction QC** (representative trace below).



Manually select the analysis region of ~150-300 bp and determine the average fragment size from the Bioanalyzer trace. This will be used as the insert size for library quantification.