

## USER GUIDE

# Chromium Next GEM Automated Single Cell 5' Reagent Kits v2



### FOR USE WITH

Chromium Next GEM Automated Single Cell 5' Kit v2, 24 rxns PN-1000290

Chromium Next GEM Automated Single Cell 5' Kit v2, 4 rxns PN-1000298

Chromium Automated Single Cell Human TCR Amplification & Library Construction Kit, 24 rxns PN-1000300

Chromium Automated Single Cell Mouse TCR Amplification & Library Construction Kit, 24 rxns PN-1000310

Chromium Automated Single Cell Human BCR Amplification & Library Construction Kit, 24 rxns PN-1000305

Chromium Automated Single Cell Mouse BCR Amplification & Library Construction Kit, 24 rxns PN-1000311

Chromium Next GEM Chip K Automated Single Cell Kit, 48 rxns PN-1000289

Chromium Next GEM Chip K Automated Single Cell Kit, 16 rxns PN-1000297

Dual Index Kit TT Set A, 96 rxns PN-1000215

## Notices

### Document Number

CG000384 • Rev D

### Legal Notices

© 2023 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: [www.10xgenomics.com/trademarks](http://www.10xgenomics.com/trademarks). 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: [www.10xgenomics.com/patents](http://www.10xgenomics.com/patents). The use of products described herein is subject to 10x Genomics Terms and Conditions of Sale, available at [www.10xgenomics.com/legal-notices](http://www.10xgenomics.com/legal-notices), or such other terms that have been agreed to in writing between 10x Genomics and user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

### Instrument & Licensed Software Updates Warranties

Updates to existing Instruments and Licensed Software may be required to enable customers to use new or existing products. In the event of an Instrument failure resulting from an update, such failed Instrument will be replaced or repaired in accordance with the 10x Limited Warranty, Assurance Plan or service agreement, only if such Instrument is covered by any of the foregoing at the time of such failure. Instruments not covered under a current 10x Limited Warranty, Assurance Plan or service agreement will not be replaced or repaired.

### Support

Email: [support@10xgenomics.com](mailto:support@10xgenomics.com)

10x Genomics

6230 Stoneridge Mall Road

Pleasanton, CA 94588 USA

## Document Revision Summary

|                        |   |
|------------------------|---|
| <b>Document Number</b> | CG000384  |
| <b>Title</b>           | Chromium Next GEM Automated Single Cell 5' Reagent Kits v2 User Guide |
| <b>Revision</b>        | Rev C to Rev D  |
| <b>Revision Date</b>   | May 2023  |

### Specific Changes:

Updated to include information regarding two different sample input volumes (pages 25, 31-34, 40-42 ). Includes an updated Cell Suspension Volume Calculator table for 10  $\mu$ l input volume and an additional table for 32  $\mu$ l input volume.

Updated Dynabead resuspension instructions (pages 23, 37).

### General Changes:

Updated for general minor consistency of language and terms throughout.

# Table of Contents

|  |           |
|--|-----------|
| Notices  | 2         |
| Document Revision Summary  | 3         |
| <b>Introduction</b>  | <b>6</b>  |
| Chromium Automated Single Cell 5' Workflows                                      | 7         |
| Additional Kits, Reagents & Equipment  | 8         |
| Recommended Thermal Cyclers  | 9         |
| Recommended Real Time qPCR System  | 9         |
| Protocol Steps & Timing  | 10        |
| Stepwise Objectives  | 11        |
| Chromium Next GEM Automated Single Cell 5' Reagent Kits v2                       | 15        |
| Chromium Next GEM Automated 5' Quick Planner Card                                | 23        |
| <b>Tips &amp; Best Practices</b>   | <b>24</b> |
| Consumables  | 25        |
| Cell Concentration   | 25        |
| Cell Preparation   | 25        |
| Reagent Handling   | 26        |
| Chromium Automated Chip Handling   | 26        |
| <b>Chromium Connect</b>  | <b>27</b> |
| Instrument Orientation   | 28        |
| Deck Orientation   | 29        |
| CSV Setup  | 31        |
| Sample Input Files   | 32        |
| Sample Input Template Files  | 33        |
| Uploading Sample Information Using a Template File                               | 34        |
| <b>Items &amp; Reagents_cDNA Amplification &amp; 5' GEX Library Construction</b> | <b>35</b> |
| Gather Items & Reagents  | 36        |
| Thaw & Prep Reagents   | 37        |
| <b>Sample Preparation Guidelines</b>   | <b>39</b> |
| Sample Preparation Guidelines  | 40        |
| 10 $\mu$ l Sample Input – Cell Suspension Volume Calculator Table                | 41        |
| 32 $\mu$ l Sample Input – Cell Suspension Volume Calculator Table                | 42        |

---

|   |                    |
|---|--------------------|
| <a href="#">Carrier Loading Guidelines</a>                                | <a href="#">43</a> |
| Carrier Loading Guidelines  | 44                 |
| <a href="#">Additional Protocol Guidelines</a>                            | <a href="#">46</a> |
| Confirm GEM Generation  | 47                 |
| cDNA Amplification Cycle Number   | 47                 |
| cDNA QC & Quantification  | 48                 |
| <a href="#">5' Gene Expression (GEX) Library Construction Guidelines</a>  | <a href="#">49</a> |
| Sample Index PCR  | 50                 |
| Post Library Construction QC  | 50                 |
| <a href="#">V(D)J Amplification &amp; Library Construction Guidelines</a> | <a href="#">51</a> |
| Deck Orientation – V(D)J Amplification & Library Construction             | 52                 |
| Gather Items & Reagents for V(D)J Amplification and Library Construction  | 54                 |
| Thaw & Prep Reagents for V(D)J Amplification & Library Construction       | 55                 |
| V(D)J Amplification QC & Quantification                                   | 56                 |
| V(D)J Library Construction QC   | 57                 |
| <a href="#">Post Library Construction Quantification &amp; Pooling</a>    | <a href="#">58</a> |
| Deck Orientation – Library Quantification                                 | 59                 |
| Gather Items & Reagents   | 60                 |
| Post Library Construction Quantification                                  | 61                 |
| Deck Orientation – Library Pooling  | 63                 |
| Gather Items & Reagents   | 64                 |
| Library Pooling   | 64                 |
| <a href="#">Sequencing</a>  | <a href="#">65</a> |
| Sequencing Libraries  | 66                 |
| Illumina Sequencer Compatibility  | 66                 |
| Sample Indices  | 66                 |
| Library Sequencing Depth & Run Parameters                                 | 67                 |
| Library Loading   | 67                 |
| Library Pooling   | 67                 |
| <a href="#">Appendix</a>  | <a href="#">68</a> |
| LabChip Traces  | 69                 |
| Oligonucleotide Sequences   | 70                 |

# Introduction

[Chromium Automated Single Cell 5' Workflow](#)

[Additional Kits, Reagents & Equipment](#)

[Recommended Thermal Cyclers](#)

[Recommended Real Time qPCR System](#)


[Protocol Steps & Timing](#)

[Stepwise Objectives](#)


[Next GEM Automated Single Cell 5' Reagent Kits v2](#)

[Next GEM Automated Single Cell 5' Quick Planner Card](#)

## Chromium Automated Single Cell 5' Workflows



\*For Automated Gene Expression flexible workflow with cDNA storage option, refer to Chromium Next GEM Single Cell 5' cDNA Kit v2 User Guide Supplement (CG000473) and Automated Library Construction User Guide (CG000474).



\*\*Unused 5' gene expression library construction reagents will be lost, if the 5' Gene Expression – cDNA Only assay is selected. A more suitable option in this case would be to use the Chromium Next GEM Single Cell 5' cDNA Kit (PN-1000425) for cDNA generation.

## Additional Kits, Reagents & Equipment

The items in the table below have been tested by 10x Genomics and are required for the Chromium Connect Automated Single Cell 5' protocol. DO NOT substitute any of the listed materials.

| Supplier   | Description   | Part Number (US) |
|--|---|------------------|
| <b>Plastics</b>  |   |                  |
| Hamilton   | CO-RE/CO-RE II Tips 50 µl Filtered Tips   | 235948           |
|  | CO-RE/CO-RE II Tips 300 µl Filtered Tips  | 235903           |
|  | 60 ml Reagent Reservoir Self-Standing   | 194051           |
|  | Hamilton PCR ComfortLid   | 814300           |
| Eppendorf  | 96-well Full Skirted Plate  | 951020460        |
|  | 96-well Semi Skirted Plate<br><i>(Blue color listed; other colors are acceptable)</i> | 951020362        |
| Thermo Fisher Scientific   | MicroAmp 8-Tube Strip, 0.2 ml   | N8010580         |
|  | MicroAmp 8-Cap Strip, clear   | 4323032          |
| <b>Kits &amp; Reagents</b>   |   |                  |
| Thermo Fisher Scientific   | Nuclease-free Water   | AM9937           |
| Millipore Sigma  | Ethanol, Pure (200 Proof, anhydrous)  | E7023-500ML      |
| Qiagen   | Qiagen Buffer EB  | 19086            |
| <b>Equipment</b>   |   |                  |
| 10x Genomics   | 10x Vortex Adapter  | 330002           |
| -  | Benchtop Vortex   | -                |
| -  | Benchtop Centrifuge   | -                |
| -  | Plate Centrifuge  | -                |
| -  | Benchtop Thermal Cycler   | -                |
| <b>Additional materials ONLY for optional assays – qPCR and pooling</b>  |   |                  |
| Bio-Rad  | 10% Tween 20  | 1662404          |
|  | 96-well PCR Plates  | HSP9665          |
| Thermo Fisher Scientific   | 2 ml-Screw-cap Tubes, NonKnurl  | 3488NK           |
|  | 0.5 ml-Screw-cap Tubes, NonKnurl  | 3472NK           |
| KAPA Biosystems  | KAPA Library Quantification Kit for Illumina Platforms                                | KK4824           |
| Qiagen   | Qiagen Buffer EB  | 19086            |
| <b>Additional materials for Chromium Connect maintenance</b><br><i>Use only indicated cleaning agents. DO NOT use bleach or organic oxidizers.</i> |   |                  |
| Thor Labs  | Lens tissues  | MC-5             |
| VWR  | Microcide SQ Broad Spectrum Disinfectant  | 25099            |
| Contec   | 70% Isopropanol<br>(alternative to VWR disinfectant))                                 | SB167030IR       |



## Additional Kits, Reagents & Equipment

| Supplier   | Description   | Part Number (US)           |
|--|---|----------------------------|
| Quantification & Quality Control   |   |                            |
| Agilent  | 2100 Bioanalyzer Laptop Bundle (discontinued)<br>(Replacement 2100 Bioanalyzer Instrument/ 2100 Expert Laptop Bundle) | G2943CA<br>G2939BA/G2953CA |
| Choose Bioanalyzer, TapeStation, LabChip, or Qubit based on availability & preference. | High Sensitivity DNA Kit  | 5067-4626                  |
|  | 4200 TapeStation  | G2991AA                    |
|  | High Sensitivity D1000 ScreenTape/Reagents  | 5067-5592/ 5067-5593       |
|  | High Sensitivity D5000 ScreenTape/Reagents  | 5067-5584/ 5067-5585       |
| Thermo Fisher Scientific   | Qubit 4.0 Fluorometer<br>Qubit dsDNA HS Assay Kit   | Q33238<br>Q32854           |
| PerkinElmer  | LabChip GX Touch HT Nucleic Acid Analyzer<br>DNA High Sensitivity Reagent Kit   | CLS137031<br>CLS760672     |

## Recommended Thermal Cyclers

Thermal cyclers for off-deck use.

| Supplier                 | Description  | Part Number   |
|--------------------------|--|---|
| Bio-Rad                  | C1000 Touch Thermal Cycler with 96-Deep Well Reaction Module | 1851197   |
| Eppendorf                | MasterCycler Pro (discontinued)                              | North America 950030010<br>International 6321 000.019 |
| Thermo Fisher Scientific | Veriti 96-Well Thermal Cycler                                | 4375786   |

## Recommended Real Time qPCR System

| Supplier | Description            | Part Number |
|----------|------------------------|-------------|
| Bio-Rad  | CFX96 Real-time System | 1855096     |

The qPCR system should be compatible with Bio-Rad 96-well PCR Plates, P/N HSP9665 and with the KAPA Library Quantification Kit dye. Refer to manufacturer's recommendation.

## Protocol Steps & Timing

|           |           | Steps   | Timing                   |
|-----------|-----------|---|--------------------------|
| 3 h       | MANUAL    | Cell Preparation (Dependent on Cell Type)<br>Gather & Load Reagents and Consumables   | ~1-1.5 h<br>~60 min      |
|           |           | Master Mix Preparation<br>Chromium Automated Controller Loading<br>GEM Generation   |                          |
| 6 h       | AUTOMATED | <b>OPTIONAL</b><br>Confirm GEM Generation (Manual, 5 min) ~45 min after starting  | ~3.5 h<br>Walk-away time |
|           |           | Post GEM RT-Cleanup – Dynabead<br>cDNA Amplification<br>cDNA Cleanup – SPRIselect   |                          |
| 9 h       | MANUAL    | cDNA QC & Quantification  | ~60 min                  |
|           | AUTOMATED | 5' Gene Expression Library Construction<br>Fragmentation, End Repair & A-tailing<br>Post Fragmentation, End Repair & A-tailing Double Sided Size Selection – SPRIselect<br>Adaptor Ligation<br>Post Ligation Cleanup – SPRIselect<br>Sample Index PCR<br>Post Sample Index PCR Double Sided Size Selection – SPRIselect | ~4.5 h<br>Walk-away time |
| 12 h plus | MANUAL    | Post Library Construction QC  | ~60 min                  |
|           | MANUAL    | V(D)J Amplification & Library Construction<br>Gather & Load Reagents and Consumables  | ~45 min                  |
|           | AUTOMATED | V(D)J Amplification 1<br>V(D)J Amplification 1 Double Sided Size Selection – SPRIselect<br>V(D)J Amplification 2<br>V(D)J Amplification 2 Double Sided Size Selection – SPRIselect  | ~3.5 h<br>Walk-away time |
|           | MANUAL    | V(D)J Amplification QC & Quantification   | ~60 min                  |
|           | AUTOMATED | Fragmentation, End Repair & A-tailing<br>Adaptor Ligation<br>Post Ligation Cleanup – SPRIselect<br>Sample Index PCR<br>Post Sample Index PCR Cleanup – SPRIselect   | ~4.5 h<br>Walk-away time |
|           | MANUAL    | Post Library Construction QC  | ~60 min                  |

**OPTIONAL**  
Library Quantification qPCR & Library Pooling

## Stepwise Objectives


Chromium Connect automates the preparation of sequencing-ready, single cell libraries from input samples with walk-away convenience. Generation of Chromium Single Cell 5' Gene Expression and V(D)J libraries on the Chromium Connect instrument includes automated Gel Beads-in-emulsion (GEM) generation, barcoding, and library preparation from single cell suspensions, along with additional functionalities for library quantification and pooling.

The Chromium Connect platform for 5' digital gene expression profiles 500-10,000 individual cells per sample. A pool of ~750,000 10x Barcodes is sampled separately to index each cell's transcriptome. It is done by partitioning thousands of cells into nanoliter-scale GEMs, where all generated cDNA share a common 10x Barcode. Libraries are generated and sequenced from the cDNA and 10x Barcodes are used to associate individual reads back to the individual partitions.

This document outlines the key automated protocol steps for generating Single Cell 5' Gene Expression and V(D)J libraries.


## Single Cell 5' Gel Beads

The Single Cell VDJ 5' Gel Beads primer enables the production of barcoded, full-length cDNA from poly-adenylated mRNA, for generating Single Cell 5' Gene Expression and V(D)J libraries.



## Automated GEM Generation & Barcoding


Automated GEM generation is done by combining barcoded Single Cell VDJ 5' Gel Beads, a Master Mix containing cells and enzymes, and Partitioning Oil onto Chromium Next GEM Automated Chip K. To achieve single cell resolution, cells are delivered at a limiting dilution, such that the majority (~90-99%) of generated GEMs contain no cell, while the remainder largely contain a single cell.



## Automated GEM Generation & Barcoding

Immediately following GEM generation, the Gel Bead is dissolved and any co-partitioned cell is lysed. Oligonucleotides containing (i) an Illumina R1 sequence (read 1 sequencing primer), (ii) a 16 nt 10x Barcode, (iii) a 10 nt unique molecular identifier (UMI), and (iv) 13 nt template switch oligo (TSO) are released and mixed with the cell lysate and a Master Mix containing reverse transcription (RT) reagents and poly(dT) RT primers. Incubation of the GEMs produces 10x Barcoded, full-length cDNA from poly-adenylated mRNA.

### Inside individual GEMs




## Automated Post GEM-RT Cleanup & cDNA Amplification

GEMs are broken and pooled after GEM-RT reaction mixtures are recovered. Silane magnetic beads are used to purify the 10x Barcoded first-strand cDNA from the post GEM-RT reaction mixture, which includes leftover biochemical reagents and primers.

10x Barcoded, full-length cDNA is amplified via PCR with primers against common 5' and 3' ends added during GEM-RT. Amplification generates sufficient material to construct multiple libraries from the same cells, e.g. both T cell and/or B cell libraries and 5' Gene Expression libraries.

### Pooled cDNA amplification




If 5' Gene Expression libraries are not desired, stop the automated protocol after cDNA amplification and proceed directly to V(D)J amplification. Unused 5' Gene Expression library construction reagents will be lost, if the Gene Expression Library Construction is not performed as the next step. A more suitable option in this case would be to use the Chromium Next GEM Single Cell 5' cDNA Kit (PN-1000425) for cDNA generation. Refer to Chromium Next GEM Single Cell 5' cDNA Kit v2 User Guide Supplement (CG000473) and Automated Library Construction User Guide (CG000474) for more details.

## Automated 5' Gene Expression Library Construction

Amplified full-length cDNA from poly-adenylated mRNA is used to generate 5' Gene Expression library. Enzymatic fragmentation and size selection are used to optimize the cDNA amplicon size prior to 5' gene expression library construction. P5, P7, i5 and i7 sample indexes, and Illumina R2 sequence (read 2 primer sequence) are added via End Repair, A-tailing, Adaptor Ligation, and Sample Index PCR. The final libraries contain the P5 and P7 priming sites used in Illumina sequencers.


Pooled amplified cDNA processed in bulk



## Automated V(D)J Amplification from cDNA

Amplified full-length cDNA from poly-adenylated mRNA is used to amplify full-length V(D)J segments (10x Barcoded) via PCR amplification with primers specific to either the TCR or BCR constant regions. If both T and B cells are expected to be present in the partitioned cell population, TCR and BCR transcripts can be amplified in separate reactions from the same amplified cDNA material.

Pooled amplified cDNA processed in bulk



## Automated V(D)J Library Construction

Enzymatic fragmentation and size selection are used to generate variable length fragments that collectively span the V(D)J segments of the amplified TCR or BCR transcripts prior to library construction.

P5, P7, i5 and i7 sample indexes, and an Illumina R2 sequence (read 2 primer sequence) are added via End Repair, A-tailing, Adaptor Ligation, and Sample Index PCR. The final libraries contain the P5 and P7 priming sites used in Illumina sequencing.

## Sequencing

Illumina-ready dual index libraries can be sequenced at the recommended depth & run parameters. Illumina sequencer compatibility, sample indices, library loading and pooling for sequencing are summarized in the Sequencing chapter.

### Chromium Single Cell V(D)J Dual Index Library



### Chromium Single Cell 5' Gene Expression Dual Index Library



[See Appendix for Oligonucleotide Sequences](#)

## Chromium Next GEM Automated Single Cell 5' Reagent Kits v2

### Chromium Next GEM Automated Single Cell 5' Kit v2, 24 rxns PN-1000290

Reagent volumes and colors are different in each of the module types.

#### Chromium Next GEM Automated Single Cell 5' Kit v2, Module 1, 24 rxns PN-1000292 (store at 4°C)

**Chromium**  
Next GEM Automated  
Single Cell 5' v2, Module 1, 24 rxns

|  | #                       |
|--|-------------------------|
| <input checked="" type="radio"/> Module 1      | 24 tube strips          |
| <input type="radio"/> Dynabeads™ MyOne™ SILANE | 6 tubes<br>(PN-2000048) |

10xGenomics.com

10x  
GENOMICS



#### Chromium Next GEM Automated Single Cell 5' Kit v2, Module 2, 24 rxns PN-1000293 (store at -20°C)

**Chromium**  
Next GEM Automated  
Single Cell 5' v2, Module 2, 24 rxns

|   | #              |
|---|----------------|
| <input checked="" type="radio"/> Module 2 | 24 tube strips |

10xGenomics.com

10x  
GENOMICS



#### Chromium Next GEM Automated Single Cell 5' Kit v2, Module 3, 24 rxns PN-1000294 (store at -20°C)

**Chromium**  
Next GEM Automated  
Single Cell 5' v2, Module 3, 24 rxns

|                                | #                       |
|--------------------------------|-------------------------|
| <input type="radio"/> Module 3 | 24 tube strips          |
| Poly-dT RT Primer              | 6 tubes<br>(PN-2000007) |

10xGenomics.com

10x  
GENOMICS




**Chromium Next GEM Automated Single Cell 5' Gel Bead Kit v2,  
24 rxns PN-1000291 (store at -80°C)**

**Chromium**  
Next GEM Automated  
Single Cell 5' Gel Bead Kit v2, 24 rxns

|                              | #             |
|------------------------------|---------------|
| Single Cell VDJ 5' Gel Beads | 3 tube strips |

10xGenomics.com

**10x**  
GENOMICS





## Chromium Next GEM Automated Single Cell 5' Kit v2, 4 rxns PN-1000298

Reagent volumes and colors are different in each of the module types.

### Chromium Next GEM Automated Single Cell 5' Kit v2, Module 1, 4 rxns PN-1000295 (store at 4°C)

**Chromium**  
Next GEM Automated  
Single Cell 5' v2, Module 1, 4 rxns

|  | #                       |
|--|-------------------------|
| <input checked="" type="radio"/> Module 1      | 4 tube strips           |
| <input type="radio"/> Dynabeads™ MyOne™ SILANE | 2 tubes<br>(PN-2000048) |

10xGenomics.com

10x  
GENOMICS



### Chromium Next GEM Automated Single Cell 5' Kit v2, Modules 2 & 3, 4 rxns PN-1000296 (store at -20°C)

**Chromium**  
Next GEM Automated  
Single Cell 5' v2, Modules 2 & 3, 4 rxns

|   | #                       |
|---|-------------------------|
| <input checked="" type="radio"/> Module 2 | 4 tube strips           |
| <input type="radio"/> Module 3            | 4 tube strips           |
| Poly-dT RT Primer                         | 2 tubes<br>(PN-2000007) |

10xGenomics.com

10x  
GENOMICS




### Chromium Next GEM Automated Single Cell 5' Gel Bead Kit v2, 4 rxns PN-1000299 (store at -80°C)

**Chromium**  
Next GEM Automated  
Single Cell 5' Gel Bead Kit v2, 4 rxns

|                              | #                |
|------------------------------|------------------|
| Single Cell VDJ 5' Gel Beads | 4 rxn-tube strip |

10xGenomics.com

10x  
GENOMICS



## Chromium Automated Single Cell Human TCR Amplification & Library Construction Kit, 24 rxns PN-1000300

Reagent volumes and colors are different in each of the module types.

### Chromium Automated Single Cell Human TCR Amplification & Library Construction, V(D)J Module 1, 24 rxns PN-1000301 (store at 4°C)

| Chromium<br>Automated Single Cell Human TCR<br>Amplification & Library Construction,<br>V(D)J Module 1, 24 rxns |                | #              |
|---|----------------|----------------|
| <input checked="" type="checkbox"/>   | V(D)J Module 1 | 24 tube strips |


10xGenomics.com **10x** GENOMICS



### Chromium Automated Single Cell Human TCR Amplification & Library Construction, V(D)J Module 2, 24 rxns PN-1000302 (store at -20°C)

| Chromium<br>Automated Single Cell Human TCR<br>Amplification & Library Construction,<br>V(D)J Module 2, 24 rxns |                              | #                       |
|---|------------------------------|-------------------------|
| <input type="checkbox"/>  | V(D)J Module 2               | 24 tube strips          |
|   | Human T Cell Primer Mix 1 v2 | 6 tubes<br>(PN-2000242) |
|   | Human T Cell Primer Mix 2 v2 | 6 tubes<br>(PN-2000246) |

10xGenomics.com **10x** GENOMICS



## Chromium Automated Single Cell Mouse TCR Amplification &amp; Library Construction Kit, 24 rxns PN-1000310

### Chromium Automated Single Cell Mouse TCR Amplification & Library Construction, V(D)J Module 1, 24 rxns PN-1000303 (store at 4°C)

**Chromium**

Automated Single Cell Mouse TCR  
Amplification & Library Construction,  
V(D)J Module 1, 24 rxns

|  | #              |
|--|----------------|
| <input checked="" type="checkbox"/> V(D)J Module 1 | 24 tube strips |

10xGenomics.com

10x  
GENOMICS

V(D)J Module 1



### Chromium Automated Single Cell Mouse TCR Amplification & Library Construction, V(D)J Module 2, 24 rxns PN-1000304 (store at -20°C)

**Chromium**

Automated Single Cell Mouse TCR  
Amplification & Library Construction,  
V(D)J Module 2, 24 rxns

|   | #                       |
|---|-------------------------|
| <input type="checkbox"/> V(D)J Module 2 | 24 tube strips          |
| Mouse T Cell Primer Mix 1 v2            | 6 tubes<br>(PN-2000256) |
| Mouse T Cell Primer Mix 2 v2            | 6 tubes<br>(PN-2000257) |

10xGenomics.com

10x  
GENOMICS

V(D)J Module 2



## Chromium Automated Single Cell Human BCR Amplification &amp; Library Construction Kit, 24 rxns PN-1000305

### Chromium Automated Single Cell Human BCR Amplification & Library Construction, V(D)J Module 1, 24 rxns PN-1000306 (store at 4°C)

**Chromium**

Automated Single Cell Human BCR Amplification & Library Construction, V(D)J Module 1, 24 rxns

|  | #              |
|--|----------------|
| <input checked="" type="checkbox"/> V(D)J Module 1 | 24 tube strips |

10xGenomics.com

10x  
GENOMICS

V(D)J Module 1



### Chromium Automated Single Cell Human BCR Amplification & Library Construction, V(D)J Module 2, 24 rxns PN-1000307 (store at -20°C)

**Chromium**

Automated Single Cell Human BCR Amplification & Library Construction, V(D)J Module 2, 24 rxns

|   | #                    |
|---|----------------------|
| <input type="checkbox"/> V(D)J Module 2 | 24 tube strips       |
| Human B Cell Primer Mix 1 v2            | 6 tubes (PN-2000254) |
| Human B Cell Primer Mix 2 v2            | 6 tubes (PN-2000255) |

10xGenomics.com

10x  
GENOMICS

V(D)J Module 2



## Chromium Automated Single Cell Mouse BCR Amplification &amp; Library Construction Kit, 24 rxns PN-1000311

### Chromium Automated Single Cell Mouse BCR Amplification & Library Construction Kit, V(D)J Module 1, 24 rxns PN-1000308 (store at 4°C)

**Chromium**

Automated Single Cell Mouse BCR  
Amplification & Library Construction,  
V(D)J Module 1, 24 rxns

|  | #              |
|--|----------------|
| <input checked="" type="checkbox"/> V(D)J Module 1 | 24 tube strips |

10xGenomics.com

10x  
GENOMICS

V(D)J Module 1



### Chromium Automated Single Cell Mouse BCR Amplification & Library Construction Kit, V(D)J Module 2, 24 rxns PN-1000309 (store at -20°C)

**Chromium**

Automated Single Cell Mouse BCR  
Amplification & Library Construction,  
V(D)J Module 2, 24 rxns

|   | #              |
|---|----------------|
| <input type="checkbox"/> V(D)J Module 2 | 24 tube strips |

|                              |                         |
|------------------------------|-------------------------|
| Mouse B Cell Primer Mix 1 v2 | 6 tubes<br>(PN-2000258) |
|------------------------------|-------------------------|

|                              |                         |
|------------------------------|-------------------------|
| Mouse B Cell Primer Mix 2 v2 | 6 tubes<br>(PN-2000259) |
|------------------------------|-------------------------|

10xGenomics.com

10x  
GENOMICS

V(D)J Module 2



### Chromium Next GEM Chip K Automated Single Cell Kit, 48 rxns PN-1000289 (store at ambient temperature)

| Chromium Partitioning Oil                               |   |         |                                    | Chromium 50% Glycerol |         |    |  |
|---|---|---------|------------------------------------|-----------------------|---------|----|--|
|   | # | PN      |                                    |                       | #       | PN |  |
| <input checked="" type="radio"/> Partitioning Oil       | 6 | 2000190 | <input type="radio"/> 50% Glycerol | 6                     | 2000109 |    |  |
| Chromium Next GEM Chip K Automated Single Cell          |   |         |                                    | #                     | PN      |    |  |
| Next GEM Chip K Automated Single Cell (gasket attached) |   |         |                                    | 6                     | 2000371 |    |  |

10xGenomics.com **10x** GENOMICS

Chip (gasket attached)



### Chromium Next GEM Chip K Automated Single Cell Kit, 16 rxns PN-1000297 (store at ambient temperature)

| Chromium Partitioning Oil                               |   |         |                                    | Chromium 50% Glycerol |         |    |  |
|---|---|---------|------------------------------------|-----------------------|---------|----|--|
|   | # | PN      |                                    |                       | #       | PN |  |
| <input checked="" type="radio"/> Partitioning Oil       | 2 | 2000190 | <input type="radio"/> 50% Glycerol | 2                     | 2000109 |    |  |
| Chromium Next GEM Chip K Automated Single Cell          |   |         |                                    | #                     | PN      |    |  |
| Next GEM Chip K Automated Single Cell (gasket attached) |   |         |                                    | 2                     | 2000371 |    |  |

10xGenomics.com **10x** GENOMICS

Chip (gasket attached)



### Dual Index Kit TT Set A, 96 rxns PN-1000215 (store at -20°C)

| Dual Index Kit TT Set A   |   |         |  |
|---------------------------|---|---------|--|
|                           | # | PN      |  |
| Dual Index Plate TT Set A | 1 | 3000431 |  |

## Chromium Next GEM Automated 5' Quick Planner Card

Gather the listed items & reagents before running the assay. Follow the touchscreen prompts for detailed information.

### Gather indicated items prior to running the assay

|  |   |
|--|---|
| <input type="checkbox"/> Set thermal cycler to 37°C and lid to 50°C  | <input type="checkbox"/> Eppendorf 96-well Semi skirted plate, 96 well – 1 per run  |
| <input type="checkbox"/> Nuclease free water – 10 ml   | <input type="checkbox"/> Eppendorf 96-well Full skirted plate, 96 well – 1 per run  |
| <input type="checkbox"/> Ethanol, Pure (200 Proof, anhydrous) – 40 ml<br><input type="checkbox"/> Combine 40 ml EtOH and 10 ml nuclease free water to prepare 80% EtOH | <input type="checkbox"/> 50 µl Black CO-RE/CO-RE II Pipette Tips, with filter <ul style="list-style-type: none"> <li>• 7-8 samples: 2 racks</li> <li>• 4-6 samples: 2 racks</li> <li>• 1-3 samples: 1 rack</li> </ul>   |
| <input type="checkbox"/> ComfortLids – 6 per run   |   |
| <input type="checkbox"/> MicroAmp 8-tube strips, 0.2 ml – 2 per run  | <input type="checkbox"/> 300 µl Black CO-RE/ CO-RE II Pipette Tips, with filter <ul style="list-style-type: none"> <li>• 7-8 samples: 4 racks</li> <li>• 4-6 samples: 3 racks</li> <li>• 1-3 samples: 2 rack</li> </ul> |
| <input type="checkbox"/> Reagent reservoirs, 60 ml – 3 per run   |   |

| 10x Reagents   | Storage                  | Preparation & Handling   |
|--|--------------------------|--|
| <input type="checkbox"/> Next GEM Chip K Automated<br>1 per run  | Room temp.               | Set aside, keep sealed. Follow the touchscreen prompts to load on deck.  |
| <input type="checkbox"/> Partitioning oil,<br>50% Glycerol<br>1 tube each per run  | Room temp.<br>(Chip box) | Keep capped. Follow the touchscreen prompts to remove the cap after cells are loaded on the deck.  |
| <input type="checkbox"/> Library Module 1 (black tube strip)<br>1 tube strip per sample  | 4°C                      | Use a thermal cycler (lid temp 50°C) to thaw for <b>30 min</b> at 37°C. Vortex at <b>15 min</b> and again at <b>30 min</b> , centrifuge at <b>300 rcf</b> for <b>1 min</b> .   |
| <input type="checkbox"/> Dynabeads MyOne Silane - 1 tube per run<br><b>DO NOT save excess</b><br>2 tubes/4rxn kit; 6 tubes/24rxn kit | 4°C<br>(Module 1 Box)    | Equilibrate to room temperature. Immediately before use, vortex (≥30 sec). Aspirate the full liquid volume with a pipette tip to verify that the beads have not settled in the bottom of the tube. Using a 200 µl pipettor (set to 150 µl), pipette mix at least 20X to fully resuspend clumps. If clumps are still present, repeat vortex ≥30 sec; pipette mix 20X or until fully resuspended. <b>DO NOT centrifuge. DO NOT remove cap until prompted on touchscreen.</b> |
| <input type="checkbox"/> Library Module 2 (grey tube strip)<br>1 tube strip per sample   | -20°C                    | Thaw at room temperature for <b>30 min</b> . Vortex, centrifuge at <b>300 rcf</b> for <b>1 min</b> .   |
| <input type="checkbox"/> Library Module 3 (white strip tube)<br>1 tube strip per sample  | -20°C                    | Thaw at 4°C or on ice. Maintain on ice until ready to load. Before loading, invert mix ( <b>DO NOT vortex</b> ), centrifuge at <b>300 rcf</b> for <b>1 min</b> .   |
| <input type="checkbox"/> Poly-dT RT Primer - 1 tube per run<br><b>DO NOT save excess</b><br>2 tubes/4rxn kit; 6 tubes/24rxn kit      | -20°C<br>(Module 3 box)  | Thaw at room temperature for <b>30 min</b> . Vortex, centrifuge briefly.   |
| <input type="checkbox"/> Dual Index TT Set A Plate (SI Plate)<br>1 plate per run, 1 well per sample                                  | -20°C                    | Thaw at room temperature for <b>30 min</b> . Vortex, centrifuge briefly.   |
| <input type="checkbox"/> Gel Beads Strip(s) – 1 well per sample  | -80°C                    | Thaw at room temperature ≥ <b>30 min</b> . Vortex <b>30 sec</b> , centrifuge <b>5 sec</b> .  |

# Tips & Best Practices







Consult the Chromium Connect User Guide (CG000180) and follow the Chromium Connect Touchscreen prompts for specifics of assay execution.

## Consumables

- Use validated and recommended emulsion-safe plastic consumables as some plastics can destabilize GEMs.

## Cell Concentration

- **Resuspend samples in PBS+ 0.04% BSA.** Total volume loaded onto the sample plate can be either 10  $\mu$ l or 32  $\mu$ l (must be the same volumes for all samples in a given run).
- **Based on cell stock concentration, do sequential stock dilutions, if needed.**
- **Use 3 independent cell counts to determine cell concentration.**
- The presence of dead cells in the suspension may also reduce the recovery rate. Consult the 10x Genomics Single Cell Protocols Cell Preparation Guide and the Guidelines for Optimal Sample Preparation flowchart (Documents CG00053 and CG000126 respectively) for more information on preparing cells.
- Refer to the 10x Genomics Support website for more information regarding cell type specific sample preparation.

| Multiplet Rate (%) | # of Cells Loaded | # of Cells Recovered |
|--------------------|-------------------|----------------------|
| ~0.4%              | ~870              | ~500                 |
| ~0.8%              | ~1,700            | ~1,000               |
| ~1.6%              | ~3,500            | ~2,000               |
| ~2.3%              | ~5,300            | ~3,000               |
| ~3.1%              | ~7,000            | ~4,000               |
| ~3.9%              | ~8,700            | ~5,000               |
| ~4.6%              | ~10,500           | ~6,000               |
| ~5.4%              | ~12,200           | ~7,000               |
| ~6.1%              | ~14,000           | ~8,000               |
| ~6.9%              | ~15,700           | ~9,000               |
| ~7.6%              | ~17,400           | ~10,000              |


## Cell Preparation

- Ensure cell counts are accurate.
- Based on cell stock concentration, do sequential stock dilutions, if needed. Based on sample input volume (10  $\mu$ l or 32  $\mu$ l), refer the applicable [Cell Suspension Volume Calculator Table](#) for optimal pipetting volumes and concentrations.
- Load cell samples when prompted on the touchscreen.
- The cDNA amplification cycle number will be based on the targeted cell recovery. The cycle number chosen for one sample will apply to all the samples in a run. Refer to [cDNA Amplification Cycle Number](#) for more information.

## Reagent Handling

- Fully thaw and thoroughly mix reagents before use.
- Resuspend Dynabeads and Poly-dT RT Primers at the end of loading.
- Ensure there are no air bubbles at the bottoms of reagent tubes.
- Follow the prompts on the touchscreen for handling Library Modules 1, 2, and 3 during setup and use.
- Follow the prompts on the touchscreen for handling V(D)J Modules 1 and 2.
- Ensure correct reagent tube barcode orientation (on tubes and racks) as prompted by the touchscreen.
- Prepare and dispense 80% ethanol off-deck to avoid spilling on consumables.
- When indicated, promptly move reagents back to the recommended storage.

Barcode Orientation



## Chromium Automated Chip Handling

- The automated chip includes a pre-installed gasket.
- Minimize exposure of reagents and chips to sources of particles and fibers, laboratory wipes, frequently opened flip-cap tubes, clothing that sheds fibers, and dusty surfaces.
- Keep chip and gasket in sealed package until prompted to load.
- After removing the chip from the sealed bag, use in  $\leq$  24 h.
- Avoid contacting the bottom surface of the chip with gloved hands and other surfaces.
- DO NOT use chips or gaskets specific to other 10x Genomics protocols.

Chip (gasket attached)




# Chromium Connect

Instrument Orientation

Deck Orientation


CSV Setup

## Instrument Orientation




Refer to the Chromium Connect Instrument User Guide (CG000180) and Quick Reference Cards (CG000256) for more information.

## Deck Orientation



T2 - Dynabeads MyOne SILANE  
 T4 - Poly-dT RT Primer  
 T5 - 50% Glycerol  
 T7 - Partitioning Oil




Stationary Carriers

Sliding Carriers

Refer to the Chromium Connect Instrument User Guide (CG000180) and Quick Reference Cards (CG000256) for more information.

| Deck Layout Reagents/Consumables<br>Chromium Next GEM Automated Single Cell 5' Gene Expression v2 Assay |                 |  |
|---|-----------------|--|
| Carrier   | Zone            | Item   |
| <b>1</b><br>Stationary  | Zone 1 (Black)  | 22°C Block, Reagent Strips, Module 1           |
|   | Zone 2 (Gray)   | 22°C Block, Reagent Strips, Module 2           |
|   | Zone 3 (White)  | 10°C Block, Reagent Strips, Module 3           |
|   | -               | Magnetic Plate                                 |
| <b>2</b><br>Stationary  | -               | Thermal Cycler                                 |
|   | -               | Full Skirted PCR Plate (within Thermal Cycler) |
|   | Zone 4          | ComfortLids                                    |
| <b>3*</b><br>Sliding<br>Deck Rails: 15-18<br>Number of Lights: 4  | Position W      | Waste Reservoirs                               |
|   | Position TS     | Tube Strips (TS-1, TS-2, TS-3 & TS-4)          |
|   | Position T2     | Dynabeads™ MyOne™ SILANE                       |
|   | Position T4     | Poly-dT RT Primer                              |
|   | Position T5     | 50% Glycerol                                   |
|   | Position T7     | Partitioning Oil                               |
|   | -               | CO-RE Paddles                                  |
|   | Position E      | Ethanol Reservoir                              |
|   | Position LID    | Lid for Ethanol Reservoir                      |
|   | Position GB-1°  | Gel Beads Primary                              |
|   | Position GB-Alt | Gel Beads Alternate                            |
| <b>4</b><br>Sliding<br>Deck Rails: 19-24<br>Number of Lights: 6   | Zone 5          | Sample Index Plate                             |
|   | Zone 6          | Pipette Tips 50 µl                             |
|   | Zone 7          | Semi Skirted Plate                             |
| <b>5</b><br>Sliding<br>Deck Rails: 25-30<br>Number of Lights: 6   | Zone 8          | Pipette Tips 300 µl                            |



## CSV Setup

Sample information can also be uploaded using a CSV file at the run setup screen. Use the folder icon to search a network file system or USB drive. Navigate to the appropriate CSV file and click "SELECT".

For 5' Gene Expression Library construction, use Chromium Connect Single Cell 5' Gene Expression Input File (CG000430) and for V(D)J Library construction, use Chromium Connect Single Cell 5' V(D)J Input File (CG000432). All the files are available on the 10x Genomics support website.

Alternatively, customer's CSV files can also be generated using the customer's LIMS system. If using a LIMS system to generate CSV files, use ChromiumConnect\_SC5-GEX\_InputSampleInfo\_Template file (CG000429) and for V(D)J Library construction, use ChromiumConnect\_SC5-VDJ\_InputSampleInfo\_Template file (CG000431).

### Run Setup Screen

The image displays two screenshots of the Chromium Connect Run Setup Screen. The left screenshot shows the 'Setup' step of the process, with a folder icon highlighted in a yellow box. The right screenshot shows the 'Select File' dialog box, which is used to upload a CSV file. The dialog box includes a file list with columns for 'Z:' and 'Name', a 'File Name' field, and a file type dropdown set to 'CSV files (\*.csv)'. The 'SELECT' button is highlighted in blue.

## Sample Input Files

Sample input files for Gene Expression and V(D)J Amplification are shown below, refer to Chromium Connect SC5'-GEX Input File (CG000430) for more info. The columns highlighted in blue are mandatory to start a run. Any missing fields/corrections can be added during sample information setup. Final selections will be recorded in the final run report CSV file.

### Gene Expression Sample Input File

#### Chromium Connect Single Cell 5' Gene Expression Input File | 10xgenomics.com

After entering assay information below, click this button to autofill the LIMS file ----

Export Data to csv

| Run Parameters           | Selection  | Notes (included in run logs) | Legend                           |
|--------------------------|------------|------------------------------|----------------------------------|
| ExperimentName           | Sample Run |                              | <b>Blue</b> Enter info manually  |
| Instruction Level        | Standard   |                              | <b>Light Blue</b> Drop-down menu |
| Sample Input Volume      | 10uL       |                              | <b>Red</b> Invalid Entry         |
| GEM Check?               | Yes        |                              |                                  |
| Feature Barcode?         | No         |                              |                                  |
| FB Library Construction? | No         |                              |                                  |
| V(D)J Amplification?     | Yes        |                              |                                  |
| qPCR Setup?              | No         |                              |                                  |
| Pooling?                 | No         |                              |                                  |
| cDNA Cycles              | 13         |                              |                                  |

| Number of Samples | 8         |                       |           |               |                 |
|-------------------|-----------|-----------------------|-----------|---------------|-----------------|
| Sample Number     | Sample ID | Sample Index (A1-H12) | CellCount | Input Type    | ExpressionLevel |
| ID1               | aaa       | A1                    | 2001-6000 | Primary Cells | High            |
| ID2               | bbb       | B1                    | 2001-6000 | Primary Cells | High            |
| ID3               | ccc       | C1                    | 2001-6000 | Primary Cells | High            |
| ID4               | ddd       | D1                    | 2001-6000 | Primary Cells | High            |
| ID5               | eee       | E1                    | 2001-6000 | Primary Cells | High            |
| ID6               | fff       | F1                    | 2001-6000 | Primary Cells | High            |
| ID7               | ggg       | G1                    | 2001-6000 | Primary Cells | High            |
| ID8               | hhh       | H1                    | 2001-6000 | Primary Cells | High            |

### V(D)J Sample Input File

#### Chromium Connect Single Cell 5' V(D)J Input File | 10xgenomics.com

After entering assay information below, click this button to autofill the LIMS file ----

Export Data to csv

| Run Parameters     | Selection  | Notes (included in run logs) | Legend                           |
|--------------------|------------|------------------------------|----------------------------------|
| ExperimentName     | Sample Run |                              | <b>Blue</b> Enter info manually  |
| Instruction Level  | Standard   |                              | <b>Light Blue</b> Drop-down menu |
| Species-Cell Type? | Human-TCR  |                              | <b>Red</b> Invalid Entry         |
| qPCR Setup?        | No         |                              |                                  |
| Pooling?           | No         |                              |                                  |

| Number of Samples | 8         |                       |           |                 |
|-------------------|-----------|-----------------------|-----------|-----------------|
| Sample Number     | Sample ID | Sample Index (A1-H12) | CellCount | ExpressionLevel |
| ID1               | aaa       | A1                    | 2001-6000 | High            |
| ID2               | bbb       | B1                    | 2001-6000 | High            |
| ID3               | ccc       | C1                    | 2001-6000 | High            |
| ID4               | ddd       | D1                    | 2001-6000 | High            |
| ID5               | eee       | E1                    | 2001-6000 | High            |
| ID6               | fff       | F1                    | 2001-6000 | High            |
| ID7               | ggg       | G1                    | 2001-6000 | High            |
| ID8               | hhh       | H1                    | 2001-6000 | High            |

≤ 32 characters/  
symbols/spaces

Rows A-D not  
accessible for  
ID5-ID8



## Sample Input Template Files

Sample input template files for Gene Expression and V(D)J Amplification are shown below, refer to Chromium Connect SC5'-GEX Input Sample Info Template (CG000429) for more info . The columns highlighted in blue are mandatory to start a run. Any missing fields/corrections can be added during sample information setup. Final selections will be recorded in the final run report CSV file.

### Gene Expression Sample Input Template File

| SAMPLE#                  | SAMPLEN   | SIINDEX | VOLUME | CellCount | InputType | Expressio | USERDEFI1 | USERDEFI2 | USERDEFI3 | USERDEFI4 | USERDEFI5 | USERDEFI6 | USERDEFI7 | USERDEFI8 |
|--------------------------|-----------|---------|--------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| ID1                      |           |         |        |           |           |           |           |           |           |           |           |           |           |           |
| ID2                      |           |         |        |           |           |           |           |           |           |           |           |           |           |           |
| ID3                      |           |         |        |           |           |           |           |           |           |           |           |           |           |           |
| ID4                      |           |         |        |           |           |           |           |           |           |           |           |           |           |           |
| ID5                      |           |         |        |           |           |           |           |           |           |           |           |           |           |           |
| ID6                      |           |         |        |           |           |           |           |           |           |           |           |           |           |           |
| ID7                      |           |         |        |           |           |           |           |           |           |           |           |           |           |           |
| ID8                      |           |         |        |           |           |           |           |           |           |           |           |           |           |           |
| RUNPARA                  | SELECTION |         |        |           |           |           |           |           |           |           |           |           |           |           |
| runName                  |           |         |        |           |           |           |           |           |           |           |           |           |           |           |
| Instruction Level        |           |         |        |           |           |           |           |           |           |           |           |           |           |           |
| Sample Input Volume      |           |         |        |           |           |           |           |           |           |           |           |           |           |           |
| GEM Check?               |           |         |        |           |           |           |           |           |           |           |           |           |           |           |
| Feature Barcode?         |           |         |        |           |           |           |           |           |           |           |           |           |           |           |
| FB Library Construction? |           |         |        |           |           |           |           |           |           |           |           |           |           |           |
| V(D)J Amplification?     |           |         |        |           |           |           |           |           |           |           |           |           |           |           |
| qPCR Setup?              |           |         |        |           |           |           |           |           |           |           |           |           |           |           |
| Pooling?                 |           |         |        |           |           |           |           |           |           |           |           |           |           |           |
| cDNA Cycles              |           |         |        |           |           |           |           |           |           |           |           |           |           |           |
| SI Cycles                |           |         |        |           |           |           |           |           |           |           |           |           |           |           |
| Notes                    |           |         |        |           |           |           |           |           |           |           |           |           |           |           |

### V(D)J Sample Input Template File

| SAMPLE#            | SAMPLEN   | SIINDEX | VOLUME | CellCount | ExpressiorCycles | USERDEFI1 | USERDEFI2 | USERDEFI3 | USERDEFI4 | USERDEFI5 | USERDEFI6 | USERDEFI7 | USERDEFI8 |
|--------------------|-----------|---------|--------|-----------|------------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| ID1                |           |         |        |           |                  |           |           |           |           |           |           |           |           |
| ID2                |           |         |        |           |                  |           |           |           |           |           |           |           |           |
| ID3                |           |         |        |           |                  |           |           |           |           |           |           |           |           |
| ID4                |           |         |        |           |                  |           |           |           |           |           |           |           |           |
| ID5                |           |         |        |           |                  |           |           |           |           |           |           |           |           |
| ID6                |           |         |        |           |                  |           |           |           |           |           |           |           |           |
| ID7                |           |         |        |           |                  |           |           |           |           |           |           |           |           |
| ID8                |           |         |        |           |                  |           |           |           |           |           |           |           |           |
| RUNPARA            | SELECTION |         |        |           |                  |           |           |           |           |           |           |           |           |
| runName            |           |         |        |           |                  |           |           |           |           |           |           |           |           |
| Instruction Level  |           |         |        |           |                  |           |           |           |           |           |           |           |           |
| Species-Cell Type? |           |         |        |           |                  |           |           |           |           |           |           |           |           |
| qPCR Setup?        |           |         |        |           |                  |           |           |           |           |           |           |           |           |
| Pooling?           |           |         |        |           |                  |           |           |           |           |           |           |           |           |
| Notes              |           |         |        |           |                  |           |           |           |           |           |           |           |           |

## Uploading Sample Information Using a Template File

The following tables provide specific guidelines on sample entry in the template file.

| Sample Parameters | Information   |
|-------------------|---|
| Sample Name       | Alphanumeric and up to 32 characters  |
| SI Index          | Location on sample index plate to be used for each sample during SI PCR   |
| Expression Level  | User defined field for tracking<br>Example:<br>High cell expression: Cell lines<br>Low cell expression: PBMCs   |
| Cell Count        | User defined field for tracking<br>(enter applicable option EXACTLY as shown below)<br>500-2000<br>2001-6000<br>6001-10000<br>DO NOT use commas.<br>Space between symbol & number required. |

Up to four user-defined fields (LIMS data) can be passed through the instrument for additional sample tracking.

| Run Parameters                       | Information  |
|--------------------------------------|--|
| Run Name                             | Alphanumeric and up to 32 characters   |
| Instruction Level                    | Standard, Advanced, Expert<br>Refer to the Chromium Connect Instrument User Guide (CG000180) for details     |
| Sample Input Volume                  | 10 µl or 32 µl   |
| Run Steps                            | GEX/cDNA only  |
| GEM Check                            | Opt-in for optional QC step: Yes/No  |
| Feature Barcode                      | Opt-in for optional assay step: Yes/No   |
| Feature Barcode Library Construction | Opt-in for optional assay step: Yes/No   |
| V(D)J Amplification                  | Opt-in for optional assay step: Yes/No   |
| qPCR Setup                           | Opt-in for optional assay step: Yes/No   |
| Pooling                              | Opt-in for optional assay step: Yes/No   |
| Species-Cell Type                    | Human-TCR/Human-BCR<br>Mouse-TCR/Mouse-BCR   |
| cDNA Cycles                          | User defined field. Refer to <a href="#">cDNA Amplification Cycle Number</a> for guidance on optimal cycles. |
| SI Cycles                            | User defined field. Refer to appropriate section in this User Guide for guidance.                            |

# Items & Reagents for cDNA Amplification & 5' GEX Library Construction

## Gather Items & Reagents

Follow prompts on the Chromium Connect touchscreen to gather the listed items and reagents for loading the Deck Carriers.

Gather the quantities specified for each of the items and reagents.

| Item   | Qty                               |
|--|-----------------------------------|
| Nuclease-free Water  | 10 ml                             |
| Ethanol, Pure (200 Proof, anhydrous)   | 40 ml                             |
| <b>Hamilton</b>  |                                   |
| ComfortLids  | 6                                 |
| 50 µl CO-RE Pipette Tips, with filter (Black, Conductive)  | 2 racks                           |
| 300 µl CO-RE Pipette Tips, with filter (Black, Conductive)   | 4 racks                           |
| Reagent Reservoir, 60 ml   | 3                                 |
| <b>Eppendorf</b>   |                                   |
| 96-well Semi Skirted Plate   | 1                                 |
| 96-well Full Skirted Plate   | 1                                 |
| <b>Thermo Fisher Scientific</b>  |                                   |
| MicroAmp 8-Tube Strip, 0.2 ml  | 2                                 |
| <b>10x Genomics</b>  |                                   |
| Chromium Next GEM Chip K Automated Single Cell Kit<br>(stored at room temperature)<br><i>Partitioning Oil</i><br><i>50% Glycerol</i><br><i>Chip K (keep chip sealed)</i> | 1                                 |
| Chromium Next GEM Automated Single Cell 5' Gel Bead Kit v2<br>(stored at -80°C)<br><i>Single Cell VDJ 5' Gel Bead v2</i>   | 1 tube/sample                     |
| Chromium Next GEM Automated Single Cell 5' Kit v2  |                                   |
| Module 1 (stored at 4°C)<br><i>Black tube strip</i><br><i>Dynabeads</i>  | 1 tube strip/sample<br>1 tube/run |
| Module 2 (stored at -20°C)<br><i>Gray tube strip</i>   | 1 tube strip/sample               |
| Module 3 (stored at -20°C)<br><i>White tube strip</i><br><i>Poly-dT RT Primer</i>  | 1 tube strip/sample<br>1 tube/run |
| Dual Index Plate TT Set A (stored at -20°C)  | 1 plate                           |

See [Additional Kits, Reagents & Equipment](#) list for performing optional assays and/or QC.

## Thaw & Prep Reagents

Follow prompts on the touchscreen to thaw and prepare reagents. Some important guidelines are highlighted below.

| ACTION            | GUIDELINES<br><i>Follow touchscreen prompts for specifics and timing</i>   |
|-------------------|--|
| Thaw Reagents     | Thaw reagents as indicated on the touchscreen. Verify no precipitate is present. Ensure that the correct thawing locations and temperatures are used. During reagent thaw load the consumables following touchscreen prompts.  |
| Prepare Ethanol   | Prepare <b>50 ml</b> 80% Ethanol in Nuclease-free water and dispense in Ethanol Reservoir when prompted.   |
| Poly-dT RT Primer | Vortex only when prompted on the touchscreen. Centrifuge briefly before loading.   |
| Dynabeads         | Equilibrate to room temperature. Immediately before use: Vortex Dynabeads for <b>≥30 sec</b> . Aspirate the full liquid volume with a pipette tip to verify that the beads have not settled in the bottom of the tube.<br>Using a <b>200 µl pipettor (set to 150 µl)</b> , pipette mix at least 20X to fully resuspend clumps. If clumps are still present, repeat vortex ≥30 sec; pipette mix 20X or until fully resuspended. <b>DO NOT</b> centrifuge. <b>DO NOT</b> remove cap until prompted on touchscreen. Confirm there are no bubbles at the bottom of the tube. |
| Library Modules   | Thaw Library Modules as prompted on the touchscreen. After reagent thaw, invert rack holding Module tube strips and vortex Library Modules 1 and 2 for <b>30 sec</b> ; verify no precipitate. Confirm there are no bubbles at the bottoms of any module tubes. Centrifuge Library Modules 1 and 2 at <b>300 rcf</b> for 1 min at 22°C. Retrieve Library Module 3 from 4°C thaw. <b>DO NOT</b> vortex. Invert-mix and centrifuge at <b>300 rcf</b> for 1 min at 22°C.   |



Resuspend Clump




Confirm that there are no bubbles at the bottoms of any module tubes, Dual Index Plate wells, or Primer tubes.

## Thaw & Prep Reagents

Follow prompts on the touchscreen to thaw and prepare reagents. Some important guidelines are highlighted below.

| ACTION                          | GUIDELINES<br><i>Follow touchscreen prompts for specifics and timing</i>   |
|---------------------------------|--|
| <p><b>Prepare Gel Beads</b></p> | <ul style="list-style-type: none"> <li>• Equilibrate the Gel Beads for <b>30 min at room temperature</b> before use.</li> <li>• Snap the tube strip holder with the Gel Bead strip into a 10x Vortex Adapter. Vortex <b>30 sec</b>.</li> <li>• Centrifuge the Gel Bead strip for <b>~5 sec</b> after removing from the holder. Confirm there are no bubbles at the bottom of the tubes and the liquid levels look even.</li> <li>• Place the Gel Bead strip back in the holder and secure the holder lid.</li> <li>• Store unused Gel Beads at <b>-80°C</b> and avoid more than 12 freeze-thaw cycles. <b>DO NOT</b> leave Gel Beads at room temperature for <b>&gt;24 h</b>.</li> <li>• Remove Gel Beads from the Deck during any of the QCs and store the holder with the unused Gel Beads at <b>-80°C</b>.</li> </ul> |




# Sample Preparation Guidelines


## Sample Preparation Guidelines



- Resuspend samples in PBS + 0.04% BSA. Total volume loaded per sample onto the sample plate is either 10  $\mu$ l or 32  $\mu$ l (must be the same volumes for all samples in a given run).
- Based on the sample input volume (10  $\mu$ l or 32  $\mu$ l), refer to the applicable [Cell Suspension Volume Calculator Table](#) for the cell suspension and buffer volumes.
- Based on cell stock concentration, do sequential stock dilutions, if needed.
- It is recommended to use 3 independent cell counts to determine cell concentration.
- The presence of dead cells in the suspension may also reduce the recovery rate. Consult the 10x Genomics Single Cell Protocols Cell Preparation Guide and the Guidelines for Optimal Sample Preparation flowchart (Documents CG00053 and CG000126, respectively) for more information on preparing cells.
- The cell load impacts PCR cycle numbers for cDNA amplification and other downstream steps in the assay. Refer to [Additional Protocol Guidelines](#) chapter for more information.
- Differences in manual and automated sample preparation are outlined below:

|                     | Manual   | Automated   |
|---------------------|--|---|
| Sample Prep         | Using 10x Genomics Demonstrated Protocols for cell prep and QC |   |
| Sample Input Volume | Up to 38.7 $\mu$ l   | 10 $\mu$ l or 32 $\mu$ l<br> Refer to the applicable Cell Suspension Volume Calculator Table |
| Sample Loading      | PCR strip tubes  | 96-well skirted plate   |
| Samples per Chip    | 1-8  | 1-8   |
| Samples Tested      | Various  | Human PBMCs, mouse PBMCs, mouse splenocytes, human melanoma   |



 Samples are loaded in column 1, starting at A1. It is not necessary to add glycerol to unused sample wells when running <8 samples.



## 10 $\mu$ l Sample Input – Cell Suspension Volume Calculator Table

(Chromium Connect Automated Single Cell 5' v2 protocol)

Volume of Cell Suspension Stock per reaction ( $\mu$ l) | Volume of PBS + 0.04% BSA ( $\mu$ l)

| Cell Stock Conc.<br>(cells/ $\mu$ l) | Targeted Cell Recovery |      |      |      |      |      |      |      |      |      |       |     |
|--------------------------------------|------------------------|------|------|------|------|------|------|------|------|------|-------|-----|
|                                      | 500                    | 1000 | 2000 | 3000 | 4000 | 5000 | 6000 | 7000 | 8000 | 9000 | 10000 |     |
| 100                                  | 8.3                    | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a   | n/a |
|                                      | 1.8                    | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a   | n/a |
| 200                                  | 4.1                    | 8.3  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a   | n/a |
|                                      | 5.9                    | 1.8  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a   | n/a |
| 300                                  | 2.8                    | 5.5  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a   | n/a |
|                                      | 7.3                    | 4.5  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a   | n/a |
| 400                                  | 2.1                    | 4.1  | 8.3  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a   | n/a |
|                                      | 7.9                    | 5.9  | 1.8  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a   | n/a |
| 500                                  | 1.7                    | 3.3  | 6.6  | 9.9  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a   | n/a |
|                                      | 8.4                    | 6.7  | 3.4  | 0.1  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a   | n/a |
| 600                                  | 1.4                    | 2.8  | 5.5  | 8.3  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a   | n/a |
|                                      | 8.6                    | 7.3  | 4.5  | 1.8  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a   | n/a |
| 700                                  | 1.2                    | 2.4  | 4.7  | 7.1  | 9.4  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a   | n/a |
|                                      | 8.8                    | 7.6  | 5.3  | 2.9  | 0.6  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a   | n/a |
| 800                                  | 1.0                    | 2.1  | 4.1  | 6.2  | 8.3  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a   | n/a |
|                                      | 9.0                    | 7.9  | 5.9  | 3.8  | 1.8  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a   | n/a |
| 900                                  | 0.9                    | 1.8  | 3.7  | 5.5  | 7.3  | 9.2  | n/a  | n/a  | n/a  | n/a  | n/a   | n/a |
|                                      | 9.1                    | 8.2  | 6.3  | 4.5  | 2.7  | 0.8  | n/a  | n/a  | n/a  | n/a  | n/a   | n/a |
| 1000                                 | 0.8                    | 1.7  | 3.3  | 5.0  | 6.6  | 8.3  | 9.9  | n/a  | n/a  | n/a  | n/a   | n/a |
|                                      | 9.2                    | 8.4  | 6.7  | 5.1  | 3.4  | 1.8  | 0.1  | n/a  | n/a  | n/a  | n/a   | n/a |
| 1100                                 | 0.8                    | 1.5  | 3.0  | 4.5  | 6.0  | 7.5  | 9.0  | n/a  | n/a  | n/a  | n/a   | n/a |
|                                      | 9.3                    | 8.5  | 7.0  | 5.5  | 4.0  | 2.5  | 1.0  | n/a  | n/a  | n/a  | n/a   | n/a |
| 1200                                 | 0.7                    | 1.4  | 2.8  | 4.1  | 5.5  | 6.9  | 8.3  | 9.6  | n/a  | n/a  | n/a   | n/a |
|                                      | 9.3                    | 8.6  | 7.3  | 5.9  | 4.5  | 3.1  | 1.8  | 0.4  | n/a  | n/a  | n/a   | n/a |
| 1300                                 | 0.6                    | 1.3  | 2.5  | 3.8  | 5.1  | 6.3  | 7.6  | 8.9  | n/a  | n/a  | n/a   | n/a |
|                                      | 9.4                    | 8.7  | 7.5  | 6.2  | 4.9  | 3.7  | 2.4  | 1.1  | n/a  | n/a  | n/a   | n/a |
| 1400                                 | 0.6                    | 1.2  | 2.4  | 3.5  | 4.7  | 5.9  | 7.1  | 8.3  | 9.4  | n/a  | n/a   | n/a |
|                                      | 9.4                    | 8.8  | 7.6  | 6.5  | 5.3  | 4.1  | 2.9  | 1.8  | 0.6  | n/a  | n/a   | n/a |
| 1500                                 | 0.6                    | 1.1  | 2.2  | 3.3  | 4.4  | 5.5  | 6.6  | 7.7  | 8.8  | 9.9  | n/a   | n/a |
|                                      | 9.5                    | 8.9  | 7.8  | 6.7  | 5.6  | 4.5  | 3.4  | 2.3  | 1.2  | 0.1  | n/a   | n/a |
| 1600                                 | 0.5                    | 1.0  | 2.1  | 3.1  | 4.1  | 5.2  | 6.2  | 7.2  | 8.3  | 9.3  | n/a   | n/a |
|                                      | 9.5                    | 9.0  | 7.9  | 6.9  | 5.9  | 4.8  | 3.8  | 2.8  | 1.8  | 0.7  | n/a   | n/a |
| 1700                                 | 0.5                    | 1.0  | 1.9  | 2.9  | 3.9  | 4.9  | 5.8  | 6.8  | 7.8  | 8.7  | 9.7   | n/a |
|                                      | 9.5                    | 9.0  | 8.1  | 7.1  | 6.1  | 5.1  | 4.2  | 3.2  | 2.2  | 1.3  | 0.3   | n/a |
| 1800                                 | 0.5                    | 0.9  | 1.8  | 2.8  | 3.7  | 4.6  | 5.5  | 6.4  | 7.3  | 8.3  | 9.2   | n/a |
|                                      | 9.5                    | 9.1  | 8.2  | 7.3  | 6.3  | 5.4  | 4.5  | 3.6  | 2.7  | 1.8  | 0.8   | n/a |
| 1900                                 | 0.4                    | 0.9  | 1.7  | 2.6  | 3.5  | 4.3  | 5.2  | 6.1  | 6.9  | 7.8  | 8.7   | n/a |
|                                      | 9.6                    | 9.1  | 8.3  | 7.4  | 6.5  | 5.7  | 4.8  | 3.9  | 3.1  | 2.2  | 1.3   | n/a |
| 2000                                 | 0.4                    | 0.8  | 1.7  | 2.5  | 3.3  | 4.1  | 5.0  | 5.8  | 6.6  | 7.4  | 8.3   | n/a |
|                                      | 9.6                    | 9.2  | 8.4  | 7.5  | 6.7  | 5.9  | 5.1  | 4.2  | 3.4  | 2.6  | 1.8   | n/a |

Grey boxes: Volumes that would exceed the allowable buffer volume in each reaction

Yellow boxes: Indicate a low transfer volume that may result in higher cell load variability

Blue boxes: Optimal range of cell stock concentration to maximize the likelihood of achieving the desired cell recovery target

## 32 $\mu$ l Sample Input – Cell Suspension Volume Calculator Table

(Chromium Connect Automated Single Cell 5' v2 protocol)

Volume of Cell Suspension Stock per reaction ( $\mu$ l) | Volume of PBS + 0.04% BSA ( $\mu$ l)

| Cell Stock Conc. (cells/ $\mu$ l) | Targeted Cell Recovery |      |      |      |      |      |      |      |      |      |       |
|-----------------------------------|------------------------|------|------|------|------|------|------|------|------|------|-------|
|                                   | 500                    | 1000 | 2000 | 3000 | 4000 | 5000 | 6000 | 7000 | 8000 | 9000 | 10000 |
| 100                               | 8.3                    | 16.5 | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a   |
|                                   | 23.8                   | 15.5 | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a   |
| 200                               | 4.1                    | 8.3  | 16.5 | 24.8 | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a   |
|                                   | 27.9                   | 23.8 | 15.5 | 7.3  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a  | n/a   |
| 300                               | 2.8                    | 5.5  | 11.0 | 16.5 | 22.0 | 27.5 | n/a  | n/a  | n/a  | n/a  | n/a   |
|                                   | 29.3                   | 26.5 | 21.0 | 15.5 | 10.0 | 4.5  | n/a  | n/a  | n/a  | n/a  | n/a   |
| 400                               | 2.1                    | 4.1  | 8.3  | 12.4 | 16.5 | 20.6 | 24.8 | 28.9 | n/a  | n/a  | n/a   |
|                                   | 29.9                   | 27.9 | 23.8 | 19.6 | 15.5 | 11.4 | 7.3  | 3.1  | n/a  | n/a  | n/a   |
| 500                               | 1.7                    | 3.3  | 6.6  | 9.9  | 13.2 | 16.5 | 19.8 | 23.1 | 26.4 | 29.7 | n/a   |
|                                   | 30.4                   | 28.7 | 25.4 | 22.1 | 18.8 | 15.5 | 12.2 | 8.9  | 5.6  | 2.3  | n/a   |
| 600                               | 1.4                    | 2.8  | 5.5  | 8.3  | 11.0 | 13.8 | 16.5 | 19.3 | 22.0 | 24.8 | 27.5  |
|                                   | 30.6                   | 29.3 | 26.5 | 23.8 | 21.0 | 18.3 | 15.5 | 12.8 | 10.0 | 7.3  | 4.5   |
| 700                               | 1.2                    | 2.4  | 4.7  | 7.1  | 9.4  | 11.8 | 14.1 | 16.5 | 18.9 | 21.2 | 23.6  |
|                                   | 30.8                   | 29.6 | 27.3 | 24.9 | 22.6 | 20.2 | 17.9 | 15.5 | 13.1 | 10.8 | 8.4   |
| 800                               | 1.0                    | 2.1  | 4.1  | 6.2  | 8.3  | 10.3 | 12.4 | 14.4 | 16.5 | 18.6 | 20.6  |
|                                   | 31.0                   | 29.9 | 27.9 | 25.8 | 23.8 | 21.7 | 19.6 | 17.6 | 15.5 | 13.4 | 11.4  |
| 900                               | 0.9                    | 1.8  | 3.7  | 5.5  | 7.3  | 9.2  | 11.0 | 12.8 | 14.7 | 16.5 | 18.3  |
|                                   | 31.1                   | 30.2 | 28.3 | 26.5 | 24.7 | 22.8 | 21.0 | 19.2 | 17.3 | 15.5 | 13.7  |
| 1000                              | 0.8                    | 1.7  | 3.3  | 5.0  | 6.6  | 8.3  | 9.9  | 11.6 | 13.2 | 14.9 | 16.5  |
|                                   | 31.2                   | 30.4 | 28.7 | 27.1 | 25.4 | 23.8 | 22.1 | 20.5 | 18.8 | 17.2 | 15.5  |
| 1100                              | 0.8                    | 1.5  | 3.0  | 4.5  | 6.0  | 7.5  | 9.0  | 10.5 | 12.0 | 13.5 | 15.0  |
|                                   | 31.3                   | 30.5 | 29.0 | 27.5 | 26.0 | 24.5 | 23.0 | 21.5 | 20.0 | 18.5 | 17.0  |
| 1200                              | 0.7                    | 1.4  | 2.8  | 4.1  | 5.5  | 6.9  | 8.3  | 9.6  | 11.0 | 12.4 | 13.8  |
|                                   | 31.3                   | 30.6 | 29.3 | 27.9 | 26.5 | 25.1 | 23.8 | 22.4 | 21.0 | 19.6 | 18.3  |
| 1300                              | 0.6                    | 1.3  | 2.5  | 3.8  | 5.1  | 6.3  | 7.6  | 8.9  | 10.2 | 11.4 | 12.7  |
|                                   | 31.4                   | 30.7 | 29.5 | 28.2 | 26.9 | 25.7 | 24.4 | 23.1 | 21.8 | 20.6 | 19.3  |
| 1400                              | 0.6                    | 1.2  | 2.4  | 3.5  | 4.7  | 5.9  | 7.1  | 8.3  | 9.4  | 10.6 | 11.8  |
|                                   | 31.4                   | 30.8 | 29.6 | 28.5 | 27.3 | 26.1 | 24.9 | 23.8 | 22.6 | 21.4 | 20.2  |
| 1500                              | 0.6                    | 1.1  | 2.2  | 3.3  | 4.4  | 5.5  | 6.6  | 7.7  | 8.8  | 9.9  | 11.0  |
|                                   | 31.5                   | 30.9 | 29.8 | 28.7 | 27.6 | 26.5 | 25.4 | 24.3 | 23.2 | 22.1 | 21.0  |
| 1600                              | 0.5                    | 1.0  | 2.1  | 3.1  | 4.1  | 5.2  | 6.2  | 7.2  | 8.3  | 9.3  | 10.3  |
|                                   | 31.5                   | 31.0 | 29.9 | 28.9 | 27.9 | 26.8 | 25.8 | 24.8 | 23.8 | 22.7 | 21.7  |
| 1700                              | 0.5                    | 1.0  | 1.9  | 2.9  | 3.9  | 4.9  | 5.8  | 6.8  | 7.8  | 8.7  | 9.7   |
|                                   | 31.5                   | 31.0 | 30.1 | 29.1 | 28.1 | 27.1 | 26.2 | 25.2 | 24.2 | 23.3 | 22.3  |
| 1800                              | 0.5                    | 0.9  | 1.8  | 2.8  | 3.7  | 4.6  | 5.5  | 6.4  | 7.3  | 8.3  | 9.2   |
|                                   | 31.5                   | 31.1 | 30.2 | 29.3 | 28.3 | 27.4 | 26.5 | 25.6 | 24.7 | 23.8 | 22.8  |
| 1900                              | 0.4                    | 0.9  | 1.7  | 2.6  | 3.5  | 4.3  | 5.2  | 6.1  | 6.9  | 7.8  | 8.7   |
|                                   | 31.6                   | 31.1 | 30.3 | 29.4 | 28.5 | 27.7 | 26.8 | 25.9 | 25.1 | 24.2 | 23.3  |
| 2000                              | 0.4                    | 0.8  | 1.7  | 2.5  | 3.3  | 4.1  | 5.0  | 5.8  | 6.6  | 7.4  | 8.3   |
|                                   | 31.6                   | 31.2 | 30.4 | 29.5 | 28.7 | 27.9 | 27.1 | 26.2 | 25.4 | 24.6 | 23.8  |

Grey boxes: Volumes that would exceed the allowable buffer volume in each reaction  
 Yellow boxes: Indicate a low transfer volume that may result in higher cell load variability  
 Blue boxes: Optimal range of cell stock concentration to maximize the likelihood of achieving the desired cell recovery target


# Carrier Loading Guidelines

## Carrier Loading Guidelines


Follow the instructions on the touchscreen to load the carriers.

### Carriers

- Handle the carriers as prompted.
- Ensure that Carriers 3, 4, and 5 are completely slid out and placed on an off-deck workspace before loading.
- Align the carriers to the corresponding Deck Rails when sliding them in or out of the deck.
- Ensure correct orientation of tube labels with barcodes to enable Barcode Scanning.




### Barcode Orientation



### Modules


- Load one tube strip/sample of each of the indicated modules in the corresponding positions on the Carrier, starting from back to front row.
- DO NOT skip any rows when loading.
- Use pinhole alignment to place module tube strips in the correct orientation (as shown on the touchscreen).



### Label Tube Strip Orientation

- The cDNA tube strip will be at Position 1 and the final library tube strip will be at Position 4 in the Tube Strip Holder.
- Label tube strip orientation for collecting cDNA and final libraries.

### Tube Strip (TS) Holder




Consult the Chromium Connect User Guide (CG000180) for more information.

## Carrier Loading Guidelines


### Load Gel Beads

- Up to 2 Gel Bead tube strips may be loaded in the primary (GB-1°) and alternate (GB-Alt) positions. One Gel Bead tube is required/sample.
- If only loading one tube strip, load in the primary position.
- Select the location of the loaded Gel Bead tube on the touchscreen.
- Examples of various Gel Bead loading combinations are illustrated below.




### Example 1

8 samples run with 1 Gel Bead tube strip loaded in GB-1° location.




### Example 2

8 samples run with 2 Gel Bead tube strips loaded in GB-1° and GB-Alt locations.



### Example 3

4 samples run with 2 Gel Bead tube strips loaded in GB-1° and GB-Alt locations.




# Additional Protocol Guidelines

Confirm GEM Generation  
cDNA Amplification Cycles  
cDNA QC & Quantification

## Confirm GEM Generation

- Instrument will pause for 5 min during GEM QC.
- Carefully remove Full Skirted Plate from the Thermal Cycler.
- Hold up the Full Skirted Plate and view the bottom of the wells in Column 3 to confirm GEM generation (shown below).
- Reload Full Skirted Plate in the Thermal Cycler.




### GEM Generation Example

Tube 1 shows normal GEM generation

Tube 2 shows a wetting failure

Tube 3 indicates a clog



## cDNA Amplification Cycle Number

- cDNA amplification cycles are determined by target cell number.
- Recommended guidelines for selecting optimal amplification cycle numbers

Recommended starting point for cycle number optimization.

| Targeted Cell Recovery | Low RNA Content Cells<br>e.g., Primary Cells<br>Total Cycles | High RNA Content Cells<br>e.g., Cell Lines<br>Total Cycles |
|------------------------|--|--|
| 500-2,000              | 16   | 14   |
| 2,001-6,000            | 14   | 12   |
| 6,001-10,000           | 13   | 11   |

- The optimal number of cycles is a trade-off between generating sufficient final mass for library construction and minimizing PCR amplification artifacts. The number of cDNA cycles should also be reduced if large numbers of cells are sampled.



Cycle number selected for one sample will apply to all the other samples in the run.

For V(D)J Target Amplification workflow, save 2µl cDNA for every target reaction.


## cDNA QC & Quantification



For V(D)J + GEX Library Construction proceed directly to GEX Library Construction first, followed by V(D)J Amplification and V(D)J Library Construction. If GEX library is not desired, proceed directly to V(D)J Amplification.

- Follow the instruction on the touchscreen for cDNA QC & quantification.
- Run sample on an Agilent Bioanalyzer High Sensitivity chip. Run 1  $\mu$ l undiluted product for input cells with low RNA content (<1 pg total RNA/cell), and 1  $\mu$ l of 1:10 diluted product for input cells with high RNA content.

Representative Trace for PBMCs




- If proceeding to 5' GEX Library Construction, determine cDNA yield for each sample. Example calculation below.
- Enter the cDNA concentration (pg/ $\mu$ l) and the calculated input volume ( $\mu$ l) on the touchscreen to proceed with GEX library construction.


### EXAMPLE CALCULATION

#### i. Select Region

Under the "Electropherogram" view choose the "Region Table". Manually select the region of ~200 – ~9000 bp



#### ii. Note Concentration [pg/ $\mu$ l]



#### iii. Calculate

Concentration: 2787.20 pg/ $\mu$ l  
Dilution Factor: 1

$$\text{cDNA Conc.} = \frac{\text{Conc. (pg/\mu l)} \times \text{Dilution Factor}}{1000 \text{ (pg/ng)}} = \frac{2787.20 \times 1}{1000} = 2.79 \text{ ng/\mu l}$$

#### Example Calculation for Carrying Forward 60 ng Sample for 5' GEX Library Construction

$$\text{Volume for 60 ng} = \frac{60 \text{ ng}}{2.79 \text{ (ng/\mu l)}} = 21.5 \mu \text{l}$$

- If the volume for 60 ng exceeds 22  $\mu$ l, carry ONLY 22  $\mu$ l sample into library construction. The sample input volume should be in the 5-22  $\mu$ l range.

$$\text{Sample volume for library construction} = 21.5 \mu \text{l}$$

If <60\* ng available, carry forward 22  $\mu$ l sample (2-60 ng) into 5' GEX Library Construction.

\*Note that the intended sample amount differs from manual protocol to account for pipetting differences in automation.



DO NOT exceed a mass of 60 ng in the 22  $\mu$ l carry forward volume.

### Alternate Quantification Methods:

- Agilent TapeStation

- PerkinElmer LabChip (See Appendix for representative traces)



# 5' Gene Expression (GEX) Library Construction Guidelines

Sample Index PCR

Post Library Construction QC

## Sample Index PCR

- The cycle numbers can be manually selected based on cDNA input.
- Recommended guidelines for selecting optimal Sample Index PCR cycle number.


| cDNA Input | Total Cycles |
|------------|--------------|
| 1-30 ng    | 16           |
| 31-60 ng   | 14           |



Cycle number selected will apply to all the samples in the run.

## Post Library Construction QC

Run sample on an Agilent Bioanalyzer High Sensitivity chip.



Determine the average fragment size from the Bioanalyzer trace. This will be used as the insert size for library quantification.

### Alternate QC Method:

- Agilent TapeStation
- PerkinElmer LabChip ([See Appendix for representative traces](#))

# V(D)J Amplification & Library Construction Guidelines

Deck Orientation for V(D)J Amplification


Gather Items & Reagents

Thaw & Prep Reagents


V(D)J Amplification

Post Library Construction QC





## Deck Orientation – V(D)J Amplification & Library Construction



P2 - Primer Mix 2  
P1 - Primer Mix 1



Refer to the Chromium Connect Instrument User Guide (CG000180) and Quick Reference Cards (CG000256) for more information.

|  |                | <b>Deck Layout Reagents/Consumables</b><br>V(D)J Amplification and Library Construction |  |
|--|----------------|---|--|
| Carrier  | Zone           | Item  |  |
|  <p>1</p>   | Zone 1 (Black) | 22°C Block, Reagent Strips, Module 1  |  |
|  | Zone 3 (White) | 10°C Block, Reagent Strips, Module 2  |  |
|  | -              | Magnetic Plate  |  |
|  <p>2</p>  | -              | Thermal Cycler  |  |
|  | -              | Full Skirted PCR Plate (within Thermal Cycler)  |  |
|  | Zone 4         | ComfortLids   |  |
|  <p>3*</p> <p>Deck Rails: 15-18<br/>Number of Lights: 4</p> <p>*Assay choices determine items loaded in Carrier 3</p> | Position W     | Waste Reservoirs  |  |
|  | Position TS    | Tube Strips (positions 1 & 4)   |  |
|  | Position P2    | Primer Mix 2  |  |
|  | Position P1    | Primer Mix 1  |  |
|  | Position CP    | CO-RE Paddles   |  |
|  | Position E     | Ethanol Reservoir   |  |
|  | Position Lid   | Lid for Ethanol Reservoir   |  |
|  <p>4</p>   | Zone 5         | Sample Index Plate  |  |
|  | Zone 6         | Pipette Tips 50 µl  |  |
|  | Zone 7         | Semi Skirted Plate  |  |
|  <p>5</p> <p>Deck Rails: 25-30<br/>Number of Lights: 6</p>  | Zone 8         | Pipette Tips 300 µl   |  |

## Gather Items & Reagents for V(D)J Amplification and Library Construction



2 µl of cDNA is needed for each V(D)J Target Amplification reaction.

Follow prompts on the Chromium Connect touchscreen to gather the listed items and reagents for loading the Deck Carriers. Gather the quantities specified for each of the items and reagents.

| Item   | Qty                 |
|--|---------------------|
| Nuclease-free Water  | 10 ml               |
| Ethanol, Pure (200 Proof, anhydrous)   | 40 ml               |
| <b>Hamilton</b>  |                     |
| ComfortLids  | 6                   |
| 50 µl CO-RE Pipette Tips, with filter (Black, Conductive)                        | 2 racks             |
| 300 µl CO-RE Pipette Tips, with filter (Black, Conductive)                       | 3 racks             |
| Reagent Reservoir, 60 ml   | 2                   |
| <b>Eppendorf</b>   |                     |
| 96-well Semi Skirted Plate   | 1                   |
| 96-well Full Skirted Plate   | 1                   |
| <b>Thermo Fisher Scientific</b>  |                     |
| MicroAmp 8-Tube Strip, 0.2 ml  | 2                   |
| <b>10x Genomics</b>  |                     |
| Chromium Automated Single Cell Human TCR Amplification & Library Construction v2 |                     |
| V(D)J Module 1 (stored at 4°C) <i>Black tube strip</i>                           | 1 tube strip/sample |
| V(D)J Module 2 (stored at -20°C) <i>White tube strip</i>                         | 1 tube strip/sample |
| <i>Human T Cell Primer Mix 1 v2</i>  | 1 tube/run          |
| <i>Human T Cell Primer Mix 2 v2</i>  | 1 tube/run          |
| Chromium Automated Single Cell Mouse TCR Amplification & Library Construction v2 |                     |
| V(D)J Module 1 (stored at 4°C) <i>Black tube strip</i>                           | 1 tube strip/sample |
| V(D)J Module 2 (stored at -20°C) <i>White tube strip</i>                         | 1 tube strip/sample |
| <i>Mouse T Cell Primer Mix 1 v2</i>  | 1 tube/run          |
| <i>Mouse T Cell Primer Mix 2 v2</i>  | 1 tube/run          |
| Chromium Automated Single Cell Human BCR Amplification & Library Construction v2 |                     |
| V(D)J Module 1 (stored at 4°C) <i>Black tube strip</i>                           | 1 tube strip/sample |
| V(D)J Module 2 (stored at -20°C) <i>White tube strip</i>                         | 1 tube strip/sample |
| <i>Human B Cell Primer Mix 1 v2</i>  | 1 tube/run          |
| <i>Human B Cell Primer Mix 2 v2</i>  | 1 tube/run          |

## Thaw & Prep Reagents for V(D)J Amplification & Library Construction

Follow prompts on the touchscreen to thaw and prepare reagents. Some important guidelines are highlighted below.

| ACTION                    | GUIDELINES<br><i>Follow touchscreen prompts for specifics and timing</i>   |
|---------------------------|--|
| Thaw Reagents             | <ul style="list-style-type: none"> <li>Thaw reagents as indicated on the touchscreen. Verify no precipitate is present.</li> <li>Ensure that the correct thawing locations and temperatures are used.</li> <li>During reagent thaw load the consumables following touchscreen prompts.</li> </ul>  |
| Prepare Ethanol           | <ul style="list-style-type: none"> <li>Prepare <b>50 ml</b> 80% Ethanol in Nuclease-free water and dispense in Ethanol Reservoir when prompted.</li> </ul>   |
| V(D)J Modules             | <ul style="list-style-type: none"> <li>Thaw V(D)J Modules as prompted on the touchscreen.</li> <li>After reagent thaw, invert rack holding Module tube strips and vortex V(D)J Modules 1 for <b>30 sec</b>; verify no precipitate.</li> <li>Confirm that there are no bubbles at the bottoms of any module tubes.</li> <li>Centrifuge V(D)J Module 1 at <b>300 rcf</b> for <b>1 min</b> at <b>22°C</b>.</li> <li>Retrieve V(D)J Module 2 from <b>4°C</b> thaw. <b>DO NOT</b> vortex. Invert-mix and centrifuge at <b>300 rcf</b> for <b>1 min</b> at <b>22°C</b>.</li> </ul> |
| Dual Index Plate TT Set A | <ul style="list-style-type: none"> <li>Vortex Dual Index Plate for 15 sec at maximum speed and centrifuge at <b>300 rcf</b> for <b>1 min</b> at <b>22°C</b>.</li> </ul>  |
| Primer Mix 1 & 2          | <ul style="list-style-type: none"> <li>Vortex and centrifuge before loading.</li> </ul>  |




Confirm that there are no bubbles at the bottoms of any module tubes, Dual Index Plate wells, or Primer Mix tubes.

## V(D)J Amplification QC & Quantification


- Follow the instruction on the touchscreen for V(D)J Amplification QC & quantification.
- Run 1  $\mu\text{l}$  sample at 1:5 dilution (Dilution Factor 5) on an Agilent Bioanalyzer High Sensitivity chip.

Samples of RNA-rich cells may require additional dilution in nuclease-free water. The number of distinct peaks may vary. Higher molecular weight product (2,000- 9,000 bp) may be present. This does not affect sequencing.

Representative Trace - PBMCs amplified for TCR



Representative Trace - PBMCs amplified for BCR




- Determine yield for each sample using the example calculation below.
- Enter the V(D)J amplified product concentration (pg/ $\mu\text{l}$ ) and the calculated input volume ( $\mu\text{l}$ ) on the touchscreen to proceed with V(D)J library construction.


### EXAMPLE CALCULATION

#### i. Select Region

Under the “Electropherogram” view choose the “Region Table”. Manually select the region of ~200 – ~9000 bp.



#### ii. Note Concentration [pg/ $\mu\text{l}$ ]



#### iii. Calculate

Concentration: 1773.07 pg/ $\mu\text{l}$   
Dilution Factor:

V(D)J Amplified Product Conc.

$$\frac{\text{Conc. (pg/\mu l)} \times \text{Dilution Factor}}{1000 \text{ (pg/ng)}} = \frac{1773.07 \times 5}{1000} = 8.9 \text{ ng/\mu l}$$

#### Example Calculation for Carrying Forward 60 ng Sample for V(D)J Library Construction

$$\text{Volume for 60 ng} = \frac{60 \text{ ng}}{8.9 \text{ (ng/\mu l)}} = 6.7 \mu\text{l}$$

- The sample input volume should be in the 5-22  $\mu\text{l}$  range.

If <60\* ng available, carry forward 22  $\mu\text{l}$  sample (2-60 ng) into V(D)J Library Construction.

\*Note that the intended sample amount differs from manual protocol to account for pipetting differences in automation.



DO NOT exceed a mass of 60 ng in the 22  $\mu\text{l}$  carry forward volume.

### Alternate Quantification Methods

- PerkinElmer LabChip (See Appendix for representative traces)
- Agilent TapeStation




## V(D)J Library Construction QC


Run sample on an Agilent Bioanalyzer High Sensitivity chip.

### Representative Trace

Representative Trace - PBMCs amplified for TCR



Representative Trace - PBMCs amplified for BCR



Determine the average fragment size from the Bioanalyzer trace. This will be used as the insert size for library quantification.

### Alternate QC Methods

- LabChip ([See Appendix for representative traces](#))
- Agilent TapeStation

# Post Library Construction Quantification & Pooling

[Deck Orientation – Library Quantification](#)


[Post Library Construction Quantification](#)

[Deck Orientation – Library Pooling](#)

[Library Pooling](#)

## Deck Orientation – Library Quantification

Library quantification using qPCR is recommended for accurate pooling and loading on sequencers. If the option is selected during gene expression run-setup, automated qPCR plate-setup can be run directly on Chromium Connect after library generation and final library QC. Alternatively, the option can be selected from the instrument home screen, at the user's convenience. Up to 8 samples can be quantified on a 96 well reaction plate, including duplicates for each sample. The minimum sample volume required is 25  $\mu\text{l}$ . Only 6  $\mu\text{l}$  of the sample will be used for qPCR plate setup.



## Gather Items & Reagents

Follow prompts on the Chromium Connect touchscreen to gather the listed items and reagents for loading the Deck Carriers for Library Quantification.

Gather the quantities specified for each of the items and reagents.

| Item   | Qty     |
|--|---------|
| <b>Hamilton</b>  |         |
| ComfortLid   | 1       |
| 50 µl CO-RE Pipette Tips, with filter (Black, Conductive)  | 2 racks |
| 300 µl CO-RE Pipette Tips, with filter (Black, Conductive) | 2 racks |
| 60-ml Reagent Reservoir                                    | 1       |
| <b>Eppendorf</b>   |         |
| 96-well Semi Skirted Plate                                 | 1       |
| <b>Thermo Fisher Scientific</b>                            |         |
| 2-ml Tube with Screw Cap                                   | 1       |
| <b>Bio-Rad</b>   |         |
| 96-well Hard-Shell Full Skirted Plate                      | 1       |
| <b>Reagent</b>   |         |
| Qiagen Buffer EB   | 50 ml   |
| Nuclease-free Water  | 1 ml    |
| 10% Tween-20   | 250 µl  |
| Libraries (in an 8-tube strip)                             | 1-8     |
| KAPA Library Quantification Kit, thawed                    |         |
| SYBR FAST Master Mix                                       | 5 ml    |
| Primer Mix   | 1 ml    |
| Standards  | 6       |

## Post Library Construction Quantification

- Prepare reagents as prompted on the touchscreen.
- Vortex and centrifuge KAPA standards and libraries before use.
- Retrieve previously prepared Master Mix + Primer Mix  
OR  
Add 1 ml Primer Mix to 5 ml SYBR FAST Master Mix.
- Prepare specified Quantification Master Mix in the 2-ml tube using the guidance below.

| # Sample | Master Mix + Primer Mix (μl) | Water (μl) | Total Vol (μl) |
|----------|------------------------------|------------|----------------|
| 8        | 1305                         | 435        | 1740           |
| 7        | 1200                         | 400        | 1600           |
| 6        | 1095                         | 365        | 1460           |
| 5        | 990                          | 330        | 1320           |
| 4        | 885                          | 295        | 1180           |
| 3        | 780                          | 260        | 1040           |
| 2        | 675                          | 225        | 900            |
| 1        | 570                          | 190        | 760            |


Volumes listed take into account volume for 6 standards

- Follow the touchscreen prompts for loading, scanning, and executing the run.
- During the run, the following steps will be executed by the instrument:
  - KAPA Master Mix transfer to the 96-well Hard Shell Full Skirted Plate (layout below)
  - Diluent transfer to dilution plate
  - Serial dilutions of libraries
  - Addition of library dilutions, KAPA Standards, and negative controls to the plate

Total reaction volume (20  $\mu$ l)=  
 16  $\mu$ l Master Mix  
 +  
 4  $\mu$ l Library Dilution/  
 KAPA Standard/  
 Negative Control (NTC)

Dilutions:

1:12,500  
 1:62,500  
 1:312,500  
 1:1,562,500




- After the run is completed, follow the unloading instructions on the touchscreen.
- Cap and store libraries at 4°C  $\leq$ 72 h or -20°C  $\leq$ 4 weeks.
- Remove Full Skirted Plate. Seal plate and centrifuge at 300 rcf for 1 min at 22°C.
- Follow the manufacturer's recommendations for qPCR-based quantification. For library quantification for sequencer clustering, determine the concentration based on average size (bp) derived from the Bioanalyzer/TapeStation trace.

| Step | Temperature                         | Run Time |
|------|-------------------------------------|----------|
| 1    | 95°C                                | 00:05:00 |
| 2    | 95°C                                | 00:00:30 |
| 3    | 60°C                                | 00:00:45 |
| 4    | Go to Step 2, 29X (Total 30 cycles) |          |

- **Resource:**  
 Use the Chromium Connect Library Quantification Worksheet (CG000157) provided on the 10x Genomics Support website for calculating library concentrations.

## Deck Orientation – Library Pooling

The libraries may be pooled on the Chromium Connect instrument and used for sequencing, taking into account the preferred cell numbers and per-cell read depth requirements for each library. Samples utilizing the same sample index should not be pooled together, or run on the same flow cell lane, as this would prevent correct sample demultiplexing. The Chromium Connect deck layout for Library Pooling is shown below.



## Gather Items & Reagents

Follow prompts on the Chromium Connect touchscreen to gather the listed items and reagents for loading the Deck Carriers for Library Pooling.

Gather the quantities specified for each of the items and reagents.

| Item   | Qty                |
|--|--------------------|
| <b>Hamilton</b>  |                    |
| 50 µl CO-RE Pipette Tips, with filter (Black, Conductive)  | 1-2 rack           |
| 300 µl CO-RE Pipette Tips, with filter (Black, Conductive) | 1-2 rack           |
| Reagent Reservoir, 60 ml                                   | 1                  |
| <b>Eppendorf</b>   |                    |
| 96-well Semi Skirted Plate                                 | 1                  |
| <b>Thermo Fisher Scientific</b>                            |                    |
| 0.5-ml Tube with Screw Cap                                 | 6                  |
| MicroAmp 8-Tube Strip, 0.2 ml                              | 1-4                |
| <b>Reagent</b>   |                    |
| Qiagen Buffer EB   | 50 ml              |
| Libraries (in up to four 8-tube strips)                    | up to 32 libraries |

## Library Pooling

- Follow the touchscreen prompts for loading, scanning, and executing the run.
- Briefly vortex and centrifuge libraries in the 8-tube strip.
- Confirm that there are no bubbles at the bottoms of any library tubes.
- Ensure a minimum **25 µl** library volume is available in the tubes.
- After run is complete, follow touchscreen prompts to unload and store the libraries.
- Unload remaining items and clean as prompted on the touchscreen.
- **Resource:** Use the Chromium Connect Library Pooling Worksheet (CG000466) provided on the 10x Genomics Support website to calculate volumes to be pooled. The calculated volumes can be input into the instrument either manually, or via the CSV file generated from this workbook.




# Sequencing

## Sequencing Libraries


Chromium Single Cell V(D)J and 5' Gene Expression Dual Index libraries comprise standard Illumina paired-end constructs which begin with P5 and end with P7. These libraries include 16 bp 10x Barcodes encoded at the start of TruSeq Read 1. Sample index sequences are incorporated as the 10 bp i5 and i7 index reads.

TruSeq Read 1 and TruSeq Read 2 are standard Illumina sequencing primer sites used in paired-end sequencing of V(D)J and 5' Gene Expression libraries. Sequencing these libraries produce a standard Illumina BCL data output folder.

### Chromium Single Cell V(D)J Dual Index Library



### Chromium Single Cell 5' Gene Expression Dual Index Library



## Illumina Sequencer Compatibility

The compatibility of the listed sequencers has been verified by 10x Genomics. Some variation in assay performance is expected based on sequencer choice. For more information about performance variation, visit the 10x Genomics Support website.

- MiSeq
- NextSeq 500/550/2000
- HiSeq 2500 (Rapid Run)
- HiSeq 3000/4000
- NovaSeq

## Sample Indices

Each well of the Dual Index Kit TT Set A (PN-1000215) contains a mix of one unique i7 and one unique i5 sample index. If multiple samples are pooled in a sequence lane, the sample index name (i.e. the Dual Index plate well ID) is needed in the sample sheet used for generating FASTQs with “cellranger mkfastq”.

If multiple libraries are pooled in a sequence lane, a separate sample index is needed with each library (see [Tips & Best Practices](#)).

## Library Sequencing Depth & Run Parameters

|                  |  |
|------------------|--|
| Sequencing Depth | Minimum 5,000 read pairs per cell for V(D)J library<br>Minimum 20,000 read pairs per cell for 5' Gene Expression library |
| Sequencing Type  | Paired-end, Dual indexing  |
| Sequencing Read  | Read 1: 26 cycles<br>i7 Index: 10 cycles<br>i5 Index: 10 cycles<br>Read 2: 90 cycles                                     |

## Library Loading

Once quantified and normalized, V(D)J and 5' Gene Expression libraries should be denatured and diluted as recommended for Illumina sequencing platforms. Refer to Illumina documentation for denaturing and diluting libraries. Refer to the 10x Genomics Support website for more information.

| Instrument   | Loading Concentration (pM) | PhiX (%) |
|--------------|----------------------------|----------|
| MiSeq        | 10                         | 1        |
| NextSeq 500  | 1.5                        | 1        |
| NovaSeq      | 150*/300                   | 1        |
| NextSeq 2000 | 650                        | 1        |

\* Use 150 pM loading concentration for Illumina XP workflow.

## Library Pooling

V(D)J and 5' Gene Expression libraries may be pooled for sequencing, taking into account the differences in depth requirements between the pooled libraries. 5' Gene Expression libraries may be sequenced using enriched library parameters, however the cost of sequencing using enriched library parameters is higher.

Refer to [Post Library Construction Quantification & Pooling](#) chapter for library pooling on the Chromium Connect instrument.

Library Pooling Examples:

| Libraries                  | Sequencing Depth (read pairs per cell) | Library Pooling Ratio |
|----------------------------|--|-----------------------|
| <b>Example 1</b>           |  |                       |
| V(D)J library              | 5,000                                  | 1                     |
| 5' Gene Expression library | 20,000                                 | 4                     |
| <b>Example 2</b>           |  |                       |
| V(D)J library              | 5,000                                  | 1                     |
| 5' Gene Expression library | 50,000                                 | 10                    |

# Appendix


LabChip Traces

Oligonucleotide Sequences


## LabChip Traces

LabChip Traces DNA High Sensitivity Reagent Kit was used.


### cDNA QC & Quantification




### Post TCR Amplification QC





### Post BCR Amplification QC



### GEX Post Library Construction QC



### Post V(D)J Library Construction QC (PBMCs amplified for TCR)    Post V(D)J Library Construction QC (PBMCs amplified for BCR)



Alternate QC Method:

[Qubit Fluorometer and Qubit dsDNA HS Assay Kit](#)

## Oligonucleotide Sequences

Protocol steps correspond to the Chromium Next GEM Automated Single Cell 5' v2 protocol.

### GEM-RT Incubation

#### Gel Bead Primer



5'-CTACACGACGCTCTCCGATCT-NNNNNNNNNNNNNNNNNNNN-NNNNNNNNNN-TTCTTATATrGrGrG-3'

#### Poly-dT RT Primer PN-200007



5'-AAGCAGTGGTATCAACGCAGAGTACTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTT-3'

#### GEM-RT Products



3'-GATGTGCTGCGAGAAGGCTAGA-N16-N10-AAAGAATATACC-cDNA\_Insert-NVTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTT-CATGAGACGCAACTATGGTGACGAA-5'

### Automated Protocol Step – cDNA Amplification


#### cDNA Primers

Forward Primer: Partial Read 1  
5'-CTACACGACGCTCTCCGATCT-3'

Reverse Primer: Non-poly(dT)  
5'-AAGCAGTGGTATCAACGCAGAG-3'

#### Amplified Products

##### Amplified cDNA from poly-adenylated mRNA



### Automated Protocol Step – GEX Adaptor Ligation (for 5' Gene Expression (GEX) Library Construction)


#### Adaptor Read 2

Read 2

5' -GATCGGAAGAGCACACGTCTGAACTCCAGTCAC-3'

3' -TCTAGCCTTCTCG-5'

#### Ligation Product



5'-GATCTACACTCTTCCCTACACGACGCTCTCCGATCT-N16-N10-TTCTTATATGGG-cDNA\_Insert-AGATCGGAAGAGCACACGTCTGAACTCCAGTCAC-3'

3'-CTAGATGTGAGAAAGGGATGTGCTGCGAGAAGGCTAGA-N16-N10-AAAGAATATACC-cDNA\_Insert-TCTAGCCTTCTCG-5'

### Automated Protocol Step – Sample Index PCR (for 5' Gene Expression (GEX) Library Construction)

#### Dual Indexing Dual Index TT Set A PN-1000215

Forward Primer: P5 Sample Partial Read 1  
Index (i5)

5'-AATGATACGGCGACCACCGAGATCT-N10-ACACTCTTCCCTACACGACGCTC-3'

Reverse Primer: P7 Sample Partial Read 2  
Index (i7)

5'-CAAGCAGAAGACGGGCATCGAGAT-N10-GTGACTGGAGTTCAAGCTGTG-3'

#### Sample Index PCR Product



5'-AATGATACGGCGACCACCGAGATCTACAC-N10-ACACTCTTCCCTACACGACGCTCTCCGATCT-N16-N10-TTCTTATATGGG-cDNA\_Insert-AGATCGGAAGAGCACACGTCTGAACTCCAGTCAC-N10-ATCTCGTATGCCGTCTTCTGCTTG-3'

3'-TTACTATGCCCTGGTGGCTCTAGATGTG-N10-TGTGAGAAAGGGATGTGCTGCGAGAAGGCTAGA-N16-N10-AAAGAATATACC-cDNA\_Insert-TCTAGCCTTCTCGTGTGCGAGACTTGAAGTCAGTG-N10-TAGACATACCGCAGAGCAGAC-5'

## Automated Protocol Step – V(D)J Amplification 1

|  |  |   |              |
|--|--|---|--------------|
| Human T Cell Mix<br>1 v2<br>PN-2000242 | Forward Primer:<br>PCR Primer<br>5'-GATCTACTACTCTTCCCTACACGACGC-3' | Reverse Outer Primers:<br>5'-TGAAGGCGTTGCACATGCA-3'<br>5'-TCAGGCAGTATCTGGAGTCATTGAG-3'  | Outer Primer |
| Human B Cell Mix<br>1 v2<br>PN-2000254 | Forward Primer:<br>PCR Primer<br>5'-GATCTACTACTCTTCCCTACACGACGC-3' | Reverse Outer Primers:<br>5'-CAGGGCACAGTCACATCCT-3'<br>5'-TGCTGGACCACGATTTGTA-3'<br>5'-GGTTTTGTTGTCGACCCAGTCT-3'<br>5'-TTGTCCACCTTGGTGTGCT-3'<br>5'-CATGACGTCCTTGAAGGCA-3'<br>5'-TGTGGGACTTCCACTG-3'<br>5'-TTCTCGTAGTCTGCTTTGCTCAG-3'   | Outer Primer |
| Mouse T Cell Mix<br>1 v2<br>PN-2000256 | Forward Primer:<br>PCR Primer<br>5'-GATCTACTACTCTTCCCTACACGACGC-3' | Reverse Outer Primers:<br>5'-CTGTTGCTCCAGGCAATGG-3'<br>5'-TGTAGGCTGAGGGTCCGT-3'   | Outer Primer |
| Mouse B Cell Mix<br>1 v2<br>PN-2000258 | Forward Primer:<br>PCR Primer<br>5'-GATCTACTACTCTTCCCTACACGACGC-3' | Reverse Outer Primers:<br>5'-TCAGCACGGGACAACTCTCT-3'<br>5'-GCAGGAGACAGACTCTTCTCCA-3'<br>5'-AACTGGCTGCTCATGGTGT-3'<br>5'-TGGTCAAGTGTGGTTGAGGT-3'<br>5'-TGGTCACTTGGCTGGTGGT-3'<br>5'-CACTTGGCAGGTGAACTGTTTCT-3'<br>5'-AACCTTCAAGGATGCTTTGGGA-3'<br>5'-GGACAGGGATCCAGAGTTCCA-3'<br>5'-AGGTGACGGTCTGACTTGGC-3'<br>5'-GCTGGACAGGGCTCATAAGTT-3'<br>5'-GGCACCTTGTCCAATCATGTTCC-3'<br>5'-ATGTCGTTACATCTGCTTGGT-3' | Outer Primer |

## Automated Protocol Step – V(D)J Amplification 2

|  |  |   |              |
|--|--|---|--------------|
| Human T Cell Mix<br>2 v2<br>PN-2000246 | Forward Primer:<br>PCR Primer<br>5'-GATCTACTACTCTTCCCTACACGACGC-3' | Reverse Inner Primers:<br>5'-AGTCTCTCAGCTGGTACAGC-3'<br>5'-TCTGATGGCTCAAACACAGC-3'  | Inner Primer |
| Human B Cell Mix<br>2 v2<br>PN-2000255 | Forward Primer:<br>PCR Primer<br>5'-GATCTACTACTCTTCCCTACACGACGC-3' | Reverse Inner Primers:<br>5'-GGGAAGTTTCTGGCGTCA-3'<br>5'-GGTGGTACCAGTTATCAAGCAT-3'<br>5'-GTGTCCAGGTCAACATCAC-3'<br>5'-TCTTGGAGACTGTAGGACAGC-3'<br>5'-CACGCTGCTCGTATCCGA-3'<br>5'-TAGCTGCTGGCCGC-3'<br>5'-GCGTTATCCACTTCCACTGT-3'  | Inner Primer |
| Mouse T Cell Mix<br>2 v2<br>PN-2000257 | Forward Primer:<br>PCR Primer<br>5'-GATCTACTACTCTTCCCTACACGACGC-3' | Reverse Inner Primers:<br>5'-AGTCAAAGTCGGTGAACAGGCA-3'<br>5'-GGCCAAGCACAGGGGTA-3'   | Inner Primer |
| Mouse B Cell Mix<br>2 v2<br>PN-2000259 | Forward Primer:<br>PCR Primer<br>5'-GATCTACTACTCTTCCCTACACGACGC-3' | Reverse Inner Primers:<br>5'-TACACACAGTGTGGCCTT-3'<br>5'-CAGGCCACTGTCACACCACT-3'<br>5'-CAGGTCACATTCATCGTCCGC-3'<br>5'-GAGGCCAGCACAGTGACCT-3'<br>5'-GCAGGGAAGTTCACAGTGCT-3'<br>5'-CTGTTTGAGATCAGTTGCCATCCT-3'<br>5'-TGCGAGGTGGCTAGGTAAGT-3'<br>5'-CCCTTGACCAGGCATCC-3'<br>5'-AGGTCACGGAGGAACCAAGTTG-3'<br>5'-GGCATCCCAGTGTACCCGA-3'<br>5'-AGAAGATCCACTTCACTTGAAC-3'<br>5'-GAAGCACAGACTGAGGCAC-3' | Inner Primer |

## V(D)J Amplified Product




5'-GATCTACACTCTTTCCCTACACGACGCTCTTCCGATCT-N16-N10-TTTCTTATATGGG-cDNA\_Insert-Inner\_Primer-3'  
 3'-CTAGATGTGAGAAAGGGATGTGCTGCGAGAAGGCTAGA-N16-N10-AAAGAATATACCC-cDNA\_Insert-Inner\_Primer-5'

## Automated Protocol Step – Adaptor Ligation (for V(D)J Library Construction)

## Adaptor (Read 2)

5'-GATCGGAAGAGCACACGTCTGAACTCCAGTCAC-3'  
 3'-TCTAGCCTTCTCG-5'

## Ligation Product



5'-GATCTACACTCTTTCCCTACACGACGCTCTTCCGATCT-N16-N10-TTTCTTATATGGG-cDNA\_Insert-AGATCGGAAGAGCACACGTCTGAACTCCAGTCAC-3'  
 3'-CTAGATGTGAGAAAGGGATGTGCTGCGAGAAGGCTAGA-N16-N10-AAAGAATATACCC-cDNA\_Insert-TCTAGCCTTCTCG-5'

## Automated Protocol Step – Sample Index PCR (for V(D)J Library Construction)

## Dual Indexing

Forward Primer:

P5 Sample Partial Read 1 Index (i5)

Reverse Primer:

P7 Sample Partial Read 2 Index (i7)

Dual Index Kit TT Set A  
 PN-1000215

5'-AATGATACGGCGACCACCGAGATCTACAC-N10-ACACTCTTTCCCTACACGACGCTCTTCCGATCT-N16-N10-TTTCTTATATGGG-Insert-AGATCGGAAGAGCACACGTCTGAACTCCAGTCAC-N10-ATCTCGTATGCGGCTCTTCTGCTTG-3'

5'-CAAGCAGAAGACGGCATACGAGAT-N10-GTGACTGGAGTTCAGACGTGT-3'

## Sample Index PCR Product



5'-AATGATACGGCGACCACCGAGATCTACAC-N10-ACACTCTTTCCCTACACGACGCTCTTCCGATCT-N16-N10-TTTCTTATATGGG-Insert-AGATCGGAAGAGCACACGTCTGAACTCCAGTCAC-N10-ATCTCGTATGCGGCTCTTCTGCTTG-3'  
 3'-TTACTATGCGGCTGGTGGCTCTAGATGTG-N10-TGTGAGAAAGGGATGTGCTGCGAGAAGGCTAGA-N16-N10-AAAGAATATACCC-Insert-TCTAGCCTTCTCGTGTGACAGACTTGAGTCACTG-N10-TAGACATACCGCAGAAGACGAAC-5'