

Visium CytAssist for FFPE Assay Protocol Planner


















Introduction

The Visium CytAssist is compatible with FFPE tissue sections placed on blank slides. This document provides workflow overviews for processing FFPE tissue and additional equipment, kits, and reagents for the following assays:

- Visium CytAssist Spatial Gene Expression
- Visium CytAssist Spatial Gene and Protein Expression

Reagents and hardware produced by 10x Genomics are not listed in this document. Unless otherwise indicated, reagent and equipment tables apply to both assays.

Workflow Overviews

	Sample Preparation	Deparaffinization & Staining / Imaging	Probe Hybridization & Ligation	Protein Labeling	Permeabilization & Barcoding	Transfer to Tube	Library Construction		
CytAssist Gene Expression for FFPE	FFPE Tissue Sections (human or mouse) on plain glass slide  STOP	Deparaffinization  -3 h	Immunofluorescence (IF) or Hematoxylin & Eosin (H&E)  >1 h STOP	Decrosslinking  -1 h	Whole Transcriptome Probe-Mediated mRNA Detection  Overnight STOP	Probe Release & Extension  -2 h	Pre-amplification, Cleanup, qPCR  -2 h STOP	SI-PCR, Cleanup & QC  -3 h STOP	
CytAssist Gene & Protein Expression for FFPE	FFPE Tissue Sections (human) on plain glass slide  STOP	Deparaffinization  -3 h	Immunofluorescence (IF) or Hematoxylin & Eosin (H&E)  >1 h STOP	Decrosslinking, DNase Treatment  -2 h	Whole Transcriptome Probe-Mediated mRNA Detection  Overnight STOP	Protein Blocking & Labeling  Overnight	Probe & Antibody Tag Release & Extension  -2 h	Pre-amplification, Cleanup, qPCR (2x)  -2 h STOP	SI-PCR, Cleanup & QC (2x)  -4.5 h STOP

*for both assays, stopping points during staining apply only to H&E Staining.

Contents

- 2 Pipette and Plpette Tip Recommendations
- 3 Sample Preparation Reagents
- 4 RNA Assessment (DV200) Reagents
- 5 Add-on Antibody Optimization
- 7 Deparaffinization, H&E Staining, & Decrosslinking Reagents
- 8 Deparaffinization, Decrosslinking, Immunofluorescence Staining & Imaging Reagents
- 9 Imaging Recommendations
- 10 Additional Reagents and Consumables
- 12 Recommended Thermal Cyclers
- 13 Data Analysis Overview

Kits, Reagents & Equipment

The items in the tables below have been validated by 10x Genomics and are highly recommended for the Visium CytAssist assay. **Substituting materials may adversely affect system performance.** This list does not include standard laboratory equipment such as water baths, centrifuges, vortex mixers, pH meters, freezers etc.

Pipette and Pipette Tip Recommendations

Supplier	Description	Part Number (US)
Rainin (pipettes)	<input type="checkbox"/> Pipet-Lite Multi Pipette L8-200XLS+	17013805
	<input type="checkbox"/> Pipet-Lite LTS Pipette L-2XLS+	17014393
	<input type="checkbox"/> Pipet-Lite LTS Pipette L-10XLS+	17014388
	<input type="checkbox"/> Pipet-Lite LTS Pipette L-20XLS+	17014392
	<input type="checkbox"/> Pipet-Lite LTS Pipette L-100XLS+	17014384
	<input type="checkbox"/> Pipet-Lite LTS Pipette L-200XLS+	17014391
	<input type="checkbox"/> Pipet-Lite LTS Pipette L-1000XLS+	17014382
Rainin (tips)	<input type="checkbox"/> Tips LTS 200UL Filter RT-L200FLR	30389240
	<input type="checkbox"/> Tips LTS 1ML Filter RT-L1000FLR	30389213
	<input type="checkbox"/> Tips LTS 20UL Filter RT-L10FLR	30389226
Wide Bore Pipette tips	<input type="checkbox"/> Tips RT LTS 200UL FLW	30389241
	<input type="checkbox"/> Tips RT LTS 1000UL FLW	30389218

Sample Preparation Reagents

For items with multiple options listed, choose option based on availability and preference.

For FFPE Tissue Sectioning & Section Placement			
Item	Alternatives/Options	Vendor	Part Number
<input type="checkbox"/> Microtome	Epredia HM 355S Automatic Microtome <i>Or any standard histology grade microtome</i>	Fisher Scientific	23-900-672
<input type="checkbox"/> Microtome blade	Epredia MX35 Premier Disposable Microtome Blades, Low Profile	Fisher Scientific	3052835
<input type="checkbox"/> Cool-Cut, <i>Optional</i>	Thermo Scientific Cool-Cut, <i>Optional</i>	Fisher Scientific	77-112-0
<input type="checkbox"/> Section transfer system (STS) <i>Optional</i>	Thermo Scientific Section Transfer System (STS), <i>Optional - If using Section Transfer System</i>	Fisher Scientific	771200
<input type="checkbox"/> Probes	Fisherbrand Fine Precision Probe	Fisher Scientific	12-000-153
<input type="checkbox"/> Forceps	Fisherbrand Curved Medium Point General Purpose Forceps	Fisher Scientific	16-100-110
<input type="checkbox"/> Water bath	Premiere Tissue Floating Bath, Lighted <i>Or any equivalent water bath</i>	Fisher Scientific	76278-820
<input type="checkbox"/> Section dryer oven	Epredia High Capacity Section Dryer	Fisher Scientific	A84600051
<input type="checkbox"/> Slide drying rack <i>If using section dryer oven</i>	StatLab TISSUE-TEK VERT 24 SLIDE RACK	Fisher Scientific	LWS2124
<input type="checkbox"/> Thermal cycler <i>Alternative to section dryer oven</i>	C1000 Touch Thermal Cycler with 96-Deep Well Reaction Module	Bio-Rad	1851197
<input type="checkbox"/> Green Marker, <i>Optional if marking area of interest</i>	Sharpie Argyle Green Permanent Marker <i>Or any equivalent hue</i>	Sharpie	1785396
<input type="checkbox"/> Brushes	Camel Hair Brushes, <i>Or any equivalent paintbrush</i>	Ted Pella	11859
<input type="checkbox"/> Blank Slides	Epredia Shandon ColorFrost Plus Slides	Fisher Scientific	6776214
	Fisherbrand Superfrost Plus Microscope Slides	Fisher Scientific	12-550-15
	Poly-Prep Slides	Millipore Sigma	P0425
	VWR Superfrost Plus Micro Slides, Premium	VWR	48311-703
Additional Materials			
<input type="checkbox"/> Razor blades		-	-
<input type="checkbox"/> Ice bucket		-	-
<input type="checkbox"/> Dry ice		-	-
<input type="checkbox"/> Ultrapure/Milli-Q Water for Water Bath, <i>from Milli-Q Integral Ultrapure Water System or equivalent</i>		-	-

This list may not include some standard laboratory equipment.

RNA Quality Assessment (DV200) Reagents

For items with multiple options listed, choose option based on availability and preference.

For RNA Quality Assessment			
<input type="checkbox"/>	Xylene	Xylene, Reagent Grade	Millipore Sigma 214736
		Xylene, Histological Grade	Millipore Sigma 534056
<input type="checkbox"/>	Ethanol	Ethyl Alcohol, Pure (200 Proof, anhydrous)	Millipore Sigma E7023-500ML
		Ethanol absolute ≥99.5%, TechniSolv, pure (Europe Only)	VWR 83813.360DP
<input type="checkbox"/>	RNA extraction	RNeasy FFPE Kit (50)	Qiagen 73504
<input type="checkbox"/>	Deparaffinization solution	Deparaffinization Solution	Qiagen 19093
<input type="checkbox"/>	RNase decontamination solution	RNaseZap RNase Decontamination Solution	Thermo Fisher Scientific AM9780
<input type="checkbox"/>	Nuclease-free water	Nuclease-free Water (not DEPC-Treated)	Thermo Fisher Scientific AM9937
<input type="checkbox"/>	Staining dishes <i>If working with hardset coverslipped slides</i>	Staining Dish, White	VWR 25608-906
<input type="checkbox"/>	Coplin jars <i>If working with hardset coverslipped slides. Alternative to staining dishes.</i>	Glass Coplin Staining Jar, Screw Cap	VWR 2100500-232
<input type="checkbox"/>	Metal Block <i>If working with hardset coverslipped slides</i>	Light Labs 96-WELL ALUMINUM BLOCK	Fisher Scientific NC9820161
<input type="checkbox"/>	1.5-ml centrifuge tubes	DNA LoBind Tubes, 1.5 ml	Eppendorf 022431021
<input type="checkbox"/>	0.2 ml PCR 8-tube strips <i>Choose either Eppendorf, USA Scientific, or Thermo Fisher Scientific PCR 8-tube strips</i>	PCR Tubes 0.2 ml 8-tube strips	Eppendorf 951010022
		TempAssure PCR 8-tube strip	USA Scientific 1402-4700
		MicroAmp 8-TubeStrip, 0.2 ml	Thermo Fisher Scientific N8020580
		MicroAmp 8-cap Strip, clear	Thermo Fisher Scientific N8010535
<input type="checkbox"/>	Nanodrop/Qubit fluorometer	Nanodrop 2000c Spectrophotometers Or any equivalent Nanodrop Alternative to Qubit Fluorometer	Thermo Fisher Scientific ND-2000C
		Qubit RNA BR Assay Kit	Thermo Fisher Scientific Q10210
		Qubit Assay Tubes	Thermo Fisher Scientific Q32856
		Qubit 4 Fluorometer	Thermo Fisher Scientific Q33238
<input type="checkbox"/>	Bioanalyzer/Tapestation <i>Choose Bioanalyzer or TapeStation based on availability & preference.</i>	RNA 6000 Pico Kit	Agilent 5067-1513
		2100 Bioanalyzer Laptop Bundle	Agilent G2953CA
		4200 TapeStation	Agilent G2991AA
		High Sensitivity RNA ScreenTape	Agilent 5067-5579
		High Sensitivity RNA ScreenTape Ladder	Agilent 5067-5581
		High Sensitivity RNA ScreenTape Sample Buffer	Agilent 5067-5580

This list may not include some standard laboratory equipment

Add-on Antibody Optimization Demonstrated Protocol

For items with multiple options listed, choose option based on availability and preference.

Item	Alternatives/Options	Vendor	Part Number
Amp Mix	AllTaq PCR Core Kit	Qiagen	203123
Xylene	Xylene, Reagent Grade	Millipore Sigma	214736
	Xylene, Histological Grade	Millipore Sigma	534056
Ethanol	Ethyl Alcohol, 200 Proof	Millipore Sigma	E7023
	Ethanol absolute ≥99.5%, TechniSolv, pure (Europe Only)	VWR	83813.360DP
Citrate Buffer	Citrate Buffer, pH 6.0, 10x, Antigen Retriever	Millipore Sigma	C9999-1000ML
Tris 1 M (Tris-HCl)	Tris 1 M, pH 8.0, RNase-free	Thermo Fisher Scientific	AM9855G
	Tris 1 M, pH 8.0	Teknova	T5080
	UltraPure 1M Tris-HCl, pH 8.0	Invitrogen	15568025
PBS	PBS - Phosphate Buffer Saline (10X) pH 7.4, RNase-free	Thermo Fisher Scientific	AM9624
10% Tween-20	Tween 20 Surfact-Amps Detergent Solution (10% solution)	Thermo Fisher Scientific	28320
qPCR Mix		KAPA Biosystems (US Only)	
	KAPA SYBR Fast qPCR Master Mix (2X)	Millipore Sigma (Europe, Asia & Canada)	KK4600
Nuclease-free water	Nuclease-free Water (not DEPC-Treated)	Thermo Fisher Scientific	AM9937
DAPI	DAPI Solution (1 mg/ml)	Thermo Fisher Scientific	62248
8M KOH	Potassium Hydroxide Solution, 8M	Millipore Sigma	P4494-50ML
20X SSC Buffer	SSC Buffer 20X Concentration	Millipore Sigma	S66391L
Oligo-conjugated Antibodies	TotalSeq-Bn Antibodies	BioLegend	-
Antibody Oligo Primers*	Forward Primer (5'-3'): GTGACTGGAGTTCAGACGTG Reverse Primer (3'-5'): GCTAGGACCGGCCTTAAAGC	-	-
Secondary Antibodies**	Donkey anti-Mouse IgG (H+L) Highly Cross-Adsorbed Secondary Antibody, Alexa Fluor 568	Thermo Fisher Scientific	A10037
	Donkey anti-Rabbit IgG (H+L) Highly Cross-Adsorbed Secondary Antibody, Alexa Fluor 568	Thermo Fisher Scientific	A10042
Autofluorescence Quencher	TrueBlack Lipofuscin Autofluorescence Quencher <i>Optional</i>	Biotium	23007
Mounting medium	SlowFade Diamond Antifade Mountant	Thermo Fisher Scientific	S36967
Glycerol	Glycerol	Acros Organics	327255000
	Glycerol Solution	Millipore Sigma	49781
Buffer EB	Qiagen Buffer EB	Qiagen	19086
Coplin jar/ staining dishes	Coplin Jar	VWR	100500-232
	Staining Dishes	VWR	25608-906
Section dryer oven	Epredia High Capacity Section Dryer <i>Or any equivalent product. Thermal cyclers may also be used for section drying.</i>	Fisher Scientific	A84600051

* Order primers at 100 micromole scale and choose standard desalting purification.

**Select secondary antibodies based on primary antibody host. Avoid secondary antibodies with fluorescent dyes with emission maxima of 520 ± 20 nm such as FITC, Alexa Fluor 488, and ATTO 488 due to background fluorescence.

Forceps	Excelta 3C-S-SE Tweezers Very Fine Point Straight SS	Excelta	3C-S-SE
1.5 ml tubes	DNA LoBind Tubes, 1.5 ml	Eppendorf	022431021
	Low DNA Binding Tubes, 1.5 ml	Sarstedt	72.706.700
PCR Plates and Sealing Film	Hard-shell PCR Plates 96-well, thin wall (pkg of 50) (or any compatible PCR plate)	Bio-Rad	HSP9665
	Microseal 'B' PCR Plate Sealing Film adhesive (Or any compatible PCR Plate sealing adhesive)	Bio-Rad	MSB1001
Additional Materials			
Razor Blades		-	-
Ultrapure/Milli-Q Water, <i>from Milli-Q Integral Ultrapure Water System or equivalent</i>		-	-

This list may not include some standard laboratory equipment.

Gene and Protein
Expression Assay Only

Deparaffinization, H&E Staining, & Decrosslinking Reagents

For items with multiple options listed, choose option based on availability and preference.

Item	Alternatives/Options	Vendor	Part Number
<input type="checkbox"/> Xylene	Xylene, Reagent Grade	Millipore Sigma	214736
	Xylene, Histological Grade	Millipore Sigma	534056
<input type="checkbox"/> Ethanol	Ethyl Alcohol, Pure (200 Proof, anhydrous)	Millipore Sigma	E7023-500ML
	Ethanol absolute ≥99.5%, TechniSolv, pure (Europe Only)	VWR	83813.360DP
<input type="checkbox"/> Eosin	Eosin Y-solution, Alcoholic	Millipore Sigma	HT110116
	Eosin Y Solution (Modified Alcoholic)	Abcam	ab246824
	Eosin Y with Phloxine 1% alcoholic solution <i>Gene Expression Assay Only</i>	VWR	10143
	Eosin Y Solution, Alcoholic, with Phloxine	Millipore Sigma	HT110332
<input type="checkbox"/> Hematoxylin	Hematoxylin Solution, Mayer's	Millipore Sigma	MHS16
	Hematoxylin Solution According to Mayer	Millipore Sigma	51275
	Hematoxylin, Mayer's	Agilent	S330930-2
<input type="checkbox"/> Bluing reagent	Bluing Reagent, Dako	Agilent	CS70230-2
	Thermo Scientific Shandon Bluing Reagent	Fisher Scientific	6769001
	Scott's Bluing Agent	Ricca Chemical Company	6697
<input type="checkbox"/> Glycerol	Glycerol Solution	Millipore Sigma	49781
	Glycerol	Acros Organics	327255000
<input type="checkbox"/> 0.1 N HCl	Hydrochloric Acid Solution, 0.1 N <i>Or any equivalent HCl</i>	Fisher Chemical	SA54-1
<input type="checkbox"/> Staining jar/dishes	Coplin Jar	VWR	100500-232
	Staining Dishes	VWR	25608-906
<input type="checkbox"/> Section dryer oven	Epredia High Capacity Section Dryer <i>Or any equivalent product. Thermal cyclers may also be used for section drying.</i>	Fisher Scientific	A84-600051
<input type="checkbox"/> Green Marker, <i>Optional, if annotating slide</i>	Sharpie Argyle Green Permanent Marker <i>Or any equivalent hue</i>	Sharpie	1785396
<input type="checkbox"/> Slide holders	Slide Holders, 24-place	VWR	25608-868
<input type="checkbox"/> Coverslips (50 x 24 mm)	Fisherbrand Cover Glasses: Rectangles	Fisher Scientific	12-544-EP
	Cover Glasses, Rectangular	VWR	16004-322
Additional Materials			
<input type="checkbox"/> 1000-ml Beakers (7)		-	-
<input type="checkbox"/> Ultrapure/Milli-Q Water for Water Bath, <i>from Milli-Q Integral Ultrapure Water System or equivalent</i>		-	-

This list may not include some standard laboratory equipment.

Deparaffinization, Decrosslinking, Immunofluorescence Staining & Imaging Reagents

For items with multiple options listed, choose option based on availability and preference.

Item	Alternatives/Options	Vendor	Part Number
<input type="checkbox"/> Xylene	Xylene, Reagent Grade	Millipore Sigma	214736
	Xylene, Histological Grade	Millipore Sigma	534056
<input type="checkbox"/> Glycerol	Glycerol Solution	Millipore Sigma	49781
	Glycerol	Acros Organics	327255000
<input type="checkbox"/> Ethanol	Ethyl Alcohol, Pure	Millipore Sigma	E7023
	Ethanol absolute ≥99.5%, TechniSolv, pure (Europe Only)	VWR	83813.360DP
<input type="checkbox"/> PBS	PBS - Phosphate Buffered Saline (10X) pH 7.4, RNase-free	Thermo Fisher Scientific	AM9624
<input type="checkbox"/> 10% Tween-20	Tween 20 Surfact-Amps Detergent Solution (10% solution)	Thermo Fisher Scientific	28320
<input type="checkbox"/> Nuclease-free water	Nuclease-free Water (not DEPC-Treated)	Thermo Fisher Scientific	AM9937
<input type="checkbox"/> DAPI	DAPI Solution (1 mg/ml)	Thermo Fisher Scientific	62248
<input type="checkbox"/> Antibodies	-	-	-
<input type="checkbox"/> RNase inhibitor	Protector RNase Inhibitor	Millipore Sigma	3335399001
<input type="checkbox"/> Mounting medium	SlowFade Diamond Antifade Mountant	Thermo Fisher Scientific	S36967
	SlowFade Gold Antifade Mountant	Thermo Fisher Scientific	S36966
<input type="checkbox"/> Section dryer oven	Epredia High Capacity Section Dryer	Fisher Scientific	A84600051
<input type="checkbox"/> Slide holders	Slide Holders, 24-place	VWR	25608-868
<input type="checkbox"/> Coverslips (50 x 24 mm)	Fisherbrand Cover Glasses: Rectangles	Fisher Scientific	12-544-EP
	Cover Glasses, Rectangular	VWR	16004-322
<input type="checkbox"/> Staining jar/dishes	Coplin Jar	VWR	100500-232
	Staining Dishes	VWR	25608-906
<input type="checkbox"/> Green Marker, <i>Optional, if annotating slide</i>	Sharpie Argyle Green Permanent Marker <i>Or any equivalent hue</i>	Sharpie	1785396
<input type="checkbox"/> Autofluorescence Quencher	TrueBlack Lipofuscin Autofluorescence Quencher <i>Optional</i>	Biotium	23007
Additional Materials			
<input type="checkbox"/> Beakers		-	-
<input type="checkbox"/> Ultrapure/Milli-Q Water for Water Bath, <i>from Milli-Q Integral Ultrapure Water System or equivalent</i>		-	-

This list may not include some standard laboratory equipment.

Imaging Recommendations

An imaging system with a scanning stage that can accurately capture and stitch together large tiled images is recommended for Visium CytAssist assays. The microscope objective should be dry and free of immersion oil. Low magnification, high Numeric Aperture (NA) microscopes, e.g., 10x/0.4NA or 20x/0.8NA, are recommended to balance high resolution and high acquisition speed. Slides are loosely mounted with glycerol and a coverslip and therefore should not be inverted for imaging. If using an inverted microscope, perform imaging through the slide rather than the coverslip and use an objective (with a correction collar) designed to image through microscope slides.

Supplier	Description	Configuration
Thermo Fisher Scientific	EVOS M7000	Inverted
Leica	Aperio Versa 8 Leica DMI8	Upright Inverted
MetaSystems	Metafer	Upright
Nikon	Nikon Eclipse Ti2	Inverted
BioTek	Cytation 7	Inverted or Upright
Keyence	Keyence BZX800	Inverted
Olympus	VS200	Upright
Zeiss	Imager.Z2	Upright

Brightfield Configuration

Color camera (3 x 8 bit, 2,424 x 2,424 pixel resolution)

White balancing functionality

2.18 µm/pixel minimum capture resolution

Exposure times 2-10 milli sec *May vary depending on imaging system*

Fluorescence Configuration (if performing IF staining)

Light source (or equivalent) with a wavelength range of 380-680 nm

Monochrome camera (14 bit, 2,424 x 2,424 pixel resolution)

DAPI filter cube (Excitation 392/23, Emission 447/60)

Additional filters - select based on fluorophore/antibodies used

2.18 µm/pixel minimum capture resolution

Exposure times 100 milli sec - 2 sec *May vary depending on imaging system*

Additional Reagents and Consumables for the Visium CytAssist Assays

The items in the table below have been validated by 10x Genomics and are highly recommended for both the Visium CytAssist Gene Expression and Visium CytAssist Gene and Protein Expression assays.

Substituting materials may adversely affect system performance. This list does not include standard laboratory equipment such as water baths, centrifuges, vortex mixers, pH meters, freezers etc.

For items with multiple options listed, choose option based on availability and preference.

Item	Description	Vendor	Part Number
Plastics			
<input type="checkbox"/> 1.5 ml tubes	DNA LoBind Tubes, 1.5 ml	Eppendorf	022431021
	Low DNA Binding Tubes, 1.5 ml	Sarstedt	72.706.700
<input type="checkbox"/> 2.0 ml tubes	DNA LoBind Tubes, 2.0 ml when processing more than two slides	Eppendorf	022431048
	Low DNA Binding Tubes, 2.0 ml	Sarstedt	72.695.700
<input type="checkbox"/> 15 ml tubes	15 ml PP Centrifuge Tubes	Corning	730791
<input type="checkbox"/> 50 ml tubes	Self-Standing Polypropylene Centrifuge Tubes (50 ml), sterile	Corning	430921
<input type="checkbox"/> 0.2 ml PCR 8-tube strips	PCR Tubes 0.2 ml 8-tube strips	Eppendorf	951010022
	TempAssure PCR 8-tube strip	USA Scientific	1402-4700
	MicroAmp 8-Tube Strip, 0.2 ml	Thermo Fisher Scientific	N8010580
	MicroAmp 8-Cap Strip, clear	Thermo Fisher Scientific	N8010535
<input type="checkbox"/> Slide mailer/tube	Simport Scientific LockMailer Tamper Evidence Slide Mailer	Thermo Fisher Scientific	22-038-399
	Self-Standing Polypropylene Centrifuge Tubes (50 ml), sterile <i>Alternative to slide mailer</i>	Corning	430921
<input type="checkbox"/> PCR plates and sealing film	Hard-shell PCR Plates 96-well, thin wall (pkg of 50) Or any compatible PCR Plate	Bio-Rad	HSP9665
	Microseal 'B' PCR Plate Sealing Film, adhesive	Bio-Rad	MSB1001
<input type="checkbox"/> Reagent reservoirs	Divided Polystyrene Reservoirs	VWR	41428-958
Kits & Reagents			
<input type="checkbox"/> Nuclease-free water	Nuclease-free Water (not DEPC-Treated)	Thermo Fisher Scientific	AM9937
<input type="checkbox"/> Tris 1 M (Tris-HCl)	Tris 1 M, pH 8.0, RNase-free	Thermo Fisher Scientific	AM9855G
	Tris 1 M, pH 8.0	Teknova	T5080
	UltraPure 1M Tris-HCl, pH 8.0	Invitrogen	15568025
<input type="checkbox"/> 10X PBS	PBS - Phosphate Buffered Saline (10X) pH 7.4, RNase-free	Thermo Fisher Scientific	AM9624
<input type="checkbox"/> Tween 20	Tween 20 Surfact-Amps Detergent Solution (10% solution)	Thermo Fisher Scientific	28320

Additional Reagents and Consumables for Visium CytAssist Assays

For items with multiple options listed, choose option based on availability and preference.

<input type="checkbox"/>	qPCR Mix	KAPA SYBR FAST qPCR Master Mix (2X)	KAPA Biosystems (US Only) Millipore Sigma (Europe, Asia, & Canada)	KK4600
<input type="checkbox"/>	SPRIselect reagent	SPRIselect Reagent Kit	Beckman Coulter	B23318
<input type="checkbox"/>	Ethanol	Ethanol, Pure (200 Proof, anhydrous)	Millipore Sigma	E7023-500ML
<input type="checkbox"/>	8M KOH	Potassium Hydroxide Solution, 8M	Millipore Sigma	P4494-50ML
<input type="checkbox"/>	20X SSC Buffer	SSC Buffer 20X Concentrate	Millipore Sigma	S66391L
<input type="checkbox"/>	Buffer EB	Qiagen Buffer EB	Qiagen	19086
<input type="checkbox"/>	Eosin	Eosin Y-solution, Alcoholic	Millipore Sigma	HT110116
		Eosin Y Solution (Modified Alcoholic)	Abcam	ab246824
		Eosin Y with Phloxine 1% alcoholic solution <i>Gene Expression Assay Only</i>	VWR	10143
		Eosin Y Solution, Alcoholic, with Phloxine	Millipore Sigma	HT110332
<input type="checkbox"/>	Ultrapure Water	Ultrapure/Milli-Q water		
<input type="checkbox"/>	Mini Centrifuge	VWR Mini Centrifuge Or any equivalent mini centrifuge	VWR	76269-064
Quantification & Quality Control				
<input type="checkbox"/>	Bioanalyzer & associated reagents	2100 Bioanalyzer Laptop Bundle (discontinued)	Agilent	G2943CA
		Replacement 2100 Bioanalyzer Instrument/ 2100 Expert Laptop Bundle	Agilent	G2939BA/ G2953CA
		High Sensitivity DNA Kit	Agilent	5067-4626
<input type="checkbox"/>	TapeStation & associated reagents	4200 TapeStation	Agilent	G2991AA
		High Sensitivity D1000: ScreenTape/ Reagents	Agilent	5067-5592/ 5067-5593
		High Sensitivity D5000: ScreenTape/ Reagents	Agilent	5067-5584/ 5067-5585
<input type="checkbox"/>	LabChip & associated reagents	LabChip GX Touch HT Nucleic Acid Analyzer	PerkinElmer	CLS137031
		DNA High Sensitivity Reagent Kit	PerkinElmer	CLS760672
<input type="checkbox"/>	Qubit & associated reagents (for determining dilution factor for TapeStation, Bioanalyzer, or LabChip)	Qubit Fluorometer	Thermo Fisher Scientific	Q33238
		Qubit Assay Tubes Reagents	Thermo Fisher Scientific	Q32856
		Qubit dsDNA HS and BR Assay Kits	Thermo Fisher Scientific	Q32854
<input type="checkbox"/>	Library Quantification Kit	KAPA Library Quantification Kit for Illumina Platforms	KAPA Biosystems	KK4824

Recommended Thermal Cyclers

Supplier	Description	Part Number
Bio-Rad	C1000 Touch Thermal Cycler with 96-Deep Well Reaction Module (discontinued)	1851197
	PTC Tempo Deepwell Thermal Cycler*	12015392
Eppendorf	MasterCycler Pro (discontinued)	North America 950030010 International 6321 000.019
	MasterCycler X50S	North America 6311000010
Thermo Fisher Scientific	Veriti 96-Well Thermal Cycler (discontinued)	4375786
Analytik Jena	Biometra TAdvanced 96 SG with 96-well block (silver, 0.2 mL) and gradient function	846-x-070-241

*If used with a single cassette, place the cassette on the side of the adapter furthest away from the thermal cycler hinge to prevent lid lock errors.

! If using thermal cyclers other than the BioRad C1000, ramp rates should be adjusted for all protocol steps as described below:


- Eppendorf MasterCycler X50s: 3°C/sec heating and 2°C/sec heating cooling
- Analytik Jena Biometra TAdvanced: 2°C/sec heating and cooling

Thermal cycler requirements if reactions are performed on a slide:

- Thermal cycler must be able to accommodate the Low Profile Thermocycler Adapter:
- Well depth: 4.5 mm
- Distance between block and heated lid: 12 mm
- Reaction block dimensions 95.5 x 73 mm

Thermal cycler requirements if reactions are performed in a tube (Elution, Step 3.0 onward):

- Uniform heating of 100 µl volumes
- Temperature-controlled lid
- 96 deep-well block or 0.2 ml block configuration




Recommended Real Time qPCR Systems

Supplier	Description	Part Number
Applied Biosystems	QuantStudio 12K Flex System	4471087
Bio-Rad	CFX96 Real-time System (discontinued)	1855096


Analysis of Visium CytAssist Data

For detailed information on analysis of Visium CytAssist data, refer to 10x Genomics [Software Support](#) site.

Data Analysis Overview – Visium CytAssist Spatial Gene Expression



Data Analysis Overview – Visium CytAssist Spatial Gene and Protein Expression



Key information that should be collected and provided to the data analyst/bioinformatician to ensure that parameters are set appropriately for Space Ranger to analyze the data:

- Sample name, library type, and either Illumina adaptor index set or index sequences
- If a microscope image of the tissue section on plain glass is provided, then the corresponding imaging type (Brightfield vs Fluorescence (color composite) vs Fluorescence (multipage TIFF) for the supplied image
- Visium slide serial number and Capture Area on the slide corresponding to each sample
- The CytAssist instrument captured Brightfield image with fiducial frame

Document Revision Summary

Document Number	CG000556
Title	Visium CytAssist for FFPE - Protocol Planner
Revision	Rev F
Revision Date	August 2023
General Changes	Updated for general minor consistency of language and terms throughout
Specific Changes	<ul style="list-style-type: none">Updated to include reagents from the Visium CytAssist Gene and Protein Expression Add-on Antibody Optimization protocol (CG000664).Added PTC Tempo Deepwell to list of recommended thermal cyclers.

© 2023 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: www.10xgenomics.com/trademarks. 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: www.10xgenomics.com/patents. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of 10x Genomics products in practicing the methods set forth herein has not been validated by 10x Genomics, and such non-validated use is NOT COVERED BY 10X GENOMICS STANDARD WARRANTY, AND 10X GENOMICS HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics terms and conditions of sale for the Chromium Controller or the Chromium Single Cell Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics that it currently or will at any time in the future offer or in any way support any application set forth herein.

Contact:**support@10xgenomics.com**

10x Genomics

6230 Stoneridge Mall Road

Pleasanton, CA 94588 USA

