

Cell Surface Protein Labeling for Chromium Fixed RNA Profiling

Getting Started

Cell surface proteins can be labeled using a specific protein binding molecule, such as an antibody conjugated to a Feature Barcode oligonucleotide. This protocol outlines cell surface protein labeling for use with Chromium Fixed RNA Profiling with Feature Barcode technology and is compatible with both TotalSeq™-C (Singleplex and Multiplex; barcode oligo capture) and TotalSeq™-B (Singleplex) antibody conjugates. It also provides guidance for enriching labeled cells using flow sorting.

Additional Guidance

Pre-read and have available Fixation of Cells & Nuclei for Chromium Fixed RNA Profiling Demonstrated Protocol (CG000478) before starting the cell surface protein labeling protocol.

Consult Cell Preparation Guide Handbook (CG00053) for Tips & Best Practices on handling cells and Guidelines on Accurate Target Cell Counts Technical Note (CG000091) for determining accurate cell counts.

Consult Thawing Dissociated Tumor Cells for Single RNA Sequencing Demonstrated Protocol (CG000233) for guidance on thawing dissociated tumor cells.

Consult Chromium Fixed RNA Profiling - Protocol Planner (CG000528) for details on workflow overview, document resources, and guidance on selecting the appropriate sample preparation and library construction protocols for different Chromium Fixed RNA workflows.

Cells carry potentially hazardous pathogens. Follow material supplier recommendations and local laboratory procedures and regulations for the safe handling, storage and disposal of biological materials.

Specific Reagents & Consumables

Vendor	Item	Part Number
For Cell Surface Protein Labeling		
BioLegend	Human TruStain FcX (Fc Receptor Blocking Solution)	422301
	TruStain FcX PLUS (anti-mouse CD16/32) Antibody	156604
	True-Stain Monocyte Blocker	426101
	TotalSeq™ Antibody-Oligonucleotide Conjugates*	-
	Cell Staining Buffer	420201
Thermo Fisher Scientific	UltraPure Bovine Serum Albumin (BSA, 50 mg/ml)	AM2616
	Fetal Bovine Serum, qualified, heat inactivated	16140071

Continued onto next page.

Millipore Sigma	Bovine Serum Albumin In DPBS (10%) <i>Alternative to Thermo Fisher product</i>	A1595
Miltenyi Biotec	MACS BSA Stock Solution <i>Alternative to Thermo Fisher product</i>	130-091-376
Corning	Phosphate-Buffered Saline, 1X without Calcium and Magnesium	21-040-CV
VWR	Fetal Bovine Serum (FBS) <i>Alternative to Thermo Fisher product</i>	97068-085

For Cell Counting

Thermo Fisher Scientific	Countess II FL Automated Cell Counter	AMQAF1000
	Countess Cell Counting Chamber Slides	C10228
	Trypan Blue Stain (0.4%)	T10282
Nexcelom	Celleca MX High-throughput Automated Cell Counter	MX-112-0127
	ViaStain AOPI Staining Solution	CS2-0106-5mL

This list may not include some standard laboratory equipment.

**Choose appropriate TotalSeq™ antibody-oligonucleotide conjugates based on the Chromium Fixed RNA Profiling workflow.*

Chromium Fixed RNA Profiling Workflows	TotalSeq™ Antibody-Oligonucleotide Conjugates
Gene & Protein Expression using Barcode Oligo Capture – Singleplex & Multiplex Workflows	TotalSeq™-C antibody-oligonucleotide conjugates
Gene & Protein Expression – Singleplex Workflow	TotalSeq™-B antibody-oligonucleotide conjugates

Tips & Best Practices

Cell Viability

- Determine sample viability before starting the cell surface protein labeling protocol.

Labeling & Wash Buffer

- Labeling or washing buffer for most sample types is PBS + 1% BSA.
- For samples containing <70% viable cells, PBS + 10% FBS can be used.
- BioLegend's Cell Staining Buffer can also be used for labeling cells. However, this buffer may not be optimal for all sample types. Cell Staining Buffer should only be used for the labeling step. PBS + 1% BSA should be used for the washing steps.
- If using lyophilized antibody panel/cocktails, please follow BioLegend's instructions for use regarding reconstitution and labeling volumes.

- The addition of True-Stain Monocyte Blocker to Human TruStain FcX or TruStain FcX PLUS (anti-mouse CD16/32) Antibody during labeling can generally help to reduce staining background in most sample types.

Centrifugation Conditions

- Centrifugation speed and time depend upon the sample type.
- Larger or fragile cell types may require slower centrifugation speeds. See Table 1 for more information about centrifugation conditions for different sample types.

Table 1. Sample type specific centrifugation conditions

Sample Type	Centrifugation Conditions
Samples containing >70% viable cells, e.g., PBMCs	400 rcf for 5 min
Samples containing <70% viable cells, e.g., tumor cells	150 rcf for 10 min

Optimal Antibody Concentration

- Consult Quality Control of Cell Surface Protein Labeling using Flow Cytometry Technical Note (CG000231) for guidance on how to titrate antibodies using flow cytometry for determining optimal antibody concentrations.
- It is highly recommended to titrate antibodies before use for optimal performance. Refer to the Manufacturer's Instructions for recommended dilutions for each antibody. A dilution of 0.5 µg per antibody for up to 2×10^6 cells is suggested as a starting point.

Sample Washing

- This protocol provides three different wash options after incubation with the antibody-oligonucleotide conjugate. See the following page for guidance on choosing the appropriate wash option.

Sample Fixation

- Consult Fixation of Cells & Nuclei for Chromium Fixed RNA Profiling Demonstrated Protocol (CG000478) for guidance on fixing single cell or nuclei suspensions following cell surface protein labeling.

Sample Batching

- For convenience, samples may be processed in smaller batches and stored long-term. If batching samples together, it is recommended that samples with different fixation times in one experiment are not mixed.
- After batching, samples can be stored at -80°C for up to 6 months with appropriate storage reagents.
- For storage recommendations that address fixation time and temperature, along with reagents used for storage, consult Fixation of Cells & Nuclei for Chromium Fixed RNA Profiling Demonstrated Protocol (CG000478), Fixation Conditions & Fixed Sample Storage and Appendix sections.

Antibody-Oligonucleotide Conjugation Guidance

- Consult Cell Surface Protein Labeling for Single Cell RNA Sequencing Protocols with Feature Barcode technology Demonstrated Protocol (CG000149) for guidance on antibody-oligonucleotide conjugation.

Table 2. TotalSeq™ Antibody Oligonucleotide Conjugates for Antibody Conjugation

Chromium RNA Fixed RNA Profiling Workflows	TotalSeq™ Antibody-Oligonucleotide Conjugates			
Gene & Protein Expression using Barcode Oligo Capture - Singleplex & Multiplex Workflows	/5AmMC12/CGGAGATGTGTATAAGAGACAG-NNNNNNNNNN-NNNNNNNNNNNNNNNNNNNNNN-CCCATATAAGAAA Nextera partial Read 2 10 nt Feature Barcode (15 nt) 9 nt Capture Sequence			
	TotalSeq™-C antibody-oligonucleotide conjugates			
Gene & Protein Expression - Singleplex Workflow	/5AmMC12/GTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT-NNNNNNNNNN-NNNNNNNNNNNNNNNNNNNNNN-GCTTTAAGGCCGGTCTAGCAA TruSeq Read 2 10 nt Feature Barcode (15 nt) 9 nt Capture Sequence 1			
	TotalSeq™-B antibody-oligonucleotide conjugates			

See Appendix for an illustrative overview of antibody-oligonucleotide conjugate capture by 10x Gel Bead primers.

Protocol Overview


Below is a general illustrative overview of the cell surface protein labeling protocol and wash options. Using the appropriate wash option is critical for obtaining high quality data and minimizing background signal. Choose a washing protocol based on the guidance outlined in the following diagram. Refer to the Appendix for supplemental data on the different wash options.

Protocol & Sample Washing Overview

Protocol Overview



Sample Washing Overview



* Refer to Fixation of Cells & Nuclei for Chromium Fixed RNA Profiling Demonstrated Protocol (CG000478).

Cell Surface Protein Labeling Protocol

This protocol was optimized using TotalSeq™-C and TotalSeq™-B antibody-oligonucleotide conjugates from BioLegend. The labeled cells can be enriched by FACS (see Appendix).



Use distinct and compatible antibody clones for FACS and cell surface protein labeling. Optimize working concentration of each antibody used.

Buffers – Preparation

For Labeling Cells

- Chilled (4°C): PBS + 1% BSA
- For samples containing <70% viable cells, chilled (4°C) PBS + 10% FBS can be used.

Prepare Antibody Mix Supernatant

Add appropriate/manufacture's recommended amount of TotalSeq™-C or TotalSeq™-B antibody-oligonucleotide conjugates to a 1.5-ml microcentrifuge tube.


- If using lyophilized antibody panel/cocktails, rehydrate lyophilized panel in the recommended volume of Cell Staining Buffer as directed from BioLegend. Follow manufacturer's instructions for use for cell labeling. Perform cell wash steps as described in this Demonstrated Protocol.
- Centrifuge the mix at **14,000 rcf** for **10 min** at **4°C**.
- Transfer the supernatant (containing Antibody Mix) to a new tube and maintain at **4°C**.

Prepare FACS Antibody Pool

- Add appropriate/manufacture's recommended amount of fluorophore antibodies to a 1.5-ml microcentrifuge tube on ice.
- Gently pipette mix and maintain at **4°C**. Avoid light exposure.

Labeling

When using custom conjugated antibodies, be sure to follow manufacturer's instructions. Refer to Cell Surface Protein Labeling for Single Cell RNA Sequencing Protocols Demonstrated Protocol (CG000149) for details. For pre-conjugated antibodies, we recommend BioLegend TotalSeq™-C or TotalSeq™-B.



Protocol Steps:

- a.** Resuspend cells in PBS + 1% BSA.

TIPS For samples containing <70% viable cells, chilled (4°C) PBS + 10% FBS can be used.

- b.** Transfer $\leq 2 \times 10^6$ cells to a new 1.5-ml microcentrifuge tube.
- c.** Centrifuge cells at **400 rcf** for **5 min** (PBMCs) or at **150 rcf** for **10 min** (dissociated tumor cells) at **4°C**. Use of a swinging-bucket rotor is recommended for higher cell recovery. Centrifugation speed and time depend upon the sample type. Larger or fragile cell types may require slower centrifugation speeds.

- d.** Remove supernatant without disturbing pellet.

- e.** Resuspend pellet in **50 µl** chilled PBS + 1% BSA. For samples containing <70% viable cells, chilled PBS + 10% FBS can be used.

- f.** Add **5 µl / 10^6 cells** Human TruStain FcX or **0.5–2 µl / 10^6 cells** TrueStain FcX PLUS (anti-mouse CD16/32 Antibody). Refer to Manufacturer's Instructions for resuspension recommendations. Gently pipette mix.

OPTIONAL: **5 µl** True-Stain Monocyte Blocker can also be added at this step in addition to TruStain FcX. See [Tips & Best Practices](#) for more information.

- g.** Incubate for **10 min** at **4°C**.
- h.** Add prepared Antibody Mix Supernatant. If performing FACS, add FACS Antibody Pool.
- i.** Add chilled PBS + 1% BSA to the cells to bring the total volume to **100 µl**. Gently pipette mix 10x (pipette set to 90 µl). For samples containing <70% viable cells, chilled PBS + 10% FBS can be used.
- j.** Incubate for **30 min** at **4°C**. If using FACS antibodies, incubate without light exposure.

- k.** Proceed to appropriate **Washing** section for **1-Wash or 2-Wash Options**


OR

Proceed directly to **Sample Fixation, Step c** of the Fixation of Cells & Nuclei for Chromium Fixed RNA Profiling Demonstrated Protocol (CG000478) (**No-Wash Option**).

TIPS

Sample fixation without washing is recommended for samples with low inputs, low viability, or time-sensitive samples where increased background is acceptable.

Washing: No Wash Option




Protocol Steps:

- a.** Proceed **immediately** to **Sample Fixation, Step c** of the Fixation of Cells & Nuclei for Chromium Fixed RNA Profiling Demonstrated Protocol (CG000478).

Step c: Add **1 ml** Fixation Buffer to the labeled cells and pipette mix 5x. Fixation Buffer preparation and fixation protocol are listed in Fixation of Cells & Nuclei for Chromium Fixed RNA Profiling Demonstrated Protocol (CG000478).

Washing: 1-Wash Option



Protocol Steps:


a. OPTIONAL: For enrichment of labeled and viable cells by FACS:

- Based on starting concentration and assuming ~50% cell loss, add an appropriate volume chilled PBS + 2% FBS (including a dead cell marker) to obtain a final cell concentration of $5\text{--}10 \times 10^6$ cells/ml.
- Proceed to **FACS (see FACS Guidance)**. After FACS, determine cell concentration and viability using an automated cell counter or a hemocytometer.

If not performing FACS enrichment:

- Add **1.4 ml** chilled PBS + 1% BSA to the labeled cells. Gently pipette mix. For samples containing <70% viable cells, chilled PBS + 10% FBS can be used.
- After resuspension, take an aliquot and determine cell concentration and viability using an automated cell counter or a hemocytometer.
- Proceed **immediately to Sample Fixation, Step a** of the Fixation of Cells & Nuclei for Chromium Fixed RNA Profiling Demonstrated Protocol (CG000478).

Washing: 2-Wash Option



Protocol Steps:

- Add **1.4 ml** chilled PBS + 1% BSA to the labeled cells. Gently pipette mix. For samples containing <70% viable cells, chilled PBS + 10% FBS can be used.
- Centrifuge cells at **400 rcf** for **5 min** (PBMCs) or at **150 rcf** for **10 min** (dissociated tumor cells) at **4°C**. Centrifugation speed and time depends upon the sample type.
- Remove the supernatant without touching the bottom of the tube to avoid dislodging the pellet.
- Resuspend the cell pellet in **1.5 ml** chilled PBS + 1% BSA and place on ice. For samples containing <70% viable cells, chilled PBS + 10% FBS can be used.
- After resuspension, take an aliquot and determine cell concentration and viability using an automated cell counter or a hemocytometer.
- Proceed **immediately** to **Sample Fixation, Step a** of the Fixation of Cells & Nuclei for Chromium Fixed RNA Profiling Demonstrated Protocol (CG000478).



Pellet
(may not
be visible)



Leaving behind excess supernatant (>30 μ l) may cause non-specific binding, which may result in increased background reads during sequencing.

Appendix

FACS Guidance

Enrich labeled cells using flow sorting prior to library generation to enable identification of rare subpopulations.


FACS Cell Collection

It is recommended to collect FACS-enriched cells in up to 20% FBS to maintain cell viability. Cells should be collected either in 20 µl volume in the collection tube/plate (96-well plate) or in 300 µl volume in a 5-ml tube. Use an optimal buffer for fragile cells to maintain a high cell viability.

Sort stream should be adjusted so that the cell-droplet falls into the collection buffer. Sorted cells should be counted and viability measured before

proceeding to the Fixation of Cells & Nuclei for Chromium Fixed RNA Profiling Demonstrated Protocol (CG000478).

Cell loss during flow sorting is common. Optimize the protocol steps accordingly.

 Once sorting is complete, proceed **immediately** to **Sample Fixation, Step a** of the Fixation of Cells & Nuclei for Chromium Fixed RNA Profiling Demonstrated Protocol (CG000478).

Post-Fixation FACS

Post-fixation samples can be sorted using FACS for advanced sample clean-up, as well as enrichment of specific populations in the Fixed RNA Profiling Assay. Refer to the 10x Genomics Support Website for more information.

Antibody-Oligonucleotide Conjugate Capture


Antibody-oligonucleotide conjugate capture by protocol-specific Gel Bead primers is illustrated below.

Illustrative Overview of Antibody-Oligonucleotide Conjugate Capture

TotalSeq™-C Antibody-oligonucleotide conjugate (Singleplex and Multiplex workflow; barcode oligo capture)

Refer to the following User Guides for more information:

- Chromium Fixed RNA Profiling Reagent Kits for Multiplexed Samples with Feature Barcode technology for Protein using Barcode Oligo Capture User Guide (CG000673).
- Chromium Fixed RNA Profiling Reagent Kits for Singleplexed Samples with Feature Barcode technology for Protein using Barcode Oligo Capture User Guide (CG000674).



TotalSeq™-B Antibody-oligonucleotide conjugate (Singleplex workflow)

Refer to the following User Guide for more information:

- Chromium Fixed RNA Profiling Reagent Kits for Singleplexed Samples with Feature Barcode technology for Protein User Guide (CG000477).




Figure 1. A. Illustrative overview of antibody-oligonucleotide conjugate capture for TotalSeq™-C (Singleplex and Multiplex; barcode oligo capture). **B.** Illustrative overview of antibody-oligonucleotide conjugate capture for TotalSeq™-B (Singleplex).


Appendix, *contd.*

Figure 2. t-SNE plots of human PBMCs comparing Feature Barcode staining for T cell marker CD3 across differing number of washes following antibody incubation (left), showing equivalent data quality between TotalSeq™-C (**A**) and TotalSeq™-B (**B**) configuration of the same antibody clone. Corresponding CD3E gene expression measured with Fixed RNA Profiling (FRP). Scale adjusted for better visualization with top 1% of cells removed. Violin plot showing the expression of the anti-CD3 antibody Feature Barcode (right).


Appendix, *contd.*

Figure 3. t-SNE plots of human PBMCs comparing Feature Barcode staining for Natural Killer and T cell marker TIGIT across differing number of washes following antibody incubation (left), showing equivalent data quality between TotalSeq™-C (**A**) and TotalSeq™-B (**B**) configuration of the same antibody clone. Corresponding TIGIT gene expression measured with Fixed RNA Profiling (FRP). Scale adjusted for better visualization with top 1% of cells removed for each feature. Violin plot showing the expression of the anti-TIGIT antibody Feature Barcode (right).

Appendix, *contd.*

Antibody	Comparison	TotalSeq™-C			TotalSeq™-B		
		Signal to Noise* (On-target cells / all other cells)			Signal to Noise* (On-target cells / all other cells)		
		No-Wash	1-Wash	2-Wash	No-Wash	1-Wash	2-Wash
anti-CD3	T cells / Other	12.44	15.72	16.99	13.09	14.89	17.58
anti-CD4	T cells / Other	9.92	10.71	9.36	7.61	7.17	8.05
anti-CD8	T cells / Other	3.2	5.81	6.32	3.11	4.94	5.93
anti-TIGIT	T & NK cells / Other	1.14	1.4	1.54	0.85	0.89	1.1
anti-CD19	B cells / Other	27.11	39.67	45.61	16.94	23.26	35.48
anti-CD15	Monocytes / Other	2.9	4.51	4.61	1.94	2.24	2.88
anti-CD14	Monocytes / Other	6.62	7.4	8.26	6.45	7.53	9.77
anti-CD11c	DCs & Monocytes / Other	13.5	12.89	13.07	14.07	12.1	14.81

*Ratio of mean expression level

Table 1. Signal to noise was generated by calculating the ratio of counts from a target group comprised of cell type(s) expected to display antibody signal (e.g. CD3 on T cells) to counts from a background group (e.g. CD3 on non-T cells).

Conclusion

Equivalent antibody Feature Barcoding data quality can be obtained regardless of format, either with TotalSeq™-C or TotalSeq™-B antibody oligonucleotide conjugates. However, only TotalSeq™-C enables sample multiplexing. Choosing an appropriate wash protocol following cell surface protein labeling is critical for experimental success in the Chromium Fixed RNA Profiling assay. The data presented in this Demonstrated Protocol show that a lower number of washes following cell labeling can be used for antibodies with distinct positive and negative populations (Figures 2A and 2B, CD3 data). However, lower number of washes cause a reduction in separation between the positive and negative populations (Figures 2A and 2B, violin plot; Table 1, lower signal to noise for CD3, CD14, CD19, etc). Feature Barcode data from antibodies with poorer separation between the positive and negative populations is adequate with the 1-Wash protocol, but the No-Wash protocol is not typically recommended due to poor signal to noise (Figures 3A and 3B, TIGIT; Table 1, TIGIT).

References

1. Chromium Fixed RNA Profiling Reagent Kits for Singleplexed Samples with Feature Barcode technology for Cell Surface Protein User Guide (CG000477)
2. Chromium Fixed RNA Profiling Reagent Kits for Multiplexed Samples User Guide (CG000527)
3. Fixation of Cells & Nuclei for Chromium Fixed RNA Profiling Demonstrated Protocol (CG000478)
4. Cell Surface Protein Labeling for Single Cell RNA Sequencing Protocols with Feature Barcode technology Demonstrated Protocol (CG000149)
5. Chromium Fixed RNA Profiling - Protocol Planner (CG000528)
6. Chromium Fixed RNA Profiling Reagent Kits for Multiplexed Samples with Feature Barcode technology for Protein using Barcode Oligo Capture User Guide (CG000673)
7. Chromium Fixed RNA Profiling Reagent Kits for Singleplexed Samples with Feature Barcode technology for Protein using Barcode Oligo Capture User Guide (CG000674)

Document Revision Summary

Document Number	CG000529
Title	Cell Surface Protein Labeling for Chromium Fixed RNA Profiling Demonstrated Protocol
Revision	Rev B
Revision Date	August 2023
Specific Changes	<p>Updated introduction for TotalSeq™-B and TotalSeq™-C antibody conjugate compatibility (page 1).</p> <p>Added TotalSeq™-C Oligonucleotide Conjugates to Specific Reagents & Consumables (page 1).</p> <p>Added TruStain FcX (anti-mouse CD16/32) Antibody to Specific Reagents & Consumables (page 1).</p> <p>Added Countess II FL Automated Cell Counter and Counting Chamber Slides to Specific Reagents & Consumables (page 2).</p> <p>Updated Tips & Best Practices with reference to custom antibody conjugation information (page 3).</p> <p>Added TruStain FcX (anti-mouse CD16/32) Antibody volume (0.5 - 2µl / 10⁶) for mouse cells to Labeling protocol (page 6).</p> <p>Changed Labeling suspension to PBS + 1% BSA (page 6).</p> <p>Updated Antibody-Oligonucleotide Conjugate Capture illustration in Appendix to include TotalSeq™-C antibody-oligo conjugate (page 10).</p> <p>Added figures to Appendix for TotalSeq™-C data (pages 11-13).</p>
General Changes	Updated for general minor consistency of language and terms throughout.

© 2023 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: www.10xgenomics.com/trademarks. 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: www.10xgenomics.com/patents. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of 10x Genomics products in practicing the methods set forth herein has not been validated by 10x Genomics, and such non-validated use is NOT COVERED BY 10X GENOMICS STANDARD WARRANTY, AND 10X GENOMICS HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics terms and conditions of sale for the Chromium Controller or the Chromium Single Cell Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics that it currently or will at any time in the future offer or in any way support any application set forth herein.

Contact:**support@10xgenomics.com**

10x Genomics
6230 Stoneridge Mall Road
Pleasanton, CA 94588 USA

