DEMONSTRATED PROTOCOL CG000410 | Rev E

Visium Spatial Gene Expression for FFPE –

Deparaffinization, Decrosslinking, Immunofluorescence Staining & Imaging

Introduction

The Visium Spatial Gene Expression for FFPE is designed to measure mRNA in tissue sections derived from formalin fixed & paraffin embedded (FFPE) tissue samples and requires a Visium Spatial slide with intact tissue sections as input. Immunostaining tissue sections with fluorescently labeled antibodies enables simultaneous protein detection. This protocol outlines deparaffinization, decrosslinking, immunofluorescence (IF) staining, and imaging of tissue for use with 10x Genomics Visium Spatial Gene Expression for FFPE assay. Deparaffinized, decrosslinked, and stained tissue sections are inputs for the downstream Visium Spatial Gene Expression for FFPE workflow.

Additional Guidance

Consult the Visium Spatial Gene Expression for FFPE - Tissue Preparation Guide (Document CG000408) for complete information on sectioning FFPE tissue blocks and placing sections on Visium Spatial slides. Ensure that tissue sections have been placed onto the appropriate slide before starting this Demonstrated Protocol. Consult the Visium Spatial Gene Expression for FFPE Imaging Guidelines (Document CG000436) to verify imaging settings before starting this Demonstrated Protocol. After completing this Demonstrated Protocol (CG000410), proceed with the Visium Spatial Gene Expression for FFPE - User Guide (CG000407).

Contents

- 2 Reagent Kits
- 3 Recommended Thermal Cyclers
- 3 Workflow Overview
- 4 Specific Reagents & Consumables
- 5 Tips & Best Practices
- 16 Protocol Steps & Timing
- 17 1. Deparaffinization, Decrosslinking & IF Staining
- 17 1.1 Preparation Buffers
- 18 Protocol Overview
- 21 1.2 Deparaffinization
- 23 1.3 Decrosslinking

- 24 1.4 Immunofluorescence Staining
- 30 2. Tissue Imaging
- 30 2.1 Imaging System Recommendations
- 31 2.2 Imaging
- 31 2.3 Coverslip Removal
- 32 Troubleshooting
- 32 Results



Reagent Kits

Visium Spatial Gene Expression for FFPE Reagent Kits

Refer to SDS for handling and disposal information

Visium Spatial Gene Expression Slide Kit, 16 rxns PN-1000185

Visium Spatial Gene Expression Slide Kit 16 rxns, PN-1000185			
store at ambient temperature or at -20°C accordi	ing to k #	rit label PN	
Visium Spatial Gene Expression Slide	4	2000233	
*Visium Slide Seals, 40-pack	1	2000284	
Visium Cassette & Gasket Assembly, 4-pack	1	2000282	
			10x

Visium Spatial Gene Expression Slide Kit, 4 rxns PN-1000188

Visium Spatial Gene Expression Slide Kit 4 rxns, PN-1000188			
store at ambient temperature or at -20°C acco	rding to k #	it label PN	
Visium Spatial Gene Expression Slide	1	2000233	
*Visium Slide Seals, 12-pack	1	2000283	
Visium Cassette & Gasket Assembly, 1-pag	ck 1	2000281	
			10x

^{*}Visium Slide Seals may come in varying dimensions and quantities in different lots. Check the number of slide seals in the kit. Additional seals may be required. Refer to page 13 (Visium Slide Seal Application & Removal) of this Demonstrated Protocol for instructions on how to resize seals or cut additional seals...


10x Genomics Accessories

Product	Part Number (Kit)	Part Number (Item)
Thermocycler Adaptor		3000380
Visium Spatial Imaging Test Slide	1000194	2000235
10x Magnetic Separator	1000194	230003
Slide Alignment Tool		3000433

Recommended Thermal Cyclers

Supplier	Description	Part Number
Bio-Rad	C1000 Touch Thermal Cycler with 96-Deep Well Reaction Module	1851197
Eppendorf	MasterCycler Pro (discontinued)	North America 950030010 International 6321 000.019
Thermo Fisher Scientific	Veriti 96-Well Thermal Cycler	4375786

Workflow Overview



Visit the 10x Genomics Support website for the most current documentation.

Specific Reagents & Consumables

For each item, a number of vendor options are listed. Choose item based on availability and preference.

Item	Alternatives/Options	Vendor	Part Number
Xylene	Xylene, Reagent Grade	Millipore Sigma	214736
	Xylene, Histological Grade	Millipore Sigma	534056
Ethanol	Ethyl Alcohol, 200 Proof	Millipore Sigma	E7023
	Ethanol absolute ≥99.5%, TechniSolv, pure (Europe Only)	VWR	83813.360DP
TE Buffer (pH 9.0)			10-0046
Or Prepare using Tris and 0.5 M	Tris Base For preparing TE buffer (pH 9.0), alternative to Genemed product	Fisher Scientific	BP152-500
EDTA	UltraPure 0.5 M EDTA, pH 8.0 For preparing TE buffer (pH 9.0), alternative to Genemed product	Thermo Fisher Scientific	15575020
	Dilute to 1.0 M HCl stock For preparing TE Buffer (pH 9.0), alternative to Genemed product	Millipore Sigma	258148
PBS	PBS - Phosphate Buffered Saline (10X) pH 7.4, RNase-free	Thermo Fisher Scientific	AM9624
BSA	Albumin, Bovine Serum, 10% Aqueous Solution, Nuclease- Free	Millipore Sigma	126615-25ML
10% Tween-20	Tween 20 Surfact-Amps Detergent Solution (10% solution)	Thermo Fisher Scientific	28320
Nuclease-free water	Nuclease-free Water (not DEPC-Treated)	Thermo Fisher Scientific	AM9937
DAPI	DAPI Solution (1 mg/ml)	Thermo Fisher Scientific	62248
Antibody	-	-	-
RNase inhibitor	Protector RNase Inhibitor	Millipore Sigma	3335399001
Mounting medium	SlowFade Diamond Antifade Mountant	Thermo Fisher Scientific	S36967
Section dryer oven	Epredia High Capacity Section Dryer	Fisher Scientific	A84600051
Slide holders	Slide Holders, 24-place	VWR	25608-868
Coverslips	Fisherbrand Cover Glasses: Rectangles	Fisher Scientific	12-544-EP
	Cover Glasses, Rectangular	VWR	16004-322
Staining jar/	Coplin Jar	VWR	100500-232
dishes	Staining Dishes	VWR	25608-906
Sealing film	Microseal 'B' PCR Plate Sealing Film, adhesive	Bio-Rad	MSB1001
Additional Mat	erials		
Beakers		-	-
Ultrapure/Milli- from Milli-Q Inte	Q Water, egral Ultrapure Water System or equivalent	-	-

Tips & Best Practices



Tips & Best Practices section includes additional guidance



Signifies critical step requiring accurate execution



Troubleshooting section includes additional guidance

General Reagent Handling


- · Thoroughly mix reagents before use.
- Use a pH meter to adjust pH as necessary during buffer preparation.

Pipette Calibration

• Follow manufacturer's calibration and maintenance schedules.

Visium Spatial Gene Expression Slide


- Visium slides include 4 Capture Areas (6.5 x 6.5 mm), each with ~5,000 unique gene expression spots.
- Each gene expression spot includes primers with a unique Spatial Barcode.
- The active surface of the slide is defined by a readable label that includes the serial number.
- The tissue sections are always placed on the active surface of the Capture Areas. For more information, consult the Visium Spatial Gene Expression for FFPE Tissue Preparation Guide (Demonstrated Protocol CG000408).
- Always store slides in a cool, dry environment. After tissue placement, store the slides at room temperature in a low moisture environment such as a desiccator.



Slide Handling

- · Always wear gloves when handling slides.
- Ensure that the active surface of a slide faces up and is never touched. The orientation of the label on the slide defines the active surface.
- The tissue sections should always be on the active surface of the slide. DO NOT touch the tissue sections on the slide.
- Minimize exposure of the slides to sources of particles and fibers.
- When immersing slides in deparaffinization solutions and water, ensure that the tissue sections are completely submerged.
- Xylene and ethanol may cause the slide label to come off. Keep the label above the surface of the liquid when immersing in xylene and ethanol.
- Keep the slide flat on the bench when adding reagents to the active surface.
- Ensure that no absorbent surface is in contact with the reagents on the slide during incubation.





Slide Incubation Guidance

Incubation at a specified temperature

Incubation using a Section Dryer Oven:

- Place the slides in a slide drying rack.
- Close the lid when incubating the slide in the oven.



Incubation using a Thermal Cycler:

- Position a Thermocycler Adaptor on a thermal cycler that is set at the incubation temperature.
- Ensure that the Thermocycler Adaptor is in contact with the thermal cycler surface uniformly.
- When incubating a slide, position the slide on the Thermocycler Adaptor with the active surface facing up.
- Ensure that the entire bottom surface of the slide is in contact with Thermocycler Adaptor.
- When incubating a slide encased in a Visium Cassette, place the assembled unit on the Thermocycler Adaptor with the wells facing up. The Visium Cassette should always be sealed when in the Thermocycler Adaptor.



Incubation at room temperature

- Place the slide/Visium Cassette on a flat, clean, non-absorbent work surface.
- Ensure that no absorbent surface is in contact with the reagents on the slide during incubation.

Visium Cassette

- The Visium Cassette encases the slide and creates leakproof wells for adding reagents.
- Place the slides in the Visium Cassette only when specified.
- The Visium Cassette includes a removable Visium Gasket.
- An Insert Clip and four tabs at the back of the Visium Cassette are used for holding the slide in the cassette, as shown.
- The removable Visium Gasket corresponds to the Capture Areas on the slides.
- The Visium Cassette may be assembled using the Slide Alignment Tool or manually. Instructions for both are provided in the following section.
- See Visium Cassette Assembly & Removal instructions for details.
- Ensure that the back of the Visium Cassette is facing the user before assembly. The active surface of the slide with tissue sections will face down such that the slide label is no longer readable.
- Practice assembly with a plain glass slide (75 x 25 x 1 mm).
- Applying excessive force to the slide may cause the slide to break.






Visium Cassette Assembly



Exercise caution when handling slide edges to prevent injury.

Slide insertion may push Visium Gasket out of alignment with slide cutouts. Adjust if necessary.

Position Visium Cassette along alignment tool ridges




Visium Cassette secured on alignment tool




Insert long edge of slide under tabs 1 & 2; ensure slide is flush



Remove Visium Cassette while pressing slide against the Visium Gasket



Push Insert Clip along the ridge & press Visium Cassette down



Position Visium Gasket to align with Visium Cassette cutouts




Press slide down until it is flush with the Visium Gasket and under tabs 3 & 4



Slides in images are representative.

Visium Cassette Removal

Position Visium Cassette along alignment tool ridges




Visium Cassette sits securely on alignment tool



Push Insert Clip along the ridge & press down



Lift slide at Visium Cassette groove



Slides in images are representative.

Manual Visium Cassette Assembly & Removal


Assembly

- a. Remove the Visium Gasket from the Visium Cassette and re-insert the Visium Gasket, ensuring that the Visium Gasket and Visium Cassette cutouts are aligned.
- **b.** Align the label on top of the slide to the top of the Visium Cassette, as shown.
- c. Insert the slide under tabs 1 and 2. Ensure that the long edge of the slide is flush with the side of the Visium Cassette.
- **d.** Press the insert clip very firmly by applying even force on the lower part of the insert clip.
- e. Place a finger in between tab 3 and the top of the cassette, and one finger between tab 4 and the bottom of the cassette. Press down on the slide evenly until the slide is under each tab and release the insert clip.

Removal

- **a.** Press the insert clip very firmly to release the slide from the cassette.
- **b.** Lift slide at Visium Cassette groove between tabs 3 and 4 until the slide can be removed.






Reagent Addition to Wells

- Place the assembled slide in the Visium Cassette flat on a clean work surface.
- · Dispense reagents along the side of the wells without touching the tissue sections and without introducing bubbles.




- Always cover the tissue section completely when adding reagents to the well. A gentle tap may help spread the reagent more evenly.
- Ensure that no bubbles are introduced in the process.



Reagent Removal from Wells

- Place the assembled slide in the Visium Cassette flat on a clean work surface.
- Slightly tilt the Visium Cassette while removing the reagent.
- Place the pipette tip to the bottom edge of the wells.
- · Remove reagents along the side of the wells without touching the tissue sections and without introducing bubbles.
- Ensure that no bubbles are introduced in the process.
- · Remove all liquid from the wells in each step. To ensure complete removal, check the bottom of the well by tilting the cassette slightly. A meniscus at the bottom of the well will indicate the presence of liquid in the well.



Visium Slide Seal Application & Removal

To generate new or resize Visium Slide Seals, use one of the provided seals (PN-2000283/2000284) as a template to cut additional seals from MicroSeal 'B' PCR Plate Sealing Film (PN-MSB1001; listed in Specific Reagents & Consumables). Contact support@10xgenomics.com if assistance is required.

Application

- Place the Visium Cassette flat on a clean work surface.
- Remove the back of the adhesive Visium Slide Seal.
- Align the Visium Slide Seal with the surface of the Visium Cassette and apply while firmly holding the Visium Cassette with one hand.
- Press on the Visium Slide Seal to ensure uniform adhesion.

Removal

- Place the Visium Cassette flat on a clean work surface.
- Pull on the Visium Slide Seal from the edge while firmly holding the Visium Cassette.
- Ensure that no liquid splashes out of the wells.



Tissue Detachment on Visium Slides



- · Monitor section adhesion on the Visium slides throughout the workflow. For more information, consult the Visium Spatial Gene Expression for FFPE - Tissue Preparation Guide (Demonstrated Protocol CG000408).
- Tissue detachment during the workflow can impact performance.


Antibody Optimization

Prior testing of the antibodies is recommended on the same tissue before performing immunofluorescence staining in combination with the Visium Spatial Gene Expression for FFPE workflow. Determination of the optimal antibody concentration is crucial for executing this protocol.

- Optimal antibody concentrations for this Demonstrated Protocol may differ from other applications. Composition of reagents and buffers may also differ from other immunofluorescence applications.
- To optimize antibody concentration, draw representative frames on the back of a 75 x 25 x 1 mm plain glass slide using the example slide layout. Place tissue sections in the frames on the front of the slide for compatibility with the Visium Cassette. Ensure that tissue sections used during optimization are similar in size to tissue sections used for Visium Spatial Gene Expression for FFPE workflow.
- Execute the Deparaffinization, Decrosslinking & Immunofluorescence Staining protocol using a range of antibody concentrations. A starting concentration of $0.01 \,\mu\text{g/µl}$ ($0.5 \,\mu\text{g/sample}$) is recommended.
- · Select the antibody concentration that results in the specific staining of desired cells, while minimizing nonspecific background staining.
- · When optimizing the antibody, ensure that stained slides can be imaged according to the imaging guidelines listed in Visium Spatial Gene Expression for FFPE Imaging Guidelines Technical Note (CG000436).



 Wash Visium Cassette and Visium Gasket after immunofluorescence staining. See Visium Cassette & Visium Gasket Cleaning for more information.




Visium Cassette and Visium Gasket Cleaning

Remove slide from Visium Cassette



Remove Visium Gasket




Spray with 70% isopropanol, then rinse with ultrapure water



Air dry



Mark the top portion of the Visium Gasket that faced the slide with an alcohol resistant marker




Rinse Visium Cassette and Gasket with ultrapure water



Spray with 70% isopropanol a second time, then rinse with ultrapure water



Re-insert the Visium Gasket such that the marked portion is at the top of the Visium Cassette and now faces the Visium Cassette.



Protocol Steps & Timing

~5.5 hours

Steps		Timing	Stop & Store
Step 1	- Deparaffinization		
1.2	Deparaffinization	60 min	
Step 2	- Decrosslinking		
1.3	Decrosslinking	60 min	
Step 3	- Immunofluorescence Staining		
1.4	Immunofluorescence Staining 60-150 min depending upon the type of immunofluorescence staining		
	performed, see below. 1.4A Fluorophore Conjugated Primary Antibodies	60 min	
	OR 1.4B Primary and Secondary Antibodies	150 min	
Step 4	- Tissue Imaging		
2.2 2.3	Imaging Coverslip Removal	variable 5 min	4°C ≤2 weeks

1. Deparaffinization, Decrosslinking & Immunofluorescence (IF) Staining

1.0 Overview

This chapter provides guidance on deparaffinization, decrosslinking, and immunofluorescence (IF) staining of Visium slides containing FFPE tissue sections that are dried overnight in a desiccator. Ensure that microscope settings have been verified and imaging programs have been created before starting this program. Consult the Visium Spatial Gene Expression for FFPE Imaging Guidelines Technical Note (CG000436) for more information. If staining using fluorophore conjugated primary antibodies, proceed directly to step 1.4A (Immunofluorescence staining-Fluorophore Conjugated Primary Antibodies) after completing step 1.3 (Decrosslinking). If staining using primary and secondary antibodies, proceed directly to step 1.4B (Immunofluorescence staining-Primary and Secondary Antibodies) after completing step 1.3 (Decrosslinking).



1.1 Preparation - Buffers


For	For Deparaffinization				
Ite	ms	Preparation & Handling			
	Xylene	Label two coplin jars as Xylene Jar 1 and 2. Dispense 30 ml xylene in each.			
	100% Ethanol	Label three coplin jars as 100% Ethanol Jar 1, 2, and 3. Dispense 30 ml 100% ethanol in each.			
	96% Ethanol	Label two coplin jars as 96% Ethanol Jar 1 and 2. Dispense 30 ml 96% ethanol in each.			
	70% Ethanol	Label one coplin jar as 70% Ethanol Jar. Dispense 30 ml 70% ethanol.			
	Milli-Q Water	Label one coplin jar as Milli-Q Water Jar. Dispense 30 ml Milli-Q water. Alternatively, use a 50-ml centrifuge tube or a beaker for each individual slide. Clean beaker with Milli-Q Water then RNase decontamination solution in between use.			



Alternatively, a slide staining dish can also be used in place of a coplin jar. Adjust the volumes of deparaffinization solutions and water, accordingly. Use xylene-resistant dishes, for immersion in xylene. Use xylene-resistant gloves or forceps for deparaffinization. Prepare fresh reagents every week. Use fresh reagents after every 20 slides or every month (whichever comes first).

For	Decrosslinking	
Ite	ms	Preparation & Handling
	1X PBS	Prepare 1X PBS using nuclease-free water.
	TE Buffer, pH 9.0 Or	
	Prepare TE Buffer, pH 9.0	 TE Buffer Preparation (Prepare fresh and maintain at room temperature) Dissolve 1.21 g Tris base in 950 ml nuclease-free water. Adjust the pH to 9.0 with 1.0 M HCl. Add 2 ml of 0.5 M EDTA and adjust the pH to 9.0 using 1.0 M HCl. Bring the volume to 1000 ml using nuclease-free water.


Protocol Overview



Remove

Remove

TE Buffer



Immunofluorescence Staining - Primary and Secondary Antibodies ~2.5 h Add 1X Blocking Buffer Add Primary Antibody (100 µl/well) & Solution Incubate at **room** (100 μ l/well) & temperature for 5 min Incubate at room temperature Decrosslinked for 60 min Remove Blocking Remove Antibody Solution slide in the Add 1X Washing Visium Cassette Buffer Buffer (100 µl/well) & Incubate at **room** temperature for **5 min** Washing Buffer Repeat 2x Add Secondary Antibody Solution (100 µl/well) & Incubate at **room** temperature for **60 min** Remove Add 1X Washing Solution Buffer (100 μ l/well) & Incubate at room temperature for **5 min** Remove Washing Buffer Repeat 2x Remove from PBS & proceed Remove from to Coverslipping Immerse 20 times in the Visium Cassette 1X PBS 1X PBS

1.2 Deparaffinization

Deparaffinization steps should be performed in a fume hood due to the hazardous nature of xylene.

Visium Spatial slide contains a readable label with a serial number on the active surface of the slide. The label may come off during deparaffinization steps. In such cases, the serial number etched on the slide can be used.

- a. Retrieve the slide with tissue sections from the desiccator after overnight drying.
- **b.** Place slides in a Section Dryer Oven and incubate uncovered at **60°C** for **2 h**. Keep the oven lid closed during incubation.





Alternatively, place a Thermocycler Adaptor on a thermal cycler set at **60°C**. Place slide on the Thermocycler Adaptor with the active surface facing up and incubate **2 h** at **60°C**. DO NOT close the thermal cycler lid.



c. Remove from the oven or thermal cycler and allow the slide to cool down to room temperature.

d. Gently immerse the slide in the xylene in Xylene Jar 1. Secure the jar cap to prevent xylene loss.




When immersing slides in xylene, ensure that the tissue sections are completely submerged and the xylene does not reach the label.

- e. Incubate for 10 min.
- **f.** During incubation, set up a thermal cycler to 37°C. Place a Thermocycler Adaptor on the thermal cycler.
- **g.** Gently immerse slide in the xylene in Xylene Jar 2 and incubate for 10 min.



When immersing slides in ethanol, ensure that the tissue sections are completely submerged and the ethanol does not reach the label.

- i. Gently immerse slide in the 100% Ethanol Jar 2 for 3 min.
- **j.** Gently immerse slide in the 100% Ethanol Jar 3 for **3 min**.
- **k.** Gently immerse slide in the 96% Ethanol Jar 1 for **3 min**.
- **1.** Gently immerse slide in the 96% Ethanol Jar 2 for **3 min**.
- m. Gently immerse slide in the 85% Ethanol Jar for 3 min.
- **n.** Gently immerse slide in the 70% Ethanol Jar for **3 min**.
- o. Gently immerse slide in the water in the tube/beaker and incubate for 20 sec.
- **p.** Let the slide air dry and proceed to Decrosslinking.



1.3 Decrosslinking

Reagent addition and removal should be done carefully. Remove reagents along the side of the wells without touching the tissue sections and without introducing bubbles.



- a. Place a Thermocycler Adaptor on a thermal cycler set at 37°C. Place slide on the Thermocycler Adaptor with the active surface facing up and incubate 3 min at **37°C**. DO NOT close the thermal cycler lid.
- **b.** Prepare the thermal cycler with the following incubation protocol and start the program.

Lid Temperature	Reaction Volume	Run Time
70°C	-	1 h
Step	Temperature	Time
Pre-equilibrate	70°C	Hold
Decrosslinking	70°C	00:60:00
Hold	22°C	Hold

c. Place the slide in the Visium Cassette. See Tips & Best Practices for assembly instructions. Practice assembly with a blank slide.



- **d.** Add **150 µl** TE Buffer (pH 9.0) along the side of the wells.
- e. Remove all TE Buffer from the wells.
- **f.** Add $100 \mu l$ TE Buffer (pH 9.0) along the side of the wells.
- g. Apply Visium Slide Seal on the Visium Cassette and place the cassette on the Thermocycler Adaptor at 70°C.
- **h.** Close the thermal cycler lid. Skip pre-equilibrate step and initiate Decrosslinking.
- i. After decrosslinking is complete, remove the cassette from the Thermocycler Adaptor and place on a flat, clean work surface.
- **j.** Remove the slide seal and using a pipette, remove all TE Buffer from the well corners.
- **k.** Add **100 μl** 1X PBS along the side of the wells to uniformly cover the tissue sections, without introducing bubbles.
- 1. Incubate for **5 min** at **room temperature**.
- m. Remove all PBS from the wells.
- **n.** Proceed **immediately** to appropriate Immunofluorescence Staining.

1.4 Immunofluorescence Staining

Choose appropriate staining protocol depending upon the antibodies used

1.4A Immunofluorescence Staining - Fluorophore Conjugated Primary Antibodies

a. Prepare 2X Blocking Buffer on ice. Pipette mix 10x and centrifuge briefly.



If using primrary and secondary antibody, proceed directly to 1.4B.

2X Blocking Buffer Add reagents in the order listed. Maintain on ice	Stock	Final	1Χ (μl)	4X+ 10% (μl)	8X + 10% (μl)
PBS	10X	2X	80.0	352.0	704.0
BSA	10%	4%	160.0	704.0	1408.0
Tween-20	10%	0.2%	8.0	35.2	70.4
Protector RNase Inhibitor	40 U/μl	2 U/μl	20.0	88.0	176.0
Nuclease-free Water	-		132.0	580.8	1161.6
Total	-		400.0	1760.0	3520.0

b. Prepare Antibody Solution on ice. Pipette mix 10x and centrifuge briefly.

Antibody Solution Add reagents in the order isted. Maintain on ice	Stock	Final	1Χ (μl)	4X+ 10% (μl)	8X + 10% (μl)
2X Blocking Buffer	2X	1X	50.0	220.0	440.0
Nuclease-free Water	Variable	Variable	Variable	Variable	Variable
Diluted DAPI Dilute 1:100 in nuclease-free water before using.	10 μg/ml	0.2 μg/ml	2.0	8.8	17.6
Antibody #1	Variable	Variable	Variable	Variable	Variable
Antibody #2 (Optional)	Variable	Variable	Variable	Variable	Variable
Antibody #3 (Optional)	Variable	Variable	Variable	Variable	Variable
Total	-	-	100.0	440.0	880.0

*Antibody dilution can change depending on the antibody, ranging from 1:50 up to 1:1000. Antibody volumes will depend on concentrations determined during antibody optimization. Add an appropriate volume of nuclease-free water based on the amount of added antibody to achieve the stated total volume. **c.** Prepare 1X Blocking Buffer on ice. Pipette mix 10x and centrifuge briefly.

1X Blocking Buffer Add reagents in the order listed. Maintain on ice	Stock	Final	1Χ (μl)	4X+ 10% (μl)	8X + 10% (μl)
2X Blocking Buffer	2X	1X	250.0	1100.0	2200.0
Nuclease-free Water	-	-	250.0	1100.0	2200.0
Total	-	-	500.0	2200.0	4400.0

- **d.** Add **100 μl** 1X Blocking Buffer along the side of the wells.
- **e.** Incubate for **5 min** at **room temperature**.
- **f.** Remove all Blocking Buffer from the wells.
- g. Add 100 µl Antibody Solution along the side of the wells. Tap gently to ensure uniform coverage.
- **h.** Incubate for **60 min** at **room temperature**.
- i. Remove Antibody Solution
- j. Add $100 \mu l$ 1X Blocking Buffer along the side of the wells.
- **k.** Incubate for **1 min** at **room temperature**.
- **l.** Remove all Blocking Buffer from the wells.
- **m. Repeat** j-l three more times for a total of four washes.



- n. Remove the slide from Visium Cassette. See Tips & Best Practices for removal instructions.
- **o.** Gently immerse slide 20 times in 1X PBS in 50-ml centrifuge tube.



p. Remove slide from the PBS and proceed immediately to Coverslipping. DO NOT let the slide dry.

1.4B Immunofluorescence Staining - Primary and Secondary Antibodies

a. Prepare 2X Blocking Buffer on ice. Pipette mix 10x and centrifuge briefly.

2X Blocking Buffer Add reagents in the order listed. Maintain on ice	Stock	Final	1Χ (μl)	4X+ 10% (μl)	8X + 10% (μl)
PBS	10X	2X	80.0	352.0	704.0
BSA	10%	4%	160.0	704.0	1408.0
Tween-20	10%	0.2%	8.0	35.2	70.4
Protector RNase Inhibitor	40 U/μl	2 U/μl	20.0	88.0	176.0
Nuclease-free Water	-		132.0	580.8	1161.6
Total	-		400.0	1760.0	3520.0

b. Prepare 1X Blocking Buffer on ice. Pipette mix 10x and centrifuge briefly.

1X Blocking Buffer Add reagents in the order listed. Maintain on ice	Stock	Final	1Χ (μl)	4X+ 10% (μl)	8X + 10% (μl)
2X Blocking Buffer	2X	1X	250.0	1100.0	2200.0
Nuclease-free Water	-	-	250.0	1100.0	2200.0
Total	-	-	500.0	2200.0	4400.0

c. Prepare Primary Antibody Solution on ice. Pipette mix 10x and centrifuge briefly.

Primary Antibody Solution Add reagents in the order listed. Maintain on ice	Stock	Final	1Χ (μl)	4X+ 10% (μl)	8X + 10% (μl)
2X Blocking Buffer	2X	1X	50.0	220.0	440.0
Antibody #1	Variable	Variable	Variable	Variable	Variable
Antibody #2 (Optional)	Variable	Variable	Variable	Variable	Variable
Antibody #3 (Optional)	Variable	Variable	Variable	Variable	Variable
Total	-	-	100.0	440.0	880.0

^{*}Antibody dilution can change depending on the antibody, ranging from 1:50 up to 1:1000. Antibody volumes will depend on concentrations determined during antibody optimization.

d. Prepare 1X Washing Buffer. Pipette mix 10x and centrifuge briefly.

1X Washing Buffer Add reagents in the order listed. Maintain on ice	Stock	Final	1Χ (μl)	4X+ 10% (μl)	8X + 10% (μl)
PBS	10X	1X	100.0	440.0	880.0
Tween-20	10%	0.1%	10.0	44.0	88.0
Nuclease-free Water	-		890.0	3916.0	7832.0
Total	-		1000.0	4400.0	8800.0

e. Prepare Secondary Antibody Solution on ice. Pipette mix 10x and centrifuge briefly.

Secondary Antibody Solution Add reagents in the order listed. Maintain on ice	Stock	Final	1Χ (μl)	4X+ 10% (μl)	8X + 10% (μl)
2X Blocking Buffer	2X	1X	50.0	220.0	440.0
Diluted DAPI Dilute 1:100 in nuclease-free water before using.	10 μg/ml	0.2 μg/ml	2.0	8.8	17.6
Antibody #1	Variable	Variable	Variable	Variable	Variable
Antibody #2 (Optional)	Variable	Variable	Variable	Variable	Variable
Antibody #3 (Optional)	Variable	Variable	Variable	Variable	Variable
Total	-	-	100.0	440.0	880.0

^{*}Antibody dilution can change depending on the antibody, ranging from 1:50 up to 1:1000. Antibody volumes will depend on concentrations determined during antibody optimization.

- f. Add 100 µl 1X Blocking Buffer along the side of the wells.
- **g.** Incubate for **5 min** at **room temperature**.
- h. Remove all Blocking Buffer from the wells.
- i. Add 100 µl Primary Antibody Solution along the side of the wells. Tap gently to ensure uniform coverage.
- **j.** Incubate for **1h** at **room temperature**.

- k. Remove Primary Antibody Solution.
- Add **100 µl** 1X Washing Buffer along the side of the wells.
- m. Incubate for 5 min at room temperature.
- **n.** Remove 1X Washing Buffer from the wells.
- o. Repeat l-n two more times for a total of three washes.
- p. Add 100 µl Secondary Antibody Solution along the side of the wells. Tap gently to ensure uniform coverage.
- **q.** Incubate for **60 min** at **room temperature**.
- Remove Secondary Antibody Solution.
- Add **100 µl** 1X Washing Buffer along the side of the wells.
- Incubate for **5 min** at **room temperature**.
- **u.** Remove 1X Washing Buffer from the wells.
- v. Repeat s-u two more times for a total of three washes.



- w. Remove the slide from Visium Cassette. See Tips & Best Practices for removal instructions.
- **x.** Gently immerse slide 20 times in 1X PBS in 50-ml centrifuge tube.
- y. Remove slide from the PBS and proceed **immediately** to Coverslipping.




DO NOT let the slide dry.

1.5 Coverslipping

- a. Place slide on a flat, clean, non-absorbent work surface. Some residual droplets may remain.
- b. Add 5 drops of SlowFade Diamond Antifade Mountant to uniformly cover all tissue sections on the slide.
- c. Apply the coverslip at an angle on one end of the slide. Slowly lower the coverslip, without introducing bubbles. Allow the mounting medium to spread and settle.
- **d.** If needed, remove any large excess of mounting media by carefully wicking away from the edge of the coverslip with a laboratory wipe. Be careful not to move coverslip and disturb the tissue.




e. Once coverslipping is complete, **immediately** proceed with imaging. DO NOT let the attached coverslip dry. DO NOT use Cytoseal or nail polish for securing the coverslip.



2. Tissue Imaging

2.0 Overview

This chapter provides guidance on imaging Visium slides containing immunofluorescent stained FFPE sections and coverslip removal.



2.1 Imaging System Recommendations

The following table shows imaging systems used by 10x Genomics in the development of this protocol. Any equivalent imaging setup can be used as an alternative.

Supplier	Model	Configuration
Thermo Fisher Scientific	EVOS M7000	Inverted
Leica	Aperio Versa 8	Upright
	Leica DMi8	Inverted
MetaSystems	Metafer	Upright
Nikon	Nikon Eclipse Ti2	Inverted
BioTek	Cytation 7	Inverted or Upright
Keyence	Keyence BZX800	Inverted

Fluorescence Recommended Configuration
Light source (or equivalent) with a wavelength range of 380–680 nm
Monochrome camera (14 bit, 2,424 x 2,424 pixel resolution)
DAPI filter cube (Excitation 392/23, Emission 447/60)
FITC filter cube (Excitation 480/40, Emission 535/50)
TRITC filter cube (Excitation 542/20, Emission 620/52)
Cy5 filter cube (Excitation 618/50, Emission 698/70)
2.18 µm/pixel minimum capture resolution
Exposure times 100 milli sec-2 sec


2.2 Imaging

- **a.** Image all Capture Areas individually at the desired magnification using fluorescence imaging settings. Ensure that fiducial frames are captured.
- **b.** Consult the Visium Spatial Gene Expression for FFPE Imaging Guidelines Technical Note (CG000436) for additional information.
- **c.** After imaging, proceed **immediately** to the Coverslip Removal.

2.3 Coverslip Removal

- a. Dispense 800 ml Milli-Q water in a beaker.
- **b.** Remove the coverslip immediately after imaging is complete.
- c. Immerse the slide sideways/horizontal in the beaker containing 800 ml Milli-Q water with the coverslipped surface fully sideways.
- **d.** Hold the slide in water until the coverslip slowly separates away from the slide. To avoid damaging the tissue sections and Capture Areas or causing tissue detachment, DO NOT move the slide up and down, shake forcibly or manually move the coverslip.





- **e.** Gently immerse 15x in the Milli-Q water to ensure all mounting medium is removed.
- f. Wipe the back of the slide with a laboratory wipe. Place on a flat, clean, nonabsorbent work surface and air dry.
- g. Incubate slide on the Thermocycler Adaptor with the thermal cycler lid open for 3 min at 37°C.
- **h.** Remove from the thermal cycler.



- Store the slide in a slide mailer with a desiccator or 50 -ml tube at 4°C for up to 2 weeks.
- j. Place the slide in the Visium Cassette. See Tips & Best Practices for assembly instructions.
- **k.** Add **100 μl** 1X PBS along the side of the wells and proceed **immediately** to Visium Spatial Gene Expression for FFPE User Guide (CG000407).

Troubleshooting

STEP	Notes
Tissue detachment	 If tissue detachment is observed during the workflow, contact support@10xgenomics.com.
2.2 Weak or no signal	 Verify that samples were not exposed to light after staining with fluorescent antibodies.
	 Verify antibody compatibility with decrosslinking conditions.
	 Verify antibody dilutions. Ensure that antibody optimization is per- formed before immunofluorescence staining.
	 Verify imaging system filter cubes and wavelength. Ensure that fluorophores and filter cubes match.
	 Protein of interest may have low expression. Test antibodies on tissues of interest before working with Visium Spatial slides.
2.2 High background	 Verify that samples did not dry out during the staining protocol. Ensure that samples always remain covered in liquid.
	 To prevent non-specific antibody binding, compare chosen anti- body with antibodies that target the same cell type. If possible, compare staining results to cells known to express higher or lower levels of the target protein.

Results

Performing the Deparaffinization, Decrosslinking, Immunofluorescence Staining & Imaging protocol will likely result in a decrease in the number of unique transcripts detected for many tissue types with no impact on the fraction of reads mapped confidently to the transcriptome, as compared to the Deparaffinization, H&E Staining, Imaging & Decrosslinking protocol (CG000409). However, this should not affect interpretation of experimental results.




Figure 1. All change in UMI counts when downsampled to the same sequencing depth of 50,000 reads per spot.

Document Revision Summary

CG000410 **Document Number**

Visium Spatial Gene Expression for FFPE - Deparaffinization, **Title**

Decrosslinking, IF Staining & Imaging

Revision Rev D to Rev E

Revision Date September 2023

Updated for general minor consistency of language and terms **General Changes**

throughout

• Updated Visium Spatial Gene Expression Slide Kit storage information

Specific Changes

LEGAL NOTICE © 2023 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: www.10xgenomics.com/trademarks, 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at:www.10xgenomics.com/patents. The use of products described herein is subject to 10x Genomics Terms and Conditions of Sale, available at www.10xgenomics.com/legal-notices, or such other terms that have been agreed to in writing between 10x Genomics and user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of 10x Genomics products in practicing the methods set forth herein has not been validated by 10x Genomics, and such non-validated use is NOT COVERED BY 10X GENOMICS STANDARD WARRANTY, AND 10X GENOMICS HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics terms and conditions of sale for the Chromium Controller or the Chromium Single Cell Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics that it currently or will at any time in the future offer or in any way support any application set forth herein.

Contact:

support@10xgenomics.com

10x Genomics 6230 Stoneridge Mall Road Pleasanton, CA 94588 USA

