

CG000705 Rev A



**USER GUIDE**

# Chromium Connect Automated Test Kit

**FOR USE WITH**

Chromium Automated Test Library and Gel Bead Kit, 24 rxns PN-1000578

## Notices

### Document Number

CG000705 • Rev A

### Legal Notices

© 2023 10x Genomics, Inc. All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: [www.10xgenomics.com/trademarks](http://www.10xgenomics.com/trademarks). 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: [www.10xgenomics.com/patents](http://www.10xgenomics.com/patents). The use of products described herein is subject to 10x Genomics Terms and Conditions of Sale, available at [www.10xgenomics.com/legal-notices](http://www.10xgenomics.com/legal-notices), or such other terms that have been agreed to in writing between 10x Genomics and user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

### Instrument & Licensed Software Updates Warranties

Updates to existing Instruments and Licensed Software may be required to enable customers to use new or existing products. In the event of an Instrument failure resulting from an update, such failed Instrument will be replaced or repaired in accordance with the 10x Limited Warranty, Assurance Plan or service agreement, only if such Instrument is covered by any of the foregoing at the time of such failure. Instruments not covered under a current 10x Limited Warranty, Assurance Plan or service agreement will not be replaced or repaired.

### Support

Email: [support@10xgenomics.com](mailto:support@10xgenomics.com)

10x Genomics

6230 Stoneridge Mall Road

Pleasanton, CA 94588 USA

**Document  
Revision  
Summary**

<b>Document Number</b>	CG000705
<b>Title</b>	Chromium Connect Automated Test Kit User Guide
<b>Revision</b>	Rev A
<b>Revision Date</b>	October 2023

---

# Table of Contents

Notices	2
Document Revision Summary	3
Objective	5
Chromium Automated Test Kit	6
Chromium Automated Test Kit	7
Additional Kits, Reagents & Equipment	8
<a href="#">Automated Test Kit Workflow</a>	<a href="#">9</a>
Instrument Orientation	10
Deck Orientation	11
Prepare Instrument	13
Initialize Run	14
Gather Items	15
Prepare Reagents	16
Load Instrument	17
Scan	18
Run	18
Unload & Check	19
<a href="#">Troubleshooting</a>	<a href="#">20</a>
Troubleshooting	21

## Objective

Chromium Connect is an automated platform for preparation of sequencing-ready libraries from input samples. Before executing an automated assay workflow on Chromium Connect using actual samples, an instrument functional test run can be performed with the Chromium Automated Test Kit. The automated functional test run uses test reagents, consumables, and a test sample to verify that the instrument hardware and software are operating as intended. After the functional test run, inspection of the final liquid volumes in plates and tubes on the instrument deck and a test run report are used as readouts to assess the run.

This document outlines the key steps for executing the Chromium Automated Test Kit workflow (see overview below).



## Chromium Automated Test Kit

### Chromium Automated Test Library and Gel Bead Kit, 24 rxns PN-1000578

#### Chromium Automated Test Library Kit Module, 24 rxns PN-1000574\*

(store at ambient temperature)

\*3 kits included per order

##### Chromium Automated Test Module, 24 rxns

	#
<input checked="" type="checkbox"/> Test Module Strips	24 tube strips
<input type="checkbox"/> Tube B	2
<input type="checkbox"/> Tube P	2

10xGenomics.com

10x  
GENOMICS



#### Chromium Automated Test Chip Kit, 48 rxns PN-1000575

(store at ambient temperature)

##### Chromium Automated Test Chip Kit, 48 rxns

	#
Test Chip	6
50% Glycerol	6
Tube O	6
Test Sample	6

10xGenomics.com

10x  
GENOMICS



## Chromium Automated Test Kit


### Chromium Automated Test Library and Gel Bead Kit, 24 rxns PN-1000578

#### Chromium Automated Test Gel Beads Kit, 24 rxns, PN-1000576 (store at ambient temperature)

##### Chromium Automated Test Gel Beads Kit, 24 rxns

	#
Test Gel Beads	3 tube strips

10xGenomics.com

10x  
GENOMICS

#### Test Index Kit, 96 rxns, PN-1000577 (store at ambient temperature)

##### Test Index Kit, 96 rxns

	#
Test Index Plate	1

10xGenomics.com

10x  
GENOMICS

## Additional Kits, Reagents & Equipment

The items in the table below have been tested by 10x Genomics and are required for this protocol. DO NOT substitute any of the listed materials.

Supplier	Description	Part Number (US)
<b>Plastics</b>		
Hamilton	CO-RE/CO-RE II Tips 50 µl Filtered Tips*	235948
	CO-RE/CO-RE II Tips 300 µl Filtered Tips*	235903
	60 ml Reagent Reservoir Self-Standing	194051
	Hamilton PCR ComfortLid	814300
<i>*CO-RE pipette tips will be phased out and replaced by new CO-RE II pipette tips (same part number as CO-RE tips) in 2022/2023. CO-RE II tips include a new sealing surface to interface with the CO-RE II stop disk. Geometry that interfaces with the current CO-RE stop disk is identical between the two tip designs and performance remains unaffected.</i>		
Eppendorf	96-well Full-Skirted Plate	951020460
	96-well Semi-Skirted Plate (Blue color listed; other colors are acceptable)	951020362
Thermo Fisher Scientific	MicroAmp 8-Tube Strip, 0.2 ml	N8010580
	MicroAmp 8-Cap Strip, clear	4323032
<b>Kits &amp; Reagents</b>		
Thermo Fisher Scientific	Nuclease-free Water	AM9937
Millipore Sigma	Ethanol, Pure (200 Proof, anhydrous)	E7023-500ML
<b>Equipment</b>		
10x Genomics	10x Vortex Adapter	330002
	Benchtop Vortex	standard lab equipment
	Benchtop Centrifuge	standard lab equipment
	Plate Centrifuge	standard lab equipment
<b>Additional materials for Chromium Connect maintenance</b>		
Use only indicated cleaning agents. DO NOT use bleach or organic oxidizers		
Thor Labs	Lens tissues	MC-5
VWR	Microcide SQ Broad Spectrum Disinfectant	25099
Contec	70% Isopropanol (alternative to VWR disinfectant)	SB167030IR


# Automated Test Kit Workflow

Instrument Orientation

Deck Orientation


Workflow Steps

## Instrument Orientation



Refer to the Chromium Connect Instrument User Guide (CG000180) and Quick Reference Cards (CG000256) for more information.

## Deck Orientation



Refer to the Chromium Connect Instrument User Guide (CG000180) and Quick Reference Cards (CG000256) for more information.

## Deck Layout (Reagents/ Consumables)



## Chromium Automated Functional Test Run

Carrier	Zone	Item
<b>1</b> Stationary	Zone 1* (Black)	22°C Block, Blue Test Module Strips
	Zone 2 *(Gray)	22°C Block, Blue Test Module Strips
	Zone 3* (White)	10°C Block, Blue Test Module Strips
	-	Magnetic Plate
*All 3 zones will be loaded with blue Test Modules		
<b>2</b> Stationary	-	Thermal Cycler
	-	Full Skirted PCR Plate (within Thermal Cycler) <i>Test Sample loaded here</i>
	Zone 4	ComfortLids
<b>3*</b> Sliding Deck Rails: 15-18 Number of Lights: 4	Position W	Waste Reservoirs
	Position B	Tube B (White cap)
	Position P	Tube P (Blue cap)
	Position O	Tube O
	-	CO-RE Paddles
	Position E	Ethanol Reservoir
	Position LID	Lid for Ethanol Reservoir
	Position GB-1°	Test Gel Beads
<b>4</b> Sliding Deck Rails: 19-24 Number of Lights: 6	Zone 5	Test Index Plate
	Zone 6	Pipette Tips 50 µl
	Zone 7	Semi Skirted Plate
<b>5</b> Sliding Deck Rails: 25-30 Number of Lights: 6	Zone 8	Pipette Tips 300 µl

## Prepare Instrument

Before executing the workflow, it is required to develop a clear understanding of instrument operation and deck orientation.



Refer to the Chromium Connect Instrument User Guide (CG000180) and Quick Reference Cards (CG000256) for more information.

- **Ensure that the Chromium Connect instrument is powered on.**  
*Ensure that the power cables are connected correctly, and the front power button is engaged, and both power switches are turned on (located at the left side of the instrument).*
- **The login or home screen should be displayed on the instrument touchscreen.**
- **Ensure that the Chromium Automated Controller located on the instrument deck is on (LED is green) and there are no errors.**  
*If the LED indicator on the Chromium Automated Controller is red, turn off the Chromium Automated Controller, wait for 1-2 min, and then turn it back on.*

## Initialize Run

- Before initializing the run, clear the instrument deck. Ensure that all consumables (tube strips, tip racks, plates, etc.), reagents, and the waste have been removed from the deck. All carriers should be placed in their respective locations on the deck and pushed back.
- Using the touchscreen, check the Chromium Automated Controller to ensure that there is no chip in the tray.
  - i. Navigate through the top right selection menu to > Support > Control Panel > Starting Control Panel.
  - ii. After that display, the gantry will move. Do not select any options until the door unlock click is heard. Select “Initialize” at the bottom right of the screen and wait for the system to initialize.
  - iii. Navigate to Chromium Automated Controller > Get Ready > Open Tray to open the tray to check for chip. Ensure that the module lids are secure before beginning operation.
  - iv. After chip check, navigate to Close Tray > Close on the bottom right of the screen to close the front panel.
- Using the touchscreen, navigate to System Verification > Functional Test.
- Follow touchscreen prompts to enter Sample Index Plate information.
- Refer to the touchscreen for an overview of the instrument and deck and follow prompts for the next steps.

## Gather Items


Follow prompts on the touchscreen to gather the listed items and reagents for loading the instrument carriers.

Item	Qty
Nuclease-free Water	10 ml
Ethanol, Pure (200 Proof, anhydrous)	40 ml
Hamilton	
ComfortLids	3
50 µl CO-RE/CO-RE II Pipette Tips, with filter (Black, Conductive)	1 rack
300 µl CO-RE/ CO-RE II Pipette Tips, with filter (Black, Conductive)	2 racks
Reagent Reservoir, 60 ml	3
Eppendorf	
96-well Semi Skirted Plate	1
96-well Full Skirted Plate	1
Thermo Fisher Scientific	
MicroAmp 8-Tube Strip, 0.2 ml	1
10x Genomics	
Automated Test Chip ( <i>keep sealed</i> )	1
Test Module ( <i>blue tube strip</i> )	24 tube strips/run
Tubes	
Tube B ( <i>white cap</i> )	1 tube/run
Tube P ( <i>blue cap</i> )	1 tube/run
Tube O	1 tube/run
Test Sample	1 tube/run
Test Index Plate	1 plate
Test Gel Bead Strip	1 tube strip/run

## Prepare Reagents

Follow prompts on the touchscreen to prepare reagents before loading the instrument deck. Some important guidelines are highlighted below.

ACTION	GUIDELINES <i>Follow touchscreen prompts for specifics and timing</i>
<b>Reagents</b>	<ul style="list-style-type: none"> <li>Centrifuge all reagent tubes to collect liquid at the bottom of the tubes. Remove all caps when loading onto the carriers.</li> </ul>
<b>Prepare Ethanol</b>	<ul style="list-style-type: none"> <li>Prepare <b>50 ml 80% Ethanol</b> in Nuclease-free water and dispense in Ethanol Reservoir when prompted.</li> </ul>
<b>Sample Index Test Plate</b>	<ul style="list-style-type: none"> <li>Vortex Sample Index Plate for <b>15 sec</b> at maximum speed and centrifuge at <b>300 rcf</b> for <b>1 min</b> at <b>22°C</b>.</li> </ul>
<b>Test Modules</b>	<ul style="list-style-type: none"> <li>Centrifuge all test modules at <b>300 rcf</b> for <b>1 min</b> at <b>22°C</b>.</li> <li>Confirm there are no bubbles at the bottoms of any module tubes.</li> </ul>
<b>Prepare Test Gel Beads</b>	<ul style="list-style-type: none"> <li>Snap the tube strip holder with the Test Gel Bead strip into a 10x Vortex Adapter. Vortex <b>30 sec</b>.</li> <li>Centrifuge the Test Gel Bead strip for <b>~5 sec</b> after removing from the holder. Confirm there are no bubbles at the bottoms of the tubes and the liquid levels look even.</li> <li>Place the Test Gel Bead strip back in the holder and secure the holder lid.</li> </ul>




## Load Instrument


Follow the instructions on the touchscreen to load all reagents and consumables on the carriers. Some key guidelines are provided below:

### Carriers

- Handle the carriers as prompted.
- Ensure that Carriers 3, 4, and 5 are completely slid out and placed on an off-deck workspace before loading.
- Align the carriers to the corresponding Deck Rails when sliding them in or out of the deck.
- Remove tube caps and ensure correct orientation of tube labels with barcodes to enable scanning (as shown on the touchscreen).
- Follow the instrument touchscreen prompts to load the Test Chip and the Test Sample.




### Barcode Orientation




### Modules

- Load test module tube strips in the corresponding positions on the Carrier, starting from back to front row. All three zones receive the same blue Test Modules for this automated functional test run.
- DO NOT skip any rows when loading.
- Use pinhole alignment to place module tube strips in the correct orientation (as shown on the touchscreen).



### Load Test Gel Beads

- Load Test Gel Bead at Position GB-1° as shown on the touchscreen.



Consult the Chromium Connect User Guide (CG000180) for more information.

---

## Scan

After loading is complete, follow prompts on the touchscreen to select the “SCAN” button to initiate scanning.



Gantry will move during scanning.

---

## Run

- After scanning is complete, follow prompts on the touchscreen to select the “PLAY” button to start the run. The run time is ~1 h 45 min.


## Unload & Check

- After run completion, unload the instrument as per prompts on the instrument touchscreen.
- While unloading, check and document (photograph) the liquid levels in the items shown below. The representative images below show the expected liquid levels after a successful functional test run.



If these items do not look similar to the representative images, refer to the [Troubleshooting](#) section for additional guidance.






- A run report can be generated and exported as per prompts on the instrument touchscreen.

Check Items	Representative Images After a Successful Functional Test Run
<b>Carrier 2</b>	
Full Skirted Plate	
<b>Carrier 3</b>	
Collection Tube Strip	
<b>Carrier 4</b>	
Semi Skirted Plate	

# Troubleshooting

## Troubleshooting

- After run completion, if the results are not as described in the [Unload & Check](#) section, check and document (photograph) the additional items in carriers 1 and 4 (representative images shown below).
- Contact [support@10xgenomics.com](mailto:support@10xgenomics.com) with all the photographed items for guidance.

Check Items	Representative Images	
<b>Carrier 1</b>		
Zone 1-Blue Tube Strip Piercing Pattern & Volumes		
Zone 2-Blue Tube Strip Piercing Pattern & Volumes		
Zone 3- Blue Tube Strip Piercing Pattern & Volumes		
<b>Carrier 3</b>		
Test Gel Beads		
Tubes		



If the functional test run was completed as part of an open 10x Genomics Support case, all above mentioned photographs and requested run logs should be submitted along with the case email.