

# Cell Surface & Intracellular Protein Labeling for GEM-X Flex Gene Expression


## Introduction

GEM-X Flex Gene Expression with Feature Barcode technology offers comprehensive, scalable solutions to measure gene and protein expression in formaldehyde fixed samples. To measure protein expression, cell surface and intracellular proteins are labeled using an antibody conjugated to a Feature Barcode oligonucleotide, followed by a fixation with formaldehyde to ensure 1) the cells are properly fixed and permeabilized for probe hybridization and 2) the antibodies firmly bind to the cells, preventing their loss during probe hybridization.

This document outlines cell surface and intracellular protein labeling protocols for use with GEM-X Flex Gene Expression. General guidance on

antibody-oligonucleotide conjugation, evaluating antibodies with flow cytometry, and best practices when staining cells is also provided here.

Both the cell surface and intracellular protein labeling protocols described in this document require a Fc receptor blocking step, after which the appropriate labeling protocol should be followed depending on the type of antibodies used. The Intracellular Protein Labeling Protocol, which includes a brief fixation before staining to allow cell permeabilization, can also be used with certain fixation compatible cell surface protein antibodies.



## Additional Guidance

Preread and have available Fixation of Cells & Nuclei for GEM-X Flex Gene Expression Demonstrated Protocol (CG000782) before starting the protocols provided in this document. Consult the Cell Preparation Guide Handbook (CG00053) for Tips & Best Practices on handling and counting cells. Consult Cell Thawing Protocols for Single Cell Assays (CG000447) for guidance on thawing dissociated tumor cells.

Consult GEM-X Flex Gene Expression - Protocol Planner (CG000790) for details on workflow overview, document resources, and guidance on selecting the appropriate sample preparation and library construction protocols for different GEM-X Flex Gene Expression workflows.

**Cells carry potentially hazardous pathogens. Follow material supplier recommendations and local laboratory procedures and regulations for the safe handling, storage, and disposal of biological materials.**

## Specific Reagents & Consumables

### For Cell Surface Protein Labeling

Vendor	Item	Part Number
BioLegend	Human TruStain FcX (Fc Receptor Blocking Solution)	422301
	TruStain FcX PLUS (anti-mouse CD16/32) Antibody	156604
	True-Stain Monocyte Blocker	426101
	TotalSeq™-C Human Universal Cocktail, V1.0	399905
	Cell Staining Buffer	420201
Thermo Fisher Scientific	UltraPure Bovine Serum Albumin (BSA, 50 mg/ml)	AM2616
	Fetal Bovine Serum, qualified, heat inactivated	16140071

Millipore Sigma	Bovine Serum Albumin In DPBS (10%) <i>Alternative to Thermo Fisher</i>	A1595
Miltenyi Biotec	MACS BSA Stock Solution <i>Alternative to Thermo Fisher</i>	130-091-376
Corning	Phosphate-Buffered Saline, 1X <i>without Calcium and Magnesium</i>	21-040-CV
VWR	Fetal Bovine Serum (FBS) <i>Alternative to Thermo Fisher</i>	97068-085

### For Cell Counting

Thermo Fisher Scientific	*Countess II FL Automated Cell Counter	AMQAF1000
	*Countess 3 FL Automated Cell Counter	A49866
	Countess Cell Counting Chamber Slides	C10228
	Trypan Blue Stain (0.4%)	T10282
	*DAPI solution, 1 mg/mL	62248
Nexcelom Bioscience	*Celleca MX High-throughput Automated Cell Counter	MX-112-0127
	*ViaStain PI Staining Solution	CS1-0109-5mL
	*ViaStain AOPI Staining Solution	CS2-0106-5mL
	Cellometer K2 Bundle w/ Matrix Software	CMT-K2-MX-150
	PD100 Counting Chambers 1 case	CHT4-PD100-003
Biotium	*NucSpot 470	40083
	<i>Dilute the stock to 1:100 and mix 1:1 with the sample. For example, add 10 µl diluted dye to 10 µl sample.</i>	

\*Choose Countess II/3, Cellaca, or equivalent fluorescent counter.  
\*Choose either AOPI, NucSpot, PI, or DAPI solution. If the sample has no debris, Trypan Blue can be used.

*This list may not include some standard laboratory equipment.*

## For Intracellular Protein Labeling

Vendor	Item	Part Number
10x Genomics	GEM-X Flex Sample Preparation v2 Kit	1000781
<i>*Kit components are listed on next page</i>		
BioLegend	Human TruStain FcX (Fc Receptor Blocking Solution)	422301
	TruStain FcX PLUS (anti-mouse CD16/32) Antibody	156604
	True-Stain Monocyte Blocker	426101
	TotalSeq™ intracellular conjugates	-
Proteintech Genomics	MultiPro Human Fixed Cell Immune Profiling Antibody cocktail	G900004
	MultiPro Human Discovery Panel	G900150
	Multipro Enhanced Blocking Reagent	G900005
Corning	Phosphate-Buffered Saline, 1X without Calcium and Magnesium	21-040-CV
Thermo Fisher Scientific	Tween 20 (10%)	28320
	NP-40 Surfact-Amps Detergent Solution (10%)	85124
	Nonidet P40 Substitute, Ultrapure, Thermo Scientific Chemicals	J19628.AP
<i>Alternative to NP-40; consult manufacturer's instructions for preparing 10% stock solution</i>		
	Formaldehyde (37% by Weight/ Molecular Biology), Fisher BioReagents	BP531-25
	Nuclease-free Water (not DEPC-Treated)	AM9937
Millipore Sigma	Albumin, Bovine Serum, Fraction V, Fatty Acid-Free, Nuclease- and Protease-Free	126609
	Dextran Sulfate Sodium Salt 8 KDa	RES2029D-A7
Invitrogen	UltraPure Salmon Sperm DNA Solution	15-632-011
	Phosphate-Buffered Saline (10X) pH 7.4, RNase-free	AM9625
Avantor	PBS 1X	K812-500ML
Roche	Protector RNase inhibitor	3335399001

## Additional Materials

Eppendorf	DNA LoBind Tubes 2.0 ml	022431048
	1.5mL Protein LoBind tubes	022431081
	1.5mL DNA LoBind tubes	022431021
	ThermoMixer C	5382000023
VWR	Vortex Mixer	10153-838
Corning	Phosphate-Buffered Saline, 1X without Calcium and Magnesium	21-040-CV
New England Biolabs	T4 Gene 32 Protein**	M0300S/ M0300L

*\*\*Optional: Only required if adding to the antibody pool + buffer to reduce background signal. For more information, see [Troubleshooting](#).*


## For Cell Counting

Thermo Fisher Scientific	*Countess II FL Automated Cell Counter	AMQAF1000
	*Countess 3 FL Automated Cell Counter	A49866
	Countess Cell Counting Chamber Slides	C10228
	Trypan Blue Stain (0.4%)	T10282
Nexcelom Bioscience	*Celleca MX High-throughput Automated Cell Counter	MX-112-0127
	ViaStain PI Staining Solution	CS1-0109-5mL
	ViaStain AOPI Staining Solution	CS2-0106-5mL
<i>Alternative to PI Staining Solution.</i>		
	Cellometer K2 Bundle w/ Matrix Software	CMT-K2-MX-150
	PD100 Counting Chambers 1 case	CHT4-PD100-003
Biotium	NucSpot 470	40083
<i>Alternative to PI Staining Solution.</i>		
<i>Dilute the stock to 1:100 and mix 1:1 with the sample. For example, add 10 µl diluted dye to 10 µl sample.</i>		

*\*Choose Countess II/3, Celleca, or equivalent fluorescent counter.*

**\*GEM-X Flex Sample Preparation v2 Kit, 48 rxns PN-1000781**

<b>GEM-X Flex Sample Preparation v2 Kit</b>		
48 rxns, PN-1000781		
<i>Shipped on dry ice</i>		
<i>Store at -20°C</i>		
	#	PN
● Conc. Fix & Perm Buffer B	2	2001301
● Conc. Quench Buffer B	6	2001300
● Enhancer	3	2000482
● Additive C	4	2001332



If following the intracellular protein labeling workflow, only 24 samples can be processed using this kit, as an extra fixation step is required before staining. If following the no wash cell surface protein workflow, only 24 samples can be processed using this kit.

## Antibody-Oligonucleotide Conjugates

- The Cell Surface Protein Labeling protocol described in this document was optimized using pre-conjugated TotalSeq™-C from BioLegend.
- The Intracellular Protein Labeling protocol was optimized using MultiPro Fixed Cell Immune Profiling Antibody Cocktail, a pool of cell surface, intracellular, & isotype control antibodies, from Proteintech Genomics & pre-conjugated TotalSeq™ intracellular conjugates from Biolegend.
- Antibodies conjugated to appropriate Feature Barcode oligonucleotide from other vendors can also be used. See the table below for antibody-oligonucleotide conjugate capture by protocol-specific Gel Bead primers.

- Both the labeling protocols are also compatible with custom conjugated antibody panels.

### Conjugation Guidance

- Follow manufacturer's instructions ([Antibody Conjugation Kit from Abcam](#) and [oYo-Link Antibody Labeling Reagents from AlphaThera](#)) for antibody-oligonucleotide conjugation.
- Antibody conjugation services from Abcam and AlphaThera can also be utilized.
- Consult Cell Surface Protein Labeling for Single Cell RNA Sequencing Protocols with Feature Barcode technology Demonstrated Protocol ([CG000149](#)) for guidance on antibody-oligonucleotide conjugation and verification of conjugation.

### Antibody-Oligonucleotide Conjugate Capture

GEM-X Flex Gene Expression Workflows	Antibody-Oligonucleotide Conjugates
<b>Gene &amp; Protein Expression using Barcode Oligo Capture* - Singleplex &amp; Multiplex Workflows</b>	/5AmMC12/CGGAGATGTGTATAAGAGACAG-N10-N15-N9-CCCATATAAGAAA Nextera partial Read 2    Feature Barcode (15 nt)    Capture Sequence TotalSeq™-C antibody-oligonucleotide conjugates
Consult the following User Guides for more information:	/5AmMC12/CGGAGATGTGTATAAGAGACAG-N15-CCCATATAAGAAA Nextera partial Read 2    Feature Barcode (15 nt)    Capture Sequence MultiPro Fixed Cell Immune Profiling Antibody Cocktail
<ul style="list-style-type: none"> <li>• GEM-X Flex Gene Expression Reagent Kits for Multiplexed Samples with Feature Barcode technology for Protein using Barcode Oligo Capture User Guide (CG000789).</li> <li>• GEM-X Flex Gene Expression Reagent Kits for Singleplexed Samples with Feature Barcode technology for Protein using Barcode Oligo Capture User Guide (CG000788).</li> </ul>	

\*Barcode Oligo Capture can be achieved by using either MultiPro Fixed Cell Immune Profiling Antibody Cocktail or TotalSeq™-C. If using both the antibodies in a single experiment, an alternative sequencing configuration will be required. Consult 10x Genomics [support website](#) for details.

## Fc Receptor Blocking

Both cell surface and intracellular labeling protocols require a Fc receptor blocking step. After this step, follow the appropriate labeling protocol based on the type of antibodies used.

### Buffers – Preparation

- Chilled (4°C): PBS + 1% BSA
- For samples containing <70% viable cells, chilled (4°C) PBS + 10% FBS can be used.
- a. Resuspend cells in PBS + 1% BSA.
- b. Transfer  $\leq 1 \times 10^6$  cells to a new 1.5-ml microcentrifuge tube.
- c. Centrifuge cells at **400 rcf** for **5 min** (PBMCs) at **4°C**. Use of a swinging-bucket rotor is recommended for higher cell recovery. Centrifugation speed and time depend upon the sample type. Larger or fragile cell types may require slower centrifugation speeds.
- d. Remove the supernatant without disturbing the pellet. Up to 30  $\mu$ l supernatant may be left behind if working with <300,000 cells.
- e. Resuspend pellet in **50  $\mu$ l** chilled PBS + 1% BSA or chilled PBS + 10% FBS (for <70% viable cells).

**TIPS** If using  $>1 \times 10^6$  cells, scale up the buffer volumes accordingly.


- f. Add **5  $\mu$ l** Human TruStain FcX or **2  $\mu$ l** TrueStain FcX PLUS (anti-mouse CD16/32 Antibody). Gently pipette mix.

**OPTIONAL:** If performing cell surface protein labeling, **5  $\mu$ l** True-Stain Monocyte Blocker can also be added at this step in addition to TruStain FcX to reduce staining background.

**TIPS** Consult manufacturer's instructions for recommended volume of Fc receptor blocking solutions.

- g. Incubate for **10 min** at **4°C**.
- h. Immediately proceed to the appropriate step:
  - Cell surface protein detection only: proceed to [Cell Surface Protein Labeling](#)
  - Intracellular protein detection only and/or using cell surface antibodies that are compatible with brief fixation: proceed to [Intracellular Protein Labeling](#)
  - Both cell surface protein and intracellular protein detection: proceed to [Cell Surface Protein Labeling](#) followed by [Intracellular Protein Labeling](#)

### Block Fc receptor



\*For samples containing <70% viable cells, PBS + 10% FBS can be used.

## Cell Surface Protein Labeling Protocol

### Cell Surface Protein Labeling Protocol Overview



Label Cells



Wash Cells  
(no wash, 1-wash, or 2-wash options)



Fix Cells\*

\* Consult Fixation of Cells & Nuclei for GEM-X Flex Gene Expression Demonstrated Protocol (CG000782).

This protocol was optimized using TotalSeq™-C antibody-oligonucleotide conjugates from BioLegend. The labeled cells can be enriched by flow sorting (see Appendix).



*Use distinct and compatible antibody clones for flow sorting and cell surface protein labeling. Optimize the working concentration of each antibody used.*

Using the appropriate wash option is critical for obtaining high-quality data and minimizing background signal. Choose a washing protocol based on the guidance outlined on page 9 in this [protocol](#).

See Appendix for supplemental data on the different wash options.

### Buffers – Preparation

#### For Labeling Cells

- Chilled (4°C): PBS + 1% BSA
- For samples containing <70% viable cells, chilled (4°C) PBS + 10% FBS can be used.

### Prepare Antibody Mix Supernatant

Add the appropriate/manufacture's recommended amount of antibody-oligonucleotide conjugates to a 1.5-ml Protein LoBind tube.

- If using lyophilized antibody panel/cocktails, rehydrate the lyophilized panel in the recommended volume of Cell Staining Buffer as directed by the manufacturer. Follow the manufacturer's instructions for use for cell labeling. After labeling, perform cell wash steps as described in this Demonstrated Protocol.
- Centrifuge the mix at **14,000 rcf for 10 min at 4°C**.
- Transfer the supernatant (containing Antibody Mix) to a new protein LoBind tube. Leave residual volume so the antibody aggregates are not carried over. Maintain at 4°C.

### Prepare Flow Sorting Antibody Pool

- Add the appropriate/manufacture's recommended amount of fluorophore-conjugated antibodies to a 1.5-ml microcentrifuge tube on ice.
- Gently pipette mix and maintain at 4°C. Avoid light exposure.

## Tips & Best Practices

### Cell Viability

- Determine sample viability before starting the cell surface protein labeling protocol.

### Labeling & Wash Buffer

- The recommended minimum cell input for cell surface labeling is 50,000 cells. Note that at this input, recovering the maximum number of cells (20,000) may not be possible.
- Labeling or washing buffer for most sample types is PBS + 1% BSA.
- For samples containing <70% viable cells, PBS + 10% FBS can be used.
- BioLegend's Cell Staining Buffer can also be used for labeling cells. However, this buffer may not be optimal for all sample types. Cell Staining Buffer should only be used for the labeling step. PBS + 1% BSA should be used for the washing steps.
- Follow the manufacturer's instructions if using lyophilized antibody panel/cocktails regarding reconstitution and labeling volumes.

### Centrifugation Conditions

- Centrifugation speed and time depend upon the sample type.
- Larger or fragile cell types may require slower centrifugation speeds.

### Optimal Antibody Concentration

- The optimal concentration is 0.1 µg to 0.5 µg.
- Consult Quality Control of Cell Surface Protein Labeling using Flow Cytometry Technical Note (CG000231) for guidance on how to titrate antibodies using flow cytometry for determining optimal antibody concentrations.

- It is highly recommended to titrate antibodies before use for optimal performance. Consult the manufacturer's instructions for recommended dilutions for each antibody. A dilution of 0.5 µg per antibody for up to 1 x 10<sup>6</sup> cells is suggested as a starting point.

### Sample Washing

- This protocol provides three wash options after incubation with the antibody-oligonucleotide conjugate. See [Cell Surface Protein Labeling: Wash Cells](#) for guidance on choosing the appropriate wash option.

### Sample Fixation

- Consult Fixation of Cells & Nuclei for GEM-X Flex Gene Expression Demonstrated Protocol (CG000782) regarding fixing single cell suspensions following cell surface protein labeling.

### Sample Storage


- After labeling, samples may be fixed in smaller batches and stored. If batching samples together, it is recommended that samples with different fixation times in one experiment are not mixed.
- Fixed samples can be stored either at -80°C for up to 12 months or at 4°C for up to 7 days with appropriate storage reagents.
- Consult Fixation of Cells & Nuclei for GEM-X Flex Gene Expression Demonstrated Protocol (CG000782) for details on storage conditions and reagents.



## Cell Surface Protein Labeling: Label Cells

### Label Cells

Prepare Antibody Mix Supernatant and Flow Sorting Antibody Pool (if performing flow sorting) as described in Buffers – Preparation.



Follow manufacturer's instructions if using custom conjugated antibodies. Consult Cell Surface Protein Labeling for Single Cell RNA Sequencing Protocols Demonstrated Protocol (CG000149) for details. The minimum recommended cell input for labeling is 50,000 cells.

- Add prepared Antibody Mix Supernatant to the Fc blocked sample. If performing flow sorting, add Flow Sorting Antibody Pool.
- Add chilled PBS + 1% BSA to the cells to bring the total volume to **100  $\mu$ l**. Gently pipette mix 10x (pipette set to 90  $\mu$ l). For samples containing <70% viable cells, chilled PBS + 10% FBS can be used.
- Incubate for **30 min** at **4°C**. If using fluorophore-conjugated antibodies, incubate without light exposure.

- Proceed to appropriate **Cell Surface Protein Labeling: Wash Cells** section for **1-Wash or 2-Wash Options**

OR

Proceed directly to **Sample Fixation, Step c** of the Fixation of Cells & Nuclei for GEM-X Flex Gene Expression Demonstrated Protocol (CG000782) (**No-Wash Option**).

#### TIPS

Sample fixation without washing is recommended for samples with low inputs, low viability, or time-sensitive samples where increased background is acceptable.

## Cell Surface Protein Labeling: Wash Cells

Using the appropriate wash option is critical for obtaining high-quality data and minimizing background signal. Choose a protocol based on the following guidance. See the Appendix for [supplemental data on the different wash options](#).

### No Wash Option

- Samples with low viability or input counts or time-sensitive samples
- Antibody background is expected to be high
- Compatible with antibodies that have clear distinction between populations


### 1-Wash Option

- When a reduced number of wash steps is desired
- Antibody background is expected to be moderate
- Compatible with most antibodies that have a clear distinction between populations
- Includes option for flow sorting enrichment

### 2-Wash Option (Recommended)

- Recommended for most sample types
- Antibody background is expected to be low
- Best quality Feature Barcode data

## No Wash Option Protocol



### a. Proceed **immediately** to:

- **Sample Fixation, Step c** of the Fixation of Cells & Nuclei for GEM-X Flex Gene Expression Demonstrated Protocol (CG000782).

**Modified step c:** Add **1 ml** Fixation Buffer B to the labeled cells and pipette mix 5x. Fixation Buffer B preparation and fixation protocol are listed in Fixation of Cells & Nuclei for GEM-X Flex Gene Expression Demonstrated Protocol (CG000782).


- After completion of Step d of the Fixation of Cells & Nuclei for GEM-X Flex Gene Expression Demonstrated Protocol (CG000782), centrifuge at 850 rcf for **5 min** at **room temperature**, remove **500 µl** of Fixation Buffer B, and proceed with step e.

OR

- Intracellular Protein Labeling

## 1-Wash Option Protocol

### 1-Wash Option



#### a. OPTIONAL: For enrichment of labeled and viable cells by flow sorting:

- Based on starting concentration and assuming ~50% cell loss, add an appropriate volume chilled PBS + 2% FBS (including a dead cell marker) to obtain a final cell concentration of  $5\text{--}10 \times 10^6$  cells/ml.
- Proceed to sorting (see Appendix for [Flow Sorting Guidance](#)). After sorting, determine cell concentration and viability using an automated cell counter or a hemocytometer.

#### If not performing flow sorting enrichment:

- Add **1.4 ml** chilled PBS + 1% BSA to the labeled cells. Gently pipette mix. For samples containing <70% viable cells, chilled PBS + 10% FBS can be used.

- b. After resuspension, take an aliquot and determine cell concentration and viability using an automated cell counter or a hemocytometer.


#### c. Proceed **immediately** to:

- **Sample Fixation, Step a** of the Fixation of Cells & Nuclei for GEM-X Flex Gene Expression Demonstrated Protocol (CG000782)

OR

- Intracellular Protein Labeling

## 2-Wash Option Protocol (Recommended)



- a. Add **1.4 ml** chilled PBS + 1% BSA to the labeled cells. Gently pipette mix. For samples containing <70% viable cells, chilled PBS + 10% FBS can be used.
- b. Centrifuge cells at **400 rcf** for **5 min** (PBMCs) at **4°C**. Centrifugation speed and time depend upon the sample type.
- c. Remove the supernatant without touching the bottom of the tube to avoid dislodging the pellet. Up to 30  $\mu$ l supernatant may be left behind if working with <300,000 cells.
- d. Resuspend the cell pellet in **1.5 ml** chilled PBS + 1% BSA and place on ice. For samples containing <70% viable cells, chilled PBS + 10% FBS can be used.
- e. After resuspension, take an aliquot and determine cell concentration and viability using an automated cell counter or a hemocytometer.
- f. Proceed **immediately** to:
  - **Sample Fixation, Step a** of the Fixation of Cells & Nuclei for GEM-X Flex Gene Expression Demonstrated Protocol (CG000782)

OR

- Intracellular Protein Labeling

**Pellet (may not be visible)**

Leaving behind excess supernatant (>30  $\mu$ l) may cause nonspecific binding, which may result in increased background reads during sequencing.

## Intracellular Protein Labeling


### Intracellular Protein Labeling Protocol Overview



**Fix Cells**  
Brief fixation




**Label Cells**



**Fix Cells**  
CG000782

This protocol was optimized using human PBMCs and the MultiPro Fixed Cell Immune Profiling Antibody Cocktail, a pool of cell surface, intracellular, and isotype control antibodies, from Proteintech Genomics and TotalSeq™ intracellular conjugates from Biolegend.

This protocol involves a short fixation before staining and can be also be used with other cell surface protein antibodies, provided they are compatible with the brief fixation step. The labeled cells can be enriched by flow sorting (see Appendix for Flow Sorting Guidance).

 Use distinct and compatible antibody clones for sorting. Optimize the working concentration of each antibody used.

### Buffers – Preparation

#### For Fixing Cells

Buffers for Fixation - Prepare fresh			
Fixation Buffer B	Stock	Final	Per Sample + 10% (µl)
<i>Maintain at room temperature</i>			
Nuclease-free Water	-	-	435
Conc. Fix & Perm Buffer B (10x Genomics PN-2001301)	10X	1X	55
<i>Thaw at room temperature. Vortex, check for precipitation, and centrifuge briefly.</i>			
Formaldehyde*	37%	4%	60
<b>Total</b>	-	-	<b>550.0</b>

\*Formaldehyde should always be used with adequate ventilation, preferably in a fume hood. Follow appropriate regulations.

#### Buffers for Fixation - Prepare fresh

Quenching Buffer B	Stock	Final	Per Sample + 10% (µl)
<i>Maintain at 4°C</i>			
Nuclease-free Water	-	-	962.5
Conc. Quench Buffer B (10x Genomics PN-2001300)	8X	1X	137.5
<i>Thaw at room temperature. Vortex and centrifuge briefly.</i>			
<b>Total</b>	-	-	<b>1,100.0</b>

#### Additional Buffers

Additive C (10x Genomics PN-2001332)  
*Thaw at room temperature. Vortex and check for precipitation. Maintain at room temperature.*

### For Labeling Cells

#### Flow Sorting Antibody Pool Preparation

- Add appropriate/manufacture's recommended amount of fluorophore-conjugated antibodies to a 1.5-ml microcentrifuge tube on ice.
- Gently pipette mix and maintain at **4°C**. Avoid light exposure.

#### Intracellular Antibody Mix Supernatant (Antibody + Buffer) Preparation

Follow the manufacturer's recommendations for resuspension and aggregate removal. General guidance is provided below.

If using a lyophilized antibody panel:

- Equilibrate the lyophilized panel vial(s) to **room temperature** for **5 min**.
- Place the lyophilized panel vial in an empty 2-ml tube and centrifuge the vial at **10,000 rcf** for **30 sec** at **room temperature**.

- a.** Prepare Antibody Buffer in a Protein LoBind tube according to the appropriate table below, depending on choice of blocker.

Buffer for Labeling (Custom Blocker) - Prepare fresh			
Antibody Buffer <i>Maintain at 4°C</i>	Stock	Final	For 1 Sample + 10% (µl)
RNase Inhibitor	40 U/µl	0.4 U/µl	1.1
Nuclease-free BSA	30%	7.55%	27.68
Tween 20	10%	0.2%	2.2
NP-40 or Nonidet P40 Substitute	10%	0.5%	5.5
Salmon Sperm DNA	10 mg/ml	0.5 mg/ml	5.5
Dextran Sulfate 8 KDa	1%	0.025%	2.75
Monocyte Blocker	-	5 µl	5.5
FcX	-	5 µl	5.5
10X PBS	10X	1X	11.0
Nuclease-free Water* <i>If using lyophilized antibody</i>	-	-	43.27*
<b>Total</b>	-	-	<b>110.0</b>

Buffer for Labeling (MultiPro Blocker) - Prepare fresh			
Antibody Buffer <i>Maintain at 4°C</i>	Stock	Final	For 1 Sample + 10% (µl)
RNase Inhibitor	40 U/µl	0.4 U/µl	1.1
Nuclease-free BSA	30%	10%	36.67
Tween 20	10%	0.2%	2.2
NP-40 or Nonidet P40 Substitute	10%	0.5%	5.5
10X PBS	10X	1X	11.0
MultiPro Enhanced Blocking Reagent	-	-	15.5
Nuclease-free Water* <i>If using lyophilized antibody</i>	-	-	38.03*
<b>Total</b>	-	-	<b>110.0</b>

\*Adjust water volume if:

- adding reconstituted antibody pool derived from non-lyophilized antibodies
- adding fluorophore-conjugated antibodies for flow sorting



Final buffer volume should not exceed 110 µl per sample.

- b.** In case of lyophilized antibody panel: Add **110 µl** Antibody Buffer to rehydrate the lyophilized antibody panel. If using non-lyophilized antibody, directly proceed to step g.
- c.** Vortex the rehydrated antibody pool for **15 sec** and centrifuge at **10,000 rcf** for **30 sec** at **room temperature**.
- d.** Incubate for **5 min** at **room temperature**.
- e.** Vortex for **15 sec** and centrifuge at **10,000 rcf** for **30 sec** at **room temperature**.
- f.** Transfer the entire volume of the reconstituted antibody pool to a new Protein LoBind tube and maintain at **4°C**.
- g.** Centrifuge at **14,000 rcf** for **10 min** at **4°C**.
- h.** Transfer **100 µl** supernatant (containing antibody + buffer) to a new Protein LoBind tube and maintain at **4°C**.

#### For Washing Cells

Buffer for Washing - Prepare fresh			
Intracellular Wash Buffer <i>Maintain at 4°C</i>	Stock	Final	For 1 Sample + 10% (µl)
RNase Inhibitor	40 U/µl	0.2 U/µl	11.0
Tween 20	10%	0.1%	22.0
Nuclease-free BSA	30%	1%	73.33
10X PBS	10X	1X	220.0
Nuclease-free Water	-	-	1873.67
<b>Total</b>	-	-	<b>2,200.0</b>

## Tips & Best Practices

### Brief Fixation & Antibody Compatibility

- The intracellular labeling protocol includes a brief fixation (20 min) step.
- Before starting the protocol, it is recommended to test the antibodies by flow cytometry to ensure compatibility with the fixation.

### Optimal Antibody Concentration

- It is highly recommended to titrate antibodies before use for optimal performance. Consult the manufacturer's instructions for recommended dilutions for each antibody.
- A dilution of 0.2 µg per antibody for up to  $1 \times 10^6$  cells is suggested as a starting point.
- Consult Quality Control of Cell Surface Protein Labeling using Flow Cytometry Technical Note (CG000231) for guidance on how to titrate antibodies using flow cytometry for determining optimal antibody concentrations.

### Antibody-Oligonucleotide Conjugation

- Custom conjugated antibodies can also be used with this protocol. See [Conjugation Guidance](#) for details.
- Consult the Demonstrated Protocol Cell Surface Protein Labeling for Single Cell RNA Sequencing Protocols with Feature Barcode technology (CG000149) for guidance on antibody-oligonucleotide conjugation.

### Antibody Quality


- The specific antibody clones used needs to be of high-quality. 10x Genomics strongly recommends validating antibody clones by flow cytometry.

### Labeling

- The recommended minimum cell input for intracellular cell surface labeling is 100,000 cells. Note that at this input, recovering the maximum number of cells (20,000) may not be possible.

### Centrifugation & Pellet Resuspension

- Use a swinging-bucket rotor for higher cell/nuclei recovery.
- When working with samples with low cell numbers (i.e. <50,000 cells), complete removal of the supernatant is not required. Up to 30 µl supernatant may be left behind to optimize cell recovery without significantly impacting assay performance.



### Sample Fixation

- Consult Fixation of Cells & Nuclei for GEM-X Flex Gene Expression Demonstrated Protocol (CG000782) regarding fixing single cell suspensions following intracellular protein labeling.

### Sample Storage After Brief Fixation


- 20 minute fixed samples can be stored prior to intracellular protein labeling at  $-80^{\circ}\text{C}$  for up to 12 months with appropriate storage reagents.
- Sample storage buffer and post-storage guidelines for brief fixation are provided in the Appendix.

### Sample Storage After Overnight Fixation

- Samples that have completed intracellular protein labeling and subsequent fixation as described in Fixation of Cells & Nuclei for GEM-X Flex Gene Expression Demonstrated Protocol (CG000782) can be stored either at  $-80^{\circ}\text{C}$  for up to 12 months or at  $4^{\circ}\text{C}$  for up to 7 days with appropriate storage reagents.

## Intracellular Protein Labeling: Fix Samples

### Fix Sample




#### a. Optional Step Cell Surface Protein Labeling

Skip this step and directly proceed to step b if not utilizing the cell surface protein labeling section in this document.

- Perform cell surface protein labeling and washing as described in the Cell Surface Protein Labeling section in this document.
- Centrifuge labeled & washed sample at 300-400 rcf for 5 min (PBMCs/cell lines) at 4°C.
- Remove the supernatant without disturbing the pellet.
- Proceed immediately to step b.

**b. Add 0.5 ml room temperature Fixation Buffer B** to the sample (Fc blocked or cell surface protein labeled) and pipette mix 5x.

**c. Incubate for 20 min at room temperature (20°C).**

 **DO NOT** agitate or mix the sample during incubation. To minimize variability for room temperature fixations, incubation at controlled 20°C temperature (e.g. with a thermomixer, heat block, or water bath) is recommended.


**d. Add 0.5 ml** room temperature Additive C to the sample in Fixation Buffer B and pipette mix 5x.

**e. Centrifuge at 2000 rcf for 5 min at room temperature.**

**f. Remove the supernatant** without disturbing the pellet. Up to 30 µl supernatant may be left behind if working with <300,000 cells.


**g. Add 1 ml chilled Quenching Buffer B** to the sample pellet and pipette mix 5x and keep on ice.

**h. Proceed to the next step (Label Samples) or store the sample** after resuspending in appropriate reagents.

 Briefly fixed samples can be stored at -80°C for up to 12 months. See Appendix for guidance on storage and post-storage processing..



## Intracellular Protein Labeling: Label Samples



- Centrifuge at **2000 rcf** for **5 min** at **4°C**.
- Remove the supernatant without disturbing the pellet.
- Add **100 µl** Intracellular Antibody Mix Supernatant (antibody + buffer) to the pellet and gently pipette mix 10x.
- Incubate for **30 min** at **4°C**. If using fluorophore-conjugated antibodies, incubate without light exposure.
- Add **1 ml** chilled Intracellular Wash Buffer to the labeled cells. Gently pipette mix.
- Centrifuge cells at **2000 rcf** for **5 min** at **4°C**.
- Remove the supernatant without touching the bottom of the tube to avoid dislodging the pellet. Up to 30 µl supernatant may be left behind if working with <300,000 cells.

### OPTIONAL: For enrichment of cells by flow sorting:

- Proceed to flow sorting. See [Post-Fixation Flow Sorting](#) for details. For collection of sorted cells, use RNase-free reagents and supplement all buffers with RNase inhibitors.
  - Flow sorting can also be performed after step i (second fixation following the Demonstrated Protocol CG000782). See [Post-Fixation Flow Sorting](#) for additional information on collection.
- Repeat** e-g one more time for a total of 2 washes.
  - Proceed **immediately** to **Sample Fixation, Step c** of the Fixation of Cells & Nuclei for GEM-X Flex Gene Expression Demonstrated Protocol (CG000782).

## Intracellular Protein Labeling – Troubleshooting

Problem	Solution
Low cell viability when beginning the protocol	Perform dead cell removal (flow sorting or bead based)
High background	<p>Remove as much residual buffer as possible during the centrifugation steps, leaving no more than 30 <math>\mu</math>l buffer</p> <p>Add single-stranded DNA binding protein (SSB) from bacteriophage T4 (T4 Gene 32 Protein) to the antibody pool + antibody buffer if the background signal is an issue. Use 8 <math>\mu</math>g SSB per <math>\mu</math>g of antibody and add <math>MgCl_2</math> at 4 mM final concentration. It is recommended to use flow cytometry to ensure that the SSB addition is not leading to a reduction in antibody signal.</p>
Poor intracellular protein labeling performance	<p>Modify fixation time ranging from 15 to 30 min; 1 h fixation may be compatible as well, but will require testing with flow cytometry. Lowering the fixation time may increase clumping while increasing the time may reduce the antibody signal.</p> <p>Modify NP-40 concentration ranging from 0.25% to 1%. Note that <math>\leq 0.25\%</math> NP-40 may not be enough for permeabilization, while <math>&gt; 0.5\%</math> may induce clumping.</p> <p>Determine optimal antibody concentration using flow cytometry</p> <p>Remove dextran sulfate in the blocking buffer, substitute with BSA</p> <p>Increase antibody incubation time to <math>&gt; 30</math> min</p>
Poor Flex assay performance	<p>Ensure that all buffers are made fresh and RNase-free</p> <p>Ensure RNase inhibitors have been added to the Antibody Buffer and the Intracellular Wash Buffer</p>
Cell clumping	<p>Perform dead cell removal (flow sorting or bead based)</p> <p>Increase fixation time to up to 1 hr (testing antibody compatibility with flow cytometry is recommended)</p> <p>Reduce NP-40 concentration to <math>&lt; 0.5\%</math>. Note that <math>\leq 0.25\%</math> NP-40 may not be enough for permeabilization. Optimization and testing with flow cytometry is recommended.</p> <p>Perform additional filtering steps, including before hybridization</p>
Reads derived from antibodies targeting cell surface antigens are a high percentage of the antibody library	<p>Titrate surface antibodies to lower the amount of antibody added</p> <p>Consult Quality Control of Cell Surface Protein Labeling using Flow Cytometry Technical Note (CG000231) for guidance on how to titrate antibodies using flow cytometry for determining optimal antibody concentrations.</p>

## Appendix

### Flow Sorting Guidance

Enrich labeled cells using flow sorting prior to library generation to enable identification of rare subpopulations.

#### Pre-Fixation Flow Sorting

It is recommended to collect enriched cells in up to 20% FBS to maintain cell viability. Cells should be collected either in 20  $\mu$ l volume in the collection tube/plate (96-well plate) or in 300  $\mu$ l volume in a 5-ml tube. Use an optimal buffer for fragile cells to maintain a high cell viability.

Sort stream should be adjusted so that the cell-droplet falls into the collection buffer. Sorted cells should be counted and viability measured before proceeding to the Fixation of Cells & Nuclei for GEM-X Flex Gene Expression Demonstrated Protocol (CG000782). Cell loss during flow sorting is common. Optimize the protocol steps accordingly.



Once sorting is complete, proceed **immediately** to **Sample Fixation, Step a** of the Fixation of Cells & Nuclei for GEM-X Flex Gene Expression Demonstrated Protocol (CG000782).

#### Post-Fixation Flow Sorting

Post fixation samples can be flow sorted for advanced sample clean-up, as well as enrichment of specific populations. For samples undergoing intracellular protein labeling, sorting can be performed either after the first wash post labeling or after second fixation using the guidance provided in Fixation of Cells & Nuclei for GEM-X Flex Gene Expression Demonstrated Protocol (CG000782). Samples should be collected in PBS + 1% nuclease-free BSA supplemented with RNase inhibitor (Protector RNase inhibitor from Sigma, PN-3335399001). A final concentration of 0.2 U/ $\mu$ l RNase inhibitor is recommended.

Consult 10x Genomics [support website](#) for more information on postfixation cell sorting.

### Overnight-Fixed Sample Storage Guidance

Following intracellular antibody labeling and fixation, follow the storage guidance in the Fixation of Cells & Nuclei for GEM-X Flex Gene Expression Demonstrated Protocol (CG000782).


#### Briefly-Fixed Sample Storage Guidance

Samples resuspended in Quenching Buffer B following brief, 20 minute fixation (see [page 16](#)) prior to intracellular labeling can be stored at **-80°C** for up to **12 months**.

- a. Add 1 U/ $\mu$ l Protector RNase inhibitor to quenched cells and pipette mix.
- b. Add 50% glycerol for a final concentration of 10%. For example: if cells are resuspended in **1 ml** Quenching Buffer B with **25  $\mu$ l** of RNase inhibitor (stock concentration 40 U/ $\mu$ l), add 256.25  $\mu$ l 50% glycerol. Pipette mix. Alternatively, to conserve RNase inhibitor, centrifuge cells at 2000 rcf, remove 500  $\mu$ l Quenching Buffer B, and add 12.5  $\mu$ l RNase inhibitor (stock concentration 40 U/ $\mu$ l), then add 128.13  $\mu$ l of 50% glycerol.
- c. Store at **-80°C** for up to **12 months**.


### Supplemental Data on Different Wash Options – Cell Surface Protein Labeling

The Cell Surface Protein Labeling protocol provides three wash options after incubation with the antibody. Using the appropriate wash option is critical for obtaining high-quality data and minimizing background signal. Supplemental data on the different wash options is provided in the following sections.



**Figure 1.** t-SNE plots of human PBMCs comparing staining for T cell marker CD3 across differing number of washes following antibody incubation (left). Corresponding CD3E gene expression measured with the GEM-X Flex Gene Expression assay. Scale adjusted for better visualization with top 1% of cells removed. Violin plot showing the expression of the anti-CD3 antibody Feature Barcode (right).

Supplemental Data on Different Wash Options – Cell Surface Protein Labeling *contd.*



**Figure 2.** t-SNE plots of human PBMCs comparing staining for Natural Killer and T cell marker TIGIT across differing number of washes following antibody incubation (left). Corresponding TIGIT gene expression measured with the GEM-X Flex Gene Expression assay. Scale adjusted for better visualization with top 1% of cells removed for each feature. Violin plot showing the expression of the anti-TIGIT antibody Feature Barcode (right).

**Supplemental Data on Different Wash Options – Cell Surface Protein Labeling** *contd.*

Antibody	Comparison	TotalSeq™-C		
		Signal to Noise* (On-target cells / all other cells)		
		No-Wash	1-Wash	2-Wash
anti-CD3	T cells / Other	12.44	15.72	16.99
anti-CD4	T cells / Other	9.92	10.71	9.36
anti-CD8	T cells / Other	3.2	5.81	6.32
anti-TIGIT	T & NK cells / Other	1.14	1.4	1.54
anti-CD19	B cells / Other	27.11	39.67	45.61
anti-CD15	Monocytes / Other	2.9	4.51	4.61
anti-CD14	Monocytes / Other	6.62	7.4	8.26
anti-CD11c	DCs & Monocytes / Other	13.5	12.89	13.07

\*Ratio of mean expression level

**Table 1.** Signal to noise was generated by calculating the ratio of counts from a target group comprised of cell type(s) expected to display antibody signal (e.g. CD3 on T cells) to counts from a background group (e.g. CD3 on non-T cells).

### Conclusion on Different Wash Options – Cell Surface Protein Labeling

Choosing an appropriate wash protocol following cell surface protein labeling is critical for experimental success in the GEM-X Flex Gene Expression assay.

The data presented in this Demonstrated Protocol show that a lower number of washes following cell labeling can be used for antibodies with distinct positive and negative populations (Figures 1A and 1B, CD3 data). However, a lower number of washes cause a reduction in separation between the positive and negative populations (Figures 1A and 1B, violin plot; Table 1, lower signal to noise for CD3, CD14, CD19, etc.). Feature Barcode data from antibodies with poorer separation between the positive and negative populations is adequate with the 1-Wash protocol, but the No-Wash protocol is not typically recommended due to poor signal to noise (Figures 2A and 2B, TIGIT; Table 1, TIGIT).

### References

1. Fixation of Cells & Nuclei for GEM-X Flex Gene Expression Demonstrated Protocol (CG000782)
2. Cell Surface Protein Labeling for Single Cell RNA Sequencing Protocols with Feature Barcode technology Demonstrated Protocol (CG000149)
3. GEM-X Flex Gene Expression - Protocol Planner (CG000790)
4. GEM-X Flex Gene Expression Reagent Kits for Singleplexed Samples with Feature Barcode technology for Cell Surface Protein User Guide (CG000788)
5. GEM-X Flex Gene Expression Reagent Kits for Multiplexed Samples with Feature Barcode technology for Protein using Barcode Oligo Capture User Guide (CG000789)

Take 1 minute to evaluate this protocol. Scan this code or [click here](#).



## Document Revision Summary

<b>Document Number</b>	CG000781
<b>Title</b>	Cell Surface & Intracellular Protein Labeling for GEM-X Flex Gene Expression Demonstrated Protocol
<b>Revision</b>	Rev A
<b>Revision Date</b>	October 2024

© 2024 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: [www.10xgenomics.com/trademarks](http://www.10xgenomics.com/trademarks). 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: [www.10xgenomics.com/patents](http://www.10xgenomics.com/patents). All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of 10x Genomics products in practicing the methods set forth herein has not been validated by 10x Genomics, and such non-validated use is NOT COVERED BY 10X GENOMICS STANDARD WARRANTY, AND 10X GENOMICS HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics terms and conditions of sale for the Chromium Controller or the Chromium Single Cell Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics that it currently or will at any time in the future offer or in any way support any application set forth herein.

**Contact:**[support@10xgenomics.com](mailto:support@10xgenomics.com)

10x Genomics

6230 Stoneridge Mall Road

Pleasanton, CA 94588 USA

