

Chromium X Series

Firmware Update Procedure

For use with:

Chromium X & Accessory Kit
12-Month Warranty, PN-1000331
24-Month Warranty, PN-1000332

Chromium Xo & Accessory Kit
12-Month Warranty, PN-1000812

Chromium iX & Accessory Kit
12-Month Warranty, PN-1000328
24-Month Warranty, PN-1000329


Chromium X Series Refurbished
PN-1000326R | PN-1000818R

Consult the Chromium X Series Instrument User Guide with Readiness Test (CG000396) for additional details. Automatic downloads of updates are enabled after firmware updates to ensure the instrument is up-to-date with the latest assay products, features, and bug fixes.

Step 1: Download and install the update


Procedure with an instrument connected to the internet and auto-download enabled (recommended):


- A pop-up on the home screen will alert the user of available updates.
- Open pop-up and click to install.
- Once the update is complete, click the “Restart Instrument” button.

 **DO NOT skip this button and use the power switch at the back of the unit; this could result in lost data.**

Procedure with an instrument connected to the internet, auto-download disabled:

- A pop-up on the home screen will alert the user of available updates.
- Open pop-up and click to download.
- Once the update has been downloaded, install the update.
- Once the update is complete, click the “Restart Instrument” button.

 **DO NOT skip this button and use the power switch at the back of the unit; this could result in lost data.**



Step 2: Reinitialize the instrument

After installing an update, run a readiness test.

- Retrieve the shipping/readiness test block with attached gasket.
See the Readiness Test section of the Instrument User Guide (CG000396) for more information on the block.
- Initiate a readiness test through the home screen by loading the block into the instrument.
- After the test is complete, follow screen instruction to remove the block from the tray.
- Store the assembled block and gasket in a plastic bag, away from sun-light and extreme temperature. Retain for future use.
- The instrument is now ready for use.

Chromium X Series

Document Revision Summary



Document Number

CG000534 | Rev D

Title

Chromium X Series Firmware Update Procedure

Revision

Rev C to Rev D

Revision Date

February 2025

Description of Changes

- References to updating the firmware now specify clicking the “Restart Instrument” button to clarify the process

Contact:

support@10xgenomics.com

10x Genomics

6230 Stoneridge Mall Road

Pleasanton, CA 94588 USA

© 2025 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: www.10xgenomics.com/trademarks. 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: www.10xgenomics.com/patents. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of 10x Genomics products in practicing the methods set forth herein has not been validated by 10x Genomics, and such non-validated use is NOT COVERED BY 10X GENOMICS STANDARD WARRANTY, AND 10X GENOMICS HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics terms and conditions of sale for the Chromium Controller or the Chromium Single Cell Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics that it currently or will at any time in the future offer or in any way support any application set forth herein.