

Quick Reference Cards | CG000730 | Rev C

# Visium Cassette S3

for Assembly & Disassembly

For use with:

Visium Slide Cassettes S3, 6.5 mm, 2pk, PN-1000847

*contains: Visium 2-port, 6.5 mm*

Visium HD Cassettes, 6.5 mm, 4 rxns, PN-1000669

*contains: Visium 2-Port Cassette S3, 6.5 mm; Visium Tissue Slide Cassette S3, 6.5 mm*

Visium 8-port Cassette S3, 4 pk, PN-1000685

Visium Tissue Slide Cassette S3, 4 pk, PN-1000684

Visium HD Cassettes, 11 mm, 2 rxns, PN-1000849

*contains: Visium 2-Port Cassette S3, 11 mm; Visium Tissue Slide Cassette S3, 11 mm*

Visium HD Tissue Slide Cassette, 11 mm, 2 rxns, PN-1001030

# Visium Cassette S3

## Document Revision Summary



### Document Number

CG000730 | Rev C

### Title

Visium Quick Reference Cards – Visium Cassette S3 for Assembly & Disassembly

### Revision

Rev B to Rev C

### Revision Date

March 2026

### Description of Changes

- Added 11 mm-sized cassette kits to the “For use with” section on the cover
- Specified 6.5 mm or 11 mm in the instructions where necessary
- Added 11 mm-sized versions of illustrations throughout
- Clarified "Install the Gasket" section on page 4

#### Contact:

[support@10xgenomics.com](mailto:support@10xgenomics.com)

10x Genomics, Inc.

6230 Stoneridge Mall Road

Pleasanton, CA 94588 USA

© 2026 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: [www.10xgenomics.com/trademarks](http://www.10xgenomics.com/trademarks). 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: [www.10xgenomics.com/patents](http://www.10xgenomics.com/patents). All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of 10x Genomics products in practicing the methods set forth herein has not been validated by 10x Genomics, and such non-validated use is NOT COVERED BY 10X GENOMICS STANDARD WARRANTY, AND 10X GENOMICS HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics terms and conditions of sale for the Chromium Controller or the Chromium Single Cell Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics that it currently or will at any time in the future offer or in any way support any application set forth herein.

# Visium Cassette S3

## Tissue Slide Cassette Assembly (1 of 3)

### Visium Tissue Slide Cassette S3, 6.5 mm

Found in: PN-1000669, PN-1000684

### Visium Tissue Slide Cassette S3, 11 mm

Found in: PN-1000849, PN-1001030

The Visium Tissue Slide Cassette S3 is a single-use item used to create wells on glass slides with tissue for reagent addition and removal.

It is distinct from the Visium 2-port Cassette S3, which encases a Visium HD Slide.

Gaskets are used to create a sealed well around tissue. They are combined with the complementary Visium Tissue Slide Cassette S3 components (cassette bottom, cassette top, translator for 6.5 mm) to create a complete cassette.

To ensure that the gasket surrounds the tissue area of interest, its position can be adjusted within the cassette during assembly (see image).

Tissue or area of interest should be centered within the gasket. Assembling against a white background may aid in tissue visualization. The tissue does not need to be smaller than the gasket; however, only tissue within the gasket will be analyzed.

Assemble slide into the cassette flat on a clean work surface. For information on tested slides and slide dimensions that are compatible with the gasket, consult the Visium HD Protocol Planner (CG000698).



Practice assembly with a blank slide (75 x 25 x 1 mm).

Take measures to reduce exposure of cassette, gasket, translator and slide to debris.

If placing the top half of the cassette on a surface, ensure the gasket faces away from the surface so it does not collect debris.

Place slides in the cassette only when specified.

Applying excessive force to the slide may cause the slide to break.




# Visium Cassette S3

## Tissue Slide Cassette Assembly (2 of 3)

### Inspection

- Ensure that the cassette, slide, gasket, and translator (for 6.5 mm) are free of debris both before and after assembly.
- In rare cases, the gasket may have raised pieces of silicone. These are considered excess and must be removed with tweezers or a pipette tip.  
*Failure to remove excess silicone may prevent adequate sealing of the cassette during heated incubation steps. Excess silicone may also lead to evaporation.*




### 6.5 mm Gasket-Translator Assembly

#### Correct assembly:

- Sides of the gasket are flush with the sides of the translator.
- Rails overlap on top of the gasket.
- Gasket slides smoothly on the tracks.

#### Incorrect assembly:


- One or both rails are to the side of the gasket rather than on top.
- There is a gap between the gasket and translator.
- There is an angle to one side of the translator.



### Gasket Installation


**6.5 mm gasket:** The opening of the translator must be oriented towards the side of the cassette, so that the gasket cannot slide out.

Ensure that the gasket does not get pushed out of the translator during cassette assembly (step 2), or when adjusting the gasket over tissue on the slide (step 4).



# Visium Cassette S3

## Tissue Slide Cassette Assembly (3 of 3)



### Visium Tissue Slide Cassette S3






Exercise caution when handling slide edges to prevent injury.  
Wear fresh gloves while assembling the cassette.  
Do not move gasket once cassette is closed, as this may damage the tissue.  
Do not touch the bottom of the gasket to avoid collecting debris.




Steps show a 6.5 mm cassette, but all instructions are universal, unless specified.

### Assembly

- Place tissue slide into lower half of cassette with tissue facing up.
- Snap pre-assembled gasket (and translator, for 6.5 mm) into place on the top half of cassette, pressing it in from the back until it clicks into place. If necessary, the gasket can be rotated 180°. Verify the gasket is still correctly placed.
- Gently place top half of cassette over bottom half. DO NOT assemble together until Step 6.






- Adjust gasket (and translator, for 6.5 mm) such that the gasket is over the tissue region of interest. Do not push down on the gasket, which could dislodge and damage the tissue (or push the gasket out of the translator, for 6.5 mm).
- Apply even pressure on top and bottom of cassette until it clicks shut. Verify that clips are completely secured over tab.
- Turn cassette over and verify that tissue area of interest is within gasket. DO NOT move gasket once cassette is closed. If necessary, open cassette and reposition gasket. *Repositioning the gasket may lead to tissue damage.*



### Removal

- Pull clip up to detach upper and lower halves of cassette.
- Open cassette by continuing to lift clip. If slide sticks to the gasket, gently hold down tissue slide near the label while lifting the top half.
- Hold the top of the slide and lift slide out from lower half.



# Visium Cassette S3

## Visium Cassette Assembly (1 of 3)

### Visium 2-port Cassette S3, 6.5 mm

Found in: PN-1000669, PN-1000847

### Visium 2-port Cassette S3, 11 mm

Found in: PN-1000849



### Visium 8-port Cassette S3

Found in: PN-1000685

The Visium 2-port Cassette S3 encases the Visium HD slide and creates leakproof wells for reagent addition. The Visium 8-port Cassette S3 is used during immunofluorescence staining optimization.

### Variations in Visium HD Cassette Features

**!** For some cassettes, the well is slightly larger towards the bottom right corner to allow space for pipetting without damaging the hydrogel. An arrow points to the correct corner for inserting the pipette.



# Visium Cassette S3

## Visium Cassette Assembly (2 of 3)


### Inspection

Cassette upper and lower halves must be from the same kit. Mixing kit parts may result in improper assembly and cassette leaking.

- Ensure that the cassette, gasket and slide are free of debris both before and after assembly.
- In rare cases, the gasket may have raised pieces of silicone. These are considered excess and must be removed with tweezers or a pipette tip (see page 4 for visual example).

*Failure to remove excess silicone may prevent adequate sealing of the cassette during heated incubation steps. Excess silicone may also lead to evaporation.*

Practice assembly with a blank slide (75 x 25 x 1 mm).

 If placing the top half of the cassette on a surface, ensure the gasket faces away from the surface so it does not collect debris.


Place slides in the cassette only when specified.

Applying excessive force to the slide may cause the slide to break.



# Visium Cassette S3

## Visium Cassette Assembly (3 of 3)






### Visium 2-port Cassette S3, Visium 8-port Cassette S3



Exercise caution when handling slide edges to prevent injury.  
 Wear fresh gloves while assembling the cassette.  
 Cassette top and bottom must be from the same kit. Mixing parts may result in improper assembly and cassette leaking.  
 Steps show a 6.5 mm 2-port cassette, but all instructions are universal.

### Assembly

- 1 Place upper and lower halves of cassettes and Visium slide on bench
- 2 Place Visium slide with active surface facing upwards into lower half of cassette.
- 3 Press slide down into grooves of the bottom half of the cassette until it sits firmly in place



- 4 Secure outer clips of top half with outer tab of bottom half
- 5 Press firmly on top and bottom of cassette until it clicks shut. Ensure that clips are completely secured over tabs.



### Removal

- 1 Pull inner clip up from inner tab to detach upper and lower halves of cassette
- 2 Open cassette by continuing to lift upper half upward
- 3 Lift slide out from lower half

