



Sample Prep User Guide | CG000505 | Rev B

# Chromium Nuclei Isolation Kit

For use with:

Chromium Nuclei Isolation Kit

16 rxns, PN-1000493

Chromium Nuclei Isolation Kit with RNase Inhibitor\*

16 rxns, PN-1000494

Nuclei Isolation Kit with RNase Inhibitor v2

16 rxns, PN-1001101

*\*Additional item required: RNase Inhibitor 40X Kit (Kit PN-1000887, Tube PN-2001488)*

# Notices

## Document Number

CG000505 | Rev B

## Legal Notices

© 2026 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: [www.10xgenomics.com/trademarks](http://www.10xgenomics.com/trademarks). 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: [www.10xgenomics.com/patents](http://www.10xgenomics.com/patents). The use of products described herein is subject to 10x Genomics Terms and Conditions of Sale, available at [www.10xgenomics.com/legal-notices](http://www.10xgenomics.com/legal-notices), or such other terms that have been agreed to in writing between 10x Genomics and user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES. Notice to Customer: The listed products and their use are the subject of United States Patent Nos. 6,159,736, and 6,294,385, European Patent No. 1115856 and related patents and patent applications licensed from the Wisconsin Alumni Research Foundation.

## Instrument & Licensed Software Updates Warranties

Updates to existing Instruments and Licensed Software may be required to enable customers to use new or existing products. In the event of an Instrument failure resulting from an update, such failed Instrument will be replaced or repaired in accordance with the 10x Limited Warranty, Assurance Plan or service agreement, only if such Instrument is covered by any of the foregoing at the time of such failure. Instruments not covered under a current 10x Limited Warranty, Assurance Plan or service agreement will not be replaced or repaired.

## Support

Email: [support@10xgenomics.com](mailto:support@10xgenomics.com)

10x Genomics

6230 Stoneridge Mall Road

Pleasanton, CA 94588 USA

# Document Revision Summary

## Document Number

CG000505 | Rev B

## Title

Chromium Nuclei Isolation Reagent Kits Sample Prep User Guide

## Revision

Rev A to Rev B

## Revision Date

April 2026

## Description of Changes

- Added information regarding Nuclei Isolation Kit with RNase Inhibitor v2 on pages 1, 9, 26, & 34, along with supported products on page 9
- Updated Lysis Buffer on pages 28 and 36
- Added additional Troubleshooting guidance on pages 53, 54
- Additional guidance on untested tissue types on page 56
- Added Cell/Nuclei Counting section in Appendix on page 59
- Updated for general minor consistency of format, language, and terms throughout

# Table of Contents

<b>Introduction</b>	<b>6</b>
Objective	7
Protocol Steps & Timing	8
Product Compatibility	9
Additional Kits, Reagents & Equipment	10
Recommended Pipette Tips	11
<b>Tips &amp; Best Practices</b>	<b>12</b>
<b>Protocol Selector</b>	<b>23</b>
<b>Chapter 1:</b>	
■ <b>Single Cell Universal Gene Expression &amp; Flex</b>	<b>25</b>
Chromium Nuclei Isolation Reagent Kits	26
Get Started	27
Buffer Preparation	28
Illustrative Overview	29
Nuclei Isolation Protocol Steps	30
<b>Chapter 2:</b>	
■ <b>Single Cell Epi Multiome ATAC + Gene Expression</b>	<b>33</b>
Chromium Nuclei Isolation Reagent Kits	34
Get Started	35
Buffer Preparation	36
Illustrative Overview	38
Nuclei Isolation Protocol Steps	39
<b>Chapter 3:</b>	
■ <b>Single Cell Epi ATAC</b>	<b>42</b>
Chromium Nuclei Isolation Reagent Kits	43
Get Started	44
Buffer Preparation	45
Illustrative Overview	47
Nuclei Isolation Protocol Steps	48
<b>Troubleshooting</b>	<b>51</b>
Troubleshooting Guide - Issues, Potential Causes & Mitigation Strategies	52

# Table of Contents

<b>Appendix</b>	<b>55</b>
Tested Tissue Types	56
Tissues with Variable Results	57
Optimization Recommendations	58
Cell/Nuclei Counting	60

# Introduction

Objective

Protocol Steps & Timing

Product Compatibility

Additional Reagents, Kits, & Equipment

Recommended Pipette Tips

## Objective


The Chromium Nuclei Isolation Kit is an all-in-one solution for the standardized isolation of nuclei from fresh and frozen tissue for use in 10x Genomics Single Cell assays. Tissue samples are homogenized with a pestle in Lysis Buffer and passed through a column. Next, debris is removed via centrifugation in Debris Removal Buffer. The isolated nuclei are then washed and resuspended and loaded directly into compatible 10x Genomics Single Cell assays.

The Chromium Nuclei Isolation Kit streamlines the nuclei isolation process into a single workflow, allowing for increased efficiency, scalability through sample batching, and reduced experimental variability using 10x Genomics pre-formulated reagents. The protocol is designed to be compatible with a wide variety of tissue types and sizes. Consult the Tested Tissue Types pages on the 10x Genomics Support website for more information.



This User Guide outlines the process for isolating nuclei from frozen tissues for use in compatible 10x Genomics Single Cell assays. Cell cultures and cell suspensions have been not tested and are not supported as inputs for the Chromium Nuclei Isolation Kit workflow. Refer to the [Product Compatibility](#) and [Protocol Selector](#) pages for additional information on choosing the appropriate nuclei isolation kit and protocol based on the intended downstream Single Cell assay

## High Level Overview of Nuclei Isolation Workflow



## Protocol Steps & Timing

The table below provides an overview of the nuclei isolation workflow steps and timing. This protocol is to be executed without any stopping points.

Steps	Timing
<b>Buffer Preparation</b>	10 min
<b>Nuclei Isolation Process</b>	
Tissue Dissociation	10 min
Nuclei Isolation & Cleanup	45 min
<b>Sample QC</b>	10 min



Keep samples on ice following QC. Proceed **immediately** to relevant 10x Genomics Single Cell protocol.

## Product Compatibility

The Chromium Nuclei Isolation Kit is compatible with 10x Genomics Single Cell assays. This table outlines the supported products for the Chromium Nuclei Isolation Kits, as well as the reagents and consumables included in each kit. Consult the [Protocol Selector](#) page of this document for additional information regarding kit and product selection.

Chromium Nuclei Isolation Kit	Kit Components	Supported Products (User Guides)*
<p>Nuclei Isolation Kit with RNase Inhibitor v2 (PN-1001101)</p> <p><i>Includes PN-1000447, PN-1000448, &amp; PN-1001099</i></p> <p>OR</p> <p>Chromium Nuclei Isolation Kit with RNase Inhibitor (PN-1000494)</p> <p><i>Includes PN-1000447, PN-1000448, PN-1000450 &amp; PN-1000449</i></p>	<ul style="list-style-type: none"> <li>Chromium Nuclei Isolation Reagents, 16 rxns (PN-1000447)</li> <li>Chromium Nuclei Isolation Consumables, 16 rxns (PN-1000448)</li> <li>Reducing Agent B (PN-1000450)</li> <li>RNase Inhibitor Kit (PN-1000449)</li> <li>RNase Inhibitor Kit v2 (PN-1001099)</li> </ul> <p><i>If using PN-1000494, additional item required: RNase Inhibitor 40X Kit (Kit PN-1000887, Tube PN-2001488)</i></p>	<ul style="list-style-type: none"> <li>Single Cell Universal 3' Gene Expression (CG000731, CG000768, CG000416, CG000315, CG000204, CG000399)</li> <li>Single Cell Universal 5' Gene Expression (CG000733, CG000770, CG000423, CG000331)</li> <li>Single Cell Epi Multiome ATAC + Gene Expression (CG001689, CG000338)</li> <li>Flex/Fixed RNA Profiling (CG000834, CG000527, CG000786, CG000787, CG000691)</li> </ul>
Chromium Nuclei Isolation Kit (PN-1000493)	<ul style="list-style-type: none"> <li>Chromium Nuclei Isolation Reagents, 16 rxns (PN-1000447)</li> <li>Chromium Nuclei Isolation Consumables, 16 rxns (PN-1000448)</li> <li>Reducing Agent B (PN-1000450)</li> </ul>	<ul style="list-style-type: none"> <li>Single Cell Epi ATAC (CG000496)</li> </ul>

*Some reagent tubes (same part number) are available in multiple kits and may be used interchangeably.*

\*Consult the 10x Genomics Support Website for the most current information

Nuclei isolation removes outer cell membranes and cytoplasmic compartments. As a result the workflow cannot be used for:

- Feature Barcode technology for Cell Surface Protein
- Feature Barcode technology for CRISPR Screening may be challenging due to low sgRNA in the nucleus and is not supported
- V(D)J capture from nucleic transcripts is not currently supported

## Additional Kits, Reagents & Equipment

The items in the table below have been validated by 10x Genomics and are highly recommended for the Nuclei Isolation protocol. Substituting materials may adversely affect system performance. This list does not include standard laboratory equipment, such as water baths, centrifuges, vortex mixers, pH meters, freezers, etc.

Item	Description	Supplier	Part Number
<b>Plastics</b>			
2-ml Tubes	DNA LoBind Tubes 2.0 ml	Eppendorf	022431048
15-ml Tubes	Corning 15 ml centrifuge tubes	Corning	CLS430791
50-ml Tubes	Corning 50 ml centrifuge tubes	Corning	CLS430829
<b>Kits &amp; Reagents</b>			
10% BSA	Bovine Serum Albumin in DPBS (10%) <i>(alternatively, use MACS BSA Stock Solution)</i>	Millipore-Sigma Miltenyi Biotec	A1595 130-091-376
1X PBS	Phosphate-Buffered Saline without Calcium & Magnesium	Corning	21-040-CV
Nuclease-free Water	Molecular Grade Nuclease-free Water	Thermo Fisher Scientific	AM9937
<b>Cell Counting</b> <i>See Appendix for a comprehensive list</i>			
Nucleic Acid Staining Fluorescent Dye	VitaStain AOPI Staining Solution <i>(alternatively, use Ethidium Homodimer-1)</i>	Nexcelom Thermo Fisher Scientific	CS2-0106-5ml E1169
Cell Counter	Cellaca MX High-throughput Automated Cell Counter <i>(alternatively, use any cell counter with fluorescent capabilities)</i>	Nexcelom	MX-112-0127
<b>Equipment</b>			
Vortex	Vortex Mixer	VWR	10153-838
Centrifuge	Refrigerated Eppendorf Centrifuge <i>(alternatively, use any equivalent centrifuge)</i>	Millipore-Sigma	5427R or 5424R

## Recommended Pipette Tips

10x Genomics recommends using only validated emulsion-safe pipette tips for all Single Cell protocols, including the Chromium Nuclei Isolation protocol. Rainin pipette tips have been extensively validated by 10x Genomics and are highly recommended for all single cell assays. If Rainin tips are unavailable, any of the listed alternate pipette tips validated by 10x Genomics may be used.

Supplier	Description	Part Number (US)
<b>Recommended Pipettes &amp; Pipette tips</b>		
Rainin (pipettes)	Pipet-Lite Multi Pipette L8-50XLS+	17013804
	Pipet-Lite Multi Pipette L8-200XLS+	17013805
	Pipet-Lite Multi Pipette L8-10XLS+	17013802
	Pipet-Lite Multi Pipette L8-20XLS+	17013803
	Pipet-Lite LTS Pipette L-2XLS+	17014393
	Pipet-Lite LTS Pipette L-10XLS+	17014388
	Pipet-Lite LTS Pipette L-20XLS+	17014392
	Pipet-Lite LTS Pipette L-100XLS+	17014384
	Pipet-Lite LTS Pipette L-200XLS+	17014391
	Pipet-Lite LTS Pipette L-1000XLS+	17014382
Rainin (pipette tips)	Tips LTS 200UL Filter RT-L200FLR	30389240
	Tips LTS 1ML Filter RT-L1000FLR	30389213
	Tips LTS 20UL Filter RT-L10FLR	30389226
<b>Alternate Recommendations</b> (If Rainin pipette tips are unavailable, any of the listed pipette tips may be used)		
Eppendorf (pipettes)	Eppendorf Research plus, 8-channel, epT.I.P.S. Box, 0.5 – 10 µL	3125000010
	Eppendorf Research plus, 8-channel, epT.I.P.S. Box, 10 – 100 µL	3125000036
	Eppendorf Research plus, 8-channel, epT.I.P.S. Box, 100 – 300 µL	3125000052
	Eppendorf Research plus, 1-channel, epT.I.P.S.® Box, 0.1 – 2.5 µL	3123000012
	Eppendorf Research plus, 1-channel, epT.I.P.S.® Box, 0.5 – 10 µL	3123000020
	Eppendorf Research plus, 1-channel, epT.I.P.S.® Box, 2 – 20 µL	3123000039
	Eppendorf Research plus, 1-channel, epT.I.P.S.® Box, 2 – 200 µL	3123000055
	Eppendorf Research plus, 1-channel, epT.I.P.S.® Box, 100 – 1000 µL	3123000063
Eppendorf (pipette tips) <i>Compatible with Eppendorf pipettes only</i>	ep Dualfilter T.I.P.S., 2-20 µL	0030078535
	ep Dualfilter T.I.P.S., 2-200 µL	0030078551
	ep Dualfilter T.I.P.S., 2-1,000 µL	0030078578
Labcon*	ZAP SLIK 20 µL Low Retention Aerosol Filter Pipet Tips for Rainin LTS	4-1143-965-008
	ZAP SLIK 200 µL Low Retention Aerosol Filter Pipet Tips for Rainin LTS	4-1144-965-008
	ZAP SLIK 1000 µL Low Retention Aerosol Filter Pipet Tips for Rainin LTS	4-1145-965-008
Biotix*	xTIP4 Racked Pipette Tips, Rainin LTS Pipette Compatible, 0.1-20uL	63300931
	xTIP4 Racked Pipette Tips, Rainin LTS Pipette Compatible, 200uL	63300001
	xTIP4 Racked Pipette Tips, Rainin LTS Pipette Compatible, 1000uL	63300003

\*Compatible with Rainin pipettes

# Tips & Best Practices



## Icons



Tips & Best Practices section includes additional guidance



Signifies critical step requiring accurate execution



Troubleshooting section includes additional guidance

## Plastics

- Use recommended plastic consumables when isolating nuclei as some plastics can introduce fibers into reagents, buffers, and solutions, leading to microfluidic failures.

## General Reagent Handling

- Fully thaw and thoroughly mix reagents before use.
- Keep all reagents on ice during setup and use. Promptly move reagents back to the recommended storage after use.
- Note that an excess of 10% of 1 reaction value is calculated for all buffer reagent volumes in this User Guide.
- Thoroughly mix samples before each step.
- If provided Lysis Reagent and Debris Removal Buffers appear cloudy or contain precipitate, warm the tubes to **40°C** and swirl until the buffers become clear again.


## Pipette Tips & Calibration

- Follow manufacturer's calibration and maintenance schedules.
- Use only recommended pipette tips.

## Experimental Setup

- Minimize exposure of reagents, columns, and tubes to sources of particles and fibers, laboratory wipes, frequently opened flip-cap tubes, clothing that sheds fibers, and dusty surfaces.
- Execute steps without pause or delay.
- Pre-chill centrifuge to **4°C** before starting protocol. Ensure centrifuge is balanced prior to loading samples.
- If using a fixed angle rotor centrifuge, orient tube caps inward to avoid breaking of caps during centrifugation steps. Alternatively, use empty balance tubes to pin the collection tube cap as indicated in the image below.

### Orient Tube Caps Facing Inward During Centrifugation



\*Caps of tubes are angled toward the balance tubes to prevent them from turning outward due to the centrifugal force.

- Label tops and sides of all tubes, as well as tops of spin columns, for clear identification of samples.
- If centrifuge model does not allow for **20 sec** spin during first centrifugation step of Nuclei Isolation Protocol, set it to **30 sec** and stop spin after **20 sec**.

## Tissue Handling & Storage

- Work quickly and minimize handling during all tissue processing steps.
- Wash tissues in a clean glass petri dish with cold PBS upon harvest and absorb excess blood using a laboratory wipe.
- Cut tissues into small pieces (i.e. the size of a rice grain) for ease of freezing and place in a cryovial.
- To flash freeze, either submerge the cryovial in liquid nitrogen or a liquid-nitrogen cooled bath (e.g. isopentane) or place the tube deep in a bucket of dry ice. Wait at least 2–3 minutes for the tissue to freeze all the way through, and transfer the tube containing the tissue to vapor phase liquid nitrogen for long-term storage.
- Tissues should be stored **long-term** in a cryovial in **liquid nitrogen** for best results. Tissues can be stored **short-term** (1–2 days) at **-80°C** if needed.
- Once removed from liquid nitrogen, tissues should be maintained at **-80°C** or on **dry ice** until use.




*Thawing tissue prior to dissociation is NOT recommended as this can result in RNA degradation.*

## Tissue Dissociation

- Weigh input tissues before proceeding with tissue dissociation. The Chromium Nuclei Isolation Kit requires input tissue masses between **3–50 mg**.




Start with sufficient material. Refer to *Nuclei Concentration for Optimal Performance* and *Nuclei Recovery* sections for guidance on nuclei recovery.




- If tissue mass is **>50 mg**, cut tissue into smaller sections on a dry ice-cooled glass petri dish and place each section into a Sample Dissociation Tube. The resulting nuclei from each isolation can be combined at the point of resuspension immediately prior to counting and loading of the downstream Single Cell assay.
- For difficult to dissociate, fibrous tissues, pre-chopping tissue into smaller pieces (**≤10 mg**) will aid in dissociation and increase nuclei recovery.
- Dissociate the tissue in the Lysis Buffer by grasping the pestle between thumb and fingers. Twist the pestle continuously within the tube. Alternate between pushing the pestle down and lifting the pestle up to allow the tissue to slide between the wall of the tube and the pestle.
- Continue dissociating the tissue until the lysate appears homogeneous and no large pieces of tissue remain.
- Use only the pestles provided in this kit for tissue dissociation steps.
- If pestles are accidentally dropped or discarded, they may be washed with 70–80% ethanol before usage.

## Tissue Dissociation

Contd.




- Note that for many sample types there may be a small amount of white or opaque particulate that remain in the tube following tissue dissociation. These debris may be left behind when transferring the lysate to the spin column in subsequent steps.



- Tissue dissociation is sample dependent. Soft tissues may only require a few pestle strokes for dissociation, while more fibrous tissue will require more pestle strokes to fully dissociate tissue. No large tissue chunks should remain after tissue dissociation and tissue homogenate should be able to easily pass through a P1000 pipette tip without clogging.

## Lysis Conditions

- Nuclei lysis should be carried out on ice using chilled reagents and tools.
- Overlysis of nuclei occurs when samples are suspended in Lysis Buffer for an extended period of time. Avoid overlysis of nuclei which can lead to leakage of nuclear content and high levels of background signal.
- Nuclei with an intact membrane will appear round and smooth, while nuclei with a compromised membrane will appear “ruffled”, an indication of blebbing.
- The Chromium Nuclei Isolation Kit lysis conditions have been shown to be compatible across a wide variety of tissue types. Follow the protocol as indicated and QC nuclei before attempting to perform additional lysis condition optimization.




## Debris & Aggregate Removal

- A high-quality nuclei suspension with minimal aggregates and debris is critical for single cell sequencing.
- For high debris tissues (e.g. brain), starting with smaller tissue input (<25 mg) will reduce levels of debris and background in final nuclei suspension.
- Aggregates can be caused by insufficient resuspension of nuclei between steps of this protocol. Vortex nuclei samples where indicated during resuspension and wash steps.
- Debris and clumping can also be a sign of overlysis. Lysis optimization may be needed if vortexing steps were followed and clumping is not resolved.



Some level of small debris is acceptable during the nuclei isolation process and is not expected to negatively impact data quality.



## Nuclei Visualization & Counting

See Appendix for additional information


- After the isolation process, visualize nuclei suspensions to determine nuclei concentration and viability, suspension quality, and nuclei sizes prior to use in 10x Genomics Single Cell protocols.
- Counting with a fluorescent nucleic acid staining dye (e.g. AOPI or Ethidium Homodimer-1) and a fluorescent capable automated counter or microscope is strongly recommended as use of Trypan Blue can lead to overestimated nuclei counts.
- Vortex samples for **3 sec** immediately prior to counting to ensure an accurate estimate of the final nuclei concentration.
- Count nuclei in replicates (2–3 reproducible counts) to ensure accuracy.

## Nuclei Concentration for Optimal Performance

- The total number of suspended nuclei used as input to 10x Genomics Single Cell protocols is determined by the nuclei recovery target. Consult the applicable 10x Genomics Single Cell protocol to determine these relationships.
- The optimal input nuclei concentration for most 10x Genomics Single Cell assays is **700–1,200 nuclei/μl**. For the Chromium Fixed RNA Profiling assay, a starting number of **1 x 10<sup>6</sup>** nuclei is recommended. Refer to the relevant User Guide for specific information.
- If possible, bring the input nuclei suspension to a concentration that is optimal for the dynamic range of counting technique used (manual or automated), allows for 2–3 reproducible counts (where the standard deviation is <25%), and requires pipetting **2.5–15 μl** nuclei suspension into the Single Cell Master Mix depending on the assay. Consult the downstream workflow user guides for more information on nuclei input.
- Higher nuclei stock concentrations will result in lower pipetting volumes that may increase nuclei input variability.
- Lower nuclei stock concentrations will result in inaccurate nuclei counts that may also increase nuclei input variability.
- Use final single nuclei suspension to estimate the number of input nuclei since nuclei are inevitably lost during washing and resuspension steps. See Appendix for Cell/Nuclei Counting details.

## Nuclei Recovery

- If using a limited sample and nuclei recovery is expected to be low, a single wash may be used in an attempt to improve recovery.
- In addition to a single wash, centrifugation time in the debris removal buffer and/or wash buffer may be extended to improve nuclei yield.
- The use of swinging-bucket rotors, which pellet nuclei to the bottom of tubes rather than the side, may also improve final nuclei recovery.



- The above graph shows that nuclei yield is impacted by starting tissue sample type. Nuclei recovery from healthy tissues typically ranges from **5,000–15,000 nuclei/mg**. Both disease state and cell type composition can impact nuclei recovery.
- Resuspend nuclei in a low initial volume (**50–200  $\mu$ l**) during final resuspension step if nuclei yield is expected to be low or is unknown.
- Do not resuspend nuclei in **<50  $\mu$ l**. Refer to specific workflow chapters for exceptions.

## Nuclei Control Sample

- The use of a quality control sample is recommended when using this protocol for the first time or testing new tissues types.
- Mouse brain tissue from BioIVT is a validated quality control sample and can be run in parallel with experimental samples. Additional validated tissues can be found in the list of Tested Tissues Types in the Appendix of this document.

## Nuclei Storage



- Freezing or cryopreservation of nuclei following isolation is NOT recommended as the freezing and thawing process can damage the nuclear membrane and lead to loss of single-nuclei behavior. Additionally, frozen nuclei have the potential to burst and increase background signal.
- Cryopreserve at the tissue stage for best results in 10x Genomics Single Cell sequencing protocols.
- Nuclei can be stored after fixation for Chromium Fixed RNA Profiling applications for up to 1 week at **4°C** or for up to 3 months at **-20°C or -80°C**.



# Protocol Selector

## Protocol Selector

The Nuclei Isolation protocol differs based on your desired downstream 10x Genomics Single Cell assay. The Chromium Nuclei Isolation Kit with RNase Inhibitor is compatible with Single Cell 3' Gene Expression, Single Cell 5' Gene Expression, Chromium Fixed RNA Profiling, and Single Cell Multiome + ATAC Gene Expression assays. The Chromium Nuclei Isolation Kit (without RNase Inhibitor) is compatible with Single Cell ATAC assays.

To ensure quality nuclei isolation preparation, proceed to the appropriate chapter outlined in the diagram below. Each chapter outlines the Reagent Kits, Get Started guidance, Buffer Preparation, and Protocol Steps for nuclei isolation of the selected Single Cell assay. The chapters are color coded with a labeled tab at the top of each page for easy navigation.



Note that following the appropriate nuclei isolation protocol is critical for the success of the nuclei isolation process as well as all downstream Single Cell applications.

### Downstream Assay

Single Cell Universal 3' and 5' Gene Expression  
and Flex

Single Cell Epi Multiome

Single Cell Epi ATAC

### Chapter and Page #s

Chapter 1:

**Nuclei Isolation Protocol for Single Cell Universal Gene Expression & Flex**  
for use with PN-1000494 or PN-1001101

Pages 25–32

Chapter 2:

**Nuclei Isolation Protocol for Epi Multiome ATAC + Gene Expression** for use with PN-1000494 or PN-1001101

Pages 33–41

Chapter 3:

**Nuclei Isolation Protocol for Single Cell Epi ATAC**  
for use with PN-1000493

Pages 42–50

---

# Nuclei Isolation Protocol:

## Single Cell Universal Gene Expression & Flex

Reagent Kits

Get Started

Buffer Preparation

Illustrative Overview

Nuclei Isolation Protocol Steps

## Reagent Kits

### Chromium Nuclei Isolation Kits

Refer to SDS for handling and disposal information

#### Nuclei Isolation Kit with RNase Inhibitor v2 PN-1001101

Includes PN-1000447,  
PN-1000448, & PN-1001099

OR

#### Chromium Nuclei Isolation Kit with RNase Inhibitor PN-1000494\*

Includes PN-1000447,  
PN-1000448, PN-1000450, &  
PN-1000449

\*Additional item required:  
RNase Inhibitor 40X Kit  
(Kit PN-1000887,  
Tube PN-2001488)




#### For use with the following 10x Genomics products:

- Universal 3' Gene Expression
- Universal 5' Gene Expression
- Chromium Fixed RNA Profiling/Flex

##### Chromium Nuclei Isolation Reagents

16 rxns, PN-1000447

Store at 4°C

	#	PN
 Lysis Reagent	4	2000558
 Surfactant A	1	2000559
 Debris Removal Reagent	4	2000560

##### Chromium Nuclei Isolation Consumables

16 rxns, PN-1000448

Ambient Temperature

	#	PN
Pestle	16	2000561
Nuclei Isolation Column	16	2000562
Collection Tube	16	2000563
Sample Dissociation Tube	16	2000564

##### Chromium Reducing Agent B

PN-1000450


Store at -20°C

	#	PN
 Reducing Agent B	1	2000087

##### Chromium RNase Inhibitor

PN-1000449




Store at -20°C

	#	PN
 RNase Inhibitor	5	2000565

##### RNase Inhibitor Kit v2

16 rxns, PN-1001099

Store at -20°C

	#	PN
 Reducing Agent B	1	2000087
 20X Nuclei Buffer	2	2000207
 RNase Inhibitor 40X	6	2001488

Some reagent tubes (same part number) are available in multiple kits and may be used interchangeably.

# Get Started

## Nuclei Isolation Protocol:

Single Cell Universal Gene Expression & Flex



If provided Lysis Reagent and Debris Removal Buffers appear cloudy or contain precipitate, warm the tubes to **40°C** and swirl until the buffers become clear again.

Action	Item	10x PN	Preparation & Handling	Storage
Place on Ice	● Lysis Reagent	2000558	Vortex, verify no precipitate, and centrifuge briefly.	4°C
	● Surfactant A	2000559	Vortex, verify no precipitate, and centrifuge briefly.	4°C
	● Debris Removal Reagent	2000560	Vortex, verify no precipitate or layering, and centrifuge briefly.	4°C
	○ Reducing Agent B	2000087	Thaw to room temperature, vortex, verify no precipitate, and centrifuge briefly.	-20°C
	● RNase Inhibitor	2000565/ 2001488	Centrifuge briefly.	-20°C
	<b>Nuclei Isolation Consumables:</b>		Pre-chill assembled Nuclei Isolation Column(s) and Collection Tube(s) on ice.	Ambient
	• Nuclei Isolation Column	2000562		
	• Collection Tube	2000563		
	<b>1X PBS</b>	—	See Buffer Preparation.	Ambient
	<b>10% BSA</b>	—	See Buffer Preparation.	4°C
Place on Dry Ice	<b>Frozen Tissue Sample</b>	—	See Tips & Best Practices.	Liquid Nitrogen (long-term) or -80°C (short-term)
	<b>Sample Dissociation Tube</b>	2000564	Pre-chill on dry ice.	Ambient
Obtain	<b>Pestles</b>	2000561	Keep on lab bench.	Ambient
	<b>Nucleic Acid Staining Fluorescent Dye</b>	—	See Tips & Best Practices.	4°C
	<b>Vortex</b>	—	See Nuclei Isolation Protocol.	—

## Buffer Preparation: Lysis Buffer & Debris Removal Buffer

Single Cell Universal Gene  
Expression & Flex

Prepare the following Lysis and Debris Removal Buffers on ice shortly before starting the Nuclei Isolation protocol. Prepare large volumes in a 15-ml or 50-ml conical tube. Vortex briefly before use.

### Lysis Buffer

Lysis Buffer (500 µl/rxn) <i>Add reagents in the order listed</i>	PN	1X+10% (µl)	4X + 10% (µl)	8X + 10% (µl)
<input checked="" type="radio"/> Lysis Reagent	2000558	536	2,144	4,288.8
<input type="radio"/> Reducing Agent B	2000087	0.5	2	4
<input checked="" type="radio"/> Surfactant A	2000559	5.5	22	44
<input checked="" type="radio"/> RNase Inhibitor	2000565/ 2001488	13.9	55.6	111.2
<b>Total</b>	-	<b>556</b>	<b>2,224</b>	<b>4,448</b>

### Debris Removal Buffer

Debris Removal Buffer (500 µl/rxn) <i>Add reagents in the order listed</i>	PN	1X+10% (µl)	4X + 10% (µl)	8X + 10% (µl)
<input checked="" type="radio"/> Debris Removal Reagent	2000560	550	2,200	4,400
<input type="radio"/> Reducing Agent B	2000087	0.5	2	4
<b>Total</b>	-	<b>550.5</b>	<b>2,202</b>	<b>4,404</b>

## Buffer Preparation: Wash and Resuspension Buffer

Single Cell Universal Gene  
Expression & Flex


Prepare the following Wash and Resuspension Buffer on ice shortly before starting the Nuclei Isolation protocol. Prepare large volumes in a 15-ml or 50-ml conical tube. Vortex briefly before use.

### Wash and Resuspension Buffer

Wash and Resuspension Buffer (3 ml/rxn) <i>Add reagents in the order listed</i>	PN	1X+10% (µl)	4X + 10% (µl)	8X + 10% (µl)
<b>1X PBS (not provided)</b>	-	<b>2,887.5</b>	11,550	23,100
<b>10% BSA (not provided)</b>	-	<b>330</b>	1,320	2,640
<input checked="" type="radio"/> RNase Inhibitor	2000565/ 2001488	82.5	330	660
<b>Total</b>	-	<b>3,300</b>	<b>13,200</b>	<b>26,400</b>

## Protocol Overview: Nuclei Isolation from Frozen Tissue: Single Cell Universal Gene Expression & Flex

### Nuclei Isolation



## Nuclei Isolation Protocol:

Single Cell Universal  
Gene Expression &  
Flex

- a. Pre-chill centrifuge to **4°C** and place reagents and tubes on ice as indicated in the [Get Started](#) guide. Label tops and sides of tubes, as well as tops of spin columns, before starting protocol.




Perform all protocol steps on ice and centrifugation steps at 4°C.

- b. Prepare Single Cell Gene Expression and/or Chromium Fixed RNA Profiling buffers according to [Buffer Preparation](#) section and place on ice.
- c. Place Sample Dissociation Tube(s) on dry ice.
- d. Obtain frozen tissue sample(s) and place **immediately** on dry ice.
- e. Transfer frozen tissue (**3–50 mg**) to pre-chilled Sample Dissociation Tube.
- f. Transfer Sample Dissociation Tubes(s) to wet ice. Add **200 µl** Lysis Buffer to Sample Dissociation Tube. Dissociate tissue with plastic pestle until homogeneous. For multiple samples, add Lysis Buffer to each tissue and then proceed to dissociate one at a time.



Perform tissue dissociation on ice. Use one pestle per sample. **DO NOT** discard pestles until nuclei isolation process is complete.

- g. Add **300 µl** Lysis Buffer. Pipette mix 10x. If pipette tip clogs with unhomogenized tissue, continue to dissociate tissue with the pestle until able to pipette mix.
- h. After mechanical dissociation is complete, incubate on ice for **10 min**.
- i. Pipette dissociated tissue into pre-chilled Nuclei Isolation Column assembled with Collection Tube using pipette set to 500 µl. Transfer all liquid from Dissociation Tube to Nuclei Isolation Column to avoid nuclei loss.
- j. Centrifuge at **16,000 rcf** for **20 sec** at **4°C**. See [Tips & Best Practices on page 14](#) for centrifuge loading guidance.



\*Caps of tubes are angled toward the balance tubes to prevent them from turning outward due to the centrifugal force.




The absence of flowthrough following centrifugation indicates a clog in the column. Consult the [Troubleshooting Guide](#) for more information.

## Nuclei Isolation Protocol:

Single Cell Universal  
Gene Expression &  
Flex

- k.** Discard column. Flowthrough in the Collection Tube will contain nuclei. Vortex **10 sec** at **3,200 rpm** or **max speed** to resuspend nuclei. Flowthrough may appear opaque or cloudy. This is normal and it is safe to proceed.
- l.** Centrifuge Collection Tube for **3 min** at **500 rcf** at **4°C**. Carefully discard supernatant using a pipette without disturbing nuclei pellet. Leave behind a small fraction (**~200 µl**) of supernatant if nuclei pellet is not apparent.



Position tubes with hinges facing in same direction within the centrifuge, which ensures that the pellet is consistently in the same place (opposite the hinge) following centrifugation.

- m.** Resuspend nuclei pellet in **500 µl** Debris Removal Buffer. Gently pipette mix at least 15x, continuing until no pellet can be visualized.
- n.** Centrifuge at **700 rcf** for **10 min** at **4°C**. Using a pipette, carefully discard supernatant without disturbing nuclei pellet. Leave behind a small fraction (**~200 µl**) of the supernatant if nuclei pellet is not visible.
- o.** Resuspend nuclei pellet in **1 ml** of Wash and Resuspension Buffer. Gently pipette mix at least 15x, continuing until no pellet can be visualized.
- p.** Centrifuge at **500 rcf** for **5 min** at **4°C**. Carefully discard supernatant using a pipette without disturbing nuclei pellet. Leave behind a small fraction (**~200 µl**) of supernatant if nuclei pellet is not apparent.
- q.** Resuspend nuclei pellet in **1 ml** of Wash and Resuspension Buffer.
- r.** Centrifuge at **500 rcf** for **5 min** at **4°C**. Carefully discard as much supernatant as possible using a pipette without disturbing nuclei pellet. Leave behind a small remaining volume if the pellet is not visible.

## Nuclei Isolation Protocol:

Single Cell Universal  
Gene Expression &  
Flex



For lower input tissue mass (<10 mg) or if low nuclei yield is expected, a single wash may be preferred to improve total nuclei yield.

- s. Resuspend nuclei pellet in **50–500 µl** Wash and Resuspension Buffer, depending on expected recovery for input tissue type and mass. Refer to [Nuclei Recovery](#) section of Tips & Best Practices for information on typical nuclei recovery. Gently pipette mix 15x using an appropriate pipette for resuspension volume.



Resuspend in a low volume if nuclei yield is expected to be low or is unknown. **DO NOT** resuspend in a volume <50 µl.

- t. Vortex nuclei for **3 sec** at **3,200 rpm or max speed** immediately prior to counting to ensure accurate nuclei count. Pulse spin the tube after vortexing to collect liquid at bottom of tube. **DO NOT** pulse spin the tube for more than 1 second to ensure that nuclei do not pellet at the bottom of the tube.
- u. Determine nuclei concentration using AOPI or Ethidium Homodimer-1 fluorescent staining dyes and dilute if necessary for target nuclei load. Follow recommendations for nuclei counting as outlined in the [Tips & Best Practices](#) and [Appendix](#) of this document. Adjust nuclei concentration as necessary for intended downstream assay.
- v. Vortex nuclei for **3 sec** at **3,200 rpm or max speed**. Pulse spin the tube after vortexing to collect liquid at bottom of tube. **DO NOT** pulse spin the tube for more than 1 second to ensure that nuclei do not pellet at the bottom of the tube.
- w. Keep samples on ice and proceed **immediately** to relevant 10x Genomics User Guide.

---

# Nuclei Isolation Protocol:

## Single Cell Epi Multiome ATAC + Gene Expression

Reagent Kits

Get Started

Buffer Preparation

Illustrative Overview

Nuclei Isolation Protocol Steps

# Reagent Kits

## Chromium Nuclei Isolation Kits

Refer to SDS for handling and disposal information

### Nuclei Isolation Kit with RNase Inhibitor v2 PN-1001101

Includes PN-1000447, PN-1000448, & PN-1001099

or

### Chromium Nuclei Isolation Kit with RNase Inhibitor PN-1000494\*

Includes PN-1000447, PN-1000448, PN-1000450, & PN-1000449

\*Additional item required: RNase Inhibitor 40X Kit (Kit PN-1000887, Tube PN-2001488)

### For use with the following 10x Genomics products:

- Single Cell Epi Multiome ATAC + Gene Expression
- Single Cell GEM-X Epi Multiome ATAC + Gene Expression

Chromium  
**Nuclei Isolation Reagents**  
16 rxns, PN-1000447  
Store at 4°C

	#	PN
<input checked="" type="radio"/> Lysis Reagent	4	2000558
<input checked="" type="radio"/> Surfactant A	1	2000559
<input checked="" type="radio"/> Debris Removal Reagent	4	2000560

Chromium  
**Nuclei Isolation Consumables**  
16 rxns, PN-1000448  
Ambient Temperature

	#	PN
<input checked="" type="radio"/> Pestle	16	2000561
<input checked="" type="radio"/> Nuclei Isolation Column	16	2000562
<input checked="" type="radio"/> Collection Tube	16	2000563
<input checked="" type="radio"/> Sample Dissociation Tube	16	2000564

Chromium  
**Reducing Agent B**  
PN-1000450  
Store at -20°C

	#	PN
<input type="radio"/> Reducing Agent B	1	2000087

Chromium  
**RNase Inhibitor**  
PN-1000449  
Store at -20°C

	#	PN
<input checked="" type="radio"/> RNase Inhibitor	5	2000565

**RNase Inhibitor Kit v2**  
16 rxns, PN-1001099  
Store at -20°C

	#	PN
<input type="radio"/> Reducing Agent B	1	2000087
<input checked="" type="radio"/> 20X Nuclei Buffer	2	2000207
<input checked="" type="radio"/> RNase Inhibitor 40X	6	2001488

Some reagent tubes (same part number) are available in multiple kits and may be used interchangeably.

## Get Started

### Nuclei Isolation Protocol:

Single Cell Epi  
Multiome ATAC +  
Gene Expression



If provided Lysis Reagent and Debris Removal Buffers appear cloudy or contain precipitate, warm the tubes to **40°C** and swirl until the buffers become clear again.

Action	Item	10x PN	Preparation & Handling	Storage
Place on Ice	● Lysis Reagent	2000558	Vortex, verify no precipitate, and centrifuge briefly.	4°C
	● Surfactant A	2000559	Vortex, verify no precipitate, and centrifuge briefly.	4°C
	● Debris Removal Reagent	2000560	Vortex, verify no precipitate or layering, and centrifuge briefly.	4°C
	○ Reducing Agent B*	2000087	Thaw to room temperature, vortex, verify no precipitate, and centrifuge briefly.	-20°C
	● RNase Inhibitor	2000565/ 2001488	Centrifuge briefly.	-20°C
	<b>Nuclei Isolation Consumables:</b>		Pre-chill assembled Nuclei Isolation Column(s) and Collection Tube(s) on ice.	Ambient
	• Nuclei Isolation Column	2000562		
	• Collection Tube	2000563		
	● 20X Nuclei Buffer*	2000207	Thaw to room temperature, vortex, verify no precipitate, and centrifuge briefly.	-20°C
		<b>Nuclease-free Water</b>	—	See Buffer Preparation.
	<b>1X PBS</b>	—	See Buffer Preparation.	Ambient
	<b>10% BSA</b>	—	See Buffer Preparation.	4°C
Place on Dry Ice	<b>Frozen Tissue Sample</b>	—	See Tips & Best Practices.	Liquid Nitrogen (long-term) or -80°C (short-term)
	<b>Sample Dissociation Tube</b>	2000564	Pre-chill on dry ice.	Ambient
Obtain	<b>Pestles</b>	2000561	Keep on lab bench.	Ambient
	<b>Nucleic Acid Staining Fluorescent Dye</b>	—	See Tips & Best Practices.	4°C
	<b>Vortex</b>	—	See Nuclei Isolation Protocol.	—

\*Included in the 10x Genomics Single Cell Epi Multiome Kits and Nuclei Isolation Kit with RNase Inhibitor v2.

## Buffer Preparation: Lysis Buffer & Debris Removal Buffer

Single Cell Epi Multiome  
ATAC + Gene Expression

Prepare the following Lysis and Debris Removal Buffers on ice shortly before starting the Nuclei Isolation protocol. Prepare large volumes in a 15-ml or 50-ml conical tube. Vortex briefly before use.

### Lysis Buffer

Lysis Buffer (500 µl/rxn) <i>Add reagents in the order listed</i>	PN	1X+10% (µl)	4X + 10% (µl)	8X + 10% (µl)
<input checked="" type="radio"/> Lysis Reagent	2000558	536	2,144	4,288.8
<input type="radio"/> Reducing Agent B	2000087	0.5	2	4
<input checked="" type="radio"/> Surfactant A	2000559	5.5	22	44
<input checked="" type="radio"/> RNase Inhibitor	2000565/ 2001488	13.9	55.6	111.2
<b>Total</b>	–	<b>556</b>	<b>2,224</b>	<b>4,448</b>

### Debris Removal Buffer

Debris Removal Buffer (500 µl/rxn) <i>Add reagents in the order listed</i>	PN	1X+10% (µl)	4X + 10% (µl)	8X + 10% (µl)
<input checked="" type="radio"/> Debris Removal Reagent	2000560	550	2,200	4,400
<input type="radio"/> Reducing Agent B	2000087	0.5	2	4
<b>Total</b>	–	<b>550.5</b>	<b>2,202</b>	<b>4,404</b>

## Buffer Preparation: Wash & Resuspension Buffers

Single Cell Epi Multiome  
ATAC + Gene Expression

Prepare the following Wash and Resuspension Buffers on ice shortly before starting the Nuclei Isolation protocol. Prepare large volumes in a 15-ml or 50-ml conical tube. Vortex briefly before use.

### Wash Buffer

Wash Buffer (2 ml/rxn) <i>Add reagents in the order listed</i>	PN	1X+10% ( $\mu$ l)	4X + 10% ( $\mu$ l)	8X + 10% ( $\mu$ l)
<b>1X PBS (not provided)</b>	-	<b>1,925</b>	7,700	15,400
<b>10% BSA (not provided)</b>	-	<b>220</b>	880	1,760
● <b>RNase Inhibitor</b>	2000565/ 2001488	<b>55</b>	220	440
<b>Total</b>	-	<b>2,200</b>	8,800	17,600

### Resuspension Buffer


Resuspension Buffer (1 ml/rxn) <i>Add reagents in the order listed</i>	PN	1X+10% ( $\mu$ l)	4X + 10% ( $\mu$ l)	8X + 10% ( $\mu$ l)
● <b>20X Nuclei Buffer*</b>	2000207	<b>55</b>	220	440
○ <b>Reducing Agent B*</b>	2000087	<b>1</b>	4	8
<b>Nuclease-free Water (not provided)</b>	-	<b>1,016</b>	4,066	8,131
● <b>RNase Inhibitor</b>	2000565/ 2001488	<b>27.5</b>	110	220
<b>Total</b>	-	<b>1,099.5</b>	4,400	8,799

\*Included in the 10x Genomics Single Cell Epi Multiome Kits and Nuclei Isolation Kit with RNase Inhibitor v2.

## Protocol Overview: Nuclei Isolation from Frozen Tissue:

### Single Cell Epi Multiome ATAC + Gene Expression


#### Nuclei Isolation




Single Cell Epi Multiome  
ATAC + GEX


## Nuclei Isolation Protocol:

Single Cell Epi  
Multiome ATAC +  
Gene Expression

- a. Pre-chill centrifuge to **4°C** and place reagents and tubes on ice as indicated in the [Get Started](#) guide. Label tops and sides of tubes, as well as tops of spin columns, before starting protocol.
-  Perform all protocol steps on ice and centrifugation steps at 4°C.
- b. Prepare Single Cell Multiome ATAC + Gene Expression buffers according to [Buffer Preparation](#) section and place on ice.
  - c. Place Sample Dissociation Tube(s) on dry ice.
  - d. Obtain frozen tissue sample(s) and place **immediately** on dry ice.
  - e. Transfer frozen tissue (**3–50 mg**) to pre-chilled Sample Dissociation Tube.
  - f. Transfer Sample Dissociation Tubes(s) to wet ice. Add **200 µl** Lysis Buffer to Sample Dissociation Tube. Dissociate tissue with plastic pestle until homogeneous. For multiple samples, add Lysis Buffer to each tissue and then proceed to dissociate one at a time.

 Perform tissue dissociation on ice. Use one pestle per sample. **DO NOT** discard pestles until nuclei isolation process is complete.

- g. Add **300 µl** Lysis Buffer. Pipette mix 10x. If pipette tip clogs with unhomogenized tissue, continue to dissociate tissue with the pestle until able to pipette mix.
- h. After mechanical dissociation is complete, incubate on ice for **10 min**.
- i. Pipette dissociated tissue into pre-chilled Nuclei Isolation Column assembled with Collection Tube using pipette set to 500 µl. Transfer all liquid from Dissociation Tube to Nuclei Isolation Column to avoid nuclei loss.
- j. Centrifuge at **16,000 rcf** for **20 sec** at **4°C**. See [Tips & Best Practices on page 14](#) for centrifuge loading guidance.



\*Caps of tubes are angled toward the balance tubes to prevent them from turning outward due to the centrifugal force.




The absence of flowthrough following centrifugation indicates a clog in the column. Consult the [Troubleshooting Guide](#) for more information.

## Nuclei Isolation Protocol:

Single Cell Epi  
Multiome ATAC +  
Gene Expression

- k. Discard column. Flowthrough in the Collection Tube will contain nuclei. Vortex **10 sec** at **3,200 rpm** or **max speed** to resuspend nuclei. Flowthrough may appear opaque or cloudy. This is normal and it is safe to proceed.
- l. Centrifuge Collection Tube for **3 min** at **500 rcf** at **4°C**. Carefully discard supernatant using a pipette without disturbing nuclei pellet. Leave behind a small fraction (**~200 µl**) of supernatant if nuclei pellet is not apparent.



**TIPS** Position tubes with hinges facing in same direction within the centrifuge, which ensures that the pellet is consistently in the same place (opposite the hinge) following centrifugation.

- m. Resuspend nuclei pellet in **500 µl** Debris Removal Buffer. Gently pipette mix at least 15x, continuing until no pellet can be visualized.
- n. Centrifuge at **700 rcf** for **10 min** at **4°C**. Carefully discard supernatant using a pipette without disturbing nuclei pellet. Leave behind a small fraction (**~200 µl**) of supernatant if nuclei pellet is not apparent.
- o. Resuspend nuclei pellet in **1 ml** of Wash Buffer.
- p. Centrifuge at **500 rcf** for **5 min** at **4°C**. Carefully discard supernatant using a pipette without disturbing nuclei pellet. Leave behind a small fraction (**~200 µl**) of supernatant if nuclei pellet is not apparent.
- q. Resuspend nuclei pellet in **1 ml** of Wash Buffer.
- r. Centrifuge at **500 rcf** for **5 min** at **4°C**. Carefully discard as much supernatant as possible using a pipette without disturbing nuclei pellet. Leave behind a small remaining volume if the pellet is not visible.

**TIPS** For lower input tissue mass (<10 mg) or if low nuclei yield is expected, a single wash may be preferred to improve total nuclei yield.

## Nuclei Isolation Protocol:

Single Cell Epi  
Multiome ATAC +  
Gene Expression

- s. Resuspend nuclei pellet in **50–500  $\mu$ l** Resuspension Buffer, depending on expected recovery for input tissue type and mass. Refer to [Nuclei Recovery](#) section of Tips & Best Practices for information on typical nuclei recovery. Gently pipette mix 15x using an appropriate pipette for resuspension volume.



*Resuspend in a low volume (~25  $\mu$ l) if nuclei yield is expected to be low or is unknown.*

- t. Vortex nuclei for **3 sec** at **3,200 rpm or max speed** immediately prior to counting to ensure accurate nuclei count. Pulse spin the tube after vortexing to collect liquid at bottom of tube. **DO NOT** pulse spin the tube for more than 1 second to ensure that nuclei do not pellet at the bottom of the tube.
- u. Determine nuclei concentration using AOPI or Ethidium Homodimer-1 fluorescent staining dyes and dilute if necessary for target nuclei load. Follow recommendations for nuclei counting as outlined in the [Tips & Best Practices](#) and [Appendix](#) of this document. Adjust nuclei concentration as necessary for intended downstream assay.
- v. Vortex nuclei for **3 sec** at **3,200 rpm or max speed**. Pulse spin the tube after vortexing to collect liquid at bottom of tube. **DO NOT** pulse spin the tube for more than 1 second to ensure that nuclei do not pellet at the bottom of the tube.
- w. Keep samples on ice and proceed **immediately** to relevant 10x Genomics User Guide.

---

# Nuclei Isolation Protocol:

## Single Cell Epi ATAC

Reagent Kits

Get Started

Buffer Preparation

Illustrative Overview

Nuclei Isolation Protocol Steps

## Chromium Nuclei Isolation Kit PN-1000493

### For use with the following 10x Genomics products:

- Single Cell ATAC

Chromium  
**Nuclei Isolation Reagents**  
16 rxns, PN-1000447  
Store at 4°C

	#	PN
<input type="radio"/> Lysis Reagent	4	2000558
<input checked="" type="radio"/> Surfactant A	1	2000559
<input type="radio"/> Debris Removal Reagent	4	2000560

Chromium  
**Nuclei Isolation Consumables**  
16 rxns, PN-1000448  
Ambient Temperature

	#	PN
Pestle	16	2000561
Nuclei Isolation Column	16	2000562
Collection Tube	16	2000563
Sample Dissociation Tube	16	2000564

Chromium  
**Reducing Agent B**  
PN-1000450  
Store at -20°C

	#	PN
<input type="radio"/> Reducing Agent B	1	2000087

Single Cell Epi  
ATAC

## Get Started

### Nuclei Isolation Protocol:

Single Cell Epi ATAC



If provided Lysis Reagent and Debris Removal Buffers appear cloudy or contain precipitate, warm the tubes to **40°C** and swirl until the buffers become clear again.

Action	Item	10x PN	Preparation & Handling	Storage
Place on Ice	● Lysis Reagent	2000558	Vortex, verify no precipitate, and centrifuge briefly.	4°C
	● Surfactant A	2000559	Vortex, verify no precipitate, and centrifuge briefly.	4°C
	● Debris Removal Reagent	2000560	Vortex, verify no precipitate or layering, and centrifuge briefly.	4°C
	○ Reducing Agent B	2000087	Thaw to room temperature, vortex, verify no precipitate, and centrifuge briefly.	-20°C
	<b>Nuclei Isolation Consumables:</b>		Pre-chill assembled Nuclei Isolation Column(s) and Collection Tube(s) on ice.	Ambient
	• Nuclei Isolation Column	2000562		
	• Collection Tube	2000563		
● 20X Nuclei Buffer*	2000207	Thaw to room temperature, vortex, verify no precipitate, and centrifuge briefly.	-20°C	
	<b>Nuclease-free Water</b>	—	See Buffer Preparation.	Ambient
	<b>1X PBS</b>	—	See Buffer Preparation.	Ambient
	<b>10% BSA</b>	—	See Buffer Preparation.	4°C
Place on Dry Ice	<b>Frozen Tissue Sample</b>	—	See Tips & Best Practices.	Liquid Nitrogen (long-term) or -80°C (short-term)
	<b>Sample Dissociation Tube</b>	2000564	Pre-chill on dry ice.	Ambient
Obtain	<b>Pestles</b>	2000561	Keep on lab bench.	Ambient
	<b>Nucleic Acid Staining Fluorescent Dye</b>	—	See Tips & Best Practices.	4°C
	<b>Vortex</b>	—	See Nuclei Isolation Protocol.	—

\*20X Nuclei Buffer is included in the 10x Genomics Single Cell ATAC v2 Reagent Kits

## Buffer Preparation: Lysis Buffer & Debris Removal Buffer

Single Cell Epi ATAC

Prepare the following Lysis and Debris Removal Buffers on ice shortly before starting the Nuclei Isolation protocol. Prepare large volumes in a 15-ml or 50-ml conical tube. Vortex briefly before use.

### Lysis Buffer

Lysis Buffer (500 µl/rxn) <i>Add reagents in the order listed</i>	PN	1X+10% (µl)	4X + 10% (µl)	8X + 10% (µl)
<input checked="" type="radio"/> Lysis Reagent	2000558	550	2,200	4,400
<input type="radio"/> Reducing Agent B	2000087	0.5	2	4
<input checked="" type="radio"/> Surfactant A	2000559	5.5	22	44
PBS	–	1,668	6,672	13,344
<b>Total</b>	–	<b>2,224</b>	<b>8,896</b>	<b>17,792</b>

### Debris Removal Buffer

Debris Removal Buffer (500 µl/rxn) <i>Add reagents in the order listed</i>	PN	1X+10% (µl)	4X + 10% (µl)	8X + 10% (µl)
<input checked="" type="radio"/> Debris Removal Reagent	2000560	550	2,200	4,400
<input type="radio"/> Reducing Agent B	2000087	0.5	2	4
<b>Total</b>	–	<b>550.5</b>	<b>2,202</b>	<b>4,404</b>

## Buffer Preparation: Wash & Resuspension Buffers

Single Cell Epi ATAC

Prepare the following Wash and Resuspension Buffers on ice shortly before starting the Nuclei Isolation protocol. Prepare large volumes in a 15-ml or 50-ml conical tube. Vortex briefly before use.

### Wash Buffer

Wash Buffer (2 ml/rxn) <i>Add reagents in the order listed</i>	PN	1X+10% ( $\mu$ l)	4X + 10% ( $\mu$ l)	8X + 10% ( $\mu$ l)
<b>1X PBS (not provided)</b>	-	<b>1,980</b>	7,920	15,840
<b>10% BSA (not provided)</b>	-	<b>220</b>	880	1,760
<b>Total</b>	-	<b>2,200</b>	8,800	17,600


### Resuspension Buffer

Resuspension Buffer (1 ml/rxn) <i>Add reagents in the order listed</i>	PN	1X+10% ( $\mu$ l)	4X + 10% ( $\mu$ l)	8X + 10% ( $\mu$ l)
<input checked="" type="radio"/> <b>20X Nuclei Buffer*</b>	2000207	<b>55</b>	220	440
<input type="radio"/> <b>Reducing Agent B</b>	2000087	<b>1</b>	4	8
<b>Nuclease-free Water (not provided)</b>	-	<b>1,044</b>	4,176	8,351
<b>Total</b>	-	<b>1,100</b>	4,400	8,799

\*20X Nuclei Buffer is included in the 10x Genomics Single Cell ATAC v2 Reagent Kits

## Protocol Overview: Nuclei Isolation from Frozen Tissue: Single Cell Epi ATAC

### Nuclei Isolation



## Nuclei Isolation Protocol:

### Single Cell Epi ATAC

- a. Pre-chill centrifuge to **4°C** and place reagents and tubes on ice as indicated in the [Get Started](#) guide. Label tops and sides of tubes, as well as tops of spin columns, before starting protocol.




Perform all protocol steps on ice and centrifugation steps at 4°C.

- b. Prepare Single Cell ATAC buffers according to [Buffer Preparation](#) section and place on ice.
- c. Place Sample Dissociation Tube(s) on dry ice.
- d. Obtain frozen tissue sample(s) and place **immediately** on dry ice.
- e. Transfer frozen tissue (**3–50 mg**) to pre-chilled Sample Dissociation Tube.
- f. Transfer Sample Dissociation Tubes(s) to wet ice. Add **200 µl** Lysis Buffer to Sample Dissociation Tube. Dissociate tissue with plastic pestle until homogeneous. For multiple samples, add Lysis Buffer to each tissue and then proceed to dissociate one at a time.



Perform tissue dissociation on ice. Use one pestle per sample. DO NOT discard pestles until nuclei isolation process is complete.

- g. Add **300 µl** Lysis Buffer. Pipette mix 10x. If pipette tip clogs with unhomogenized tissue, continue to dissociate tissue with the pestle until able to pipette mix.
- h. After mechanical dissociation is complete, incubate on ice for **10 min**.
- i. Pipette dissociated tissue into pre-chilled Nuclei Isolation Column assembled with Collection Tube using pipette set to 500 µl. Transfer all liquid from Dissociation Tube to Nuclei Isolation Column to avoid nuclei loss.
- j. Centrifuge at **16,000 rcf** for **20 secs** at **4°C**. See [Tips & Best Practices](#) section on page 14 for centrifuge loading guidance.



\*Caps of tubes are angled toward the balance tubes to prevent them from turning outward due to the centrifugal force.




The absence of flowthrough following centrifugation indicates a clog in the column. Consult the [Troubleshooting Guide](#) for more information.

## Nuclei Isolation Protocol:

### Single Cell Epi ATAC

- k.** Discard column. Flowthrough in the Collection Tube will contain nuclei. Vortex **10 sec** at **3,200 rpm** or **max speed** to resuspend nuclei. Flowthrough may appear opaque or cloudy. This is normal and it is safe to proceed.
- l.** Centrifuge Collection Tube for **3 min** at **500 rcf** at **4°C**. Carefully discard supernatant using a pipette without disturbing nuclei pellet. Leave behind a small fraction (**~200 µl**) of supernatant if nuclei pellet is not apparent.



Position tubes with hinges facing in same direction within the centrifuge, which ensures that the pellet is consistently in the same place (opposite the hinge) following centrifugation.

- m.** Resuspend nuclei pellet in **500 µl** Debris Removal Buffer. Gently pipette mix at least 15x, continuing until no pellet can be visualized.
- n.** Centrifuge at **700 rcf** for **10 min** at **4°C**. Carefully discard supernatant using a pipette without disturbing nuclei pellet. Leave behind a small fraction (**~200 µl**) of supernatant if nuclei pellet is not apparent.
- o.** Resuspend nuclei pellet in **1 ml** of Wash Buffer.
- p.** Centrifuge at **500 rcf** for **5 min** at **4°C**. Carefully discard supernatant using a pipette without disturbing nuclei pellet. Leave behind a small fraction (**~200 µl**) of supernatant if nuclei pellet is not apparent.
- q.** Resuspend nuclei pellet in **1 ml** of Wash Buffer.
- r.** Centrifuge at **500 rcf** for **5 min** at **4°C**. Carefully discard as much supernatant as possible using a pipette without disturbing nuclei pellet. Leave behind a small remaining volume if the pellet is not visible.



For lower input tissue mass (<10 mg) or if low nuclei yield is expected, a single wash may be preferred to improve total nuclei yield.

## Nuclei Isolation Protocol:

### Single Cell Epi ATAC

- s. Resuspend nuclei pellet in **50–500  $\mu$ l** Resuspension Buffer, depending on expected recovery for input tissue type and mass. Refer to [Nuclei Recovery](#) section of Tips & Best Practices for information on typical nuclei recovery. Gently pipette mix 15x using an appropriate pipette for resuspension volume.





*Resuspend in a low volume if nuclei yield is expected to be low or is unknown. DO NOT resuspend in a volume <50  $\mu$ l.*

- t. Vortex nuclei for **3 sec at 3,200 rpm or max speed** immediately prior to counting to ensure accurate nuclei count. Pulse spin the tube after vortexing to collect liquid at bottom of tube. DO NOT pulse spin the tube for more than 1 second to ensure that nuclei do not pellet at the bottom of the tube.
- u. Determine nuclei concentration using AOPI or Ethidium Homodimer-1 fluorescent staining dyes and dilute if necessary for target nuclei load. Follow recommendations for nuclei counting as outlined in the [Tips & Best Practices](#) and [Appendix](#) of this document. Adjust nuclei concentration as necessary for intended downstream assay.
- v. Vortex nuclei for **3 sec at 3,200 rpm or max speed**. Pulse spin the tube after vortexing to collect liquid at bottom of tube. DO NOT pulse spin the tube for more than 1 second to ensure that nuclei do not pellet at the bottom of the tube.
- w. Keep samples on ice and proceed **immediately** to relevant 10x Genomics User Guide.

# Troubleshooting





## Troubleshooting Guide

Issue	Potential Causes	Mitigation Strategies
<b>Tissue not dissociated</b> (A) Dissociated tissue and (B) Tissue that is not fully dissociated		
 <p data-bbox="170 672 267 703"><b>A</b> Normal</p> <p data-bbox="365 672 397 703"><b>B</b></p>	<ul data-bbox="592 430 998 703" style="list-style-type: none"> <li>• Tissue mass outside range of recommended tissue sizes.</li> <li>• Lysis time insufficient for dissociation.</li> <li>• Tissue type not compatible with Chromium Nuclei Isolation Kit.</li> <li>• Insufficient dissociation with pestle.</li> </ul>	<ul data-bbox="1047 430 1469 651" style="list-style-type: none"> <li>• Lyse tissue for a longer period of time.</li> <li>• Increase number of pestle strokes for complete tissue dissociation.</li> <li>• Ensure tissue mass is within recommended tissue size range.</li> </ul>
<b>Clogged nuclei isolation column</b> (A) Column with successful flowthrough and (B) Clogged column		
 <p data-bbox="170 1344 267 1375"><b>A</b> Normal</p> <p data-bbox="365 1344 397 1375"><b>B</b></p>	<ul data-bbox="592 829 998 1102" style="list-style-type: none"> <li>• Tissue mass outside range of recommended tissue sizes.</li> <li>• Tissue not fully dissociated.</li> <li>• Lysis time not sufficient for dissociation of tissue.</li> <li>• Tissue type not compatible with Chromium Nuclei Isolation Kit.</li> </ul>	<ul data-bbox="1047 829 1469 1512" style="list-style-type: none"> <li>• Lyse tissue for a longer period of time.</li> <li>• Increase number of pestle strokes for complete tissue dissociation.</li> <li>• Ensure tissue mass is within recommended tissue size range.</li> <li>• Ensure tissue is not on list of Incompatible Tissue Types.</li> <li>• Use a wide-bore pipette to transfer tissues with excessive fibrous debris to limit tip clogging and maximize transfer of nuclei to the Nuclei Isolation Column.</li> <li>• In the event of a column clog, transfer liquid remaining on top of the Nuclei Isolation Column to a new column leaving behind any debris. Place into same Collection Tube and repeat centrifugation step.</li> </ul>







If a column clog occurs, contact 10x Genomics Technical Support at [support@10xgenomics.com](mailto:support@10xgenomics.com) for further assistance.

## Troubleshooting Guide

Issue	Potential Causes	Mitigation Strategies
<b>Low nuclei recovery</b> (A) <i>Acceptable nuclei recovery</i> and (B) <i>Low nuclei recovery</i>		
 <p><b>A</b> <i>Normal</i></p> <p><b>B</b></p>	<ul style="list-style-type: none"> <li>• Improper freezing and storage of input tissue.</li> <li>• Not enough starting tissue material.</li> <li>• Reagents and samples not kept on ice during Nuclei Isolation process.</li> <li>• Inaccurate counting of nuclei.</li> <li>• Tissue not fully dissociated.</li> <li>• Low cellularity (e.g., high levels of extracellular matrix or lower cell density).</li> <li>• Pellet lost or discarded during wash steps.</li> </ul>	<ul style="list-style-type: none"> <li>• Store tissue samples in liquid nitrogen for best results and at -80°C for short-term needs.</li> <li>• Minimize handling steps during tissue processing.</li> <li>• Mince the tissue before dissociation.</li> <li>• Perform entire Nuclei Isolation protocol on ice with chilled reagents and tools.</li> <li>• Reduce number of wash steps during Nuclei Isolation protocol.</li> <li>• Increase centrifugation times to 10 minutes during wash steps.</li> <li>• Use DNA-binding dye and a fluorescent microscope to count nuclei.</li> <li>• Use a swinging-bucket rotor centrifuge.</li> <li>• Resuspend nuclei pellet in 25 µl Resuspension Buffer and instead of vortexing, use a narrow bore pipette tip to mix ~10x to increase the final nuclei concentration.</li> <li>• Perform two or three parallel preparations from the same sample, then pool the isolated nuclei prior to the final wash step.</li> </ul>
<b>Poor nuclei quality</b> (A) <i>Good nuclei quality</i> and (B) <i>poor nuclei quality with blebbing</i>		
 <p><b>A</b> <i>Normal</i></p> <p><b>B</b></p>	<ul style="list-style-type: none"> <li>• Improper freezing and storage of input tissue.</li> <li>• Overlysis of nuclei.</li> <li>• Reagents and samples not kept on ice during Nuclei Isolation process.</li> <li>• Resuspension of nuclei is too harsh.</li> <li>• Excess delays between nuclei isolation and counting/loading of chip.</li> </ul>	<ul style="list-style-type: none"> <li>• Store tissue samples in liquid nitrogen for best results and at -80°C for short-term needs.</li> <li>• Minimize handling steps during tissue processing.</li> <li>• Perform entire Nuclei Isolation protocol on ice with chilled reagents and tools.</li> <li>• Minimize delays between start of lysis and loading of chip.</li> <li>• Reduce the lysis buffer strength to 0.25x (the default for standalone ATAC assays) and check the suspension under the microscope following lysis to determine the nuclei quality.</li> </ul>

## Troubleshooting Guide

Issue	Potential Causes	Mitigation Strategies
<p><b>Debris and aggregates</b> (A) Acceptable and (B) unacceptable amounts of debris and aggregates</p>	<ul style="list-style-type: none"> <li>• Overlysis of nuclei.</li> <li>• Cell lysis not performed on ice.</li> <li>• Resuspension between steps not sufficient.</li> <li>• Transfer of lysate to Debris Removal step without spin column filtration.</li> </ul>	<ul style="list-style-type: none"> <li>• Assess lysis efficacy via microscopy after incubation, if necessary.</li> <li>• Optimize lysis time for new tissues types or sizes.</li> <li>• Perform entire Nuclei Isolation protocol on ice.</li> <li>• Ensure sufficient vortex time and speed for resuspension.</li> <li>• If aggregates &gt;50 µm are observed, the nuclei suspension may optionally be passed through a 40 µm FlowMi filter.*</li> </ul>
<p><b>Low cDNA/Poor sequencing data</b> (A) Normal barcode rank plot and (B) abnormal barcode rank plot</p>	<ul style="list-style-type: none"> <li>• Improper freezing and storage of input tissue.</li> <li>• Samples not vortexed after resuspension steps.</li> <li>• Incorrect amount of RNase inhibitor added to buffers.</li> <li>• Samples degraded during isolation process.</li> </ul>	<ul style="list-style-type: none"> <li>• Store tissue samples in liquid nitrogen for best results and at -80°C for short-term needs.</li> <li>• Vortex samples where indicated during resuspension and wash steps.</li> <li>• Add the recommended amount of RNase inhibitor to buffers according to compatible assay, as outlined in the Buffer Preparation sections of this document. Test using smaller tissue fragments and shorter lysis times.</li> <li>• Perform entire Nuclei Isolation protocol on ice.</li> <li>• Increase cDNA amplification by 1–2 cycles.</li> </ul>



\*Note that this additional filtration may cause a significant reduction in nuclei recovery.

# Appendix

Tested Tissue Types

Tissues with Variable Results

Optimization Recommendations

Cell/Nuclei Counting

## Tested Tissue Types



Human tumors are highly complex tissues. Results can vary from sample to sample and within a single tumor biopsy.

The Chromium Nuclei Isolation Kit was optimized using human and mouse samples and is expected to be compatible with most mammalian tissues. Optimization may be needed for more challenging or untested sample types. The following tissues have been successfully tested using the Chromium Nuclei Isolation Kit within the supported mass range of the kit (**3–50 mg**). Refer to the 10x Genomics Website for a complete and updated list of Tested Tissue Types.\*

For tissues that have not been tested, perform a pilot experiment before scaling up. If a tissue of interest is not in the tested list, consider whether it shares characteristics with tissues that have been tested. For example, fibrous or difficult-to-dissociate tissues may behave similarly to mouse heart, which is fibrous and has been successfully tested for nuclei isolation. Similarly, tissues with high RNase activity may be comparable to pancreas or spleen, both of which have been successfully processed for nuclei isolation and downstream gene expression data. These examples can help guide expectations for performance with the Nuclei Isolation Kit and provide general recommendations based on tissue type.

Organism	Tissue (Healthy/Tumor)
Mouse	Kidney (Healthy)
	Liver (Healthy)
	Lung (Healthy)
	Brain (Healthy)
	Heart (Healthy)
	Small Intestine (Healthy)
	Eye (Healthy)
	Skeletal Muscle (Healthy)
	Spinal Cord (Healthy)
	Bladder (Healthy)
	Ovary (Healthy)
	Colon (Healthy)
	Adipose (Healthy)
	Stomach (Healthy)
Testis (Healthy)	
Human	Jejunum (Healthy)
	Duodenum (Healthy)
	Ileum (Healthy)
	Testis (Healthy)
	Breast (Tumor)
	Prostate (Tumor)
	Melanoma (Tumor)

## Tested Tissue Types

Contd.



Human tumors are highly complex tissues. Results can vary from sample to sample and within a single tumor biopsy.

Organism	Tissue (Healthy/Tumor)
Human	Ovarian (Tumor)
	Colorectal (Tumor)
	Pancreas (Tumor)
	Kidney (Tumor)
	Lung (Tumor)

\*Note that tissues not listed may still be compatible with the Chromium Nuclei Isolation Kit. This list only summarizes tissues validated specifically by 10x Genomics.

## Tissues with Variable Results

The following tissues may require further optimization for use with the Chromium Nuclei Isolation Kit. Note that tissues with high RNase content are sensitive to collection and storage conditions. Refer to [Tip & Best Practices](#) for guidance on tissue handling and storage. Cut tissues with lower yields into small pieces ( $\leq 10$  mg) before dissociation in Lysis Buffer. See [Tissue Dissociation](#) section of Tips & Best Practices for more information.

Organism	Tissue (Healthy/Tumor)	Additional Notes
Mouse	Pancreas (Healthy)	High RNase tissue
	Spleen (Healthy)	High RNase tissue
	Skin (Healthy)	Lower yield tissue
	Tongue (Healthy)	Lower yield tissue
Human	Spleen (Healthy)	High RNase tissue

## Incompatible Tissue Types

The following tissues have been tested by 10x Genomics and are NOT recommended for use with the Chromium Nuclei Isolation Kit:

- Cell suspensions (i.e. cultured cells, PBMCs)
- Plants
- Insects
- Calcified tissue (i.e. bone)
- FFPE tissue

## Optimization Recommendations

The Chromium Nuclei Isolation Kit protocol was validated without modification using the indicated Tested Tissue Types. If expected performance is not achieved for tissue(s) of interest using the recommended protocol, the following optimization of some protocol steps may improve performance based on the unique properties of target tissue(s).

### **Lysis time:**

- Perform a lysis timeline to determine appropriate lysis incubation time for new tissue types.

### **Lysis buffer strength:**

- If nuclei quality is poor, buffer detergent strength can be decreased for a gentler lysis by diluting with 1X PBS.

### **Sample cleanup steps:**

- Additional washes may further reduce small debris and ambient RNA in the sample.

## Cell/Nuclei Counting

Vendor	Item	Part Number
<b>For Cell/Counting</b>		
Nexcelom Biosciences	*ViaStain PI Staining Solution	CS1-0109-5mL
	*ViaStain AOPI Staining Solution	CS2-0106-5mL
	**Cellaca MX High-throughput Automated Cell Counter	MX-112-0127
	**Cellometer K2 Fluorescent Cell Counter	CMT-K2-MX-150
	PD100 Counting Chambers 1 case	CHT4-PD100-003
Biotium	*NucSpot 470	40083
	<i>Dilute the stock to 1:100 and mix 1:1 with the sample. For example, add 10 <math>\mu</math>l diluted dye to 10 <math>\mu</math>l sample.</i>	
Thermo Fisher Scientific	**Countess II FL Automated Cell Counter	AMAQAF1000
	<i>Discontinued</i>	
	Countess Automated Cell Counting Chamber Slides	C10228
	**Countess 3 FL Automated Cell Counter	AMQAF2000
	Trypan Blue Stain (0.4%)	T10282
	*DAPI solution, 1 mg/mL	62248
<i>*Choose either AOPI, NucSpot, PI, or DAPI solution. If the sample has no debris, Trypan Blue can be used.</i>		
<i>**Choose Countess II/3, Cellaca, Cellometer, or equivalent fluorescent counter.</i>		

## Cell/Nuclei Counting



Where recommended, accurate counting is essential for workflow execution. A fluorescent dye is strongly recommended. Trypan blue may be used in debris free samples only if fluorescent dye and counter are not available.

- Accurate sample counting is critical for achieving desired cell/nuclei recovery. Table below shows the combination of counters and dyes tested for counting nuclei.
- It is strongly recommended that the sample be stained with a fluorescent dye such as PI staining solution and counted using an automated fluorescent cell counter or hemocytometer.
- The use of fluorescent dye during counting enables accurate quantification even in the presence of debris.
- Automated fluorescent cell counters are strongly recommended when counting nuclei.
- Ensure that the counter laser/filter setup is compatible with the fluorescent dye used.
- Ensure cell/nuclei are well-focused under brightfield before switching to the fluorescent channel for counting.
- Increase exposure time to help adjust signal to noise during counting.
- Perform visual inspection to confirm that the counting number is accurate. For example, after obtaining the counting number, switch between the brightfield and fluorescent channels to make sure the counts include minimal debris and the most cell/nuclei.

Counter Type	Fluorescent Dye	Counting Comparison
<b>Cellaca</b> Range: $1 \times 10^5$ – $1 \times 10^7$ cells/ml Automated exclusion of debris from cell count	<ul style="list-style-type: none"> <li>• Propidium Iodide</li> <li>• NucSpot 470</li> <li>• DAPI</li> </ul>	Comparable counting results at both counting steps for all three dyes
<b>Countess II FL/Countess 3 FL</b> Range: $1 \times 10^4$ – $1 \times 10^7$ cells/ml (optimal $1 \times 10^5$ – $4 \times 10^6$ ) Manual debris exclusion from cell count post-image capture, using gates on the instrument program	<ul style="list-style-type: none"> <li>• Propidium Iodide</li> <li>• NucSpot 470</li> <li>• DAPI</li> </ul>	Comparable counting results at both counting steps for the three dyes
<b>Cellometer K2</b> Range: $1 \times 10^5$ – $1 \times 10^7$ cells/ml Debris exclusion from cell count by adjusting instrument program settings before image capture	<ul style="list-style-type: none"> <li>• Propidium Iodide</li> <li>• NucSpot 470</li> </ul>	Comparable counting results at both counting steps for the two dyes  Propidium Iodide stained nuclei require longer exposure compared to NucSpot 470 but can still be relatively dimmer

## Cell/Nuclei Counting

*Contd.*

### Counting using PI Staining Solution

This protocol provides instructions for counting samples using PI staining solution and the Cellaca Counter to enable accurate quantification even in the presence of subcellular debris. The optimal cell/nuclei concentration for the Cellaca Counter is 100-10,000 cells/ $\mu$ l. Refer to manufacturer's instructions for details on operations.

- Add **25  $\mu$ l** PI Staining Solution into Mixing Row of Cellaca plate
- Gently mix the sample. If the sample is too concentrated, a 1:1 dilution in PBS can also be prepared. For example, add 15  $\mu$ l nuclei suspension to 15  $\mu$ l PBS. Ensure that this dilution factor is accounted for during counting. For example, because a 1:1 dilution was performed, the final cell concentration should be multiplied by two.
- Add **25  $\mu$ l** sample to Mixing Row of plate containing PI Staining Solution. Gently pipette mix 8x.
- Transfer stained sample to Loading Row of Cellaca plate.

Samples stained with PI staining solution can also be counted using the Countess II/3 or K2 Automated Cell Counter. Refer to manufacturer's instructions for details.

### Counting using Trypan Blue (Only for Debris-free Samples)



Debris-free samples (cells or nuclei suspensions) can also be counted using trypan blue. This protocol provides instructions for counting sample using trypan blue and a hemocytometer or Countess II Automated Cell Counter.

- Mix **1 part** 0.4% trypan blue and **1 part** sample.
- Transfer **10  $\mu$ l** sample to a Countess II Cell Counting Slide chamber or a hemocytometer.
- Insert the slide into the Countess II Cell Counter and determine the cell concentration. Or if using a hemocytometer, count by placing the hemocytometer under the microscope.
- The majority of nuclei suspensions will be stained with trypan blue stain and appear nonviable.

## Cell/Nuclei Counting

### Representative Examples

**Brightfield & Fluorescent Images from Automated Cell Counter - Samples Stained with Propidium Iodide**



Example of a Cell Not Counted

**Brightfield image from Automated Cell Counter - Samples Stained with Trypan Blue**

Fresh PBMCs

