

SAFETY DATA SHEETS

Visium Spatial Gene Expression for FFPE Reagent Kits v2

Visium FFPE Reagents v2, Small PN-1000436

REAGENTS	PN
Amp Mix B	2000567
Extension Enzyme	2000389
Extension Buffer	2000409
RNase Enzyme	3000593
RNase Buffer B	2000551
Tissue Removal Enzyme	3000387
Tissue Removal Buffer B	2000543
Tissue Removal Buffer Enhancer	2000557
Decrosslinking Buffer	2000566
TS Primer Mix B	2000537
Block and Stain Buffer	2000554

SECTION 1: Identification

1.1. Identification

Product form : Mixtures
Trade name : Amp Mix B
Product code : 2000567

1.2. Recommended use and restrictions on use

Use of the substance/mixture : Reagent

1.3. Supplier

Manufacturer/Supplier:
10x Genomics
6230 Stoneridge Mall Road Pleasanton,
CA 94588-3260
T:+1 925 401 7300
E: info@10xgenomics.com

1.4. Emergency telephone number

Emergency number : +1 925 401 7300

SECTION 2: Hazard(s) identification

2.1. Classification of the substance or mixture

GHS US classification

Acute toxicity (oral) Category 4	Harmful if swallowed
Skin corrosion/irritation Category 2	Causes skin irritation
Specific target organ toxicity (single exposure) Category 1	Causes damage to organs

2.2. GHS Label elements, including precautionary statements

GHS US labeling

Hazard pictograms (GHS US) :



Signal word (GHS US) :

Danger

Hazard statements (GHS US) :

Harmful if swallowed
Causes skin irritation
Causes damage to organs

Precautionary statements (GHS US) :

Do not breathe mist, spray, vapors.
Wash hands thoroughly after handling.
Do not eat, drink or smoke when using this product.
Wear protective gloves.
If swallowed: Call a poison center or doctor if you feel unwell.
If on skin: Wash with plenty of water.
If exposed: Call a poison center/doctor.
Rinse mouth.
If skin irritation occurs: Get medical advice/attention.
Take off contaminated clothing and wash it before reuse.
Store locked up.
Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation.

Amp Mix B

Safety Data Sheet

according to US HazCom 2012

2.3. Other hazards which do not result in classification

No additional information available

2.4. Unknown acute toxicity (GHS US)

Not applicable

SECTION 3: Composition/Information on ingredients

3.1. Substances

Not applicable

3.2. Mixtures

Name	Product identifier	%	GHS US classification
Methanaminium, N,N,N-trimethyl-, chloride	CAS-No.: 75-57-0	2 – 6	Acute Tox. 2 (Oral), H300 Acute Tox. 3 (Dermal), H311 Skin Irrit. 2, H315 STOT SE 1, H370
sulphuric acid	CAS-No.: 7664-93-9	1 – 3	Skin Corr. 1A, H314

*Chemical name, CAS number and/or exact concentration have been withheld as a trade secret

SECTION 4: First-aid measures

4.1. Description of first aid measures

First-aid measures after inhalation	: Remove person to fresh air and keep comfortable for breathing.
First-aid measures after skin contact	: Rinse immediately with plenty of water for 15 minutes. Take off contaminated clothing. If skin irritation occurs: Get medical advice/attention.
First-aid measures after eye contact	: Rinse eyes with water as a precaution. Remove contact lenses, if present and easy to do. Continue rinsing. Seek medical attention if ill effect or irritation develops.
First-aid measures after ingestion	: Rinse mouth. Do NOT induce vomiting. Call a poison center/doctor/physician if you feel unwell.

4.2. Most important symptoms and effects (acute and delayed)

Symptoms/effects	: May cause damage to organs.
Symptoms/effects after skin contact	: Irritation.
Symptoms/effects after ingestion	: Harmful if swallowed.

4.3. Immediate medical attention and special treatment, if necessary

Treat symptomatically.

SECTION 5: Fire-fighting measures

5.1. Suitable (and unsuitable) extinguishing media

Suitable extinguishing media	: Use extinguishing media appropriate for surrounding fire.
Unsuitable extinguishing media	: None known.

5.2. Specific hazards arising from the chemical

Fire hazard	: The product is not flammable. Does not sustain combustion. In the event of fire, may decompose : Carbon oxides (CO, CO ₂). Sulphur oxides.
Explosion hazard	: No hazard identified.

Amp Mix B

Safety Data Sheet

according to US HazCom 2012

Hazardous decomposition products in case of fire : Thermal decomposition can lead to the release of irritating gases and vapors.

5.3. Special protective equipment and precautions for fire-fighters

Firefighting instructions : Fight fire with normal precautions from a reasonable distance.
Protection during firefighting : Do not attempt to take action without suitable protective equipment.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

General measures : Evacuate unnecessary personnel.

6.1.1. For non-emergency personnel

Protective equipment : Wear personal protective equipment. For further information refer to section 8: "Exposure controls/personal protection".
Emergency procedures : Ventilate spillage area. Do not breathe mist, spray, vapors. Avoid contact with skin.

6.1.2. For emergency responders

Protective equipment : Wear recommended personal protective equipment. For further information refer to section 8: "Exposure controls/personal protection".
Emergency procedures : Stop leak if safe to do so. Evacuate unnecessary personnel. Ventilate area.

6.2. Environmental precautions

Avoid release to the environment.

6.3. Methods and material for containment and cleaning up

For containment : Soak up with inert absorbent material (for example sand, sawdust, a universal binder, silica gel).
Methods for cleaning up : Collect spillage. Store away from other materials.
Other information : Dispose in a safe manner in accordance with local/national regulations.

6.4. Reference to other sections

For further information refer to section 8: "Exposure controls/personal protection". For disposal of residues refer to section 13 : "Disposal considerations".

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : Ensure good ventilation of the work station. Do not breathe mist, spray, vapors. Avoid contact with skin. Wear personal protective equipment.
Hygiene measures : Wash hands and other exposed areas with mild soap and water before eating, drinking or smoking and when leaving work. Take off contaminated clothing and wash it before reuse. Do not eat, drink or smoke when using this product. Always wash hands after handling the product. Handle in accordance with good industrial hygiene and safety practice.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Keep only in the original container in a cool well ventilated place. Keep container closed when not in use. Keep cool. Store in a well-ventilated place. Store locked up.
Incompatible materials : Strong oxidizing agents. Strong acids. Strong bases. Strong reducing agents.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Amp Mix B

Safety Data Sheet

according to US HazCom 2012

Amp Mix B	
No additional information available	
sulphuric acid (7664-93-9)	
USA - ACGIH - Occupational Exposure Limits	
Local name	Sulfuric acid
ACGIH OEL TWA	0.2 mg/m ³ (T - Thoracic particulate matter)
Remark (ACGIH)	TLV® Basis: Pulm func. Notations: A2 (Suspected Human Carcinogen. Classification refers to sulfuric acid contained in strong inorganic acid mists)
Regulatory reference	ACGIH 2022
Methanaminium, N,N,N-trimethyl-, chloride (75-57-0)	
No additional information available	

8.2. Appropriate engineering controls

Appropriate engineering controls : Ensure good ventilation of the work station.
Environmental exposure controls : Avoid release to the environment.

8.3. Individual protection measures/Personal protective equipment

Hand protection:

Wear suitable gloves resistant to chemical penetration. Please follow the instructions related to the permeability and the penetration time provided by the manufacturer

Eye protection:

If there is a risk of liquid being splashed : Safety glasses

Respiratory protection:

In case of insufficient ventilation, wear suitable respiratory equipment. Appropriate dust or mist respirator should be used if airborne particles are generated when handling this material

Other information:

Do not eat, drink or smoke during use.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state	: Liquid
Appearance	: Clear.
Color	: Colourless
Odor	: odorless
Odor threshold	: No data available
pH	: No data available
Melting point	: No data available
Freezing point	: No data available
Boiling point	: No data available
Flash point	: No data available
Relative evaporation rate (butyl acetate=1)	: No data available
Flammability (solid, gas)	: Not applicable.
Vapor pressure	: No data available
Relative vapor density at 20 °C	: No data available
Relative density	: No data available

Amp Mix B

Safety Data Sheet

according to US HazCom 2012

Solubility	: No data available
Partition coefficient n-octanol/water (Log Pow)	: No data available
Auto-ignition temperature	: No data available
Decomposition temperature	: No data available
Viscosity, kinematic	: No data available
Viscosity, dynamic	: No data available
Explosion limits	: No data available
Explosive properties	: No data available
Oxidizing properties	: No data available

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

The product is non-reactive under normal conditions of use, storage and transport.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

No dangerous reactions known under normal conditions of use.

10.4. Conditions to avoid

None under recommended storage and handling conditions (see section 7).

10.5. Incompatible materials

Strong oxidizing agents. Strong acids. Strong bases. Strong reducing agents.

10.6. Hazardous decomposition products

Under normal conditions of storage and use, hazardous decomposition products should not be produced. In the event of fire, may decompose : Carbon oxides (CO, CO₂). Sulphur oxides.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity (oral)	: Harmful if swallowed.
Acute toxicity (dermal)	: Not classified (Based on available data, the classification criteria are not met)
Acute toxicity (inhalation)	: Not classified (Based on available data, the classification criteria are not met)

Amp Mix B

ATE US (oral)	833 mg/kg body weight
---------------	-----------------------

Methanaminium, N,N,N-trimethyl-, chloride (75-57-0)

LD50 oral rat	50 mg/kg
LD50 dermal rabbit	200 – 500 mg/kg

Skin corrosion/irritation	: Causes skin irritation.
Serious eye damage/irritation	: Not classified (Based on available data, the classification criteria are not met)
Respiratory or skin sensitization	: Not classified (Based on available data, the classification criteria are not met)
Germ cell mutagenicity	: Not classified (Based on available data, the classification criteria are not met)

Amp Mix B

Safety Data Sheet

according to US HazCom 2012

Carcinogenicity	: Not classified (Based on available data, the classification criteria are not met)
Reproductive toxicity	: Not classified (Based on available data, the classification criteria are not met)
STOT-single exposure	: Causes damage to organs.

Methanaminium, N,N,N-trimethyl-, chloride (75-57-0)

STOT-single exposure	Causes damage to organs.
STOT-repeated exposure	: Not classified (Based on available data, the classification criteria are not met)
Aspiration hazard	: Not classified (Based on available data, the classification criteria are not met)
Viscosity, kinematic	: No data available
Symptoms/effects	: May cause damage to organs.
Symptoms/effects after skin contact	: Irritation.
Symptoms/effects after ingestion	: Harmful if swallowed.
Other information	: Likely routes of exposure: ingestion, inhalation, skin and eye.

SECTION 12: Ecological information

12.1. Toxicity

Methanaminium, N,N,N-trimethyl-, chloride (75-57-0)

LC50 - Fish [1]	431 – 495 mg/l (Exposure time: 96 h - Species: Pimephales promelas [flow-through])
-----------------	--

12.2. Persistence and degradability

No additional information available

12.3. Bioaccumulative potential

Methanaminium, N,N,N-trimethyl-, chloride (75-57-0)

Partition coefficient n-octanol/water (Log Pow)	< -1.6 (at 20 °C (at pH 7))
---	-----------------------------

12.4. Mobility in soil

No additional information available

12.5. Other adverse effects

Other information : Avoid release to the environment.

SECTION 13: Disposal considerations

13.1. Disposal methods

Waste treatment methods	: Dispose of contents/container in accordance with licensed collector's sorting instructions.
Product/Packaging disposal recommendations	: Dispose in a safe manner in accordance with local/national regulations.
Ecology - waste materials	: Avoid release to the environment.

SECTION 14: Transport information

In accordance with DOT / TDG / IMDG / IATA

DOT	TDG	IMDG	IATA
14.1. UN number			
Not regulated for transport			

Amp Mix B

Safety Data Sheet

according to US HazCom 2012

DOT	TDG	IMDG	IATA
14.2. Proper Shipping Name			
Not applicable	Not applicable	Not applicable	Not applicable
14.3. Transport hazard class(es)			
Not applicable	Not applicable	Not applicable	Not applicable
14.4. Packing group			
Not applicable	Not applicable	Not applicable	Not applicable
14.5. Environmental hazards			
Not applicable	Not applicable	Not applicable	Not applicable
No supplementary information available			

SECTION 15: Regulatory information

15.1. US Federal regulations

sulphuric acid (7664-93-9)

CERCLA RQ	1000 lb
RQ (Reportable quantity, section 304 of EPA's List of Lists)	1000 lb
SARA Section 302 Threshold Planning Quantity (TPQ)	1000 lb

15.2. International regulations

CANADA

sulphuric acid (7664-93-9)

Listed on the Canadian DSL (Domestic Substances List)

Methanaminium, N,N,N-trimethyl-, chloride (75-57-0)

Listed on the Canadian DSL (Domestic Substances List)

EU-Regulations

Methanaminium, N,N,N-trimethyl-, chloride (75-57-0)

Listed on the EEC inventory EINECS (European Inventory of Existing Commercial Chemical Substances)

National regulations

sulphuric acid (7664-93-9)

Listed as carcinogen on NTP (National Toxicology Program)
Listed on INSQ (Mexican National Inventory of Chemical Substances)

Amp Mix B

Safety Data Sheet

according to US HazCom 2012

Methanaminium, N,N,N-trimethyl-, chloride (75-57-0)

Listed introduction on Australian Industrial Chemicals Introduction Scheme (AICIS Inventory)
Listed on PICCS (Philippines Inventory of Chemicals and Chemical Substances)
Listed on the Japanese ENCS (Existing New Chemical Substances) inventory
Listed on KECL/KECI (Korean Existing Chemicals Inventory)
Listed on IECSC (Inventory of Existing Chemical Substances Produced or Imported in China)
Listed on NZIoC (New Zealand Inventory of Chemicals)
Listed on the Japanese ISHL (Industrial Safety and Health Law)
Listed on INSQ (Mexican National Inventory of Chemical Substances)
Listed on the TCSI (Taiwan Chemical Substance Inventory)
Listed on the NCI (Vietnam - National Chemical Inventory)

15.3. US State regulations

Component	State or local regulations
sulphuric acid(7664-93-9)	U.S. - Massachusetts - Right To Know List; U.S. - New Jersey - Right to Know Hazardous Substance List; U.S. - New York City - Right to Know Hazardous Substances List; U.S. - Pennsylvania - RTK (Right to Know) List

SECTION 16: Other information

according to US HazCom 2012

Revision date : 6 June 2022
Other information : None.

Safety Data Sheet (SDS), USA

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product.

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 04/20/2021

Reviewed on 04/20/2021

* 1 Identification

- **Product Identifier**
- **Trade Name:** *Barcoding Enzyme
Second Strand Enzyme
Extension enzyme
Extension enzyme*
- **Product Number:**
2000139
2000125
2000183
2000218
2000389
2000427
- **Relevant identified uses of the substance or mixture and uses advised against:**
- **Product Description:** Not relevant information.
- **Details of the Supplier of the Safety Data Sheet:**
- **Manufacturer/Supplier:**
10x Genomics
6230 Stoneridge Mall Road
Pleasanton, CA 94588 USA
1 925 401 7300
<https://www.10xgenomics.com/>
- **Emergency telephone number:** 1 925 401 7300

2 Hazard(s) Identification

- **Classification of the substance or mixture:**
The product does not need classification according to OSHA HazCom Standard 29 CFR paragraph (d) of §1910.1200(g) and GHS Rev 03.
- **Label elements:**
- **Hazard pictograms:** Non-Regulated Material
- **Signal word:** Non-Regulated Material
- **Hazard statements:** Non-Regulated Material
- **Precautionary statements:**
P264 Wash thoroughly after handling.
P270 Do not eat, drink or smoke when using this product.
- **Unknown acute toxicity:**
This value refers to knowledge of known, established toxicological or ecotoxicological values.
0 % of the mixture consists of component(s) of unknown toxicity.
- **Classification system:** NFPA/HMIS Definitions: 0-Least, 1-Slight, 2-Moderate, 3-High, 4-Extreme
- **NFPA ratings (scale 0 - 4)**



- **HMIS-ratings (scale 0 - 4)**
- | | | |
|------------|---|---------------------|
| HEALTH | 0 | Health = 0 |
| FIRE | 0 | Fire = 0 |
| REACTIVITY | 0 | Physical Hazard = 0 |
- **Hazard(s) not otherwise classified (HNOC):** None known

(Contd. on page 2)

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 04/20/2021

Reviewed on 04/20/2021

**Trade Name: *Barcoding Enzyme
Second Strand Enzyme
Extension enzyme
Extension enzyme***

3 Composition/Information on Ingredients

· **Non-hazardous components:**

Non-Hazardous Components	90-99%
--------------------------	--------

· **Chemical characterization: Substance**

· **Description:** Mixture of substances listed below with non-hazardous additions.

· **Dangerous Components:** Non-Regulated Material

· **Additional information:**

The exact percentages of the ingredients of this mixture are considered to be proprietary and are withheld in accordance with the provisions of paragraph (i) of §1910.1200 of 29 CFR 1910.1200 Trade Secrets.

4 First-Aid Measures

· **Description of first aid measures**

· **General information:** No special measures required.

· **After inhalation:**

If breathed in, move person into fresh air. If not breathing, give artificial respiration. In case of discomfort seek medical attention.

· **After skin contact:**

Wash off with soap and plenty of water for at least 15 minutes. In case of discomfort seek medical attention.

· **After eye contact:**

Flush eyes thoroughly with water for 15 minutes. Remove contact lenses after the initial 1-2 minutes and continue flushing. In case of discomfort seek medical attention.

· **After swallowing:**

Do NOT induce vomiting. Never give anything by mouth to an unconscious person. Rinse mouth with water. In case of discomfort seek medical attention.

· **Information for doctor**

· **Most important symptoms and effects, both acute and delayed:** No further relevant information available.

· **Indication of any immediate medical attention and special treatment needed:**

No further relevant information available.

5 Fire-Fighting Measures

· **Extinguishing media**

· **Suitable extinguishing agents:**

CO₂, extinguishing powder or water spray. Fight larger fires with water spray.

Use fire fighting measures that suit the environment.

· **For safety reasons unsuitable extinguishing agents:** No further relevant information.

· **Special hazards arising from the substance or mixture:** No further relevant information available.

· **Advice for firefighters**

· **Special protective equipment for firefighters:**

As in any fire, wear self-contained breathing apparatus pressure-demand (NIOSH approved or equivalent) and full protective gear to prevent contact with skin and eyes.

6 Accidental Release Measures

· **Personal precautions, protective equipment and emergency procedures:**

Avoid contact with skin, eyes and clothing.

· **Methods and material for containment and cleaning up:**

Absorb with liquid-binding material (i.e. sand, diatomite, acid binders, universal binders, sawdust).

(Contd. on page 3)

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 04/20/2021

Reviewed on 04/20/2021

**Trade Name: *Barcoding Enzyme
Second Strand Enzyme
Extension enzyme
Extension enzyme***

Dispose of the collected material according to regulations.

- **Reference to other sections:**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.
- **Protective Action Criteria for Chemicals**

· PAC-1:
None of the ingredients are listed.
· PAC-2:
None of the ingredients are listed.
· PAC-3:
None of the ingredients are listed.

7 Handling and Storage

- **Handling**
- **Precautions for safe handling:**
Avoid contact with skin, eyes and clothing
No special measures required.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** None.
- **Specific end use(s):** No further relevant information available.

8 Exposure Controls/Personal Protection

- **Additional information about design of technical systems:** No further data; see section 7.
- **Control parameters:**
- **Components with occupational exposure limits:**
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:** The lists that were valid during the creation of this SDS were used as basis.
- **Exposure controls:**
- **Personal protective equipment**
- **General protective and hygienic measures:**
The usual precautionary measures for handling chemicals should be followed.
- **Breathing equipment:** Not required.
- **Protection of hands:**



Protective gloves

- **Material of gloves:** Not applicable.
- **Penetration time of glove material:** Not applicable.

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 04/20/2021

Reviewed on 04/20/2021

**Trade Name: *Barcoding Enzyme
Second Strand Enzyme
Extension enzyme
Extension enzyme***

· **Eye protection:**



Safety glasses

· **Body protection:** Not required.

9 Physical and Chemical Properties

· **Information on basic physical and chemical properties**

· **General Information**

· **Appearance:**

· **Form:**

Liquid

· **Color:**

Clear

· **Odor:**

Odorless

· **Odor threshold:**

Not determined.

· **pH-value:**

Not applicable.

· **Change in condition**

· **Melting point/Melting range:**

Not determined.

· **Flash point:**

None

· **Flammability (solid, gaseous):**

Not applicable.

· **Ignition temperature:**

Not applicable

· **Decomposition temperature:**

Not determined.

· **Auto igniting:**

Product is not self-igniting.

· **Danger of explosion:**

Product does not present an explosion hazard.

· **Explosion limits:**

· **Lower:**

Not determined.

· **Upper:**

Not determined.

· **Vapor pressure:**

Not determined.

· **Density @ 20 °C (68 °F):**

1.1983 g/cm³ (9.9998 lbs/gal)

· **Relative density:**

Not determined.

· **Vapor density:**

Not determined.

· **Evaporation rate:**

Not determined.

· **Solubility in / Miscibility with:**

· **Water:**

Not miscible or difficult to mix.

· **Partition coefficient (n-octanol/water):** Not determined.

· **Viscosity:**

· **Dynamic:**

Not determined.

· **Kinematic:**

Not determined.

· **Solvent content:**

· **VOC content:**

0.00 %

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 04/20/2021

Reviewed on 04/20/2021

Trade Name: *Barcoding Enzyme
Second Strand Enzyme
Extension enzyme
Extension enzyme*

Solids content: 0.0 %

· **Other information:** No further relevant information available.

10 Stability and Reactivity

- **Reactivity:** No further relevant information available.
- **Chemical stability:** Stable under normal conditions.
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions:** No dangerous reactions known.
- **Conditions to avoid:** Do not store next to heat source, in direct sunlight, or elevated storage temperature.
- **Incompatible materials:** Strong acids, strong bases
- **Hazardous decomposition products:** No dangerous decomposition products known.

11 Toxicological Information

- **Information on toxicological effects:**
- **Acute toxicity:**
- **LD/LC50 values that are relevant for classification:** No data available.
- **Primary irritant effect:**
- **On the skin:** No irritating effect.
- **On the eye:** No irritating effect.
- **Additional toxicological information:**
The product is not subject to classification according to internally approved calculation methods for preparations.
When used and handled according to specifications, the product does not have any harmful effects according to our experience and the information provided to us.
- **Carcinogenic categories:**

· **IARC (International Agency for Research on Cancer):**

None of the ingredients are listed.

· **NTP (National Toxicology Program):**

None of the ingredients are listed.

· **OSHA-Ca (Occupational Safety & Health Administration):**

None of the ingredients are listed.

12 Ecological Information

- **Toxicity:**
- **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability:** No further relevant information available.
- **Behavior in environmental systems:**
- **Bioaccumulative potential:** No further relevant information available.
- **Mobility in soil:** No further relevant information available.
- **Additional ecological information:**
- **General notes:**
Do not allow undiluted product or product that has not been neutralized to reach ground water, water course or sewage system.
- **Results of PBT and vPvB assessment:**
- **PBT:** Not applicable.

(Contd. on page 6)

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 04/20/2021

Reviewed on 04/20/2021

Trade Name: *Barcoding Enzyme
Second Strand Enzyme
Extension enzyme
Extension enzyme*

- **vPvB:** Not applicable.
- **Other adverse effects:** No further relevant information available.

13 Disposal Considerations

- **Waste treatment methods**
- **Recommendation:**
Smaller quantities can be disposed of with household waste.
Observe all federal, state and local environmental regulations when disposing of this material.
- **Uncleaned packaging**
- **Recommendation:** Disposal must be made according to official regulations.

14 Transport Information

- **UN-Number:**
- **DOT, ADR/ADN, ADN, IMDG, IATA** Non-Regulated Material
- **UN proper shipping name:**
- **DOT, ADR/ADN, ADN, IMDG, IATA** Non-Regulated Material
- **Transport hazard class(es):**
- **DOT, ADR/ADN, ADN, IMDG, IATA**
- **Class:** Non-Regulated Material
- **Packing group:**
- **DOT, ADR/ADN, IMDG, IATA** Non-Regulated Material
- **Environmental hazards:** Not applicable.
- **Special precautions for user:** Not applicable.
- **Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code:** Not applicable.
- **UN "Model Regulation":** Non-Regulated Material

15 Regulatory Information

- **Safety, health and environmental regulations/legislation specific for the substance or mixture:**
No further relevant information available.
- **SARA (Superfund Amendments and Reauthorization):**

· **Section 355 (extremely hazardous substances):**

None of the ingredients are listed.

· **Section 313 (Specific toxic chemical listings):**

None of the ingredients are listed.

· **TSCA (Toxic Substances Control Act):**

None of the ingredients are listed.

· **California Proposition 65:**

· **Chemicals known to cause cancer:**

None of the ingredients are listed.

· **Chemicals known to cause reproductive toxicity for females:**

None of the ingredients are listed.

· **Chemicals known to cause reproductive toxicity for males:**

None of the ingredients are listed.

(Contd. on page 7)

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 04/20/2021

Reviewed on 04/20/2021

**Trade Name: *Barcoding Enzyme
Second Strand Enzyme
Extension enzyme
Extension enzyme***

<ul style="list-style-type: none"> · Chemicals known to cause developmental toxicity:
None of the ingredients are listed.
<ul style="list-style-type: none"> · New Jersey Right-to-Know List:
None of the ingredients are listed.
<ul style="list-style-type: none"> · New Jersey Special Hazardous Substance List:
None of the ingredients are listed.
<ul style="list-style-type: none"> · Pennsylvania Right-to-Know List:
None of the ingredients are listed.
<ul style="list-style-type: none"> · Pennsylvania Special Hazardous Substance List:
None of the ingredients are listed.
<ul style="list-style-type: none"> · Carcinogenic categories:
<ul style="list-style-type: none"> · EPA (Environmental Protection Agency):
None of the ingredients are listed.
<ul style="list-style-type: none"> · TLV (Threshold Limit Value established by ACGIH):
None of the ingredients are listed.
<ul style="list-style-type: none"> · NIOSH-Ca (National Institute for Occupational Safety and Health):
None of the ingredients are listed.
<ul style="list-style-type: none"> · GHS label elements Non-Regulated Material · Hazard pictograms: Non-Regulated Material · Signal word: Non-Regulated Material · Hazard statements: Non-Regulated Material · Precautionary statements:
P264 Wash thoroughly after handling.
P270 Do not eat, drink or smoke when using this product.
<ul style="list-style-type: none"> · National regulations:
None of the ingredients are listed.
<ul style="list-style-type: none"> · Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

16 Other Information

The information and recommendations in this safety data sheet are, to the best of our knowledge, accurate as of the date of issue. Nothing herein shall be deemed to create warranty, expressed or implied, and shall not establish a legally valid contractual relationship. It is the responsibility of the user to determine applicability of this information and the suitability of the material or product for any particular purpose.

- **Contact:**
- **Date of last revision/ revision number:** 04/20/2021 / 16
- **Abbreviations and acronyms:**
 ADR: The European Agreement concerning the International Carriage of Dangerous Goods by Road
 ADN: The European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways
 IMDG: International Maritime Code for Dangerous Goods
 DOT: US Department of Transportation
 IATA: International Air Transport Association
 EINECS: European Inventory of Existing Commercial Chemical Substances
 ELINCS: European List of Notified Chemical Substances
 CAS: Chemical Abstracts Service (division of the American Chemical Society)
 NFPA: National Fire Protection Association (USA)
 HMIS: Hazardous Materials Identification System (USA)

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 04/20/2021

Reviewed on 04/20/2021

**Trade Name: *Barcoding Enzyme
Second Strand Enzyme
Extension enzyme
Extension enzyme***

VOC: Volatile Organic Compounds (USA, EU)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety and Health

OSHA: Occupational Safety & Health Administration

TLV: Threshold Limit Value

PEL: Permissible Exposure Limit

REL: Recommended Exposure Limit

*** *Data compared to the previous version altered.***

SDS created by MSDS Authoring Services www.msdsauthoring.com +1-877-204-9106

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 05/05/2021

Reviewed on 05/05/2021

*** 1 Identification**


- **Product Identifier**
- **Trade Name: *Barcoding Reagent
Second Strand Reagent
Extension buffer***
- **Product Number:**
2000124
2000219
2000220
2000409
2000408
- **Relevant identified uses of the substance or mixture and uses advised against:**
- **Product Description:** Not relevant information.
- **Details of the Supplier of the Safety Data Sheet:**
- **Manufacturer/Supplier:**
10x Genomics
6230 Stoneridge Mall Road
Pleasanton, CA 94588 USA
1 925 401 7300
<https://www.10xgenomics.com/>
- **Emergency telephone number:** 1 925 401 7300

2 Hazard(s) Identification

- **Classification of the substance or mixture:**
Eye Irrit. 2B H320 Causes eye irritation.
- **Label elements:**
- **Hazard pictograms:** Non-Regulated Material
- **Signal word:** Warning
- **Hazard statements:**
H320 Causes eye irritation.
- **Precautionary statements:**
P264 Wash thoroughly after handling.
P305+P351+P338 If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337+P313 If eye irritation persists: Get medical advice/attention.
- **Unknown acute toxicity:**
This value refers to knowledge of known, established toxicological or ecotoxicological values.
0 % of the mixture consists of component(s) of unknown toxicity.
- **Classification system:** NFPA/HMIS Definitions: 0-Least, 1-Slight, 2-Moderate, 3-High, 4-Extreme
- **NFPA ratings (scale 0 - 4)**



- **HMIS-ratings (scale 0 - 4)**



- **Hazard(s) not otherwise classified (HNOC):** None known

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 05/05/2021

Reviewed on 05/05/2021

**Trade Name: *Barcoding Reagent
Second Strand Reagent
Extension buffer***

3 Composition/Information on Ingredients

- **Chemical characterization: Substance**
- **Description:** Mixture of substances listed below with non-hazardous additions.

- **Dangerous Components:**

CAS: 56-81-5	Glycerine	5-10%
RTECS: MA 8050000		

- **Additional information:**

The exact percentages of the ingredients of this mixture are considered to be proprietary and are withheld in accordance with the provisions of paragraph (i) of §1910.1200 of 29 CFR 1910.1200 Trade Secrets.

4 First-Aid Measures

- **Description of first aid measures**
- **After inhalation:** Generally the product does not irritate with inhalation.
- **After skin contact:**
Wash with soap and water.
If skin irritation occurs, consult a doctor.
- **After eye contact:**
Rinse opened eye for at least 15 minutes under running water. If symptoms persist, consult a doctor.
If easy to do so, remove contact lenses if worn.
- **After swallowing:**
Call a doctor if irritation or symptoms develop.
Do not induce vomiting.
- **Information for doctor**
- **Most important symptoms and effects, both acute and delayed:** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed:**
No further relevant information available.

5 Fire-Fighting Measures

- **Extinguishing media**
- **Suitable extinguishing agents:**
CO₂, extinguishing powder or water spray. Fight larger fires with water spray.
Use fire fighting measures that suit the environment.
- **For safety reasons unsuitable extinguishing agents:** No further relevant information.
- **Special hazards arising from the substance or mixture:** No further relevant information available.
- **Advice for firefighters**
- **Special protective equipment for firefighters:**
As in any fire, wear self-contained breathing apparatus pressure-demand (NIOSH approved or equivalent) and full protective gear to prevent contact with skin and eyes.

6 Accidental Release Measures

- **Personal precautions, protective equipment and emergency procedures:**
Avoid contact with skin, eyes and clothing.
- **Methods and material for containment and cleaning up:**
Absorb with liquid-binding material (i.e. sand, diatomite, acid binders, universal binders, sawdust).
Ensure adequate ventilation.
Dispose of the collected material according to regulations.

(Contd. on page 3)

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 05/05/2021

Reviewed on 05/05/2021

**Trade Name: *Barcoding Reagent
Second Strand Reagent
Extension buffer***

- **Reference to other sections:**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.

- **Protective Action Criteria for Chemicals**

· PAC-1:		
56-81-5	Glycerine	45 mg/m ³
	Trade Secret	6.9 mg/m ³
	Trade Secret	12 mg/m ³
	Trade Secret	11 mg/m ³
· PAC-2:		
56-81-5	Glycerine	180 mg/m ³
	Trade Secret	76 mg/m ³
	Trade Secret	130 mg/m ³
	Trade Secret	120 mg/m ³
· PAC-3:		
56-81-5	Glycerine	1,100 mg/m ³
	Trade Secret	460 mg/m ³
	Trade Secret	790 mg/m ³
	Trade Secret	550 mg/m ³

7 Handling and Storage

- **Handling**
- **Precautions for safe handling:**
Avoid contact with skin, eyes and clothing
Ensure good ventilation/exhaustion at the workplace.
Prevent formation of aerosols.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** Keep receptacle tightly sealed.
- **Specific end use(s):** No further relevant information available.

8 Exposure Controls/Personal Protection

- **Additional information about design of technical systems:** No further data; see section 7.
- **Control parameters:**

· Components with occupational exposure limits:	
56-81-5 Glycerine	
PEL	Long-term value: 15* 5** mg/m ³ mist; *total dust **respirable fraction
TLV	TLV withdrawn-insufficient data human occup. exp.

(Contd. on page 4)

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 05/05/2021

Reviewed on 05/05/2021

**Trade Name: *Barcoding Reagent
Second Strand Reagent
Extension buffer***

- **Additional information:** The lists that were valid during the creation of this SDS were used as basis.
- **Exposure controls:**
- **Personal protective equipment**
- **General protective and hygienic measures:**
The usual precautionary measures for handling chemicals should be followed.
Keep away from foodstuffs, beverages and feed.
Immediately remove all soiled and contaminated clothing and wash before reuse.
Wash hands before breaks and at the end of work.
Avoid contact with the eyes.
Avoid contact with the eyes and skin.
- **Breathing equipment:** Not required.
- **Protection of hands:**



Protective gloves

- **Material of gloves:** Not applicable.
- **Penetration time of glove material:** Not applicable.
- **Eye protection:**



Safety glasses

- **Body protection:** Not required.

9 Physical and Chemical Properties

· **Information on basic physical and chemical properties**

· **General Information**

· **Appearance:**

· **Form:** Liquid

· **Color:** Clear

· **Odor:** Odorless

· **Odor threshold:** Not determined.

· **pH-value:** Not applicable.

· **Change in condition**

· **Melting point/Melting range:** Not determined.

· **Flash point:** None

· **Flammability (solid, gaseous):** Not applicable.

· **Ignition temperature:** Not applicable

· **Decomposition temperature:** Not determined.

· **Auto igniting:** Product is not self-igniting.

· **Danger of explosion:** Product does not present an explosion hazard.

· **Explosion limits:**

· **Lower:** Not determined.

(Contd. on page 5)

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 05/05/2021

Reviewed on 05/05/2021

**Trade Name: *Barcoding Reagent
Second Strand Reagent
Extension buffer***

- Upper:** Not determined.
- Vapor pressure:** Not determined.
- Density:**
 - Relative density:** Not determined.
 - Vapor density:** Not determined.
 - Evaporation rate:** Not determined.
- Solubility in / Miscibility with:**
 - Water:** Not miscible or difficult to mix.
- Partition coefficient (n-octanol/water):** Not determined.
- Viscosity:**
 - Dynamic:** Not determined.
 - Kinematic:** Not determined.
- Solvent content:**
 - VOC content:** 0.00 %
 - Solids content:** 0.8 %
- Other information:** No further relevant information available.

10 Stability and Reactivity

- Reactivity:** No further relevant information available.
- Chemical stability:** Stable under normal conditions.
- Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- Possibility of hazardous reactions:** No dangerous reactions known.
- Conditions to avoid:** Strong acids, strong bases.
- Incompatible materials:** Strong acids, strong bases
- Hazardous decomposition products:** No dangerous decomposition products known.

11 Toxicological Information

- Information on toxicological effects:**
- Acute toxicity:**

LD/LC50 values that are relevant for classification:		
56-81-5 Glycerine		
Oral	LD50	12,600 mg/kg (Rat)
Dermal	LD50	>10,000 mg/kg (Rabbit)
Inhalative	LC50/96 hours	>50,000 mg/l (Trout)

- Primary irritant effect:**
- On the skin:** No irritating effect.
- On the eye:** Irritating effect.
- Additional toxicological information:**
The product shows the following dangers according to internally approved calculation methods for preparations:
Irritant

(Contd. on page 6)

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 05/05/2021

Reviewed on 05/05/2021

**Trade Name: *Barcoding Reagent
Second Strand Reagent
Extension buffer***

· **Carcinogenic categories:**

· **IARC (International Agency for Research on Cancer):**

None of the ingredients are listed.

· **NTP (National Toxicology Program):**

None of the ingredients are listed.

· **OSHA-Ca (Occupational Safety & Health Administration):**

None of the ingredients are listed.

12 Ecological Information

· **Toxicity:**

· **Aquatic toxicity:** No further relevant information available.

· **Persistence and degradability:** No further relevant information available.

· **Behavior in environmental systems:**

· **Bioaccumulative potential:** No further relevant information available.

· **Mobility in soil:** No further relevant information available.

· **Additional ecological information:**

· **General notes:**

Do not allow undiluted product or product that has not been neutralized to reach ground water, water course or sewage system.

· **Results of PBT and vPvB assessment:**

· **PBT:** Not applicable.

· **vPvB:** Not applicable.

· **Other adverse effects:** No further relevant information available.

13 Disposal Considerations

· **Waste treatment methods**

· **Recommendation:**

Observe all federal, state and local environmental regulations when disposing of this material.

· **Uncleaned packaging**

· **Recommendation:** Disposal must be made according to official regulations.

14 Transport Information

· **UN-Number:**

· **DOT, ADR/ADN, ADN, IMDG, IATA** Non-Regulated Material

· **UN proper shipping name:**

· **DOT, ADR/ADN, ADN, IMDG, IATA** Non-Regulated Material

· **Transport hazard class(es):**

· **DOT, ADR/ADN, ADN, IMDG, IATA**

· **Class:** Non-Regulated Material

· **Packing group:**

· **DOT, ADR/ADN, IMDG, IATA** Non-Regulated Material

· **Environmental hazards:** Not applicable.

· **Special precautions for user:** Not applicable.

· **Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code:** Not applicable.

(Contd. on page 7)

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 05/05/2021

Reviewed on 05/05/2021

**Trade Name: *Barcoding Reagent
Second Strand Reagent
Extension buffer***

· **UN "Model Regulation":** Non-Regulated Material

15 Regulatory Information

· **Safety, health and environmental regulations/legislation specific for the substance or mixture:**
No further relevant information available.

· **SARA (Superfund Amendments and Reauthorization):**

· **Section 355 (extremely hazardous substances):**

None of the ingredients are listed.

· **Section 313 (Specific toxic chemical listings):**

None of the ingredients are listed.

· **TSCA (Toxic Substances Control Act):**

56-81-5	Glycerine
---------	-----------

	Trade Secret
--	--------------

	Trade Secret
--	--------------

7447-40-7	Potassium chloride
-----------	--------------------

12125-02-9	Ammonium Chloride
------------	-------------------

	Trade Secret
--	--------------

· **California Proposition 65:**

· **Chemicals known to cause cancer:**

None of the ingredients are listed.

· **Chemicals known to cause reproductive toxicity for females:**

None of the ingredients are listed.

· **Chemicals known to cause reproductive toxicity for males:**

None of the ingredients are listed.

· **Chemicals known to cause developmental toxicity:**

None of the ingredients are listed.

· **New Jersey Right-to-Know List:**

56-81-5	Glycerine
---------	-----------

12125-02-9	Ammonium Chloride
------------	-------------------

· **New Jersey Special Hazardous Substance List:**

None of the ingredients are listed.

· **Pennsylvania Right-to-Know List:**

56-81-5	Glycerine
---------	-----------

12125-02-9	Ammonium Chloride
------------	-------------------

· **Pennsylvania Special Hazardous Substance List:**

12125-02-9	Ammonium Chloride	E
------------	-------------------	---

· **Carcinogenic categories:**

· **EPA (Environmental Protection Agency):**

None of the ingredients are listed.

· **TLV (Threshold Limit Value established by ACGIH):**

None of the ingredients are listed.

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 05/05/2021

Reviewed on 05/05/2021

**Trade Name: *Barcoding Reagent
Second Strand Reagent
Extension buffer***

· **NIOSH-Ca (National Institute for Occupational Safety and Health):**

None of the ingredients are listed.

· **GHS label elements**

The product is classified and labeled according to the Globally Harmonized System (GHS).

· **Hazard pictograms:** Non-Regulated Material

· **Signal word:** Warning

· **Hazard statements:**

H320 Causes eye irritation.

· **Precautionary statements:**

P264 Wash thoroughly after handling.

P305+P351+P338 If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337+P313 If eye irritation persists: Get medical advice/attention.

· **National regulations:**

None of the ingredients are listed.

· **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

16 Other Information

The information and recommendations in this safety data sheet are, to the best of our knowledge, accurate as of the date of issue. Nothing herein shall be deemed to create warranty, expressed or implied, and shall not establish a legally valid contractual relationship. It is the responsibility of the user to determine applicability of this information and the suitability of the material or product for any particular purpose.

· **Contact:**

· **Date of last revision/ revision number:** 05/05/2021 / 16

· **Abbreviations and acronyms:**

ADR: The European Agreement concerning the International Carriage of Dangerous Goods by Road

ADN: The European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

VOC: Volatile Organic Compounds (USA, EU)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety and Health

OSHA: Occupational Safety & Health Administration

TLV: Threshold Limit Value

PEL: Permissible Exposure Limit

REL: Recommended Exposure Limit

Eye Irrit. 2B: Serious eye damage/eye irritation – Category 2B

· *** Data compared to the previous version altered.**

SDS created by MSDS Authoring Services www.msdsauthoring.com +1-877-204-9106

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 04/20/2021

Reviewed on 04/20/2021

1 Identification


- **Product Identifier**
- **Trade Name: RNase Enzyme**
RNase Enzyme
- **Product Number:**
3000605
3000593
- **Relevant identified uses of the substance or mixture and uses advised against:**
- **Product Description:** No further relevant information available.
- **Details of the Supplier of the Safety Data Sheet:**
- **Manufacturer/Supplier:**
10x Genomics
6230 Stoneridge Mall Road
Pleasanton, CA 94588 USA
1 925 401 7300
<https://www.10xgenomics.com/>
- **Emergency telephone number:** 1 925 401 7300

2 Hazard(s) Identification

- **Classification of the substance or mixture:**
The product does not need classification according to OSHA HazCom Standard 29 CFR paragraph (d) of §1910.1200(g) and GHS Rev 03.
- **Label elements:**
- **Hazard pictograms:** Non-Regulated Material
- **Signal word:** Non-Regulated Material
- **Hazard statements:** Non-Regulated Material
- **Unknown acute toxicity:**
This value refers to knowledge of known, established toxicological or ecotoxicological values.
0 % of the mixture consists of component(s) of unknown toxicity.
- **Classification system:** NFPA/HMIS Definitions: 0-Least, 1-Slight, 2-Moderate, 3-High, 4-Extreme
- **NFPA ratings (scale 0 - 4)**



- **HMIS-ratings (scale 0 - 4)**



- **Hazard(s) not otherwise classified (HNOC):** None known

3 Composition/Information on Ingredients

- **Chemical characterization: Substance**
- **Description:** Mixture of substances listed below with non-hazardous additions.

- **Dangerous Components:**

Trade Secret

25-50%

(Contd. on page 2)

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 04/20/2021

Reviewed on 04/20/2021

Trade Name: RNase Enzyme
RNase Enzyme

· **Additional information:**

Trade secret made in accordance with paragraph (i) of §1910.1200 of 29 CFR 1910.1200, the OSHA Hazard Communication Standard and U.S. Code of Federal Regulations.

The exact percentages of the ingredients of this mixture are considered to be proprietary and are withheld in accordance with the provisions of paragraph (i) of §1910.1200 of 29 CFR 1910.1200 Trade Secrets.

4 First-Aid Measures

· **Description of first aid measures**

· **General information:** If symptoms persist, call a physician.

· **After inhalation:** Supply fresh air; consult doctor in case of complaints.

· **After skin contact:** If skin irritation occurs, consult a doctor.

· **After eye contact:**

If eye irritation occurs, consult a doctor.

Rinse opened eye for several minutes under running water.

· **After swallowing:** If swallowed and symptoms occur, consult a doctor.

· **Information for doctor**

· **Most important symptoms and effects, both acute and delayed:** No further relevant information available.

· **Indication of any immediate medical attention and special treatment needed:**

No further relevant information available.

5 Fire-Fighting Measures

· **Extinguishing media**

· **Suitable extinguishing agents:** Use fire fighting measures that suit the environment.

· **For safety reasons unsuitable extinguishing agents:** No further relevant information.

· **Special hazards arising from the substance or mixture:** No further relevant information available.

· **Advice for firefighters**

· **Special protective equipment for firefighters:**

As in any fire, wear self-contained breathing apparatus pressure-demand (NIOSH approved or equivalent) and full protective gear to prevent contact with skin and eyes.

6 Accidental Release Measures

· **Personal precautions, protective equipment and emergency procedures:** Not required.

· **Environmental precautions:** Dilute with plenty of water.

· **Methods and material for containment and cleaning up:**

Absorb with liquid-binding material (i.e. sand, diatomite, acid binders, universal binders, sawdust).

Dispose of the collected material according to regulations.

· **Reference to other sections:**

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

· **Protective Action Criteria for Chemicals**

· **PAC-1:**

	Trade Secret	45 mg/m ³
	Trade Secret	12 mg/m ³
7786-30-3	Magnesium Chloride	11 mg/m ³

· **PAC-2:**

	Trade Secret	180 mg/m ³
	Trade Secret	130 mg/m ³

(Contd. on page 3)

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 04/20/2021

Reviewed on 04/20/2021

Trade Name: RNase Enzyme
RNase Enzyme

7786-30-3	Magnesium Chloride	120 mg/m ³
· PAC-3:		
	Trade Secret	1,100 mg/m ³
	Trade Secret	790 mg/m ³
7786-30-3	Magnesium Chloride	550 mg/m ³

7 Handling and Storage

- **Handling**
- **Precautions for safe handling:** No special precautions are necessary if used correctly.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** None.
- **Specific end use(s):** No further relevant information available.

8 Exposure Controls/Personal Protection

- **Additional information about design of technical systems:** No further data; see section 7.
- **Control parameters:**

· Components with occupational exposure limits:	
Trade Secret	
PEL	Long-term value: 15* 5** mg/m ³ mist; *total dust **respirable fraction
TLV	TLV withdrawn-insufficient data human occup. exp.

- **Additional information:** The lists that were valid during the creation of this SDS were used as basis.
- **Exposure controls:**
- **Personal protective equipment**
- **General protective and hygienic measures:**
The usual precautionary measures for handling chemicals should be followed.
Wash hands before breaks and at the end of work.
- **Breathing equipment:** Not required.
- **Protection of hands:** Not required.
- **Material of gloves:** Not applicable.
- **Penetration time of glove material:** Not applicable.
- **Eye protection:**



Goggles recommended during refilling.

- **Limitation and supervision of exposure into the environment:** None

(Contd. on page 4)

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 04/20/2021

Reviewed on 04/20/2021

Trade Name: RNase Enzyme
RNase Enzyme

9 Physical and Chemical Properties

· **Information on basic physical and chemical properties**

· **General Information**

· **Appearance:**

· **Form:** Liquid
 · **Color:** White
 · **Odor:** Slight
 · **Odor threshold:** Not determined.

· **pH-value:** Not determined.

· **Change in condition**

· **Melting point/Melting range:** Not determined.

· **Flash point:** ≥190 °C (≥374 °F)

· **Flammability (solid, gaseous):** Not applicable.

· **Ignition temperature:** ≥400 °C (≥752 °F)

· **Decomposition temperature:** Not determined.

· **Auto igniting:** Product is not self-igniting.

· **Danger of explosion:** Product does not present an explosion hazard.

· **Explosion limits:**

· **Lower:** Not determined.
 · **Upper:** Not determined.

· **Vapor pressure @ 20 °C (68 °F):** <0.1 hPa

· **Density @ 20 °C (68 °F):** 1.2679 g/cm³ (10.5806 lbs/gal)

· **Relative density:** Not determined.

· **Vapor density:** Not determined.

· **Evaporation rate:** Not determined.

· **Solubility in / Miscibility with:**

· **Water:** Fully miscible.

· **Partition coefficient (n-octanol/water):** Not determined.

· **Viscosity:**

· **Dynamic:** Not determined.
 · **Kinematic:** Not determined.

· **Solvent content:**

· **Organic solvents:** 50.0 %
 · **VOC content:** 0.00 %

· **Solids content:** 1.1 %

· **Other information:** No further relevant information available.

10 Stability and Reactivity

· **Reactivity:** No further relevant information available.

· **Chemical stability:** Product is stable under normal conditions.

· **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.

(Contd. on page 5)

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 04/20/2021

Reviewed on 04/20/2021

Trade Name: RNase Enzyme
RNase Enzyme

- **Possibility of hazardous reactions:** No dangerous reactions known.
- **Conditions to avoid:** No further relevant information available.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** No dangerous decomposition products known.

11 Toxicological Information

- **Information on toxicological effects:**
- **Acute toxicity:**

· **LD/LC50 values that are relevant for classification:**

Trade Secret		
Oral	LD50	12,600 mg/kg (Rat)
Dermal	LD50	>10,000 mg/kg (Rabbit)
Inhalative	LC50/96 hours	>50,000 mg/l (Trout)

- **Primary irritant effect:**
- **On the skin:** No irritating effect.
- **On the eye:** No irritating effect.
- **Additional toxicological information:**
- **Carcinogenic categories:**

· **IARC (International Agency for Research on Cancer):**

None of the ingredients are listed.

· **NTP (National Toxicology Program):**

None of the ingredients are listed.

· **OSHA-Ca (Occupational Safety & Health Administration):**

None of the ingredients are listed.

12 Ecological Information

- **Toxicity:**
- **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability:** No further relevant information available.
- **Behavior in environmental systems:**
- **Bioaccumulative potential:** No further relevant information available.
- **Mobility in soil:** No further relevant information available.
- **Additional ecological information:**
- **General notes:**
Do not allow undiluted product or product that has not been neutralized to reach ground water, water course or sewage system.
- **Results of PBT and vPvB assessment:**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **Other adverse effects:** No further relevant information available.

13 Disposal Considerations

- **Waste treatment methods**
- **Recommendation:**
Observe all federal, state and local environmental regulations when disposing of this material.

(Contd. on page 6)

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 04/20/2021

Reviewed on 04/20/2021

Trade Name: RNase Enzyme
RNase Enzyme

- **Uncleaned packaging**
- **Recommendation:** Disposal must be made according to official regulations.
- **Recommended cleansing agent:** Water, if necessary with cleansing agents.

14 Transport Information

- **UN-Number:** Non-Regulated Material
- **DOT, ADR/ADN, IMDG, IATA** Non-Regulated Material
- **UN proper shipping name:** Non-Regulated Material
- **DOT, ADR/ADN, IMDG, IATA** Non-Regulated Material
- **Transport hazard class(es):**
- **DOT, ADR/ADN, ADN, IMDG, IATA** Non-Regulated Material
- **Class:** Non-Regulated Material
- **Packing group:**
- **DOT, ADR/ADN, IMDG, IATA** Non-Regulated Material
- **Environmental hazards:** Not applicable.
- **Special precautions for user:** Not applicable.
- **Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code:** Not applicable.
- **UN "Model Regulation":** Non-Regulated Material

15 Regulatory Information

- **Safety, health and environmental regulations/legislation specific for the substance or mixture:**
 No further relevant information available.
- **SARA (Superfund Amendments and Reauthorization):**

Section 355 (extremely hazardous substances):	None of the ingredients are listed.
Section 313 (Specific toxic chemical listings):	None of the ingredients are listed.
TSCA (Toxic Substances Control Act):	All components have the value ACTIVE.
Hazardous Air Pollutants	None of the ingredients are listed.
California Proposition 65:	
Chemicals known to cause cancer:	None of the ingredients are listed.
Chemicals known to cause reproductive toxicity for females:	None of the ingredients are listed.
Chemicals known to cause reproductive toxicity for males:	None of the ingredients are listed.
Chemicals known to cause developmental toxicity:	None of the ingredients are listed.
New Jersey Right-to-Know List:	Trade Secret

(Contd. on page 7)

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 04/20/2021

Reviewed on 04/20/2021

Trade Name: RNase Enzyme
RNase Enzyme

<ul style="list-style-type: none"> · New Jersey Special Hazardous Substance List:
None of the ingredients are listed.
<ul style="list-style-type: none"> · Pennsylvania Right-to-Know List:
Trade Secret
<ul style="list-style-type: none"> · Pennsylvania Special Hazardous Substance List:
None of the ingredients are listed.
<ul style="list-style-type: none"> · Carcinogenic categories:
<ul style="list-style-type: none"> · EPA (Environmental Protection Agency):
None of the ingredients are listed.
<ul style="list-style-type: none"> · TLV (Threshold Limit Value established by ACGIH):
None of the ingredients are listed.
<ul style="list-style-type: none"> · NIOSH-Ca (National Institute for Occupational Safety and Health):
None of the ingredients are listed.
<ul style="list-style-type: none"> · GHS label elements Non-Regulated Material · Hazard pictograms: Non-Regulated Material · Signal word: Non-Regulated Material · Hazard statements: Non-Regulated Material
<ul style="list-style-type: none"> · National regulations:
None of the ingredients are listed.
<ul style="list-style-type: none"> · Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

16 Other Information

The information and recommendations in this safety data sheet are, to the best of our knowledge, accurate as of the date of issue. Nothing herein shall be deemed to create warranty, expressed or implied, and shall not establish a legally valid contractual relationship. It is the responsibility of the user to determine applicability of this information and the suitability of the material or product for any particular purpose.

- **Contact:**
- **Date of last revision/ revision number:** 04/20/2021 / -
- **Abbreviations and acronyms:**
 ADR: The European Agreement concerning the International Carriage of Dangerous Goods by Road
 ADN: The European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways
 IMDG: International Maritime Code for Dangerous Goods
 DOT: US Department of Transportation
 IATA: International Air Transport Association
 EINECS: European Inventory of Existing Commercial Chemical Substances
 ELINCS: European List of Notified Chemical Substances
 CAS: Chemical Abstracts Service (division of the American Chemical Society)
 NFPA: National Fire Protection Association (USA)
 HMIS: Hazardous Materials Identification System (USA)
 VOC: Volatile Organic Compounds (USA, EU)
 LC50: Lethal concentration, 50 percent
 LD50: Lethal dose, 50 percent
 PBT: Persistent, Bioaccumulative and Toxic
 vPvB: very Persistent and very Bioaccumulative
 NIOSH: National Institute for Occupational Safety and Health
 OSHA: Occupational Safety & Health Administration
 TLV: Threshold Limit Value
 PEL: Permissible Exposure Limit
 REL: Recommended Exposure Limit
- *** Data compared to the previous version altered.**

SECTION 1: Identification

1.1. Identification

Product form : Mixtures
Trade name : RNase Buffer
Product code B : 2000551

1.2. Recommended use and restrictions on use

Use of the substance/mixture : Reagent

1.3. Supplier

Manufacturer/Supplier:
10x Genomics
6230 Stoneridge Mall Road Pleasanton,
CA 94588-3260
T: +1 925 401 7300
E: info@10xgenomics.com

1.4. Emergency telephone number

Emergency number : +1 925 401 7300

SECTION 2: Hazard(s) identification

2.1. Classification of the substance or mixture

GHS US classification

Not classified

2.2. GHS Label elements, including precautionary statements

GHS US labeling

No labeling applicable

2.3. Other hazards which do not result in classification

No additional information available

2.4. Unknown acute toxicity (GHS US)

Not applicable

SECTION 3: Composition/Information on ingredients

3.1. Substances

Not applicable

3.2. Mixtures

This mixture does not contain any substances to be mentioned according to the criteria of section 3.2 of HazCom 2012

RNase Buffer B

Safety Data Sheet

according to US HazCom 2012

SECTION 4: First-aid measures

4.1. Description of first aid measures

First-aid measures general	: If you feel unwell, seek medical advice (show the label where possible).
First-aid measures after inhalation	: Remove person to fresh air and keep comfortable for breathing. If experiencing respiratory symptoms: Call a poison center or a doctor.
First-aid measures after skin contact	: Wash with water and soap as a precaution. If skin irritation occurs: Get medical advice/attention.
First-aid measures after eye contact	: Rinse eyes with water as a precaution. Seek medical attention if ill effect or irritation develops.
First-aid measures after ingestion	: Rinse mouth out with water. Do not induce vomiting. If you feel unwell, seek medical advice.

4.2. Most important symptoms and effects (acute and delayed)

Symptoms/effects	: No adverse effects expected.
------------------	--------------------------------

4.3. Immediate medical attention and special treatment, if necessary

Treat symptomatically.

SECTION 5: Fire-fighting measures

5.1. Suitable (and unsuitable) extinguishing media

Suitable extinguishing media	: Use extinguishing media appropriate for surrounding fire.
Unsuitable extinguishing media	: None known.

5.2. Specific hazards arising from the chemical

Fire hazard	: The product is not flammable. Does not sustain combustion. In the event of fire, may decompose: Carbon oxides (CO, CO ₂). Potassium oxides. Chlorides.
Explosion hazard	: No hazard identified.
Hazardous decomposition products in case of fire	: Thermal decomposition can lead to the release of irritating gases and vapors.

5.3. Special protective equipment and precautions for fire-fighters

Firefighting instructions	: Fight fire with normal precautions from a reasonable distance.
Protection during firefighting	: Do not attempt to take action without suitable protective equipment.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

General measures	: No special requirements.
------------------	----------------------------

6.1.1. For non-emergency personnel

Protective equipment	: No special protection required.
Emergency procedures	: No additional risk management measures required.

6.1.2. For emergency responders

Protective equipment	: No special protection required.
Emergency procedures	: No additional risk management measures required.

6.2. Environmental precautions

No special environmental precautions required.

6.3. Methods and material for containment and cleaning up

For containment	: Absorb spillage to prevent material-damage.
Methods for cleaning up	: Wipe up with absorbent material (for example cloth).

RNase Buffer B

Safety Data Sheet

according to US HazCom 2012

Other information : Dispose in a safe manner in accordance with local/national regulations.

6.4. Reference to other sections

For further information refer to section 8: "Exposure controls/personal protection". For disposal of residues refer to section 13 : "Disposal considerations".

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : No special handling advices are necessary.
Hygiene measures : Use good personal hygiene practices.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store in a well-ventilated place. Keep cool.
Incompatible materials : Oxidizing materials.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

RNase Buffer B

No additional information available

8.2. Appropriate engineering controls

Appropriate engineering controls : No special requirements.

8.3. Individual protection measures/Personal protective equipment

Hand protection:

In case of repeated or prolonged exposure : Wear suitable gloves resistant to chemical penetration. Please follow the instructions related to the permeability and the penetration time provided by the manufacturer

Eye protection:

If there is a risk of eye contact: Safety glasses with side shields

Other information:

Do not eat, drink or smoke when using this product.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state : Liquid
Appearance : Clear.
Color : Colourless
Odor : odorless
Odor threshold : No data available
pH : No data available
Melting point : No data available
Freezing point : No data available
Boiling point : No data available
Flash point : No data available

RNase Buffer B

Safety Data Sheet

according to US HazCom 2012

Relative evaporation rate (butyl acetate=1)	: No data available
Flammability (solid, gas)	: Not applicable.
Vapor pressure	: No data available
Relative vapor density at 20 °C	: No data available
Relative density	: No data available
Solubility	: No data available
Partition coefficient n-octanol/water (Log Pow)	: No data available
Auto-ignition temperature	: No data available
Decomposition temperature	: No data available
Viscosity, kinematic	: No data available
Viscosity, dynamic	: No data available
Explosion limits	: No data available
Explosive properties	: No data available
Oxidizing properties	: No data available

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

The product is non-reactive under normal conditions of use, storage and transport.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

Hazardous polymerization will not occur.

10.4. Conditions to avoid

None under recommended storage and handling conditions (see section 7).

10.5. Incompatible materials

Oxidizing materials.

10.6. Hazardous decomposition products

Under normal conditions of storage and use, hazardous decomposition products should not be produced. Thermal decomposition may produce : Carbon oxides (CO, CO₂). Potassium oxides. Chlorides.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity (oral)	: Not classified (Based on available data, the classification criteria are not met)
Acute toxicity (dermal)	: Not classified (Based on available data, the classification criteria are not met)
Acute toxicity (inhalation)	: Not classified (Based on available data, the classification criteria are not met)
Skin corrosion/irritation	: Not classified (Based on available data, the classification criteria are not met)
Serious eye damage/irritation	: Not classified (Based on available data, the classification criteria are not met)
Respiratory or skin sensitization	: Not classified (Based on available data, the classification criteria are not met)
Germ cell mutagenicity	: Not classified (Based on available data, the classification criteria are not met)
Carcinogenicity	: Not classified (Based on available data, the classification criteria are not met)
Reproductive toxicity	: Not classified (Based on available data, the classification criteria are not met)

RNase Buffer B

Safety Data Sheet

according to US HazCom 2012

STOT-single exposure	: Not classified (Based on available data, the classification criteria are not met)
STOT-repeated exposure	: Not classified (Based on available data, the classification criteria are not met)
Aspiration hazard	: Not classified (Based on available data, the classification criteria are not met)
Viscosity, kinematic	: No data available
Symptoms/effects	: No adverse effects expected.
Other information	: Likely routes of exposure: ingestion, inhalation, skin and eye.

SECTION 12: Ecological information

12.1. Toxicity

Ecology - general : This material has not been tested for environmental effects.

12.2. Persistence and degradability

No additional information available

12.3. Bioaccumulative potential

No additional information available

12.4. Mobility in soil

No additional information available

12.5. Other adverse effects

No additional information available

SECTION 13: Disposal considerations

13.1. Disposal methods

Product/Packaging disposal recommendations : Dispose of in accordance with relevant local regulations.

SECTION 14: Transport information

In accordance with DOT / TDG / IMDG / IATA

DOT	TDG	IMDG	IATA
14.1. UN number			
Not regulated for transport			
14.2. Proper Shipping Name			
Not applicable	Not applicable	Not applicable	Not applicable
14.3. Transport hazard class(es)			
Not applicable	Not applicable	Not applicable	Not applicable
14.4. Packing group			
Not applicable	Not applicable	Not applicable	Not applicable
14.5. Environmental hazards			
Not applicable	Not applicable	Not applicable	Not applicable
No supplementary information available			

RNase Buffer B

Safety Data Sheet

according to US HazCom 2012

SECTION 15: Regulatory information

15.1. US Federal regulations

No additional information available

15.2. International regulations

CANADA

No additional information available

EU-Regulations

No additional information available

National regulations

No additional information available

15.3. US State regulations

No additional information available

SECTION 16: Other information

according to US HazCom 2012

Revision date : 27 May 2022

Other information : None.

Safety Data Sheet (SDS), USA

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product.

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 05/13/2021

Reviewed on 05/13/2021

1 Identification

- **Product Identifier**
- **Trade Name: Tissue Removal Enzyme**
Perm Enzyme B
Perm Enzyme B
- **Product Number:**
3000387
3000602
3000553
- **Relevant identified uses of the substance or mixture and uses advised against:**
- **Product Description:** No further relevant information available.
- **Details of the Supplier of the Safety Data Sheet:**
- **Manufacturer/Supplier:**
10x Genomics
6230 Stoneridge Mall Road
Pleasanton, CA 94588 USA
1 925 401 7300
<https://www.10xgenomics.com/>
- **Emergency telephone number:** 1 925 401 7300

2 Hazard(s) Identification

- **Classification of the substance or mixture:**



STOT SE 3 H335 May cause respiratory irritation.

- **Label elements:**
- **Hazard pictograms:**



- **Signal word:** Warning
- **Hazard-determining components of labeling:**
Trade Secret
- **Hazard statements:**
H335 May cause respiratory irritation.
- **Precautionary statements:**
P261 Avoid breathing dust/fume/gas/mist/vapors/spray.
P271 Use only outdoors or in a well-ventilated area.
P304+P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.
P312 Call a poison center/doctor if you feel unwell.
P403+P233 Store in a well-ventilated place. Keep container tightly closed.
P405 Store locked up.
P501 Dispose of contents/container in accordance with local/regional/national/international regulations.
- **Unknown acute toxicity:**
This value refers to knowledge of known, established toxicological or ecotoxicological values.
0 % of the mixture consists of component(s) of unknown toxicity.
- **Classification system:** NFPA/HMIS Definitions: 0-Least, 1-Slight, 2-Moderate, 3-High, 4-Extreme

(Contd. on page 2)

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 05/13/2021


Reviewed on 05/13/2021

Trade Name: Tissue Removal Enzyme
Perm Enzyme B
Perm Enzyme B

· **NFPA ratings (scale 0 - 4)**



· **HMIS-ratings (scale 0 - 4)**



· **Hazard(s) not otherwise classified (HNOC):** None known

3 Composition/Information on Ingredients

- **Chemical characterization: Substance**
- **Description:** Non-Regulated Material

· **Dangerous Components:**

Trade Secret	⚠ Resp. Sens. 1, H334; ⚠ Skin Irrit. 2, H315; Eye Irrit. 2A, H319; STOT SE 3, H335	0-12%
--------------	--	-------

· **Additional information:**

Trade secret made in accordance with paragraph (i) of §1910.1200 of 29 CFR 1910.1200, the OSHA Hazard Communication Standard and U.S. Code of Federal Regulations.
 The exact percentages of the ingredients of this mixture are considered to be proprietary and are withheld in accordance with the provisions of paragraph (i) of §1910.1200 of 29 CFR 1910.1200 Trade Secrets.

4 First-Aid Measures

- **Description of first aid measures**
- **After inhalation:** Supply fresh air; consult doctor in case of complaints.
- **After skin contact:** If skin irritation occurs, consult a doctor.
- **After eye contact:**
 If eye irritation occurs, consult a doctor.
 Rinse opened eye for several minutes under running water.
- **After swallowing:** If swallowed and symptoms occur, consult a doctor.
- **Information for doctor**
- **Most important symptoms and effects, both acute and delayed:** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed:**
 No further relevant information available.

5 Fire-Fighting Measures

- **Extinguishing media**
- **Suitable extinguishing agents:** Use fire fighting measures that suit the environment.
- **For safety reasons unsuitable extinguishing agents:** No further relevant information.
- **Special hazards arising from the substance or mixture:** No further relevant information available.
- **Advice for firefighters**
- **Special protective equipment for firefighters:**
 As in any fire, wear self-contained breathing apparatus pressure-demand (NIOSH approved or equivalent) and full protective gear to prevent contact with skin and eyes.

(Contd. on page 3)

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 05/13/2021

Reviewed on 05/13/2021

Trade Name: Tissue Removal Enzyme
Perm Enzyme B
Perm Enzyme B

6 Accidental Release Measures

- **Personal precautions, protective equipment and emergency procedures:** Not required.
- **Environmental precautions:** No special measures required.
- **Methods and material for containment and cleaning up:**
 Absorb with liquid-binding material (i.e. sand, diatomite, acid binders, universal binders, sawdust).
 Dispose of the collected material according to regulations.
- **Reference to other sections:**
 See Section 7 for information on safe handling.
 See Section 8 for information on personal protection equipment.
 See Section 13 for disposal information.
- **Protective Action Criteria for Chemicals**

· **PAC-1:**

None of the ingredients are listed.

· **PAC-2:**

None of the ingredients are listed.

· **PAC-3:**

None of the ingredients are listed.

7 Handling and Storage

- **Handling**
- **Precautions for safe handling:** No special precautions are necessary if used correctly.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** None.
- **Specific end use(s):** No further relevant information available.

8 Exposure Controls/Personal Protection

- **Additional information about design of technical systems:** No further data; see section 7.
- **Control parameters:**
- **Components with occupational exposure limits:**
 The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:** The lists that were valid during the creation of this SDS were used as basis.
- **Exposure controls:**
- **Personal protective equipment**
- **General protective and hygienic measures:**
 The usual precautionary measures for handling chemicals should be followed.
 Wash hands before breaks and at the end of work.
- **Breathing equipment:** Not required.
- **Protection of hands:** Not required.
- **Material of gloves:** Not applicable.
- **Penetration time of glove material:** Not applicable.

(Contd. on page 4)

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 05/13/2021

Reviewed on 05/13/2021

Trade Name: Tissue Removal Enzyme
Perm Enzyme B
Perm Enzyme B

· **Eye protection:**



Goggles recommended during refilling.

· **Limitation and supervision of exposure into the environment:** None

9 Physical and Chemical Properties

· **Information on basic physical and chemical properties**

· **General Information**

· **Appearance:**

Form: Liquid
Color: Clear
Odor: Odorless
Odor threshold: Not determined.

· **pH-value @ 20 °C (68 °F):** 7.5

· **Change in condition**

Melting point/Melting range: Not determined.

· **Flash point:** None

· **Flammability (solid, gaseous):** Not applicable.

· **Ignition temperature:** Not applicable

· **Decomposition temperature:** Not determined.

· **Auto igniting:** Product is not self-igniting.

· **Danger of explosion:** Product does not present an explosion hazard.

· **Explosion limits:**

Lower: Not determined.
Upper: Not determined.

· **Vapor pressure:** Not determined.

· **Density @ 20 °C (68 °F):** 1.1 g/cm³ (9.1795 lbs/gal)

· **Relative density:** Not determined.

· **Vapor density:** Not determined.

· **Evaporation rate:** Not determined.

· **Solubility in / Miscibility with:**

Water: Not miscible or difficult to mix.

· **Partition coefficient (n-octanol/water):** Not determined.

· **Viscosity:**

Dynamic: Not determined.
Kinematic: Not determined.

· **Solvent content:**

VOC content: 0.00 %

Solids content: 0.0 %

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 05/13/2021

Reviewed on 05/13/2021

Trade Name: Tissue Removal Enzyme
Perm Enzyme B
Perm Enzyme B

· **Other information:** No further relevant information available.

10 Stability and Reactivity

- **Reactivity:** No further relevant information available.
- **Chemical stability:**
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions:** No dangerous reactions known.
- **Conditions to avoid:** No further relevant information available.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** No dangerous decomposition products known.

11 Toxicological Information

- **Information on toxicological effects:**
- **Acute toxicity:**
- **LD/LC50 values that are relevant for classification:** No data available.
- **Primary irritant effect:**
- **On the skin:** No irritating effect.
- **On the eye:** No irritating effect.
- **Additional toxicological information:**
- **Carcinogenic categories:**

· **IARC (International Agency for Research on Cancer):**

None of the ingredients are listed.

· **NTP (National Toxicology Program):**

None of the ingredients are listed.

· **OSHA-Ca (Occupational Safety & Health Administration):**

None of the ingredients are listed.

12 Ecological Information

- **Toxicity:**
- **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability:** No further relevant information available.
- **Behavior in environmental systems:**
- **Bioaccumulative potential:** No further relevant information available.
- **Mobility in soil:** No further relevant information available.
- **Additional ecological information:**
- **General notes:** Generally not hazardous for water.
- **Results of PBT and vPvB assessment:**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **Other adverse effects:** No further relevant information available.

13 Disposal Considerations

- **Waste treatment methods**
- **Recommendation:**
 Observe all federal, state and local environmental regulations when disposing of this material.

(Contd. on page 6)

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 05/13/2021

Reviewed on 05/13/2021

Trade Name: Tissue Removal Enzyme
Perm Enzyme B
Perm Enzyme B

- **Uncleaned packaging**
- **Recommendation:** Disposal must be made according to official regulations.

14 Transport Information

- **UN-Number:** Non-Regulated Material
- **DOT, ADR/ADN, IMDG, IATA** Non-Regulated Material
- **UN proper shipping name:** Non-Regulated Material
- **DOT, ADR/ADN, IMDG, IATA** Non-Regulated Material
- **Transport hazard class(es):**
- **DOT, ADR/ADN, ADN, IMDG, IATA** Non-Regulated Material
- **Class:** Non-Regulated Material
- **Packing group:**
- **DOT, ADR/ADN, IMDG, IATA** Non-Regulated Material
- **Environmental hazards:** Not applicable.
- **Special precautions for user:** Not applicable.
- **Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code:** Not applicable.
- **UN "Model Regulation":** Non-Regulated Material

15 Regulatory Information

- **Safety, health and environmental regulations/legislation specific for the substance or mixture:**
 No further relevant information available.
- **SARA (Superfund Amendments and Reauthorization):**

Section 355 (extremely hazardous substances):	None of the ingredients are listed.
Section 313 (Specific toxic chemical listings):	None of the ingredients are listed.
TSCA (Toxic Substances Control Act):	None of the ingredients are listed.
Hazardous Air Pollutants	None of the ingredients are listed.
California Proposition 65:	
Chemicals known to cause cancer:	None of the ingredients are listed.
Chemicals known to cause reproductive toxicity for females:	None of the ingredients are listed.
Chemicals known to cause reproductive toxicity for males:	None of the ingredients are listed.
Chemicals known to cause developmental toxicity:	None of the ingredients are listed.
New Jersey Right-to-Know List:	None of the ingredients are listed.

(Contd. on page 7)

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 05/13/2021

Reviewed on 05/13/2021

Trade Name: Tissue Removal Enzyme
Perm Enzyme B
Perm Enzyme B

· New Jersey Special Hazardous Substance List:
None of the ingredients are listed.
· Pennsylvania Right-to-Know List:
None of the ingredients are listed.
· Pennsylvania Special Hazardous Substance List:
None of the ingredients are listed.
· Carcinogenic categories:
· EPA (Environmental Protection Agency):
None of the ingredients are listed.
· TLV (Threshold Limit Value established by ACGIH):
None of the ingredients are listed.
· NIOSH-Ca (National Institute for Occupational Safety and Health):
None of the ingredients are listed.

- **GHS label elements**
The product is classified and labeled according to the Globally Harmonized System (GHS).
- **Hazard pictograms:**



- **Signal word:** Warning
- **Hazard-determining components of labeling:**
Trade Secret
- **Hazard statements:**
H335 May cause respiratory irritation.
- **Precautionary statements:**
P261 Avoid breathing dust/fume/gas/mist/vapors/spray.
P271 Use only outdoors or in a well-ventilated area.
P304+P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.
P312 Call a poison center/doctor if you feel unwell.
P403+P233 Store in a well-ventilated place. Keep container tightly closed.
P405 Store locked up.
P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

· National regulations:
None of the ingredients are listed.

- **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

16 Other Information

The information and recommendations in this safety data sheet are, to the best of our knowledge, accurate as of the date of issue. Nothing herein shall be deemed to create warranty, expressed or implied, and shall not establish a legally valid contractual relationship. It is the responsibility of the user to determine applicability of this information and the suitability of the material or product for any particular purpose.

- **Contact:**
- **Date of last revision/ revision number:** 05/13/2021 / 3

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 05/13/2021

Reviewed on 05/13/2021

Trade Name: Tissue Removal Enzyme
Perm Enzyme B
Perm Enzyme B

· Abbreviations and acronyms:

ADR: The European Agreement concerning the International Carriage of Dangerous Goods by Road
ADN: The European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways
IMDG: International Maritime Code for Dangerous Goods
DOT: US Department of Transportation
IATA: International Air Transport Association
EINECS: European Inventory of Existing Commercial Chemical Substances
ELINCS: European List of Notified Chemical Substances
CAS: Chemical Abstracts Service (division of the American Chemical Society)
NFPA: National Fire Protection Association (USA)
HMIS: Hazardous Materials Identification System (USA)
VOC: Volatile Organic Compounds (USA, EU)
LC50: Lethal concentration, 50 percent
LD50: Lethal dose, 50 percent
PBT: Persistent, Bioaccumulative and Toxic
vPvB: very Persistent and very Bioaccumulative
NIOSH: National Institute for Occupational Safety and Health
OSHA: Occupational Safety & Health Administration
TLV: Threshold Limit Value
PEL: Permissible Exposure Limit
REL: Recommended Exposure Limit
Skin Irrit. 2: Skin corrosion/irritation – Category 2
Eye Irrit. 2A: Serious eye damage/eye irritation – Category 2A
Resp. Sens. 1: Respiratory sensitisation – Category 1
STOT SE 3: Specific target organ toxicity (single exposure) – Category 3

· * Data compared to the previous version altered.SDS created by MSDS Authoring Services www.msdsauthoring.com +1-877-204-9106

SECTION 1: Identification

1.1. Identification

Product form : Mixtures
Trade name : Tissue Removal Buffer B
Product code : 2000543

1.2. Recommended use and restrictions on use

Use of the substance/mixture : Reagent

1.3. Supplier

Manufacturer/Supplier:
10x Genomics
6230 Stoneridge Mall Road Pleasanton,
CA 94588-3260
T:+1 925 401 7300
E: info@10xgenomics.com

1.4. Emergency telephone number

Emergency number : +1 925 401 7300

SECTION 2: Hazard(s) identification

2.1. Classification of the substance or mixture

GHS US classification

Acute toxicity (inhalation:dust,mist) Category 4
Serious eye damage/eye irritation Category 2

Harmful if inhaled
Causes serious eye irritation

2.2. GHS Label elements, including precautionary statements

GHS US labeling

Hazard pictograms (GHS US) :



Signal word (GHS US) :

Warning

Hazard statements (GHS US) :

Causes serious eye irritation
Harmful if inhaled

Precautionary statements (GHS US) :

Avoid breathing mist, spray.
Wash hands thoroughly after handling.
Use only outdoors or in a well-ventilated area.
Wear eye protection.
If inhaled: Remove person to fresh air and keep comfortable for breathing.
IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
Call a poison center or doctor if you feel unwell.
If eye irritation persists: Get medical advice/attention.

2.3. Other hazards which do not result in classification

No additional information available

Tissue Removal Buffer B

Safety Data Sheet

according to US HazCom 2012

2.4. Unknown acute toxicity (GHS US)

Not applicable

SECTION 3: Composition/Information on ingredients

3.1. Substances

Not applicable

3.2. Mixtures

Name	Product identifier	%	GHS US classification
SODIUM LAUROYL SARCOSINATE	CAS-No.: 137-16-6	3 – 10	Acute Tox. 2 (Inhalation), H330 Skin Irrit. 2, H315 Eye Dam. 1, H318

*Chemical name, CAS number and/or exact concentration have been withheld as a trade secret

SECTION 4: First-aid measures

4.1. Description of first aid measures

First-aid measures general	: Call a poison center/doctor/physician if you feel unwell.
First-aid measures after inhalation	: Remove person to fresh air and keep comfortable for breathing. Call a poison center/doctor/physician if you feel unwell.
First-aid measures after skin contact	: Seek medical attention if ill effect or irritation develops. Wash with water and soap as a precaution.
First-aid measures after eye contact	: Rinse immediately with plenty of water for 15 minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Obtain medical attention if pain, blinking or redness persists.
First-aid measures after ingestion	: Rinse mouth. Do NOT induce vomiting. Call a poison center/doctor/physician if you feel unwell.

4.2. Most important symptoms and effects (acute and delayed)

Symptoms/effects after inhalation	: Harmful if inhaled.
Symptoms/effects after eye contact	: Eye irritation.

4.3. Immediate medical attention and special treatment, if necessary

Treat symptomatically.

SECTION 5: Fire-fighting measures

5.1. Suitable (and unsuitable) extinguishing media

Suitable extinguishing media	: Use extinguishing media appropriate for surrounding fire.
Unsuitable extinguishing media	: None known.

5.2. Specific hazards arising from the chemical

Fire hazard	: The product is not flammable. Does not sustain combustion. In the event of fire, may decompose: Carbon oxides (CO, CO ₂). Nitrogen oxides. Chlorides.
Explosion hazard	: No hazard identified.
Hazardous decomposition products in case of fire	: Thermal decomposition can lead to the release of irritating gases and vapors.

5.3. Special protective equipment and precautions for fire-fighters

Firefighting instructions	: Fight fire with normal precautions from a reasonable distance.
Protection during firefighting	: Do not attempt to take action without suitable protective equipment.

Tissue Removal Buffer B

Safety Data Sheet

according to US HazCom 2012

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

General measures : Avoid breathing mist, spray. Avoid contact with eyes.

6.1.1. For non-emergency personnel

Protective equipment : Wear personal protective equipment. For further information refer to section 8: "Exposure controls/personal protection".

Emergency procedures : Evacuate area. Ventilate spillage area.

6.1.2. For emergency responders

Protective equipment : Wear recommended personal protective equipment. For further information refer to section 8: "Exposure controls/personal protection".

Emergency procedures : Stop leak if safe to do so. Evacuate unnecessary personnel. Ventilate area.

6.2. Environmental precautions

Avoid release to the environment.

6.3. Methods and material for containment and cleaning up

For containment : Soak up with inert absorbent material (for example sand, sawdust, a universal binder, silica gel).

Methods for cleaning up : Collect spillage. Store away from other materials.

Other information : Dispose in a safe manner in accordance with local/national regulations.

6.4. Reference to other sections

For further information refer to section 8: "Exposure controls/personal protection". For disposal of residues refer to section 13: "Disposal considerations".

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : Ensure adequate air ventilation. Use only outdoors or in a well-ventilated area. Avoid breathing mist, spray. Avoid contact with eyes. Wear personal protective equipment.

Hygiene measures : Wash hands and other exposed areas with mild soap and water before eating, drinking or smoking and when leaving work. Do not eat, drink or smoke when using this product. Always wash hands after handling the product. Handle in accordance with good industrial hygiene and safety practice.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Keep only in the original container in a cool well ventilated place. Keep container closed when not in use. Store in a well-ventilated place. Keep cool.

Incompatible materials : Oxidizing materials.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Tissue Removal Buffer B

No additional information available

SODIUM LAUROYL SARCOSINATE (137-16-6)

No additional information available

8.2. Appropriate engineering controls

Appropriate engineering controls : Ensure good ventilation of the work station.

Environmental exposure controls : Avoid release to the environment.

Tissue Removal Buffer B

Safety Data Sheet

according to US HazCom 2012

8.3. Individual protection measures/Personal protective equipment

Hand protection:

In case of repeated or prolonged exposure: Wear suitable gloves resistant to chemical penetration. Please follow the instructions related to the permeability and the penetration time provided by the manufacturer

Eye protection:

If there is a risk of eye contact: Safety glasses.

Respiratory protection:

In case of insufficient ventilation, wear suitable respiratory equipment. Appropriate dust or mist respirator should be used if airborne particles are generated when handling this material.

Other information:

Do not eat, drink or smoke during use.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state	: Liquid
Appearance	: Clear.
Color	: Colourless
Odor	: odorless
Odor threshold	: No data available
pH	: No data available
Melting point	: No data available
Freezing point	: No data available
Boiling point	: No data available
Flash point	: No data available
Relative evaporation rate (butyl acetate=1)	: No data available
Flammability (solid, gas)	: Not applicable.
Vapor pressure	: No data available
Relative vapor density at 20 °C	: No data available
Relative density	: No data available
Solubility	: No data available
Partition coefficient n-octanol/water (Log Pow)	: No data available
Auto-ignition temperature	: No data available
Decomposition temperature	: No data available
Viscosity, kinematic	: No data available
Viscosity, dynamic	: No data available
Explosion limits	: No data available
Explosive properties	: No data available
Oxidizing properties	: No data available

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

The product is non-reactive under normal conditions of use, storage and transport.

Tissue Removal Buffer B

Safety Data Sheet

according to US HazCom 2012

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

No dangerous reactions known under normal conditions of use.

10.4. Conditions to avoid

None under recommended storage and handling conditions (see section 7).

10.5. Incompatible materials

Oxidizing materials.

10.6. Hazardous decomposition products

Under normal conditions of storage and use, hazardous decomposition products should not be produced. In the event of fire, may decompose: Carbon oxides (CO, CO₂). Nitrogen oxides. Chlorides

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity (oral) : Not classified (Based on available data, the classification criteria are not met)
Acute toxicity (dermal) : Not classified (Based on available data, the classification criteria are not met)
Acute toxicity (inhalation) : Harmful if inhaled.

Tissue Removal Buffer B

ATE US (dust, mist)	1.5 mg/l/4h
---------------------	-------------

SODIUM LAUROYL SARCOSINATE (137-16-6)

ATE US (gases)	100 ppmV/4h
----------------	-------------

ATE US (vapors)	0.5 mg/l/4h
-----------------	-------------

ATE US (dust, mist)	0.05 mg/l/4h
---------------------	--------------

Skin corrosion/irritation : Not classified (Based on available data, the classification criteria are not met)
Serious eye damage/irritation : Causes serious eye irritation.
Respiratory or skin sensitization : Not classified (Based on available data, the classification criteria are not met)
Germ cell mutagenicity : Not classified (Based on available data, the classification criteria are not met)
Carcinogenicity : Not classified (Based on available data, the classification criteria are not met)
Reproductive toxicity : Not classified (Based on available data, the classification criteria are not met)
STOT-single exposure : Not classified (Based on available data, the classification criteria are not met)
STOT-repeated exposure : Not classified (Based on available data, the classification criteria are not met)
Aspiration hazard : Not classified (Based on available data, the classification criteria are not met)
Viscosity, kinematic : No data available
Symptoms/effects after inhalation : Harmful if inhaled.
Symptoms/effects after eye contact : Eye irritation.
Other information : Likely routes of exposure: ingestion, inhalation, skin and eye.

SECTION 12: Ecological information

12.1. Toxicity

No additional information available

Tissue Removal Buffer B

Safety Data Sheet

according to US HazCom 2012

12.2. Persistence and degradability

No additional information available

12.3. Bioaccumulative potential

No additional information available

12.4. Mobility in soil

No additional information available

12.5. Other adverse effects

Other information : Avoid release to the environment.

SECTION 13: Disposal considerations

13.1. Disposal methods

Waste treatment methods : Dispose of contents/container in accordance with licensed collector's sorting instructions.
Product/Packaging disposal recommendations : Dispose in a safe manner in accordance with local/national regulations.
Ecology - waste materials : Avoid release to the environment.

SECTION 14: Transport information

In accordance with DOT / TDG / IMDG / IATA

DOT	TDG	IMDG	IATA
14.1. UN number			
Not regulated for transport			
14.2. Proper Shipping Name			
Not applicable	Not applicable	Not applicable	Not applicable
14.3. Transport hazard class(es)			
Not applicable	Not applicable	Not applicable	Not applicable
14.4. Packing group			
Not applicable	Not applicable	Not applicable	Not applicable
14.5. Environmental hazards			
Not applicable	Not applicable	Not applicable	Not applicable
No supplementary information available			

SECTION 15: Regulatory information

15.1. US Federal regulations

No additional information available

Tissue Removal Buffer B

Safety Data Sheet

according to US HazCom 2012

15.2. International regulations

CANADA

SODIUM LAUROYL SARCOSINATE (137-16-6)

Listed on the Canadian DSL (Domestic Substances List)

EU-Regulations

No additional information available

National regulations

SODIUM LAUROYL SARCOSINATE (137-16-6)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

15.3. US State regulations

No additional information available

SECTION 16: Other information

according to US HazCom 2012

Revision date : 6 June 2022

Other information : None.

Safety Data Sheet (SDS), USA

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product.

SECTION 1: Identification

1.1. Identification

Product form : Mixtures
Trade name : Tissue Removal Buffer Enhancer
Product code : 2000557

1.2. Recommended use and restrictions on use

Use of the substance/mixture : Reagent

1.3. Supplier

Manufacturer/Supplier
10x Genomics
6230 Stoneridge Mall Road Pleasanton,
CA 94588-3260
T: +1 925 401 7300
E: info@10xgenomics.com

1.4. Emergency telephone number

Emergency number : +1 925 401 7300

SECTION 2: Hazard(s) identification

2.1. Classification of the substance or mixture

GHS US classification

Not classified

2.2. GHS Label elements, including precautionary statements

GHS US labeling

No labeling applicable

2.3. Other hazards which do not result in classification

No additional information available

2.4. Unknown acute toxicity (GHS US)

Not applicable

SECTION 3: Composition/Information on ingredients

3.1. Substances

Not applicable

3.2. Mixtures

This mixture does not contain any substances to be mentioned according to the criteria of section 3.2 of HazCom 2012

Tissue Removal Buffer Enhancer

Safety Data Sheet

according to US HazCom 2012

SECTION 4: First-aid measures

4.1. Description of first aid measures

First-aid measures general	: If you feel unwell, seek medical advice (show the label where possible).
First-aid measures after inhalation	: Remove person to fresh air and keep comfortable for breathing. If experiencing respiratory symptoms: Call a poison center or a doctor.
First-aid measures after skin contact	: Wash with water and soap as a precaution. If skin irritation occurs: Get medical advice/attention.
First-aid measures after eye contact	: Rinse eyes with water as a precaution. Seek medical attention if ill effect or irritation develops.
First-aid measures after ingestion	: Rinse mouth out with water. Do not induce vomiting. If you feel unwell, seek medical advice.

4.2. Most important symptoms and effects (acute and delayed)

Symptoms/effects	: No adverse effects expected.
------------------	--------------------------------

4.3. Immediate medical attention and special treatment, if necessary

Treat symptomatically.

SECTION 5: Fire-fighting measures

5.1. Suitable (and unsuitable) extinguishing media

Suitable extinguishing media	: Use extinguishing media appropriate for surrounding fire.
Unsuitable extinguishing media	: None known.

5.2. Specific hazards arising from the chemical

Fire hazard	: The product is not flammable. In the event of fire, may decompose : Carbon oxides (CO, CO ₂), Nitrogen oxides. Ammonia.
Explosion hazard	: No hazard identified.
Hazardous decomposition products in case of fire	: Thermal decomposition can lead to the release of irritating gases and vapors.

5.3. Special protective equipment and precautions for fire-fighters

Firefighting instructions	: Fight fire with normal precautions from a reasonable distance.
Protection during firefighting	: Do not attempt to take action without suitable protective equipment.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

General measures	: No special requirements.
------------------	----------------------------

6.1.1. For non-emergency personnel

Protective equipment	: No special protection required.
Emergency procedures	: No additional risk management measures required.

6.1.2. For emergency responders

Protective equipment	: No special protection required.
Emergency procedures	: No additional risk management measures required.

6.2. Environmental precautions

No special environmental precautions required.

6.3. Methods and material for containment and cleaning up

For containment	: Absorb spillage to prevent material-damage.
Methods for cleaning up	: Wipe up with absorbent material (for example cloth).

Tissue Removal Buffer Enhancer

Safety Data Sheet

according to US HazCom 2012

Other information : Dispose in a safe manner in accordance with local/national regulations.

6.4. Reference to other sections

For further information refer to section 8: "Exposure controls/personal protection". For disposal of residues refer to section 13 : "Disposal considerations".

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : No special handling advices are necessary.
Hygiene measures : Use good personal hygiene practices.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store in a well-ventilated place. Keep cool.
Incompatible materials : Acids. Strong reducing agents. oxidizing agents (peroxides, chromates, dichromates). hypochlorites. Nitrates. Fluorine.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Tissue Removal Buffer Enhancer

No additional information available

8.2. Appropriate engineering controls

Appropriate engineering controls : No special requirements.

8.3. Individual protection measures/Personal protective equipment

Hand protection:

In case of repeated or prolonged exposure : Wear suitable gloves resistant to chemical penetration. Please follow the instructions related to the permeability and the penetration time provided by the manufacturer

Eye protection:

If there is a risk of eye contact: Safety glasses with side shields

Other information:

Do not eat, drink or smoke when using this product.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state : Liquid
Appearance : Clear.
Color : Colourless
Odor : odorless
Odor threshold : No data available
pH : No data available
Melting point : No data available
Freezing point : No data available
Boiling point : No data available

Tissue Removal Buffer Enhancer

Safety Data Sheet

according to US HazCom 2012

Flash point	: No data available
Relative evaporation rate (butyl acetate=1)	: No data available
Flammability (solid, gas)	: Not applicable.
Vapor pressure	: No data available
Relative vapor density at 20 °C	: No data available
Relative density	: No data available
Solubility	: No data available
Partition coefficient n-octanol/water (Log Pow)	: No data available
Auto-ignition temperature	: No data available
Decomposition temperature	: No data available
Viscosity, kinematic	: No data available
Viscosity, dynamic	: No data available
Explosion limits	: No data available
Explosive properties	: No data available
Oxidizing properties	: No data available

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

The product is non-reactive under normal conditions of use, storage and transport.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

Hazardous reactions may occur if contact with incompatible material. An exothermic reaction may occur. Toxic fumes may be released.

10.4. Conditions to avoid

None under recommended storage and handling conditions (see section 7).

10.5. Incompatible materials

Acids. Strong reducing agents. oxidizing agents (peroxides, chromates, dichromates). hypochlorites. Nitrates. Fluorine.

10.6. Hazardous decomposition products

Under normal conditions of storage and use, hazardous decomposition products should not be produced. Thermal decomposition may produce : Carbon oxides (CO, CO₂). Nitrogen oxides. Ammonia.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity (oral)	: Not classified (Based on available data, the classification criteria are not met)
Acute toxicity (dermal)	: Not classified (Based on available data, the classification criteria are not met)
Acute toxicity (inhalation)	: Not classified (Based on available data, the classification criteria are not met)
Skin corrosion/irritation	: Not classified (Based on available data, the classification criteria are not met)
Serious eye damage/irritation	: Not classified (Based on available data, the classification criteria are not met)
Respiratory or skin sensitization	: Not classified (Based on available data, the classification criteria are not met)
Germ cell mutagenicity	: Not classified (Based on available data, the classification criteria are not met)
Carcinogenicity	: Not classified (Based on available data, the classification criteria are not met)

Tissue Removal Buffer Enhancer

Safety Data Sheet

according to US HazCom 2012

Reproductive toxicity	: Not classified (Based on available data, the classification criteria are not met)
STOT-single exposure	: Not classified (Based on available data, the classification criteria are not met)
STOT-repeated exposure	: Not classified (Based on available data, the classification criteria are not met)
Aspiration hazard	: Not classified (Based on available data, the classification criteria are not met)
Viscosity, kinematic	: No data available
Symptoms/effects	: No adverse effects expected.
Other information	: Likely routes of exposure: ingestion, inhalation, skin and eye.

SECTION 12: Ecological information

12.1. Toxicity

Ecology - general : This material has not been tested for environmental effects.

12.2. Persistence and degradability

No additional information available

12.3. Bioaccumulative potential

No additional information available

12.4. Mobility in soil

No additional information available

12.5. Other adverse effects

No additional information available

SECTION 13: Disposal considerations

13.1. Disposal methods

Product/Packaging disposal recommendations : Dispose of in accordance with relevant local regulations.

SECTION 14: Transport information

In accordance with DOT / TDG / IMDG / IATA

DOT	TDG	IMDG	IATA
14.1. UN number			
Not regulated for transport			
14.2. Proper Shipping Name			
Not applicable	Not applicable	Not applicable	Not applicable
14.3. Transport hazard class(es)			
Not applicable	Not applicable	Not applicable	Not applicable
14.4. Packing group			
Not applicable	Not applicable	Not applicable	Not applicable
14.5. Environmental hazards			
Not applicable	Not applicable	Not applicable	Not applicable

Tissue Removal Buffer Enhancer

Safety Data Sheet

according to US HazCom 2012

DOT	TDG	IMDG	IATA
No supplementary information available			

SECTION 15: Regulatory information

15.1. US Federal regulations

No additional information available

15.2. International regulations

CANADA

No additional information available

EU-Regulations

No additional information available

National regulations

No additional information available

15.3. US State regulations

No additional information available

SECTION 16: Other information

according to US HazCom 2012

Revision date : 27 May 2022

Other information : None.

Safety Data Sheet (SDS), USA

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product.

SECTION 1: Identification

1.1. Identification

Product form : Mixtures
Trade name : Decrosslinking Buffer
Product code : 2000566

1.2. Recommended use and restrictions on use

Use of the substance/mixture : Reagent

1.3. Supplier

Manufacturer/Supplier
10x Genomics
6230 Stoneridge Mall Road Pleasanton,
CA 94588-3260
T: +1 925 401 7300
E: info@10xgenomics.com

1.4. Emergency telephone number

Emergency number : +1 925 401 7300

SECTION 2: Hazard(s) identification

2.1. Classification of the substance or mixture

GHS US classification

Not classified

2.2. GHS Label elements, including precautionary statements

GHS US labeling

No labeling applicable

2.3. Other hazards which do not result in classification

Other hazards which do not result in classification : Material is hygroscopic.

2.4. Unknown acute toxicity (GHS US)

Not applicable

SECTION 3: Composition/Information on ingredients

3.1. Substances

Not applicable

3.2. Mixtures

This mixture does not contain any substances to be mentioned according to the criteria of section 3.2 of HazCom 2012

Decrosslinking Buffer

Safety Data Sheet

according to US HazCom 2012

SECTION 4: First-aid measures

4.1. Description of first aid measures

First-aid measures general	: If you feel unwell, seek medical advice (show the label where possible).
First-aid measures after inhalation	: Remove person to fresh air and keep comfortable for breathing. If experiencing respiratory symptoms: Call a poison center or a doctor.
First-aid measures after skin contact	: Wash with water and soap as a precaution. If skin irritation occurs: Get medical advice/attention.
First-aid measures after eye contact	: Rinse eyes with water as a precaution. Seek medical attention if ill effect or irritation develops.
First-aid measures after ingestion	: Rinse mouth out with water. Do not induce vomiting. If you feel unwell, seek medical advice.

4.2. Most important symptoms and effects (acute and delayed)

Symptoms/effects : No adverse effects expected.

4.3. Immediate medical attention and special treatment, if necessary

Treat symptomatically.

SECTION 5: Fire-fighting measures

5.1. Suitable (and unsuitable) extinguishing media

Suitable extinguishing media	: Use extinguishing media appropriate for surrounding fire.
Unsuitable extinguishing media	: None known.

5.2. Specific hazards arising from the chemical

Fire hazard	: The product is not flammable. Does not sustain combustion. In the event of fire, may decompose : Carbon oxides (CO, CO ₂). Sodium oxides (NaOx).
Explosion hazard	: No hazard identified.
Hazardous decomposition products in case of fire	: Thermal decomposition can lead to the release of irritating gases and vapors.

5.3. Special protective equipment and precautions for fire-fighters

Firefighting instructions	: Fight fire with normal precautions from a reasonable distance.
Protection during firefighting	: Do not attempt to take action without suitable protective equipment.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

General measures : No special requirements.

6.1.1. For non-emergency personnel

Protective equipment	: No special protection required.
Emergency procedures	: No additional risk management measures required.

6.1.2. For emergency responders

Protective equipment	: No special protection required.
Emergency procedures	: No additional risk management measures required.

6.2. Environmental precautions

No special environmental precautions required.

6.3. Methods and material for containment and cleaning up

For containment	: Absorb spillage to prevent material-damage.
Methods for cleaning up	: Wipe up with absorbent material (for example cloth).

Decrosslinking Buffer

Safety Data Sheet

according to US HazCom 2012

Other information : Dispose in a safe manner in accordance with local/national regulations.

6.4. Reference to other sections

For further information refer to section 8: "Exposure controls/personal protection". For disposal of residues refer to section 13 : "Disposal considerations".

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : No special handling advices are necessary.
Hygiene measures : Use good personal hygiene practices.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store in a well-ventilated place. Keep cool.
Incompatible materials : Oxidizing materials. Reducing agents. Strong bases.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Decrosslinking Buffer

No additional information available

8.2. Appropriate engineering controls

Appropriate engineering controls : No special requirements.

8.3. Individual protection measures/Personal protective equipment

Hand protection:

In case of repeated or prolonged exposure : Wear suitable gloves resistant to chemical penetration. Please follow the instructions related to the permeability and the penetration time provided by the manufacturer

Eye protection:

If there is a risk of eye contact: Safety glasses with side shields

Other information:

Do not eat, drink or smoke when using this product.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state : Liquid
Appearance : Clear.
Color : Colourless
Odor : odorless
Odor threshold : No data available
pH : 6 – 6.5
Melting point : No data available
Freezing point : No data available
Boiling point : No data available
Flash point : No data available

Decrosslinking Buffer

Safety Data Sheet

according to US HazCom 2012

Relative evaporation rate (butyl acetate=1)	: No data available
Flammability (solid, gas)	: Not applicable.
Vapor pressure	: No data available
Relative vapor density at 20 °C	: No data available
Relative density	: No data available
Solubility	: No data available
Partition coefficient n-octanol/water (Log Pow)	: No data available
Auto-ignition temperature	: No data available
Decomposition temperature	: No data available
Viscosity, kinematic	: No data available
Viscosity, dynamic	: No data available
Explosion limits	: No data available
Explosive properties	: No data available
Oxidizing properties	: No data available

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

The product is non-reactive under normal conditions of use, storage and transport.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

Hazardous polymerization will not occur.

10.4. Conditions to avoid

None under recommended storage and handling conditions (see section 7).

10.5. Incompatible materials

Oxidizing materials. Reducing agents. Strong bases.

10.6. Hazardous decomposition products

Under normal conditions of storage and use, hazardous decomposition products should not be produced. Thermal decomposition may produce : Carbon oxides (CO, CO₂). Sodium oxides (NaOx).

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity (oral)	: Not classified (Based on available data, the classification criteria are not met)
Acute toxicity (dermal)	: Not classified (Based on available data, the classification criteria are not met)
Acute toxicity (inhalation)	: Not classified (Based on available data, the classification criteria are not met)
Skin corrosion/irritation	: Not classified (Based on available data, the classification criteria are not met) pH: 6 – 6.5
Serious eye damage/irritation	: Not classified (Based on available data, the classification criteria are not met) pH: 6 – 6.5
Respiratory or skin sensitization	: Not classified (Based on available data, the classification criteria are not met)
Germ cell mutagenicity	: Not classified (Based on available data, the classification criteria are not met)
Carcinogenicity	: Not classified (Based on available data, the classification criteria are not met)

Decrosslinking Buffer

Safety Data Sheet

according to US HazCom 2012

Reproductive toxicity	: Not classified (Based on available data, the classification criteria are not met)
STOT-single exposure	: Not classified (Based on available data, the classification criteria are not met)
STOT-repeated exposure	: Not classified (Based on available data, the classification criteria are not met)
Aspiration hazard	: Not classified (Based on available data, the classification criteria are not met)
Viscosity, kinematic	: No data available
Symptoms/effects	: No adverse effects expected.
Other information	: Likely routes of exposure: ingestion, inhalation, skin and eye.

SECTION 12: Ecological information

12.1. Toxicity

Ecology - general : This material has not been tested for environmental effects.

12.2. Persistence and degradability

No additional information available

12.3. Bioaccumulative potential

No additional information available

12.4. Mobility in soil

No additional information available

12.5. Other adverse effects

No additional information available

SECTION 13: Disposal considerations

13.1. Disposal methods

Product/Packaging disposal recommendations : Dispose of in accordance with relevant local regulations.

SECTION 14: Transport information

In accordance with DOT / TDG / IMDG / IATA

DOT	TDG	IMDG	IATA
14.1. UN number			
Not regulated for transport			
14.2. Proper Shipping Name			
Not applicable	Not applicable	Not applicable	Not applicable
14.3. Transport hazard class(es)			
Not applicable	Not applicable	Not applicable	Not applicable
14.4. Packing group			
Not applicable	Not applicable	Not applicable	Not applicable
14.5. Environmental hazards			
Not applicable	Not applicable	Not applicable	Not applicable
No supplementary information available			

Decrosslinking Buffer

Safety Data Sheet

according to US HazCom 2012

SECTION 15: Regulatory information

15.1. US Federal regulations

No additional information available

15.2. International regulations

CANADA

No additional information available

EU-Regulations

No additional information available

National regulations

No additional information available

15.3. US State regulations

No additional information available

SECTION 16: Other information

according to US HazCom 2012

Revision date : 27 May 2022

Other information : None.

Safety Data Sheet (SDS), USA

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product.

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 12/10/2021

Reviewed on 12/10/2021

* 1 Identification

· **Product Identifier**

· **Trade Name:** *cDNA Primers, Feature cDNA Primers 1, Feature cDNA Primers 2, Feature cDNA Primers 3, Second Strand Primer, Dual Index Plate TS Set A, Human WT Probes - RHS, Human WT Probes - LHS, Mouse WT Probes - RHS, Mouse WT Probes - LHS Human WT Probes - RHS, Human WT Probes - LHS, Mouse WT Probes - RHS, Mouse WT Probes - LHS, TS Primer Mix A, Pre-Amp Primers B, Feature Pre-Amp Primers, Human WTA Probes BC001, Human WTA Probes BC002, Human WTA Probes BC003, Human WTA Probes BC004, Human WTA Probes BC005, Human WTA Probes BC006, Human WTA Probes BC007, Human WTA Probes BC008, Human WTA Probes BC009, Human WTA Probes BC010, Human WTA Probes BC011, Human WTA Probes BC012, Human WTA Probes BC013, Human WTA Probes BC014, Human WTA Probes BC015, Human WTA Probes BC016, Antibody Blocker, Mouse WTA Probes BC001, Mouse WTA Probes BC002, Mouse WTA Probes BC003, Mouse WTA Probes BC004, Mouse WTA Probes BC005, Mouse WTA Probes BC006, Mouse WTA Probes BC007, Mouse WTA Probes BC008, Mouse WTA Probes BC009, Mouse WTA Probes BC010, Mouse WTA Probes BC011, Mouse WTA Probes BC012, Mouse WTA Probes BC013, Mouse WTA Probes BC014, Mouse WTA Probes BC015, Mouse WTA Probes BC016, Human WT Probes v2, RHS, Human WT Probes v2, LHS, TS Primer Mix B, NT Primer Mix A*

· **Product Number:**

2000089, 2000096, 2000097, 2000289, 2000217,
3000511, 2000453, 2000454, 2000457, 2000458,
2000449, 2000450, 2000455, 2000456, 2000447,
2000529, 2000515, 2000495, 2000496, 2000497,
2000498, 2000499, 2000500, 2000501, 2000502,
2000503, 2000504, 2000505, 2000506, 2000507,
2000508, 2000509, 2000510, 2000569, 2000703,
2000704, 2000705, 2000706, 2000707, 2000708,
2000709, 2000710, 2000711, 2000712, 2000713,
2000714, 2000715, 2000716, 2000717, 2000718,
2000657, 2000658, 2000537, 2000540

· **Relevant identified uses of the substance or mixture and uses advised against:**

· **Product Description:** No further relevant information available.

· **Details of the Supplier of the Safety Data Sheet:**

· **Manufacturer/Supplier:**

10x Genomics
6230 Stoneridge Mall Road
Pleasanton, CA 94588 USA
1 925 401 7300
<https://www.10xgenomics.com/>

· **Emergency telephone number:** 1 925 401 7300

2 Hazard(s) Identification

· **Classification of the substance or mixture:**

The product does not need classification according to OSHA HazCom Standard 29 CFR paragraph (d) of §1910.1200(g) and GHS Rev 03.

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 12/10/2021

Reviewed on 12/10/2021

Trade Name: *cDNA Primers, Feature cDNA Primers 1, Feature cDNA Primers 2 ,Feature cDNA Primers 3, Second Strand Primer, Dual Index Plate TS Set A, Human WT Probes - RHS, Human WT Probes - LHS, Mouse WT Probes - RHS, Mouse WT Probes - LHS Human WT Probes - RHS, Human WT Probes - LHS, Mouse WT Probes - RHS, Mouse WT Probes - LHS, TS Primer Mix A, Pre-Amp Primers B, Feature Pre-Amp Primers, Human WTA Probes BC001, Human WTA Probes BC002, Human WTA Probes BC003, Human WTA Probes BC004, Human WTA Probes BC005, Human WTA Probes BC006, Human WTA Probes BC007, Human WTA Probes BC008, Human WTA Probes BC009, Human WTA Probes BC010, Human WTA Probes BC011, Human WTA Probes BC012, Human WTA Probes BC013, Human WTA Probes BC014, Human WTA Probes BC015, Human WTA Probes BC016, Antibody Blocker, Mouse WTA Probes BC001, Mouse WTA Probes BC002, Mouse WTA Probes BC003, Mouse WTA Probes BC004, Mouse WTA Probes BC005, Mouse WTA Probes BC006, Mouse WTA Probes BC007, Mouse WTA Probes BC008, Mouse WTA Probes BC009, Mouse WTA Probes BC010, Mouse WTA Probes BC011, Mouse WTA Probes BC012, Mouse WTA Probes BC013, Mouse WTA Probes BC014, Mouse WTA Probes BC015, Mouse WTA Probes BC016, Human WT Probes v2, RHS, Human WT Probes v2, LHS*

Label elements:

Hazard pictograms: Non-Regulated Material

Signal word: Non-Regulated Material

Hazard statements: Non-Regulated Material


Unknown acute toxicity:

This value refers to knowledge of known, established toxicological or ecotoxicological values.


0 % of the mixture consists of component(s) of unknown toxicity.

Classification system: NFPA/HMIS Definitions: 0-Least, 1-Slight, 2-Moderate, 3-High, 4-Extreme

NFPA ratings (scale 0 - 4)



HMIS-ratings (scale 0 - 4)



Hazard(s) not otherwise classified (HNOC): None known

3 Composition/Information on Ingredients

Chemical characterization: Substance

Description: Non-Regulated Material

Dangerous Components: Non-Regulated Material

4 First-Aid Measures

Description of first aid measures

General information: No special measures required.

After inhalation:

In case of unconsciousness place patient stably in the side position for transportation.

Supply fresh air; consult doctor in case of complaints.

After skin contact:

Wash with soap and water.

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 12/10/2021 Trade

Reviewed on 12/10/2021

Name: *cDNA Primers, Feature cDNA Primers 1, Feature cDNA Primers 2, Feature cDNA Primers 3, Second Strand Primer, Dual Index Plate TS Set A, Human WT Probes - RHS, Human WT Probes - LHS, Mouse WT Probes - RHS, Mouse WT Probes - LHS Human WT Probes - RHS, Human WT Probes - LHS, Mouse WT Probes - RHS, Mouse WT Probes - LHS, TS Primer Mix A, Pre-Amp Primers B, Feature Pre-Amp Primers, Human WTA Probes BC001, Human WTA Probes BC002, Human WTA Probes BC003, Human WTA Probes BC004, Human WTA Probes BC005, Human WTA Probes BC006, Human WTA Probes BC007, Human WTA Probes BC008, Human WTA Probes BC009, Human WTA Probes BC010, Human WTA Probes BC011, Human WTA Probes BC012, Human WTA Probes BC013, Human WTA Probes BC014, Human WTA Probes BC015, Human WTA Probes BC016, Antibody Blocker, Mouse WTA Probes BC001, Mouse WTA Probes BC002, Mouse WTA Probes BC003, Mouse WTA Probes BC004, Mouse WTA Probes BC005, Mouse WTA Probes BC006, Mouse WTA Probes BC007, Mouse WTA Probes BC008, Mouse WTA Probes BC009, Mouse WTA Probes BC010, Mouse WTA Probes BC011, Mouse WTA Probes BC012, Mouse WTA Probes BC013, Mouse WTA Probes BC014, Mouse WTA Probes BC015, Mouse WTA Probes BC016, Human WT Probes v2, RHS, Human WT Probes v2, LHS*

If skin irritation occurs, consult a doctor.

· **After eye contact:**

If eye irritation occurs, consult a doctor.

Rinse opened eye for several minutes under running water.

· **After swallowing:**

Do not induce vomiting without medical advice.

If swallowed and symptoms occur, consult a doctor.

· **Information for doctor**

· **Most important symptoms and effects, both acute and delayed:** No further relevant information available.

· **Indication of any immediate medical attention and special treatment needed:**

No further relevant information available.

5 Fire-Fighting Measures

· **Extinguishing media**

· **Suitable extinguishing agents:** Use fire fighting measures that suit the environment.

· **For safety reasons unsuitable extinguishing agents:** No further relevant information.

· **Special hazards arising from the substance or mixture:** No further relevant information available.

· **Advice for firefighters**

· **Special protective equipment for firefighters:**

As in any fire, wear self-contained breathing apparatus pressure-demand (NIOSH approved or equivalent) and full protective gear to prevent contact with skin and eyes.

6 Accidental Release Measures

· **Personal precautions, protective equipment and emergency procedures:**

Avoid contact with skin, eyes and clothing.

· **Environmental precautions:** Dilute with plenty of water.

· **Methods and material for containment and cleaning up:**

Absorb with liquid-binding material (i.e. sand, diatomite, acid binders, universal binders, sawdust).

Dispose of the collected material according to regulations.

· **Reference to other sections:**

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

(Contd. on page 4)

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 12/10/2021

Reviewed on 12/10/2021

Trade Name: *cDNA Primers, Feature cDNA Primers 1, Feature cDNA Primers 2, Feature cDNA Primers 3, Second Strand Primer, Dual Index Plate TS Set A, Human WT Probes - RHS, Human WT Probes - LHS, Mouse WT Probes - RHS, Mouse WT Probes - LHS Human WT Probes - RHS, Human WT Probes - LHS, Mouse WT Probes - RHS, Mouse WT Probes - LHS, TS Primer Mix A, Pre-Amp Primers B, Feature Pre-Amp Primers, Human WTA Probes BC001, Human WTA Probes BC002, Human WTA Probes BC003, Human WTA Probes BC004, Human WTA Probes BC005, Human WTA Probes BC006, Human WTA Probes BC007, Human WTA Probes BC008, Human WTA Probes BC009, Human WTA Probes BC010, Human WTA Probes BC011, Human WTA Probes BC012, Human WTA Probes BC013, Human WTA Probes BC014, Human WTA Probes BC015, Human WTA Probes BC016, Antibody Blocker, Mouse WTA Probes BC001, Mouse WTA Probes BC002, Mouse WTA Probes BC003, Mouse WTA Probes BC004, Mouse WTA Probes BC005, Mouse WTA Probes BC006, Mouse WTA Probes BC007, Mouse WTA Probes BC008, Mouse WTA Probes BC009, Mouse WTA Probes BC010, Mouse WTA Probes BC011, Mouse WTA Probes BC012, Mouse WTA Probes BC013, Mouse WTA Probes BC014, Mouse WTA Probes BC015, Mouse WTA Probes BC016, Human WT Probes v2, RHS, Human WT Probes v2, LHS*

· **PAC-1:**

None of the ingredients are listed.

· **PAC-2:**

None of the ingredients are listed.

· **PAC-3:**

None of the ingredients are listed.

7 Handling and Storage

· **Handling**

· **Precautions for safe handling:** Avoid breathing fume/gas/mist/vapors/spray.

· **Information about protection against explosions and fires:** No special measures required.

· **Conditions for safe storage, including any incompatibilities**

· **Storage**

· **Requirements to be met by storerooms and receptacles:** No special requirements.

· **Information about storage in one common storage facility:** Not required.

· **Further information about storage conditions:** None.

· **Specific end use(s):** No further relevant information available.

8 Exposure Controls/Personal Protection

· **Additional information about design of technical systems:** No further data; see section 7.

· **Control parameters:**

· **Components with occupational exposure limits:**

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

· **Additional information:** The lists that were valid during the creation of this SDS were used as basis.

· **Exposure controls:**

· **Personal protective equipment**

· **General protective and hygienic measures:**

The usual precautionary measures for handling chemicals should be followed.

· **Breathing equipment:** Not required.

· **Protection of hands:** Not required.

· **Material of gloves:** Not applicable.

(Contd. on page 5)

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 12/10/2021

Reviewed on 12/10/2021

Trade Name: *cDNA Primers, Feature cDNA Primers 1, Feature cDNA Primers 2, Feature cDNA Primers 3, Second Strand Primer, Dual Index Plate TS Set A, Human WT Probes - RHS, Human WT Probes - LHS, Mouse WT Probes - RHS, Mouse WT Probes - LHS Human WT Probes - RHS, Human WT Probes - LHS, Mouse WT Probes - RHS, Mouse WT Probes - LHS, TS Primer Mix A, Pre-Amp Primers B, Feature Pre-Amp Primers, Human WTA Probes BC001, Human WTA Probes BC002, Human WTA Probes BC003, Human WTA Probes BC004, Human WTA Probes BC005, Human WTA Probes BC006, Human WTA Probes BC007, Human WTA Probes BC008, Human WTA Probes BC009, Human WTA Probes BC010, Human WTA Probes BC011, Human WTA Probes BC012, Human WTA Probes BC013, Human WTA Probes BC014, Human WTA Probes BC015, Human WTA Probes BC016, Antibody Blocker, Mouse WTA Probes BC001, Mouse WTA Probes BC002, Mouse WTA Probes BC003, Mouse WTA Probes BC004, Mouse WTA Probes BC005, Mouse WTA Probes BC006, Mouse WTA Probes BC007, Mouse WTA Probes BC008, Mouse WTA Probes BC009, Mouse WTA Probes BC010, Mouse WTA Probes BC011, Mouse WTA Probes BC012, Mouse WTA Probes BC013, Mouse WTA Probes BC014, Mouse WTA Probes BC015, Mouse WTA Probes BC016, Human WT Probes v2, RHS, Human WT Probes v2, LHS*

- **Penetration time of glove material:** Not applicable.
- **Eye protection:**



Safety glasses

9 Physical and Chemical Properties

· **Information on basic physical and chemical properties**

· **General Information**

· **Appearance:**

· **Form:** Liquid

· **Color:** Clear, colorless

· **Odor:** Odorless

· **Odor threshold:** Not determined.

· **pH-value:** Not applicable.

· **Change in condition**

· **Melting point/Melting range:** Not determined.

· **Flash point:** None

· **Flammability (solid, gaseous):** Not applicable.

· **Ignition temperature:** Not applicable

· **Decomposition temperature:** Not determined.

· **Auto igniting:** Product is not self-igniting.

· **Danger of explosion:** Product does not present an explosion hazard.

· **Explosion limits:**

· **Lower:** Not determined.

· **Upper:** Not determined.

· **Vapor pressure:** Not determined.

· **Density @ 20 °C (68 °F):** \ g/cm³

· **Relative density:** Not determined.

(Contd. on page 6)

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 12/10/2021

Reviewed on 12/10/2021

Trade Name: *cDNA Primers, Feature cDNA Primers 1, Feature cDNA Primers 2, Feature cDNA Primers 3, Second Strand Primer, Dual Index Plate TS Set A, Human WT Probes - RHS, Human WT Probes - LHS, Mouse WT Probes - RHS, Mouse WT Probes - LHS Human WT Probes - RHS, Human WT Probes - LHS, Mouse WT Probes - RHS, Mouse WT Probes - LHS, TS Primer Mix A, Pre-Amp Primers B, Feature Pre-Amp Primers, Human WTA Probes BC001, Human WTA Probes BC002, Human WTA Probes BC003, Human WTA Probes BC004, Human WTA Probes BC005, Human WTA Probes BC006, Human WTA Probes BC007, Human WTA Probes BC008, Human WTA Probes BC009, Human WTA Probes BC010, Human WTA Probes BC011, Human WTA Probes BC012, Human WTA Probes BC013, Human WTA Probes BC014, Human WTA Probes BC015, Human WTA Probes BC016, Antibody Blocker, Mouse WTA Probes BC001, Mouse WTA Probes BC002, Mouse WTA Probes BC003, Mouse WTA Probes BC004, Mouse WTA Probes BC005, Mouse WTA Probes BC006, Mouse WTA Probes BC007, Mouse WTA Probes BC008, Mouse WTA Probes BC009, Mouse WTA Probes BC010, Mouse WTA Probes BC011, Mouse WTA Probes BC012, Mouse WTA Probes BC013, Mouse WTA Probes BC014, Mouse WTA Probes BC015, Mouse WTA Probes BC016, Human WT Probes v2, RHS, Human WT Probes v2, LHS*

- **Vapor density:** Not determined.
- **Evaporation rate:** Not determined.
- **Solubility in / Miscibility with:**
 - Water:** Fully miscible.
- **Partition coefficient (n-octanol/water):** Not determined.
- **Viscosity:**
 - Dynamic:** Not determined.
 - Kinematic:** Not determined.
- **Solvent content:**
 - VOC content:** 0.00 %
- **Other information:** No further relevant information available.

10 Stability and Reactivity

- **Reactivity:** No further relevant information available.
- **Chemical stability:** Stable under normal conditions.
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions:** No dangerous reactions known.
- **Conditions to avoid:** No further relevant information available.
- **Incompatible materials:** Strong acids, strong bases
- **Hazardous decomposition products:** No dangerous decomposition products known.

11 Toxicological Information

- **Information on toxicological effects:**
- **Acute toxicity:**
- **LD/LC50 values that are relevant for classification:** No data available.
- **Primary irritant effect:**
- **On the skin:** No irritating effect.
- **On the eye:** No irritating effect.
- **Additional toxicological information:**
The product is not subject to classification according to internally approved calculation methods for preparations.

(Contd. on page 7)

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 12/10/2021

Reviewed on 12/10/2021

Trade Name: *cDNA Primers, Feature cDNA Primers 1, Feature cDNA Primers 2, Feature cDNA Primers 3, Second Strand Primer, Dual Index Plate TS Set A, Human WT Probes - RHS, Human WT Probes - LHS, Mouse WT Probes - RHS, Mouse WT Probes - LHS Human WT Probes - RHS, Human WT Probes - LHS, Mouse WT Probes - RHS, Mouse WT Probes - LHS, TS Primer Mix A, Pre-Amp Primers B, Feature Pre-Amp Primers, Human WTA Probes BC001, Human WTA Probes BC002, Human WTA Probes BC003, Human WTA Probes BC004, Human WTA Probes BC005, Human WTA Probes BC006, Human WTA Probes BC007, Human WTA Probes BC008, Human WTA Probes BC009, Human WTA Probes BC010, Human WTA Probes BC011, Human WTA Probes BC012, Human WTA Probes BC013, Human WTA Probes BC014, Human WTA Probes BC015, Human WTA Probes BC016, Antibody Blocker, Mouse WTA Probes BC001, Mouse WTA Probes BC002, Mouse WTA Probes BC003, Mouse WTA Probes BC004, Mouse WTA Probes BC005, Mouse WTA Probes BC006, Mouse WTA Probes BC007, Mouse WTA Probes BC008, Mouse WTA Probes BC009, Mouse WTA Probes BC010, Mouse WTA Probes BC011, Mouse WTA Probes BC012, Mouse WTA Probes BC013, Mouse WTA Probes BC014, Mouse WTA Probes BC015, Mouse WTA Probes BC016, Human WT Probes v2, RHS, Human WT Probes v2, LHS*

When used and handled according to specifications, the product does not have any harmful effects according to our experience and the information provided to us.

· **Carcinogenic categories:**

· IARC (International Agency for Research on Cancer):
None of the ingredients are listed.
· NTP (National Toxicology Program):
None of the ingredients are listed.
· OSHA-Ca (Occupational Safety & Health Administration):
None of the ingredients are listed.

12 Ecological Information

- **Toxicity:**
- **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability:** No further relevant information available.
- **Behavior in environmental systems:**
- **Bioaccumulative potential:** No further relevant information available.
- **Mobility in soil:** No further relevant information available.
- **Additional ecological information:**
- **General notes:**
Do not allow undiluted product or product that has not been neutralized to reach ground water, water course or sewage system.
- **Results of PBT and vPvB assessment:**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **Other adverse effects:** No further relevant information available.

13 Disposal Considerations

- **Waste treatment methods**
- **Recommendation:**
Smaller quantities can be disposed of with household waste.
Observe all federal, state and local environmental regulations when disposing of this material.
- **Uncleaned packaging**
- **Recommendation:** Disposal must be made according to official regulations.

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 12/10/2021

Reviewed on 12/10/2021

Trade Name: *cDNA Primers, Feature cDNA Primers 1, Feature cDNA Primers 2, Feature cDNA Primers 3, Second Strand Primer, Dual Index Plate TS Set A, Human WT Probes - RHS, Human WT Probes - LHS, Mouse WT Probes - RHS, Mouse WT Probes - LHS Human WT Probes - RHS, Human WT Probes - LHS, Mouse WT Probes - RHS, Mouse WT Probes - LHS, TS Primer Mix A, Pre-Amp Primers B, Feature Pre-Amp Primers, Human WTA Probes BC001, Human WTA Probes BC002, Human WTA Probes BC003, Human WTA Probes BC004, Human WTA Probes BC005, Human WTA Probes BC006, Human WTA Probes BC007, Human WTA Probes BC008, Human WTA Probes BC009, Human WTA Probes BC010, Human WTA Probes BC011, Human WTA Probes BC012, Human WTA Probes BC013, Human WTA Probes BC014, Human WTA Probes BC015, Human WTA Probes BC016, Antibody Blocker, Mouse WTA Probes BC001, Mouse WTA Probes BC002, Mouse WTA Probes BC003, Mouse WTA Probes BC004, Mouse WTA Probes BC005, Mouse WTA Probes BC006, Mouse WTA Probes BC007, Mouse WTA Probes BC008, Mouse WTA Probes BC009, Mouse WTA Probes BC010, Mouse WTA Probes BC011, Mouse WTA Probes BC012, Mouse WTA Probes BC013, Mouse WTA Probes BC014, Mouse WTA Probes BC015, Mouse WTA Probes BC016, Human WT Probes v2, RHS, Human WT Probes v2, LHS*

- **Recommended cleansing agent:** Water, if necessary with cleansing agents.

14 Transport Information

- **UN-Number:**
- **DOT, ADR/ADN, ADN, IMDG, IATA** Non-Regulated Material
- **UN proper shipping name:**
- **DOT, ADR/ADN, ADN, IMDG, IATA** Non-Regulated Material
- **Transport hazard class(es):**
- **DOT, ADR/ADN, ADN, IMDG, IATA**
- **Class:** Non-Regulated Material
- **Packing group:**
- **DOT, ADR/ADN, IMDG, IATA** Non-Regulated Material
- **Environmental hazards:** Not applicable.
- **Special precautions for user:** Not applicable.
- **Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code:** Not applicable.
- **UN "Model Regulation":** Non-Regulated Material

15 Regulatory Information

- **Safety, health and environmental regulations/legislation specific for the substance or mixture:**
No further relevant information available.
- **SARA (Superfund Amendments and Reauthorization):**

· **Section 355 (extremely hazardous substances):**
None of the ingredients are listed.

· **Section 313 (Specific toxic chemical listings):**
None of the ingredients are listed.

- **TSCA (Toxic Substances Control Act):**
- **California Proposition 65:**

· **Chemicals known to cause cancer:**
None of the ingredients are listed.

· **Chemicals known to cause reproductive toxicity for females:**
None of the ingredients are listed.

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 12/10/2021

Reviewed on 12/10/2021

Trade Name: *cDNA Primers, Feature cDNA Primers 1, Feature cDNA Primers 2, Feature cDNA Primers 3, Second Strand Primer, Dual Index Plate TS Set A, Human WT Probes - RHS, Human WT Probes - LHS, Mouse WT Probes - RHS, Mouse WT Probes - LHS Human WT Probes - RHS, Human WT Probes - LHS, Mouse WT Probes - RHS, Mouse WT Probes - LHS, TS Primer Mix A, Pre-Amp Primers B, Feature Pre-Amp Primers, Human WTA Probes BC001, Human WTA Probes BC002, Human WTA Probes BC003, Human WTA Probes BC004, Human WTA Probes BC005, Human WTA Probes BC006, Human WTA Probes BC007, Human WTA Probes BC008, Human WTA Probes BC009, Human WTA Probes BC010, Human WTA Probes BC011, Human WTA Probes BC012, Human WTA Probes BC013, Human WTA Probes BC014, Human WTA Probes BC015, Human WTA Probes BC016, Antibody Blocker, Mouse WTA Probes BC001, Mouse WTA Probes BC002, Mouse WTA Probes BC003, Mouse WTA Probes BC004, Mouse WTA Probes BC005, Mouse WTA Probes BC006, Mouse WTA Probes BC007, Mouse WTA Probes BC008, Mouse WTA Probes BC009, Mouse WTA Probes BC010, Mouse WTA Probes BC011, Mouse WTA Probes BC012, Mouse WTA Probes BC013, Mouse WTA Probes BC014, Mouse WTA Probes BC015, Mouse WTA Probes BC016, Human WT Probes v2, RHS, Human WT Probes v2, LHS*

<ul style="list-style-type: none"> · Chemicals known to cause reproductive toxicity for males:
<p>None of the ingredients are listed.</p>
<ul style="list-style-type: none"> · Chemicals known to cause developmental toxicity:
<p>None of the ingredients are listed.</p>
<ul style="list-style-type: none"> · New Jersey Right-to-Know List:
<p>None of the ingredients are listed.</p>
<ul style="list-style-type: none"> · New Jersey Special Hazardous Substance List:
<p>None of the ingredients are listed.</p>
<ul style="list-style-type: none"> · Pennsylvania Right-to-Know List:
<p>None of the ingredients are listed.</p>
<ul style="list-style-type: none"> · Pennsylvania Special Hazardous Substance List:
<p>None of the ingredients are listed.</p>
<ul style="list-style-type: none"> · Carcinogenic categories:
<ul style="list-style-type: none"> · EPA (Environmental Protection Agency):
<p>None of the ingredients are listed.</p>
<ul style="list-style-type: none"> · TLV (Threshold Limit Value established by ACGIH):
<p>None of the ingredients are listed.</p>
<ul style="list-style-type: none"> · NIOSH-Ca (National Institute for Occupational Safety and Health):
<p>None of the ingredients are listed.</p>
<ul style="list-style-type: none"> · GHS label elements Non-Regulated Material · Hazard pictograms: Non-Regulated Material · Signal word: Non-Regulated Material · Hazard statements: Non-Regulated Material
<ul style="list-style-type: none"> · National regulations:
<p>The product is not subject to be labelled according with the prevailing version of the regulations on hazardous substances.</p>
<ul style="list-style-type: none"> · Chemical safety assessment: A Chemical Safety Assessment has not been carried out.
<div style="border: 1px solid black; height: 20px; width: 100%;"></div>

Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 12/10/2021

Reviewed on 12/10/2021

Trade Name: *cDNA Primers, Feature cDNA Primers 1, Feature cDNA Primers 2, Feature cDNA Primers 3, Second Strand Primer, Dual Index Plate TS Set A, Human WT Probes - RHS, Human WT Probes - LHS, Mouse WT Probes - RHS, Mouse WT Probes - LHS Human WT Probes - RHS, Human WT Probes - LHS, Mouse WT Probes - RHS, Mouse WT Probes - LHS, TS Primer Mix A, Pre-Amp Primers B, Feature Pre-Amp Primers, Human WTA Probes BC001, Human WTA Probes BC002, Human WTA Probes BC003, Human WTA Probes BC004, Human WTA Probes BC005, Human WTA Probes BC006, Human WTA Probes BC007, Human WTA Probes BC008, Human WTA Probes BC009, Human WTA Probes BC010, Human WTA Probes BC011, Human WTA Probes BC012, Human WTA Probes BC013, Human WTA Probes BC014, Human WTA Probes BC015, Human WTA Probes BC016, Antibody Blocker, Mouse WTA Probes BC001, Mouse WTA Probes BC002, Mouse WTA Probes BC003, Mouse WTA Probes BC004, Mouse WTA Probes BC005, Mouse WTA Probes BC006, Mouse WTA Probes BC007, Mouse WTA Probes BC008, Mouse WTA Probes BC009, Mouse WTA Probes BC010, Mouse WTA Probes BC011, Mouse WTA Probes BC012, Mouse WTA Probes BC013, Mouse WTA Probes BC014, Mouse WTA Probes BC015, Mouse WTA Probes BC016, Human WT Probes v2, RHS, Human WT Probes v2, LHS*

16 Other Information

The information and recommendations in this safety data sheet are, to the best of our knowledge, accurate as of the date of issue. Nothing herein shall be deemed to create warranty, expressed or implied, and shall not establish a legally valid contractual relationship. It is the responsibility of the user to determine applicability of this information and the suitability of the material or product for any particular purpose.

· **Contact:**

· **Date of last revision/ revision number:** 12/10/2021 / 2

· **Abbreviations and acronyms:**

ADR: The European Agreement concerning the International Carriage of Dangerous Goods by Road

ADN: The European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

VOC: Volatile Organic Compounds (USA, EU)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety and Health

OSHA: Occupational Safety & Health Administration

TLV: Threshold Limit Value

PEL: Permissible Exposure Limit

REL: Recommended Exposure Limit

· *** Data compared to the previous version altered.**

SDS created by MSDS Authoring Services www.msdsauthoring.com +1-877-204-9106

SECTION 1: Identification

1.1. Identification

Product form : Mixtures
Trade name : Block and Stain Buffer
Product code : 2000554

1.2. Recommended use and restrictions on use

Use of the substance/mixture : Reagent

1.3. Supplier

Manufacturer/Supplier:
10x Genomics
6230 Stoneridge Mall Road Pleasanton,
CA 94588-3260
T: +1 925 401 7300
E: info@10xgenomics.com

1.4. Emergency telephone number

Emergency number : +1 925 401 7300

SECTION 2: Hazard(s) identification

2.1. Classification of the substance or mixture

GHS US classification

Not classified

2.2. GHS Label elements, including precautionary statements

GHS US labeling

No labeling applicable

2.3. Other hazards which do not result in classification

No additional information available

2.4. Unknown acute toxicity (GHS US)

Not applicable

SECTION 3: Composition/Information on ingredients

3.1. Substances

Not applicable

3.2. Mixtures

This mixture does not contain any substances to be mentioned according to the criteria of section 3.2 of HazCom 2012

Block and Stain Buffer

Safety Data Sheet

according to US HazCom 2012

SECTION 4: First-aid measures

4.1. Description of first aid measures

First-aid measures general	: If you feel unwell, seek medical advice (show the label where possible).
First-aid measures after inhalation	: Remove person to fresh air and keep comfortable for breathing. If experiencing respiratory symptoms: Call a poison center or a doctor.
First-aid measures after skin contact	: Wash with water and soap as a precaution. If skin irritation occurs: Get medical advice/attention.
First-aid measures after eye contact	: Rinse eyes with water as a precaution. Seek medical attention if ill effect or irritation develops.
First-aid measures after ingestion	: Rinse mouth out with water. Do not induce vomiting. If you feel unwell, seek medical advice.

4.2. Most important symptoms and effects (acute and delayed)

Symptoms/effects	: No adverse effects expected.
------------------	--------------------------------

4.3. Immediate medical attention and special treatment, if necessary

Treat symptomatically.

SECTION 5: Fire-fighting measures

5.1. Suitable (and unsuitable) extinguishing media

Suitable extinguishing media	: Use extinguishing media appropriate for surrounding fire.
Unsuitable extinguishing media	: None known.

5.2. Specific hazards arising from the chemical

Fire hazard	: The product is not flammable. Does not sustain combustion. In the event of fire, may decompose : Carbon oxides (CO, CO ₂). Sulphur oxides. Phosphorus oxides.
Explosion hazard	: No hazard identified.
Hazardous decomposition products in case of fire	: Thermal decomposition can lead to the release of irritating gases and vapors.

5.3. Special protective equipment and precautions for fire-fighters

Firefighting instructions	: Fight fire with normal precautions from a reasonable distance.
Protection during firefighting	: Do not attempt to take action without suitable protective equipment.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

General measures	: No special requirements.
------------------	----------------------------

6.1.1. For non-emergency personnel

Protective equipment	: No special protection required.
Emergency procedures	: No additional risk management measures required.

6.1.2. For emergency responders

Protective equipment	: No special protection required.
Emergency procedures	: No additional risk management measures required.

6.2. Environmental precautions

No special environmental precautions required.

6.3. Methods and material for containment and cleaning up

For containment	: Absorb spillage to prevent material-damage.
Methods for cleaning up	: Wipe up with absorbent material (for example cloth).

Block and Stain Buffer

Safety Data Sheet

according to US HazCom 2012

Other information : Dispose in a safe manner in accordance with local/national regulations.

6.4. Reference to other sections

For further information refer to section 8: "Exposure controls/personal protection". For disposal of residues refer to section 13 : "Disposal considerations".

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : No special handling advices are necessary.
Hygiene measures : Use good personal hygiene practices.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store in a well-ventilated place. Keep cool.
Incompatible materials : Acids. Oxidizing materials.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Block and Stain Buffer

No additional information available

8.2. Appropriate engineering controls

Appropriate engineering controls : No special requirements.

8.3. Individual protection measures/Personal protective equipment

Hand protection:

In case of repeated or prolonged exposure : Wear suitable gloves resistant to chemical penetration. Please follow the instructions related to the permeability and the penetration time provided by the manufacturer

Eye protection:

If there is a risk of eye contact: Safety glasses with side shields

Other information:

Do not eat, drink or smoke when using this product.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state : Liquid
Appearance : Viscous.
Color : Yellow
Odor : odorless
Odor threshold : No data available
pH : No data available
Melting point : No data available
Freezing point : No data available
Boiling point : No data available
Flash point : No data available

Block and Stain Buffer

Safety Data Sheet

according to US HazCom 2012

Relative evaporation rate (butyl acetate=1)	: No data available
Flammability (solid, gas)	: Not applicable.
Vapor pressure	: No data available
Relative vapor density at 20 °C	: No data available
Relative density	: No data available
Solubility	: No data available
Partition coefficient n-octanol/water (Log Pow)	: No data available
Auto-ignition temperature	: No data available
Decomposition temperature	: No data available
Viscosity, kinematic	: No data available
Viscosity, dynamic	: No data available
Explosion limits	: No data available
Explosive properties	: No data available
Oxidizing properties	: No data available

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

The product is non-reactive under normal conditions of use, storage and transport.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

Hazardous polymerization will not occur.

10.4. Conditions to avoid

None under recommended storage and handling conditions (see section 7).

10.5. Incompatible materials

Acids. Oxidizing materials.

10.6. Hazardous decomposition products

Under normal conditions of storage and use, hazardous decomposition products should not be produced. Thermal decomposition may produce : Carbon oxides (CO, CO₂). Sulphur oxides. Phosphorus oxides.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity (oral)	: Not classified (Based on available data, the classification criteria are not met)
Acute toxicity (dermal)	: Not classified (Based on available data, the classification criteria are not met)
Acute toxicity (inhalation)	: Not classified (Based on available data, the classification criteria are not met)
Skin corrosion/irritation	: Not classified (Based on available data, the classification criteria are not met)
Serious eye damage/irritation	: Not classified (Based on available data, the classification criteria are not met)
Respiratory or skin sensitization	: Not classified (Based on available data, the classification criteria are not met)
Germ cell mutagenicity	: Not classified (Based on available data, the classification criteria are not met)
Carcinogenicity	: Not classified (Based on available data, the classification criteria are not met)
Reproductive toxicity	: Not classified (Based on available data, the classification criteria are not met)

Block and Stain Buffer

Safety Data Sheet

according to US HazCom 2012

STOT-single exposure	: Not classified (Based on available data, the classification criteria are not met)
STOT-repeated exposure	: Not classified (Based on available data, the classification criteria are not met)
Aspiration hazard	: Not classified (Based on available data, the classification criteria are not met)
Viscosity, kinematic	: No data available
Symptoms/effects	: No adverse effects expected.
Other information	: Likely routes of exposure: ingestion, inhalation, skin and eye.

SECTION 12: Ecological information

12.1. Toxicity

Ecology - general : This material has not been tested for environmental effects.

12.2. Persistence and degradability

No additional information available

12.3. Bioaccumulative potential

No additional information available

12.4. Mobility in soil

No additional information available

12.5. Other adverse effects

No additional information available

SECTION 13: Disposal considerations

13.1. Disposal methods

Product/Packaging disposal recommendations : Dispose of in accordance with relevant local regulations.

SECTION 14: Transport information

In accordance with DOT / TDG / IMDG / IATA

DOT	TDG	IMDG	IATA
14.1. UN number			
Not regulated for transport			
14.2. Proper Shipping Name			
Not applicable	Not applicable	Not applicable	Not applicable
14.3. Transport hazard class(es)			
Not applicable	Not applicable	Not applicable	Not applicable
14.4. Packing group			
Not applicable	Not applicable	Not applicable	Not applicable
14.5. Environmental hazards			
Not applicable	Not applicable	Not applicable	Not applicable
No supplementary information available			

Block and Stain Buffer

Safety Data Sheet

according to US HazCom 2012

SECTION 15: Regulatory information

15.1. US Federal regulations

No additional information available

15.2. International regulations

CANADA

No additional information available

EU-Regulations

No additional information available

National regulations

No additional information available

15.3. US State regulations

No additional information available

SECTION 16: Other information

according to US HazCom 2012

Revision date : 27 May 2022

Other information : None.

Safety Data Sheet (SDS), USA

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product.