# SAFETY DATA SHEETS

# Chromium Next GEM Automated Single Cell 5' Reagent Kits v2

Automated Feature Barcode Library Construction Kit, Module 2, 24 rxns PN-1000453

REAGENTS	PN
Feature Barcode Library Construction Module 2 Reagent 1	2000589





OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 10/26/2021 Reviewed on 01/01/2021

## 1 Identification

- · Product Identifier
- · Trade Name: Feature Barcode Library Construction Module 2 Reagent 1
- · Product Number: 2000589
- · Relevant identified uses of the substance or mixture and uses advised against:
- · Product Description: No Further Information
- · Application of the substance / the mixture: Industry Specific Application
- · Details of the Supplier of the Safety Data Sheet:
- Manufacturer/Supplier:

10x Genomics

6230 Stoneridge Mall Road Pleasanton, CA 94588 USA

1 925 401 7300

https://www.10xgenomics.com/

Emergency telephone number: 1 925 401 7300

## 2 Hazard(s) Identification

· Classification of the substance or mixture:



Skin Irrit. 2 H315 Causes skin irritation.

Eye Irrit. 2B H320 Causes eye irritation.

- · Label elements:
- · Hazard pictograms:



- · Signal word: Warning
- · Hazard statements:

H315+H320 Causes skin and eye irritation.

· Precautionary statements:

P264 Wash thoroughly after handling.

P280 Wear protective gloves.

P302+P352 If on skin: Wash with plenty of water.

P321 Specific treatment (see supplementary first aid instructions on this Safety Data Sheet).

P305+P351+P338 If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if

present and easy to do. Continue rinsing.

P332+P313 If skin irritation occurs: Get medical advice/attention.
P362+P364 Take off contaminated clothing and wash it before reuse.
P337+P313 If eye irritation persists: Get medical advice/attention.

· Unknown acute toxicity:

This value refers to knowledge of known, established toxicological or ecotoxicological values.

0 % of the mixture consists of component(s) of unknown toxicity.

· Classification system: NFPA/HMIS Definitions: 0-Least, 1-Slight, 2-Moderate, 3-High, 4-Extreme

(Contd. on page 2)



OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 10/26/2021 Reviewed on 01/01/2021

## Trade Name: Feature Barcode Library Construction Module 2 Reagent 1

· NFPA ratings (scale 0 - 4)



· HMIS-ratings (scale 0 - 4)



· Hazard(s) not otherwise classified (HNOC): None known

## 3 Composition/Information on Ingredients

- · Chemical characterization: Substance
- · Description: Mixture of substances listed below with non-hazardous additions.

· Dangerous Components:		
CAS: 56-81-5 RTECS: MA 8050000	Glycerine	0-12%
	Trade Secret	2-12%
	♦ Skin Irrit. 2, H315; Eye Irrit. 2A, H319	
	Trade Secret	≤2.5%
	♦ Skin Corr. 1A, H314	
	Specific concentration limits: Skin Corr. 1A; H314: C ≥ 15 %	
	Skin Irrit. 2; H315: 5 % ≤ C < 15 %	
	Eye Irrit. 2; H319: 5 % ≤ C < 15 %	

## Additional information:

The exact percentages of the ingredients of this mixture are considered to be proprietary and are withheld in accordance with the provisions of paragraph (i) of §1910.1200 of 29 CFR 1910.1200 Trade Secrets.

## 4 First-Aid Measures

- · Description of first aid measures
- · After inhalation: In case of unconsciousness place patient stably in the side position for transportation.
- · After skin contact:

Wash with soap and water.

If skin irritation occurs, consult a doctor.

· After eye contact:

Rinse opened eye for several minutes under running water.

If easy to do so, remove contact lenses if worn.

If eye irritation occurs, consult a doctor.

- · After swallowing: Do not induce vomiting without medical advice.
- · Information for doctor
- · Most important symptoms and effects, both acute and delayed: No further relevant information available.
- Indication of any immediate medical attention and special treatment needed:

No further relevant information available.

## 5 Fire-Fighting Measures

- · Extinguishing media
- · Suitable extinguishing agents: Use fire fighting measures that suit the environment.
- For safety reasons unsuitable extinguishing agents: No further relevant information.

(Contd. on page 3)



OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 10/26/2021 Reviewed on 01/01/2021

## Trade Name: Feature Barcode Library Construction Module 2 Reagent 1

- · Special hazards arising from the substance or mixture: No further relevant information available.
- Advice for firefighters
- Special protective equipment for firefighters:

As in any fire, wear self-contained breathing apparatus pressure-demand (NIOSH approved or equivalent) and full protective gear to prevent contact with skin and eyes.

## 6 Accidental Release Measures

## · Personal precautions, protective equipment and emergency procedures:

Do not inhale gases / fumes /aerosols

Avoid contact with skin, eyes and clothing.

- Environmental precautions: Dilute with plenty of water.
- Methods and material for containment and cleaning up:

Absorb with liquid-binding material (i.e. sand, diatomite, acid binders, universal binders, sawdust).

Ensure adequate ventilation.

Dispose of the collected material according to regulations.

· Reference to other sections:

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

PAC-1:		
56-81-5 Gly	cerine	45 mg/m³
Tra	de Secret	43 mg/m³
Tra	nde Secret	0.20 mg/m³
· PAC-2:		
56-81-5 Gly	cerine	180 mg/m³
Tra	nde Secret	470 mg/m³
Tra	de Secret	8.7 mg/m³
PAC-3:		
56-81-5 Gly	cerine	1,100 mg/m³
Tra	nde Secret	5,400 mg/m³
Tra	nde Secret	160 mg/m³

## 7 Handling and Storage

- · Handling
- · Precautions for safe handling:

Avoid release into the environment.

Ensure good ventilation/exhaustion at the workplace.

Prevent formation of aerosols.

- · Information about protection against explosions and fires: No special measures required.
- · Conditions for safe storage, including any incompatibilities
- Storage
- · Requirements to be met by storerooms and receptacles: Store at room temperature.
- Information about storage in one common storage facility: Not required.
- · Further information about storage conditions: Keep receptacle tightly sealed.
- · Specific end use(s): No further relevant information available.

(Contd. on page 4)



OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 10/26/2021 Reviewed on 01/01/2021

Trade Name: Feature Barcode Library Construction Module 2 Reagent 1

## 8 Exposure Controls/Personal Protection

- Additional information about design of technical systems: No further data; see section 7.
- · Control parameters:
- · Components with occupational exposure limits:

The following constituents are the only constituents of the product which have a PEL, TLV or other recommended exposure limit.

At this time, the remaining constituents have no known exposure limits.

	e anno, and romaniming contentations make the factor of postal community
<b>56-8</b>	1-5 Glycerine
PEL	Long-term value: 15* 5** mg/m³ mist; *total dust **respirable fraction
TLV	TLV withdrawn-insufficient data human occup. exp.
Trad	e Secret
PEL	Long-term value: 1 mg/m³
REL	Long-term value: 1 mg/m³
TLV	Long-term value: 0.2* mg/m³ *as thoracic fraction

- · Additional information: The lists that were valid during the creation of this SDS were used as basis.
- · Exposure controls:
- · Personal protective equipment
- · General protective and hygienic measures:

The usual precautionary measures for handling chemicals should be followed.

Keep away from foodstuffs, beverages and feed.

Immediately remove all soiled and contaminated clothing and wash before reuse.

Wash hands before breaks and at the end of work.

Avoid contact with the eyes and skin.

- · Breathing equipment: Not required.
- Protection of hands:



Protective gloves

## · Material of gloves:

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material cannot be calculated in advance and has therefore to be checked prior to the application.

#### · Penetration time of glove material:

The exact break-through time has to be determined and observed by the manufacturer of the protective gloves.

· Eye protection:



Tightly sealed goggles

· Limitation and supervision of exposure into the environment: None

(Contd. on page 5)



OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 10/26/2021 Reviewed on 01/01/2021

Trade Name: Feature Barcode Library Construction Module 2 Reagent 1

## 9 Physical and Chemical Properties

· Information on basic physical and chemical properties

· General Information

· Appearance:

Form: Liquid

Color:Clear, colorlessOdor:OdorlessOdor threshold:Not determined.pH-value:Not applicable.

· Change in condition

**Melting point/Melting range:** Not determined.

· Flash point: None

Flammability (solid, gaseous): Not applicable.
 Ignition temperature: ≥400 °C (≥752 °F)
 Decomposition temperature: Not determined.

· **Auto igniting:** Product is not self-igniting.

• Danger of explosion: Product does not present an explosion hazard.

· Explosion limits:

Lower: Not determined.Vapor pressure: Not determined.

· Density:

Relative density:Not determined.Vapor density:Not determined.Evaporation rate:Not determined.

· Solubility in / Miscibility with:

Water: Fully miscible.Partition coefficient (n-octanol/water): Not determined.

· Viscosity:

**Dynamic:** Not determined. **Kinematic:** Not determined.

· Other information: No further relevant information available.

## 10 Stability and Reactivity

- · Reactivity: No further relevant information available.
- · Chemical stability: Stable under normal conditions.
- · Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
- · Possibility of hazardous reactions: No dangerous reactions known.
- · Conditions to avoid: Contact with moisture, acids and bases.
- · Incompatible materials:

Strong oxidizing agents, strong reducing agents

Strong acids, strong bases

(Contd. on page 6)



OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 10/26/2021 Reviewed on 01/01/2021

## Trade Name: Feature Barcode Library Construction Module 2 Reagent 1

· Hazardous decomposition products: No dangerous decomposition products known.

## \*11 Toxicological Information

- · Information on toxicological effects:
- · Acute toxicity:

· LD/LC50 values that are relevant for classification:			
56-81-5 G	56-81-5 Glycerine		
Oral	LD50	12600 mg/kg (Rat)	
Dermal	LD50	>10000 mg/kg (Rabbit)	
Inhalative	LC50/96 hours	>50000 mg/l (Trout)	
Trade Sec	Trade Secret		
Oral	LD50	3310 mg/kg (Rat)	
Dermal	LD50	>2000 mg/kg (Rat)	
Trade Sec	Trade Secret		
Oral	LD50 Oral	2140 ml/kg (Rat)	
Inhalative	LC50/96 hours		
		Sulfuric acid is harmful by all routes of entry.	

- · Primary irritant effect:
- · On the skin: Irritant to skin and mucous membranes.
- · On the eye: Irritating effect.
- · Additional toxicological information:

The product shows the following dangers according to internally approved calculation methods for preparations:

Irritant

· Carcinogenic categories:

· IARC (International Agency for Research on Cancer):		
Trade Secret	7664-93-9	
· NTP (National Toxicology Program):		
None of the ingredients are listed.		
· OSHA-Ca (Occupational Safety & Health Administration):		
None of the ingredients are listed.		

## 12 Ecological Information

· Toxicity:

· Aquatic toxicity:		
Trade Secret		
EC50	12 mg/l (Algae - Selenastrum capricornutum) (96 HOURS)	
	0.148 mg/l (Daphnia) (48 hours)	
Trade Secret		
EC50 22 mg/l (Daphnia)		

- · Persistence and degradability: No further relevant information available.
- · Behavior in environmental systems:
- · Bioaccumulative potential: No further relevant information available.
- · Mobility in soil: No further relevant information available.

(Contd. on page 7)



OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 10/26/2021 Reviewed on 01/01/2021

## Trade Name: Feature Barcode Library Construction Module 2 Reagent 1

- · Additional ecological information:
- · General notes:

Do not allow undiluted product or product that has not been neutralized to reach ground water, water course or sewage system.

- · Results of PBT and vPvB assessment:
- · PBT: Not applicable.
- · vPvB: Not applicable.
- · Other adverse effects: No further relevant information available.

## 13 Disposal Considerations

- · Waste treatment methods
- · Recommendation:

Observe all federal, state and local environmental regulations when disposing of this material.

- · Uncleaned packaging
- Recommendation: Disposal must be made according to official regulations.
- Recommended cleansing agent: Water, if necessary with cleansing agents.

## 14 Transport Information

· UN-Number:

· DOT, ADR/ADN, ADN, IMDG, IATA Non-Regulated Material

· UN proper shipping name:

· DOT, ADR/ADN, ADN, IMDG, IATA Non-Regulated Material

· Transport hazard class(es):

· DOT, ADR/ADN, ADN, IMDG, IATA

· Class: Non-Regulated Material

· Packing group:

· DOT, ADR/ADN, IMDG, IATA Non-Regulated Material

Environmental hazards: Not applicable.Special precautions for user: Not applicable.

Transport in bulk according to Annex II of

MARPOL73/78 and the IBC Code: Not applicable.

UN "Model Regulation": Non-Regulated Material

## 15 Regulatory Information

- · Safety, health and environmental regulations/legislation specific for the substance or mixture: No further relevant information available.
- · SARA (Superfund Amendments and Reauthorization):

· Section 355	(extremel	v hazardous	substances	):
---------------	-----------	-------------	------------	----

**Trade Secret** 

· Section 313 (Specific toxic chemical listings):

Trade Secret

Trade Secret

· TSCA (Toxic Substances Control Act):

All components have the value ACTIVE.

· Hazardous Air Pollutants

None of the ingredients are listed.

(Contd. on page 8)



OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 10/26/2021 Reviewed on 01/01/2021

## Trade Name: Feature Barcode Library Construction Module 2 Reagent 1

· California Proposition 65:

· Chemicals known to cause cancer:	
None of the ingredients are listed.	
· Chemicals known to cause reproductive toxicity for females:	
None of the ingredients are listed.	
· Chemicals known to cause reproductive toxicity for males:	
None of the ingredients are listed.	
· Chemicals known to cause developmental toxicity:	
None of the ingredients are listed.	
· New Jersey Right-to-Know List:	
56-81-5 Glycerine	
Trade Secret	
· New Jersey Special Hazardous Substance List:	
Trade Secret	CA, CO, R
· Pennsylvania Right-to-Know List:	
56-81-5 Glycerine	
Trade Secret	
Pennsylvania Special Hazardous Substance List:	
Trade Secret	E

· Carcinogenic categories:	
· EPA (Environmental Protection Agency):	
None of the ingredients are listed.	
· TLV (Threshold Limit Value established by ACGIH):	
Trade Secret	A2
· NIOSH-Ca (National Institute for Occupational Safety and Health):	
None of the ingredients are listed.	

· GHS label elements

The product is classified and labeled according to the Globally Harmonized System (GHS).

· Hazard pictograms:



· Signal word: Warning

· Hazard statements:

H315+H320 Causes skin and eye irritation.

· Precautionary statements:

P264 Wash thoroughly after handling.

P280 Wear protective gloves.

P302+P352 If on skin: Wash with plenty of water.

P321 Specific treatment (see supplementary first aid instructions on this Safety Data Sheet). P305+P351+P338 If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if

present and easy to do. Continue rinsing.

P332+P313 If skin irritation occurs: Get medical advice/attention.

P362+P364 Take off contaminated clothing and wash it before reuse.



OSHA HazCom Standard 29 CFR 1910.1200(g) revised in 2012 and GHS Rev 03.

Issue date 10/26/2021 Reviewed on 01/01/2021

## Trade Name: Feature Barcode Library Construction Module 2 Reagent 1

P337+P313 If eye irritation persists: Get medical advice/attention.

## · National regulations:

The product is not subject to be labelled according with the prevailing version of the regulations on hazardous substances.

· Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

## 16 Other Information

The information and recommendations in this safety data sheet are, to the best of our knowledge, accurate as of the date of issue. Nothing herein shall be deemed to create warranty, expressed or implied, and shall not establish a legally valid contractual relationship. It is the responsibility of the user to determine applicability of this information and the suitability of the material or product for any particular purpose.

- · Contact:
- · Date of last revision/ revision number: 10/26/2021 / 3
- · Abbreviations and acronyms:

ADR: The European Agreement concerning the International Carriage of Dangerous Goods by Road

ADN: The European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

vPvB. very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety and Health

OSHA: Occupational Safety & Health Administration

TLV: Threshold Limit Value

PEL: Permissible Exposure Limit

REL: Recommended Exposure Limit

Skin Corr. 1A: Skin corrosion/irritation - Category 1A

Skin Irrit. 2: Skin corrosion/irritation – Category 2

Eye Irrit. 2A: Serious eye damage/eye irritation – Category 2A

Eye Irrit. 2B: Serious eye damage/eye irritation – Category 2B

\* \* Data compared to the previous version altered.

SDS created by MSDS Authoring Services www.msdsauthoring.com +1-877-204-9106