




Specifications | CG001694 | Rev B

Atera Instrument

For use with:

Atera Instrument, PN-1001128



Notices

Document Number

CG001694 | Rev B

Legal Notices

© 2026 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: www.10xgenomics.com/trademarks. 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: www.10xgenomics.com/patents. The use of products described herein is subject to 10x Genomics Terms and Conditions of Sale, available at www.10xgenomics.com/legal-notices, or such other terms that have been agreed to in writing between 10x Genomics and user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

Instrument & Licensed Software Updates Warranties

Updates to existing Instruments and Licensed Software may be required to enable customers to use new or existing products.

Support

Email: support@10xgenomics.com

10x Genomics, Inc.

6230 Stoneridge Mall Road

Pleasanton, CA, 94588

Document Revision Summary

Document Number

CG001694

Title

Atera Instrument Specifications

Revision

Rev B

Revision Date

May 08, 2026

Description of Changes

- Updated electrical specifications, heat output and clarified language regarding storage capacity and Ethernet link speed in [Instrument Specifications on page 4](#)

Instrument Specifications

Parameter	Atera Instrument Specifications		
Weight			
Instrument	1,350 lb/614 kg		
Instrument Live Load	500 kg/m ²		
Dimensions	Width	Depth	Height
Instrument	52.9"/134.5 cm	35.9"/91.3cm*	63.6"/161.5 cm
UPS (APC SMX2000LV for low voltage or SMX2200HV for high voltage configurations)	7"/17.8 cm	19"/48.3 cm	17"/43.2 cm
Storage Requirements			
Humidity	<60% for long-term storage (>1 week) <80% for transportation and short-term storage (<1 week)		
Temperature	13-25°C (55-77°F) for long-term storage (>1 week) 10-40°C (50-104°F) for transportation and short-term storage (<1 week)		
Operating Specifications			
Electrical Specifications	100-240 VAC, 50-60 Hz, 1,600 W maximum		
Pollution Degree	2 (Indoor Use Only)		
Operating Temperature	19-25°C (66-77°F) Use in a typical indoor laboratory environment. Extreme temperature conditions will result in run failure or possible damage to the instrument		
Humidity	30-80% Relative Humidity, noncondensing		
Altitude	Altitude up to 2,000 m (1.2 miles) above sea level		
Environmental Vibration Guidelines	ISO Office (or better) during idle ISO Operating Theater (or better) during run No bumps or shocks adjacent to the instrument during a run		
Heat Output	~1,600 W (5,460 BTU/h)		
Power Cable Length	6 ft/183 cm		
Sound Output	<60 dB at 1 meter from instrument		
Output Data Storage Capacity	30 TB NVMe		
Ethernet Link Speed	Up to 10 Gbps		
Operating System	Ubuntu 24.04 LTS		
*Instrument can temporarily be configured to pass through a standard 35"/89 cm door			

